

NON-CEMENTED MODULAR STEM BIOMEC - BM

* For applications where a metal or alloy comes into contact with each other and there is no joint, provider is given special attention to the design, surface finishing, surface treatment and metallurgical conditions.

Commercial Presentation of the Product - Package Features: Non-Cemented Modular Stem Biomec - BM is provided in two types of sterilization, Ethylene Oxide and Gamma Radiation. For each type of sterilization the packages are composed of:


Type of Sterilization	Package Description
Ethylene Oxide	PVC blister or PET (Polyethylene Terephthalate, Bed, Tyvek, Carton Box, and PVC shrink film)
Gamma Radiation	Non-toxic PET blister, Bed, Tyvek, Carton Box, and PVC shrink film

For Ethylene Oxide Sterilization is used the PVC blister package (Polyvinyl chloride) and the PET package (Polyethylene Terephthalate) for products sterilized by gamma radiation. PET (Polyethylene Terephthalate) and PVC (Polyvinyl chloride) are raw materials indicated for the manufacture of thermofomed blisters and beds (Vacuum

- Forming in the Drug, Food, Toys, Electronic Components, Cosmetics, and Automobile Industries. Both have as the main features: good impact strength; great dimensional stability; non-toxic; allowing contact with food and hospital products; good modulus; good brittleness; and self-extinguishing. They are preferred by mechanical properties, high transparency, virtually no defects, good hot fluidity, and by not presenting danger to human health.


Non-Cemented Modular Stem Biomec - BM after receiving a double blister packaging are packed in cardboard box, these boxes are customized with company logo and wrapped with PVC shrink film, containing product labels and labels stating Sterile Product (Red label), distributed unitarily.

Traceability Labels (Label) as per NBRSIO15223



Information Contained on the Cardboard Box (External Package): Cardboard boxes (external package) also have sayings of special storage, conservation and / OR handling conditions of the product with symbols. The table below describes the information contained on this cardboard box.

Special Storage, Conservation and / OR handling conditions contained on the external package of the product:



BIO MECANICA

Indication, purpose or intended use of the medical product, as indicated by the Manufacturer:

Intended Use: The use of Non-Cemented Modular Stem Biomec - BM is mainly for the following medical indications: Primary hip Osteoarthritis, Osteoarthritis Secondary to A. Rheumatic diseases (rheumatoid Arthritis, Ankylosing Spondylitis, Acetabular Protrusion of Rheumatic Origin); B. Sequels to Hip Disease in Childhood (Hip Developmental Dysplasia, Leg-Gauche-Perthes Disease, Epiphysiolysis of the Femur Head, Non-Rheumatic Acetabular Protrusion); C. Osteonecrosis of the Femur Head; D. Trauma (sequo to fractures of the femur neck and head, and acetabular fracture); E. And Osteometabolic disorders (Paget's Disease); III. Hip fractures (Neck, Head and Neck). Due to the excellent risk and cost / effectiveness ratio of hip arthroplasty, it may be subject of the total hip arthroplasty procedure, observed the contraindications and warnings contained in this document.

Recommendations on Orthopedic Prostheses: To achieve better results in hip arthroplasty, due following recommendations: 1- The use of orthopedic prostheses must be performed by surgeons qualified and trained on such procedure. It is essential a careful preoperative planning, including the aid of product transparency or "template". 2- Every care should be taken in preparing the receiver bone to get the perfect fit to the implant, avoiding undesirable radiolucency and undesirable emergence of microcomponents. Surgical instruments and tools available to assist in the surgery are recommended for the orthopedic prosthesis. Intraoperative instruments and tests used to those designed specifically for the product. Variation in design and sizes of surgical instruments and stents may compromise critical measure required for an accurate implantation. 4- In order to protect the doctor and patient from future issues, the hospital must take responsibility for noting on the patient protocol the code and batch number of the implants used. 5- During the surgery, the collection of implants, the commanding instruments are intact and in full. 6- Never mix implants from different manufacturers in the same patient. There may be differences in sizes and tolerances in the fittings of these implants, causing joint dislocations, early loosening or flaws. 7- The orthopedic prostheses are classified as "Single Use Products", that is, cannot be reused. 8- In all cases must be followed preclinical surgical practices in the procedure of the surgery. 9- If the product is found out of or correct positioning or showing any nonconformity, it is the responsibility of the surgeon to take appropriate corrective actions.

Evaluations of the Product Implanted: after the implantation, intraoperatively the professional in charge must perform radiological control to check the correct positioning of the product. The professional in charge must perform, being of his responsibility, clinical and radiological evaluations after the surgical procedure at a frequency stipulated by him to check the implant status and the evolution of bone healing. If the product is found out of or correct positioning or showing any nonconformity, it is the responsibility of the surgeon to take appropriate corrective actions.

Useful Information to Avoid Risks Arising from Implantation: To decrease risks arising from implantation, the following must be strictly observed:

contraindications, instruction for use, and all information contained in the product "Instruction for Use".

Risk Inherent to Implantation: The Family of Non-Cemented Modular Stem Biomec - BM is manufactured with materials of recognized biomedical use. Being ASTM F75 Standard Specification for Cobalt-28 Chromium-Molybdenum Alloy Castings and Casting Alloy for Surgical Implants.

Contamination Risk: There are risks of biological contamination and transmission of viral diseases, such as HIV and Hepatitis, because the components of the Family of Non-Cemented Modular Stem Biomec - BM comes into contact with tissue and body fluids. Explanted products must be treated as highly contaminating.

Sterilization: This product is provided in sterile form. According to Customer's request, it can be supplied with the following Sterilization Processes:

Process	Validity
Sterilization by Ethene Oxide	5 years
Sterilization by Gamma Ray 25 kgye	5 years

Note: Differentiation of sterilization type will be defined by the product code. ETO sterilized products have codes different from those products sterilized by gamma radiation. If the product has the sterilization expired / end / or damaged packaging, it must be returned to the responsible supplier or directly to Biomecanica.

Holding and Transport: Care for this Medical Product: storage location of the medical product must be clean and dry, and must be kept at room temperature, for storage or transport, as well as its physical and chemical integrity. The product should be stored and transported in a cool and dry place, at room temperature (-35 °C), and relative humidity of 30% to 70%. The effects of vibration, shock, corrosion, defective seating during handling and transport, inappropriate stacking during storage, all should be avoided. The product should not be stored on the floor or on shelves, as well as the capacity and willingness of the patient to follow the instructions are the two most important aspects for successful arthroplasty. The patient must be warned that failure to follow the postoperative instructions may lead to breakdown or migration of the orthopedic prosthesis, requiring new surgery for revision or removal. 12- Every effort must be made in this sense of using compatible materials and components that do not cause any reaction or migration of the prosthesis, either in the joint and the recent history of this practice. 9- The surgeon must replace the replaced joint against excessive "stress". 18- In general, orthopedic prostheses are supplied in sterile condition, in double blister, external rigid cardboard box, properly identified with labels, with all relevant legal information regarding the product, which guarantees full identification and traceability thereof. The packaging shall be intact at Receipt (Do not use the product if the packaging is damaged). Check the Sterilization Validity (Do not use the product after the sterilization validity).

Disclaimer: Biomecanica, as the Device Manufacturer, does not practice medicine nor recommends any surgical technique for use in a given patient. The surgeon should be informed about the content and type of product to be used and the techniques for implanting the prosthesis in a given patient. Biomecanica is not responsible for choosing the proper surgical technique to be used in the patient, who is responsible for choosing the method, type and size of products to be used.

Contraindication: The Only absolute contraindication to total hip arthroplasty is the presence of aseptic infection at the aseptic site or systemically. In such condition, the implantation of metal material may perpetuate an existing infection or stimulate the infection from other sites in systemic infection. Patients with remote history of hip joint infection may have recurrence of infection, thus Clinical, Laboratory and Imaging tests are mandatory for ruling out latent infections.

Precautions, restrictions, warnings, special care, and clarifications on the use of the medical product, as well as its storage and transportation:

Handling of Sterile Material: When handling aseptically sterilized material, one must observe some rules in order to keep it sterile: it is essential to wash hands with soap and water before handling sterilized material; use materials with undamaged, dry, spotless, identified packaging (type of materials and sterilization date); work facing the material; Handle the material from the waist up; Avoid coughing, sneezing or speaking over the exposed material; Do not make moves over a sterile field; Check validity and packaging suitability; Work in a clean, calm and environment without dust; Use gloves, clean clothing, mask, cap, apron, etc.; It is better to be hand dressed over the clothing; The Nursing Technique advised for handling sterilized materials follows: Open it, starting from the opposite end to the handle; protect the exposed material with the surrounding sterile field; Touch your hands only on the outside of the package; Do not keep as sterile material a previously open package.

Intraoperative Fractures: For immediate mechanical stability it is necessary coupling

by contact pressures, which means an implant size slightly greater than the cavity created interoperatorily according to the surgical technique recommended by the manufacturer. Failure to observe this aspect may lead to impacting an eventual fracture and cause fracture to the proximal femur. Early Load: Despite an appropriate initial fixation (immediate), it is advisable during the osteointegration period to protect the load support in order to not cause micro-movements on the implant-bone interface. The use of gait support is recommended for 4-6 weeks postoperatively. Implant Cementation: Due to plasma spraying with titanium implants, the implants are indicated for use without bone cement. Cementing these implants may predispose them to corrosion on gaps in environments without oxygen exchange (differential aeration cell) and thus cause fail.

Precautions and Warnings:

Product Disposal: Implants making up the Non-Cemented Modular Stem Biomec - BM Family and explanted from patients must be properly disposed by the hospital. It is the responsibility of the hospital the complete de-characterization of the implant, preventing it to be reused. It is the responsibility of the hospital the method used to de-cement the load support in order to not cause micro-movements on the implant-bone interface. The use of gait support is recommended for 4-6 weeks postoperatively. Implant Cementation: Due to plasma spraying with titanium implants, the implants are indicated for use without bone cement. Cementing these implants may predispose them to corrosion on gaps in environments without oxygen exchange (differential aeration cell) and thus cause fail.


After Sales (Customer Complaint): In case of need to do anything complaint on Non-Cemented Modular Stem Biomec - BM regarding some adverse effect affecting the user safety, such as product not working, damage to metal implantable component, severe problems or death related to these components, the surgeon in charge must communicate these adverse events to the competent Health Authority and to Biomecanica by the email sac@biomecanica.com.br or call to 0xx14 2104 7792. In case of being reprocessed.

REUSE THE PRODUCT - PROHIBITED AFTER EXPLANTED. DO NOT REUSE THE PRODUCT: The manufacturer can result in a single patient only. Although it may appear to be undamaged, previous tensions can create imperfections that may reduce the implant success, improper selection of implant can cause unsatisfactory tensile on the implant and result in subsequent fracture. The surgeon must be familiar with and have sufficient knowledge, including the pre- and post-operative period, surgical technique adopted, precautions, and potential risks. When handling the implant, it should be considered that the implants may contain some materials that may be crack susceptible such as damages, connectors, tensions and may be sites of crack nucleation and decrease corrosion resistance, and can lead to implant fracture. Check that the product is undamaged. Check the sterilization date. Aspirically open of package after making sure the size and the chosen one. Sterilization Validity: 5 years (indicated on the internal and external packaging of the product). To be considered valid to be treated by the Complementary Health System, one made available to the Supplier Control (History Distribution Record - RHD), one made available to the surgeon in charge (main), and the last one to Health Insurance, if applicable. We inform that mandatory the registration number of the implant used, name of the implant used, implant code, batch number, and the Product Registration Number at the ANVISA. The information is described in the traceability labels accompanying the product, and on the external labeling. Other information must also be considered essential, such as surgery date, name of the patient who received the implant, name of the surgeon, patient weight, patient age, and other information requested in the Patient Record must also be completed.

Traceability Labels (Label) as per NBRSIO15223

Information Contained on the Cardboard Box (External Package): Cardboard boxes (external package) also have sayings of special storage, conservation and / OR handling conditions of the product with symbols. The table below describes the information contained on this cardboard box.

Special Storage, Conservation and / OR handling conditions contained on the external package of the product:



Indication, purpose or intended use of the medical product, as indicated by the Manufacturer:

Intended Use: The use of Non-Cemented Modular Stem Biomec - BM is mainly for the following medical indications: Primary hip Osteoarthritis, Osteoarthritis Secondary to A. Rheumatic diseases (rheumatoid Arthritis, Ankylosing Spondylitis, Acetabular Protrusion of Rheumatic Origin); B. Sequels to Hip Disease in Childhood (Hip Developmental Dysplasia, Leg-Gauche-Perthes Disease, Epiphysiolysis of the Femur Head, Non-Rheumatic Acetabular Protrusion); C. Osteonecrosis of the Femur Head; D. Trauma (sequo to fractures of the femur neck and head, and acetabular fracture); E. And Osteometabolic disorders (Paget's Disease); III. Hip fractures (Neck, Head and Neck). Due to the excellent risk and cost / effectiveness ratio of hip arthroplasty, it may be subject of the total hip arthroplasty procedure, observed the contraindications and warnings contained in this document.

Recommendations on Orthopedic Prostheses: To achieve better results in hip arthroplasty, due following recommendations: 1- The use of orthopedic prostheses must be performed by surgeons qualified and trained on such procedure. It is essential a careful preoperative planning, including the aid of product transparency or "template". 2- Every care should be taken in preparing the receiver bone to get the perfect fit to the implant, avoiding undesirable radiolucency and undesirable emergence of microcomponents. Surgical instruments and tools available to assist in the surgery are recommended for the orthopedic prosthesis. Intraoperative instruments and tests used to those designed specifically for the product. Variation in design and sizes of surgical instruments and stents may compromise critical measure required for an accurate implantation. 4- In order to protect the doctor and patient from future issues, the hospital must take responsibility for noting on the patient protocol the code and batch number of the implants used. 5- During the surgery, the collection of implants, the commanding instruments are intact and in full. 6- Never mix implants from different manufacturers in the same patient. There may be differences in sizes and tolerances in the fittings of these implants, causing joint dislocations, early loosening or flaws. 7- The orthopedic prostheses are classified as "Single Use Products", that is, cannot be reused. 8- In all cases must be followed preclinical surgical practices in the procedure of the surgery. 9- If the product is found out of or correct positioning or showing any nonconformity, it is the responsibility of the surgeon to take appropriate corrective actions.

Evaluations of the Product Implanted: after the implantation, intraoperatively the professional in charge must perform radiological control to check the correct positioning of the product. The professional in charge must perform, being of his responsibility, clinical and radiological evaluations after the surgical procedure at a frequency stipulated by him to check the implant status and the evolution of bone healing. If the product is found out of or correct positioning or showing any nonconformity, it is the responsibility of the surgeon to take appropriate corrective actions.

Useful Information to Avoid Risks Arising from Implantation: To decrease risks arising from implantation, the following must be strictly observed:

contraindications, instruction for use, and all information contained in the product "Instruction for Use".

Risk Inherent to Implantation: The Family of Non-Cemented Modular Stem Biomec - BM is manufactured with materials of recognized biomedical use. Being ASTM F75 Standard Specification for Cobalt-28 Chromium-Molybdenum Alloy Castings and Casting Alloy for Surgical Implants.

Contamination Risk: There are risks of biological contamination and transmission of viral diseases, such as HIV and Hepatitis, because the components of the Family of Non-Cemented Modular Stem Biomec - BM comes into contact with tissue and body fluids. Explanted products must be treated as highly contaminating.

Sterilization: This product is provided in sterile form. According to Customer's request, it can be supplied with the following Sterilization Processes:

Process	Validity
Sterilization by Ethene Oxide	5 years
Sterilization by Gamma Ray 25 kgye	5 years

Note: Differentiation of sterilization type will be defined by the product code. ETO sterilized products have codes different from those products sterilized by gamma radiation. If the product has the sterilization expired / end / or damaged packaging, it must be returned to the responsible supplier or directly to Biomecanica.

Holding and Transport: Care for this Medical Product: storage location of the medical product must be clean and dry, and must be kept at room temperature, for storage or transport, as well as its physical and chemical integrity. The product should be stored and transported in a cool and dry place, at room temperature (-35 °C), and relative humidity of 30% to 70%. The effects of vibration, shock, corrosion, defective seating during handling and transport, inappropriate stacking during storage, all should be avoided. The product should not be stored on the floor or on shelves, as well as the capacity and willingness of the patient to follow the instructions are the two most important aspects for successful arthroplasty. The patient must be warned that failure to follow the postoperative instructions may lead to breakdown or migration of the orthopedic prosthesis, requiring new surgery for revision or removal. 12- Every effort must be made in this sense of using compatible materials and components that do not cause any reaction or migration of the prosthesis, either in the joint and the recent history of this practice. 9- The surgeon must replace the replaced joint against excessive "stress". 18- In general, orthopedic prostheses are supplied in sterile condition, in double blister, external rigid cardboard box, properly identified with labels, with all relevant legal information regarding the product, which guarantees full identification and traceability thereof. The packaging shall be intact at Receipt (Do not use the product if the packaging is damaged). Check the Sterilization Validity (Do not use the product after the sterilization validity).

Disclaimer: Biomecanica, as the Device Manufacturer, does not practice medicine nor recommends any surgical technique for use in a given patient. The surgeon should be informed about the content and type of product to be used and the techniques for implanting the prosthesis in a given patient. Biomecanica is not responsible for choosing the proper surgical technique to be used in the patient, who is responsible for choosing the method, type and size of products to be used.

Contraindication: The Only absolute contraindication to total hip arthroplasty is the presence of aseptic infection at the aseptic site or systemically. In such condition, the implantation of metal material may perpetuate an existing infection or stimulate the infection from other sites in systemic infection. Patients with remote history of hip joint infection may have recurrence of infection, thus Clinical, Laboratory and Imaging tests are mandatory for ruling out latent infections.

Precautions, restrictions, warnings, special care, and clarifications on the use of the medical product, as well as its storage and transportation:

Handling of Sterile Material: When handling aseptically sterilized material, one must observe some rules in order to keep it sterile: it is essential to wash hands with soap and water before handling sterilized material; use materials with undamaged, dry, spotless, identified packaging (type of materials and sterilization date); work facing the material; Handle the material from the waist up; Avoid coughing, sneezing or speaking over the exposed material; Do not make moves over a sterile field; Check validity and packaging suitability; Work in a clean, calm and environment without dust; Use gloves, clean clothing, mask, cap, apron, etc.; It is better to be hand dressed over the clothing; The Nursing Technique advised for handling sterilized materials follows: Open it, starting from the opposite end to the handle; protect the exposed material with the surrounding sterile field; Touch your hands only on the outside of