# MINI TWISTER – ÂNCORA DE SUTURA

## PORTEGUÊS

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# MINI TWISTER – SUTURA ANCHOR

## ENGLISH

### Quantity and Description of the Traceability Labels that Comes with the Product:

- 6 adhesive labels of identification and traceability of the product;

### Description of the Components, Quantity and Raw Material:

- 1 unit of titanium alloy Anchor Screw Ti-6AL-4V; 1 unit of White / Blue Suture Thread of polyethylene fiber UHMWPE; - 1 unit of Stainless Steel Anchor Screw Applicator AISI 304, Polypropylene and Silicone.

### Product Presentation:

The product is packed in a blister-type double package, which is manufactured from non-toxic PET (Polietileno Tereftalato), allowing a contact with hospital products. The two blister packages are sealed with a surgical grade paper, without dye and free of holes, through the heat sealing process. The double sealed blister package is also placed in an outer cardboard package, covered with security seal. The implants are sterilized in ETO with 5 years validity.

### Traceability:

Along with the product, there are 6 labels containing the necessary data for traceability of the product (Manufacturer's name, trade name of the product and code, technical product name and code, batch number, registration number at ANVISA, quantity, manufacturing date, expiration date, manufacturer's address, technical responsible name and the words: "sterile product," "do not reprocess," "single-use product, do not reuse," "before using read the label instructions" and "warnings / precautions / special care / symbology: see instructions for use").

A label must be mandatorily fixed to the patient's clinical record, another one in the report delivered to the patient, other in the tax documentation generating the collection, at Alfi, in case of a patient treated by SUS (Brazilian Unified Health System), or in the sales invoice, in case of a patient cared by the supplementary health system, another one made available to the supplier's control (historical record of distribution - RHD, another one made available to the control of the responsible surgeon (principal) and the last one to the health insurance, if any, in order to keep the traceability cycle of the product.

Only the Aplicador para Parafuso Âncora 2,0/ 2,7, present in Mini Twister - Suture Anchor, receives the laser marking, which is performed in the stainless steel manufactured part, according to the sequence below:

Logo - Batch - Product code

The other components of Mini Twister - Suture Anchor do not receive marking due to their anatomy. Even though, the traceability of all components of Mini Twister - Suture Anchor is assured through the batch number present in the traceability labels of the product, which can be identified by tracing the batch number of each component.

### Product Description

Mini Twister - Suture Anchor consists on a suture medical device, not active, composed of three components, two non-absorbable implants and an instrumental (applicator), assembled and ready for use.


Figura 1: Mini Twister - Suture Anchor

The single component of Mini Twister - Suture Anchor showing the variation is the Anchor Screw, and may have a diameter of 2.0 mm or 2.7 mm. In the presence of the Anchor Screw 2.0 mm, the code of Mini Twister - Suture Anchor becomes 3128-20-00; in the presence of the Anchor Screw 2.7 mm, the code of Mini Twister - Suture Anchor becomes 3128-27-00. Both the Anchor Screw 0,2,0, as the Anchor Screw 0,2,7, are implants, not absorbable and not active, manufactured from titanium alloy 6-aluminum 4-vanadium (Ti6Al4V), according to the ASTM F136 standard.

Both the Anchor Screw 0,2,0 and the Anchor Screw 0,2,7 have a pointed end, which assist the introduction of the screw in the bone wall. The threaded part causes the screw to be introduced in a safe and practical way, providing strength and stability for screw fixation in the bone.

Both the Anchor Screw 0,2,0 and the Anchor Screw 0,2,7 have a polished finishing.



Figura 2: Anchor Screw

Table 1: Dimensional Comparison between the Anchor Screw Ø 2,0 and Anchor Screw Ø 2,7

	Thread Diameter	Thread Length	Pointed End Length	Fitting Area Length of the Thread and the Applicator	Total length
Anchor Screw Ø 2,0	2.0 mm	4.5 mm	1.5 mm	3.0 mm	9 mm
Anchor Screw Ø 2,7	2.7 mm	6.0 mm	2.8 mm	3.2 mm	12 mm

The White / Blue Suture Thread is a non-absorbable implant of synthetic origin for surgical suture, composed of a fiber called dyneema Purity with a high mechanical resistance based on ultra-high molecular weight polyethylene (UHMWPE); multifilament braided; and not active.

The White / Blue Suture Thread has a length of 900 mm and an average diameter of 0.585 mm, and can vary from 0.470 mm to 0.705 mm.

Figure 3: Suture Thread

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Open Aseptically the Package	Abrir Asépticamente a Embalagem
It is Recommended a Temperature of 21 °C (+-4 °C), Environmental Humidity in a Closed, Vented Located and Protected against Weatherproof.	Recomendasse Temperatura de 21°C (+-4°C), Umidade Ambiente em Local Fechado, Ventilado e Protegido de Intempéries.
Do not store directly on the floor (minimum height = 20 cm) and nor in very high places, close to lamps, which could cause dryness of the package or damage to the label.	Não armazenar diretamente no chão (altura mínima = 20cm) e nem em locais muito altos, próximos a lâmpadas, o que poderia ocasionar ressecamento da embalagem ou dano no rótulo.
Do not store in locations containing contaminant substances such as, for example, cleaning materials, insecticides, pesticides, etc.	Não armazenar em locais que contenham substâncias contaminantes como, por exemplo, materiais de limpeza, inseticidas, pesticidas, etc.
Do not use if the package is damaged	Keep away from sunlight
Fragile, handle with care	keep dry
STERILE EO Esterilizado en Gas ETO ESTÉRIL	Proteger de la exposición solar directa
"PRODUCTO DE UN SÓLO USO - PROHIBIDO UTILIZAR NUEVAMENTE"	Mantener seco

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### Indications, Contraindications, Adverse Effects and Warnings

#### Indications

Mini Twister - Suture Anchor is indicated for the tendon-bone repair of the flexor tendons of the