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PRODUCT DESCRIPTION

Orthopedic Cement

Acrilex Orthopedic Cement consists of Methyl methacrylate: a mixture of Poly Methyl methacrylate, Methyl methacrylate - Styrene-copolymer. It is a cement packed in a blister with two sterile components.

Component 1 = Liquid Component: A sterile blister with an ampoule containing 20ml of a colorless, flammable, liquid monomer with the following composition:

Methyl Methacrylate Formulation	98,215% w/w
NN Dimethyl-p-toluidine	0,816% w/w
Ethyl Alcohol	0,945% w/w
Ascorbic Acid	0,022% w/w
Hydroquinone	0,002% w/w

* NN Dimethyl-p-toluidine is added to promote cold curing of the finished liquid.

**Hydroquinone is added to prevent premature polymerization, which can occur under certain conditions such as high temperatures, exposure to light.

Component 2 = Powder Component: A 40g sterile blister. Its main component is Polymethyl methacrylate, Barium Sulfate, and Benzoyl

Peroxide. Benzoyl peroxide starts polymerization when mixing the powder and liquid component.

Polymethyl Methacrylate Formulation	87,50% w/w
Barium Sulfate	10,00% w/w
Benzoyl Peroxide	2,5% w/w

PRESENTATION

Packaging containing:

LIQUID COMPONENT: A sterile blister with an ampoule containing 20ml of a colorless liquid monomer.

POWDER COMPONENT: A 40g sterile blister.

INDICATIONS

Orthopedic cement is indicated for setting the prosthesis to the live bone in joint arthroplasties.

Recommended when reconstruction is necessary due to osteoarthritis, rheumatoid arthritis, traumatic arthritis, vascular necrosis, osteoporosis due to sickle cell anemia, collagen and severe joint destruction secondary to trauma or other conditions, and in previous arthroplasty revisions.

Orthopedic Cement is indicated for use in children only in cases of limb preservation or where no other procedure seems likely to be successful. Orthopedic cement must be used with the appropriate prosthesis.

INDICATION FOR USE (Preparation and Application)

For the preparation, the product must be at room temperature (23°C). If the product has been stored at a different temperature, it takes 3 to 4 hours to achieve room temperature. Orthopedic Cements are very sensitive to heat, and any increase in the room temperature or the cement components above 23°C may reduce preparation, handling, and polymerization times. On the other hand, lower temperatures will increase the time. The packaging must be opened by a circulating nurse, and the sterile content must be transferred over a sterile surface. Break the ampoule containing the liquid and pour ALL the liquid content into the powder in the blister. ADD THE LIQUID TO THE POWDER, NOT THE POWDER TO THE LIQUID. Containers and spatulas should be inert materials such as glass, ceramics, or non-reactive plastics.

DIGITAL APPLICATION

Orthopedic Cement is mixed with a dry, sterilized spatula for approximately 1 1/4 minutes. It is very important to minimize the entry of air while mixing, as the presence of air impairs the effectiveness of the cement and can impede its mechanical performance. The time this occurs after the mixing depends on the type of cement and the temperature. If the Orthopedic Cement was taken by hand, very early, it may initially stick onto the gloves, but as soon as it is handled, any adherent cement will be quickly absorbed by the cement ball, leaving the gloves clean. The handling process should continue from 1/4 to 1 minute when the cement is ready for application to the bone. The digital application of the cement should last from 1 to 1/4 minutes (after that time the Orthopedic Cement

becomes too viscous to insert it in the narrow cavity). The prostheses must be inserted and held firmly in position until setting the Orthopedic Cement.

Excessive Orthopedic Cement must be removed before the cement has been completely hardened.

FINAL HARDENING TIME

Normal Viscosity: "10 minutes.

IMPORTANT NOTES

1 - Before applying Orthopedic Cement to the bone, the cavity must be thoroughly sanitized and dried to prevent blood and foreign bodies from mixing with the cement, which may prevent the expected application for the cement.

2 - During the digital application of Orthopedic Cement, it is important to use a thin plastic probe remove blood, air and fat. The probe must be removed before setting the prosthesis.

3 - Orthopedic Cement is a relatively fast-setting bone cement, minimizing the waiting period for polymerization during the surgery. Orthopedic Cements do not have adhesive properties; however, they have a precise mechanical blocking with the irregular surfaces of the bone as those of the prostheses.

CONTRAINDICATIONS

Orthopedic cement is contraindicated in case of hypersensitivity to any of the components of the cement or in the presence of infections.

INTERACTION WITH OTHER DRUGS

Unknown.

WARNINGS

For the safe and effective use of Orthopedic Cement, the surgeon must have received specific preparation for the use of bone cement and be aware of its handling characteristics, its use limitations and its appropriate application.

The evidence from clinical investigations indicates the need for thorough submission to good aseptic surgical techniques.

Properly cemented prostheses are stable and long-lasting. However, there is the possibility of loosening and fracture of the cement or prosthesis, or both, caused by trauma, inadequate cementation technique, mechanical failure of the materials, or latent infection. Therefore, routine long-term follow-up of all patients is recommended.

Following the operation, the patient should be advised if there has been an intercurrent infection, and immediately seek medical help to reduce the risk of infection caused by the implant.

PRECAUTIONS

- 1 - The use of Orthopedic Cement requires collaboration and consultation between surgeon and anesthetist.
- 2 - The anesthetist must be informed, during the operation, of the moment to set the cement.
- 3 - The monomer can cause hypotension, and such low blood pressure after deliberate or intentional hypotension can produce cardiac arrhythmia or ischemic myocardium.
- 4 - The hypotensive effects of methyl methacrylate can be enhanced if the patient suffers from hypovolemia.
- 5 - The setting of a foreign body in the tissues increases the normal risk of infection after surgical interventions.
- 6 - Regarding the setting of cement, it is important to maintain the position of the prosthetic component until the completion of the polymerization process, to allow an appropriate setting.
- 7 - The routine use of an appropriate prophylactic antibiotic may be desirable in the clinician's prescription.

USER PRECAUTIONS

Orthopedic Cement liquid and powder components are pre-measured to give optimal results and are produced for use as single dose. Always add the entire liquid to the powder when mixing the Orthopedic Cement. TO PREVENT ANY POSSIBLE CONTAMINATION OF THE CEMENT WITH GLASS FRAGMENTS DO NOT BREAK THE AMPOULE OVER THE CONTAINER.

The liquid is volatile and flammable. The operating room must be properly ventilated. The liquid component and its vapor must not be exposed to flames. Continuous inhalation of liquid vapor can have a soporific effect. Excessive exposure to concentrated vapors can also cause irritation of the respiratory tract and eyes, which should be avoided.

If the liquid component has direct contact with the eyes, wash them with plenty of water. The liquid component is a powerful solvent that can cause contact dermatitis in sensitive people. The possibility of hypersensitive reactions is reduced by wearing a second pair of gloves and strictly following the preparation instructions. The monomer must NOT come into contact with rubber. Concentrated vapors from the liquid component may have an adverse reaction with soft contact lenses.

PHARMACEUTICAL PRECAUTIONS

Store below 25°C and protect from light. The content is sterile, unless the packaging is opened or damaged. Re-sterilization of any of the components must not be done.

If the liquid hardens in the ampoules, it should be noted that the Liquid Component also polymerizes itself, i.e., without adding the Powder Component in some special cases, namely: (1) in case of lack of stabilizer; (2) in constant exposure to sunlight, directly or indirectly; and (3) in contact or exposure to heat (remember that the product is flammable).

ADVERSE EFFECTS

The most severe and frequent adverse effects, which can occur with the use of bone cement, are as follows. Surgeons must be aware of such reactions, and prepared to treat them if applicable. The most severe ad-

verse reactions reported with its use, but not directly reported to acrylic bone cement, are:

SEVERE:

- 1) Myocardial infarction,
- 2) Cerebrovascular accident,
- 3) Cardiac arrest,
- 4) Sudden death, and
- 5) Pulmonary embolism.

MORE FREQUENTLY:

- 1) Transient drop in blood pressure,
- 2) Thrombophlebitis,
- 3) Hemorrhage and hematoma,
- 4) Loss or displacement of the prosthesis,
- 5) Superficial or deep infection of the lesion,
- 6) Bursitis in the trochanter, and
- 7) Short-term cardiac conduction irregularity.

OTHER REACTIONS REPORTED:

- 1) New heterotopic bone, and
- 2) Trochanter separation

USE IN PREGNANCY AND BREASTFEEDING

There is no evidence of safety in the use of bone cement for associated breastfeeding procedures.

STERILIZATION

Product sterilized by Ethylene Oxide - ETO - Validity 3 years.

This product must not be re-sterilized. In case the surgical packaging is damaged, it must be disposed of.

NOTE

This product must be handled by qualified persons only.

RECOMMENDATIONS

We recommend filling in a product quality certificate and returning it to Biomecanica.

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