FIOS E PINOS PARA SÍNTESE ÓSSEA

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THREADS AND PINS FOR ENGLISH

PRODUCT DESCRIPTION

Metallic devices built in stainless steel with several diameters and lengths in cylindrical form, flat or with central thread, or with thread in one of the extremities, with a tip of very sharp penetration and of good resistance. The pins possess a chamfer that allows a better adaptation in the drill, facilitating the placement.

PRESENTATION FORM

Identified with non-sterile label, conditioned in plastic sack. Packing with 1 unit. When accompanied with Sterile External Fixators for OsteoSynthesis, they come conditioned in double blister stamped with surgical paper and sterilized.

INDICATIONS

The device for osteosynthesis used as auxiliary element for the bony fixation in the treatment of the fractures, to be used separately, or with External Fixator.

CONTRAINDICATIONS

Does not exist any contraindication related to the use of Threads and Pins for Osteo Synthesis.

WARNINGS

Do not reuse the product.

Orthopedic implants can be used in a single patient only. Apparently, it can seem not damaged, previous tensions can create imperfections that can reduce the success of the implants.

The responsible surgeon, in agreement with the surgical technique, should establish the dimension.

Due to your presentation form and characteristics, this product does not present collateral effect; however, it is suggested to monitor the patient's re-establishment through the X Ray, according to International Norms of Safety. Do not use magnetic resonance due to interferences.

PRECAUTIONS

As they are temporary implants, in the aid of bone fractures treatment, they should accompany of care: in transport, manipulation, and sterilization. Conserve at airy place, dry and far away from the action of bad weather.

The surgeon should be aware of the possibility of the patient's physical development, could have revision need of the implants.

Do not use the product if is damaged.

ADVERSE EFFECTS

Mechanical loosening can be the result of defective fixation or infection hides. Sensibility reactions to metal in patients were rarely informed.

Implantation of strange material in tissues results in histological reactions that involve several macrophages sizes and fibroblasts (reaction type strange body). The chemical importance of this effect is uncertain, as well as similar changes can happen as a precursor, or during the healing process

STERILIZATION

STERILE: When accompanied with External Fixator, it comes sterilized by Oxide of Ethylene - E.T.O. - Validity 3years.

In the eventuality of sterilization is expired or surgical paper damaged, this product must be returned to manufacturer.

NOT STERILE: When supplied not sterile, product should be sterilized by autoclave. recommended the sterilization through steam in Autoclave in the hospital. (ISO 11134: 1994 *Sterilization of health care products Requirements for validation and*

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routine control Industrial moist heat sterilization).

AUTOCLAVE

It is equipment for sterilization through steam saturated under pressure.

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FOR BETTER EXPLANATION CONSULT THE MANUAL OF INSTRUCTIONS THAT ACCOMPANIES EACH AUTOCLAVE.

RECOMMENDATIONS

Qualified people should only manipulate this product.

Do not use implant components from different manufacturers. Therefore, we recommended that the products have the same origin.

All implants are of unique use. When retired from patient, it should be clean and discarded.

PRODUCT DISCARD

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way threads and pins for osteosynthesis explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

CONTAMINATION RISKS

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

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