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Blunt Hinged Surgical Instruments - BM

RG.: 80128580116

**ENGLISH
EUA**

Table 1 - Blunt Hinged Surgical Instrument Listing - BM

Code	Description	Level	Raw Material	Indication for Use
5124-00-000	Two-Pole Assembler	Hip	Stainless Steel AISI 304 (NBR 13911)	
6177-00-000	Reduction Tweezers 130 mm	Trauma	Stainless Steel AISI 420 (NBR 13911)	
6178-00-000	Reduction Tweezers 140 mm	Trauma	Stainless Steel AISI 420 (NBR 13911)	
6042-00-000	Reduction Tweezers 4.5	Trauma	Stainless Steel AISI 420 (NBR 13911)	
6043-00-000	Self-Centering Tweezers 4.5	Trauma	Stainless Steel AISI 420 (NBR 13911)	
5115-00-000	Acetabular Positioner CP3	Hip	Stainless Steel AISI 304 (NBR 13911)	
	Handle		Celeron	
5354-00-000	Pin Extractor Pliers	Knee	Stainless Steel AISI 420 (NBR 13911)	
5384-00-000	Clamp Patellary Component	Knee	Stainless Steel AISI 420 (NBR 13911)	
5355-02-000	BKS III - Tibial Plateau Introducing Component	Knee	Stainless Steel AISI 420B (NBR 13911)	
5352-00-000	BKS III - Tibial Plateau Removing Component	Knee	Stainless Steel AISI 420B (NBR 13911)	
8518-00-000	Rod Modeler	Column	Stainless Steel AISI 420 (NBR 13911)	
8517-00-000	'ST' Hook Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
8564-00-000	Plate Modeling Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
8556-00-000	Cervical Plate Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
8522-00-000	Rod Holder Tweezers Ø 5.8 mm	Column	Stainless Steel AISI 420 (NBR 13911)	
8622-00-000	Compression Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
6206-00-000	Rod Extractor	Trauma	Stainless Steel AISI 304 (NBR 13911)	



5136-00-000	Acetabular Screw Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	
5604-00-000	Insert Extractor Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	
8584-00-000	Distractor	Column	Stainless Steel AISI 304 (NBR 13911)	
5111-00-000	Ratchet Wrench	Hip	Stainless Steel AISI 304 (NBR 13911)	
5386-00-000	BKS III - Tibial Plateau Introducing Component	Knee	Stainless Steel AISI 304 (NBR 13911)	
5169-01-000	Extractor Tips (Moore Prost.)			
6041-00-000	Self-Centering Tweezers 3.5	Trauma	Stainless Steel AISI 304 (NBR 13911)	
6040-00-000	Reduction Tweezers	Trauma	Stainless Steel AISI 304 (NBR 13911)	
6066-00-000	Angular Plate Goniometer	Trauma	Stainless Steel AISI 304 (NBR 13911)	
6067-04-000	Positioner 70° 90° 110°	Trauma	Stainless Steel AISI 304 (NBR 13911)	
6067-03-000	Positioner 40° 50° 90°			
6067-02-000	Positioner 30° 70° 80°			
6067-01-000	Positioner 20° 60° 100°			
6273-00-000	Bioloking Extractor - Proximal Bioloking	Trauma	Stainless Steel AISI 304 (NBR 13911)	
8515-00-000	Rod Tweezers Ø 5.8 mm	Column	Stainless Steel AISI 420 (NBR 13911)	
8516-00-000	Rod Tweezers Ø 4.0 mm	Column	Stainless Steel AISI 420 (NBR 13911)	
8519-00-000	Distractor Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
8520-00-000	Compression Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
6233-00-000	Head Extractor	Trauma	Stainless Steel AISI 304 (NBR 13911)	
5309-00-000	Rod Extractor	Knee	Stainless Steel AISI 304 (NBR 13911)	
6241-00-000	Hook (Kuntscher)	Trauma	Stainless Steel AISI 304 (NBR 13911)	
5503-00-000	Rod Impactor / Extractor BR3	Hip	Stainless Steel AISI 304 (NBR 13911)	
6321-00-000	Rod Extractor SP2	Trauma	Stainless Steel AISI 304 (NBR 13911)	
8520-00-000	Compression Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	
8689-00-000	Thoracolumbar Plate Tweezers	Column	Stainless Steel AISI 304 (NBR 13911)	
8694-00-000	Vertebral Body Retractor	Column	Stainless Steel AISI 304 (NBR 13911)	
5124-10-000	Constricted Acetabular / Insert Assembler	Hip	Stainless Steel AISI 304 (NBR 13911) / Aluminum (NBR ISO209:10) / Polyacetal	
8691-00-000	Kelly Tweezers (Cylindrical Cage)	Column	Stainless Steel AISI 420 (NBR 13911)	

8681-00-000	Graft Impactor Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Bone Graft Compactor Tweezers
6500-01-000	Self-Centering Tweezers (small)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Bone tweezers for osteosynthesis plate positioning
6500-01-001	Reduction Tweezers (small)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-01-002	Reduction Tweezers with Tip (small)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction.
6500-01-021	Plate Modeling Wrench	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument to implant anatomic modeling.
6500-02-000	Self-Centering Tweezers (large)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Bone tweezers for osteosynthesis plate positioning
6500-02-001	Reduction Tweezers (large)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction.
6500-02-002	Reduction Tweezers with Tip (large)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction.
6500-03-007	Multi-angle Guide 90°~150°	Trauma	Stainless Steel AISI 304 (NBR 13911)	Drilling Hinged Guide.
6500-03-032	Reduction Wrench	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-04-008	Retractor	Trauma	Stainless Steel AISI 304 (NBR 13911)	Soft Part Retractor.
6500-04-010	Reduction Tweezers with Tip	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction.
6500-04-011	Reduction Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-012	Card Wrench with Spring	Hip	Stainless Steel AISI 420, 304 and 301 (NBR 13911)	Tightening/Loosening Wrench
6500-05-013	Reduction Tweezers with Round Tips	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-014	Reduction Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-015	Small Fragment Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-016	Large Fragment Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-018	Tweezers with Triple Spherical Tip	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-021	Pelvis Retractor	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction.
6500-05-024	Straight Reduction Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-05-025	Oblique Reduction Tweezers	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument used to bone reduction
6500-09-009	Multi-Hole Parallel Guide	Trauma	Stainless Steel AISI 304 (NBR 13911)	Bone Drilling Hinged Guide

6500-11-008	Cardan Wrench sw8.0	Trauma	Stainless Steel AISI 440C, 304 and 301 (NBR13911)	Tightening/Loosening Wrench
6500-12-006	Card Wrench with Spring	Trauma	Stainless Steel AISI 440C, 304 and 301 (NBR13911)	Tightening/Loosening Wrench
6500-17-012	Distractor Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Vertebral Distractor Tweezers.
6500-17-013	Compression Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Vertebral Compression Tweezers.
6500-17-014	Bar Modeler	Column	Stainless Steel AISI 420 (NBR 13911)	Instrument to implant anatomic modeling
6500-17-016	Bar Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Column Bar Tweezers
6500-17-018	Bar Holder Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Column Bar Tweezers
6500-17-019	Positioner Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Hook and Screw Positioning Tweezers
6500-17-020	Screw Break Tweezers	Column	Stainless Steel AISI 420 (NBR 13911)	Screw Break Tweezers.
6500-17-023	Vertebral Distractor/Compressor	Column	Stainless Steel AISI 304 (NBR 13911)	Instrument to Vertebra Compression and Distraction
5612-00-001	Card Wrench with Spring	Trauma	Stainless Steel AISI 304, 420B and 301 (NBR 13911)	Tightening/Loosening Wrench with Articulation
	Handle		Polypropylene and Polychloroprene	
5443-00-000	Headless Nail Extractor	Knee	Stainless Steel 420 (NBR 13911)	Extract headless nails
5445-00-000	Fracture Junction Tweezers	Knee	Stainless Steel 420 and 304 (NBR 13911)	Joint fractures
5448-00-000	Cardan Wrench 3.5	Knee	Stainless Steel 440C and 304 (NBR 13911)	Tighten screws
	Handle		Polypropylene and Polychloroprene	
4001-00-002	Extractor Pliers	Trauma	Stainless Steel AISI 420B (NBR 13911)	Hold tightly flexible rods with pliers tip. Ratchet closing for easy handling. Hammer guide use enables easy rod extraction.
4001-00-003	Flat Pliers	Trauma	Stainless Steel AISI 420B (NBR 13911)	Hold tightly flexible rods with pliers tip. Ratchet closing for easy handling. Hammer guide use enables easy rod extraction.
4000-00-009	Bending Pliers	Maxillofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Bend and model plates.
4000-00-010	Bending Pliers 2.0	Maxillofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Bend and model plates.
4000-00-011	Pliers with Bending Tooth 1.5	Maxillofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Bend and model plates.
4000-00-012	Pliers with Bending Tooth 2.0	Maxillofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Bend and model plates.
4000-00-027	Plate Bending Tweezers	Maxillofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Bend and model plates.
4000-00-028	Fracture Reduction Tweezers	Maxillofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Fracture Reduction.

4000-00-029	Plate Support Tweezers.	Maxilofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Keep the plate in pre-set position during surgery.
4000-00-031	Plate Bending Pliers	Maxilofacial Skull	Stainless Steel AISI 420B (NBR 13911)	Bend and model plates.
4002-00-014	Impactor - Humeral Prosthesis Extractor	Trauma	Stainless Steel AISI 304 (NBR 13911)	For prosthesis impacting and extracting handle

Traceability Label Symbols (Label) as per NBR ISO 15223

Manufacturing Date 	Valid Date 	Batch No. 	Manufacturer Data 	Code 	Non-Sterile 
--	--	---	---	--	---

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

Manufacturing materials: Most surgical instruments are made of stainless steels, as provided by NBR 13911:1997 or ASTM F899:07; that characterize their chemical composition, mechanical properties and final behavior of stainless steel, to produce surgical instruments, and in some instruments handles are made of polypropylene and polychloroprene.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Blunt Hinged Surgical Instruments - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Blunt Hinged Surgical Instruments - BM are reused in another patient, it is required that it is cleaned, sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as

provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

Cleaning: Appropriate care with surgical instruments starts with appropriate cleaning. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

Instrument Sterilization: Instruments are supplied non-sterile. They have to be sterilized before use.

We recommend autoclave steam sterilization at the hospital (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a piece of equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

Table 2 - Autoclave physical sterilization parameters

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
-------	-------------	-------------------------	-----------------------

Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

For further clarification, refer to instruction manual supplied with autoclave.

Another sterilization method: that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

Storage: Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls. Keep the in ventilated, dry place, away from light and severe weather action. As they are transported, inadequate bumps and stacking must be avoided. Store and transport in dry and ventilated place, with room temperature (max. 35°C) and relative humidity around 30% to 70%. Do not store directly on the floor (minimum height = 20 cm) or in very high places, close to lamps, as the packaging might get dried or the label might be damaged. Do not store in places where contaminant substances are stored, including, for instance, cleaning materials, insecticides, pesticides, etc.

Identification: Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile

Product' ;Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

Consumer Service: In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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