# KIT DE VERTEBROPLASTIA/ PROXIMAL DE FÊMUR

**Português**

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**English**

**PRODUCT DESCRIPTION:**

The Vertebroplasty/Proximal Femur Kit is composed of only five instrumental support, used to insert, by means of compression, biomaterial – bone liquid cement (Polimetilmetacrilato – PMMA) – non integrant of this product and not purpose of this register – into the weakened vertebra of the vertebral column or into the interior of the femoral head.

**PRESENTATION OF THE PRODUCT:**

All components of the Vertebroplasty/Proximal Femur Kit are commercialized in the same packaging, ready for use, being that those components are not commercialized individually. The Kit is packed in a Double Pack blister type that is made of PET (polyethylene terephthalate) non-toxic allowing contact with food and hospital products. Both blister packs are sealed with Tyvek and without holes. It is sealed on the blister by the hot-packing process. The sealed Double Blister Pack is also put in an external cardboard pack.

**TRACEABILITY:**

Only the injector receives laser marking, being engraved the company logo, the lot number, and the indication "1turn = 1ml." The other components of the kit receive no marking due to their anatomy. Still traceability of all kit components is guaranteed by the Kit lot number, present in the product traceability label which tracking is possible to identify the lot number of each component.

In the records the following information is essential: used implant name; date of surgery;

**KIT COMPONENTS DESCRIPTION, QUANTITY AND RAW MATERIAL:**

- 1 unit of Syringe without Needle-10ml (code: 102-00-010) made of Polypropylene,
- 1 unit of Extensorfix-20cm (code: 102-05-000) made of PVC (Cloreto de Polivinila );
- 1 unit of BM 8-15 Needle (code: 8570-40-151) made of stainless steel 304 Type and Policarbonate;
- 1 unit of the Injector for the Vertebroplasty (code: 8574-00-000) made of Polypropylene;
- 1 unit of the Funnel (code: 8575-00-000) made of PVC (Cloreto de Polivinila).

**No. 1 Component – Syringe without Needle-10ml (Code: 102-00-010)**

Syringe without needle,APIROGENICA, non-toxic, with capacity of 10ml, made of a polypropylene cylinder and stick (a semi-crystalline resin, produced through the polymerization of the propeno monomer, using a stereo specific catalyzer) and by a synthetic rubber piston. This component has the finality to store and inject the biomaterial (non integrant of this product and not the purpose of this register) in the weakened vertebra of the vertebral column.

**No. 2 Component - Extensorfix-20cm (Code: 102-05-000)**

PVC Pipe – Cloreto de Polivinila (termplastic polymer which contains, in weight, 57% of chlorine, Derived from chloride of sodium and 43% of eteno derived from petroleum), aprotogênico, non-toxic with 20 cm of length. This component has the finality of extending the syringe without needle-10ml, contributing to the injection of the biomaterial (non integrant of this product and not purpose of this register) into the weakened vertebra.

**No. 3 Component - BM 8-15 Needle (Code: 8570-40-151)**

Product made of an stainless steel needle 304 type and policarbonate (long polymer chain, formed by functional groups united by carbonate groups (-O-(C=O)-O) ), a punch also made of stainless steel 304 type and policarbonate, a string also made of stainless steel 304 type a lid also made of policarbonate of 6 cm length. The BM 8-15 needle has capacity of 8G for 150mm and has the function of introducing the biomaterial (non integrant of this product and not purpose of this register) into the weakened vertebra of the vertebr a column.

**No. 4 Component – Vertebroplasty Injector (Code: 8574-00-000)**

Component made of a polypropylene manufactured injector body (a semi crystalline resin, produced through the polymerization of the propene monomer, using a stereo specific catalyst), an injector fuse and a propylene manufactured maniple. The injector body has a total length of 135,0mm (+-1,0) and a diameter of 35,0mm (+-0,5), the injector fuse has the length of 140,0mm (+0,5) and diameter of 15,0mm (-0,2), and the maniple has 70mm by 24mm, being non-toxic e non-pyrogenic. The finality of this component is to help handle the syringe without neddel-10ml, this one fitting the injector.

**No. 5 Component - Funnel (Code: 8575-00-000)**

Funnel, made of PVC – Cloreto de Polivinila (termplastic polymer which contains, in weight, 57% of chlorine, Derived from chloride sodium chloride, and 43% of eteno, derived from petroleum), showing a diameter superior of 40,0mm and diameter inferior to 13,9mm. This component has the purpose of helping place the biomaterial (non integrant of this product and not purpose of this register) to the Syringe Without Neddle-10ml.

**Ancillary Components**

The Vertebroplasty/Proximal Femur Kit has the Acrilex Orthopedic Cement as an ancillary component, registered at ANVISA nº 10171110070, not being purpose of this register, and must be acquired separately. It is allowed the use of other kind of biomaterial as ancillary component for the Vertebrospstasy Kit, in substitution to the Acrilex Orthopedic Cement, since it can be injected in the vertebra (vertebroplasty) or into the interior of the femoral head being regularly registered at ANVISA for this purpose.

**Accessories**

The Kit does not have accessories.

**INDICATIONS**

The Vertebroplasty/Proximal Femur Kit is appropriate to minimal invasive surgeries, helping in the stability to fractured vertebras, reducing the pain and making better the activities of the weakened individuals. It must be used together with Acrilex Orthopedic Cement – ANVISA register nº 10171110070 (not part of this product and not present in this register) being allowed the use of other kind of biomaterial, in substitution to the Acrilex Orthopedic Cement, since it can be injected in the vertebra (vertebroplasty) and that it is regularly registered in ANVISA to this purpose.

The Vertebroplasty, which consists in the injection of bone liquid cement (Polimetilmetacrilato – PMMA) sterilized on the fractured vertebra bodies in the vertebral column, initially it was used for severe cervical pain treatment associated to hemangiome and posterior in the osteoporosis and osteolical neoplasty treatments. The Osteoporotical vertebral injury is the result of the weakened trabecular bone collapse or implosion. The Vertebroplasty is a method which presents good treatment results on fractures by osteoporosis, there is immediate improvement on symptoms in most cases, reducing the pain and improving patient´s life quality, besides avoiding complications originated from prolonged time in bed immobilized.

The purpose of vertebroplasty is to increase resistance and reduce pain symptoms, and not the correction of the deformity.

The ideal patient to be treated by vertebroplasty is that one which the conventional medication treatment has failed and it demonstrates to have a relatively new vertebral body fracture, preferable with six months duration, it must be observed:

- Back pain from vertebral compressive fractures caused by osteoporosis and refract of the medical treatment.
- Pain from vertebral fracture or severe osteolises from a tumor;
- Pain from vertebral fracture associated to osteonecrosis;
- Instable compressive fracture with edge movement;
- Multiple compressive fractures with collapse that can result in pulmonar or gastro-intestinal compromise;
- Severe Traumatic fractures in normal bones with lack of union among fragments;
- Internal Cistical Change

The same kit can be used in the proximal femur osteosynthesis, with bone cement injection into the femoral head via cannulated screws holes of intramedullary nails lock. It should be used in conjunction with orthopedic cement Acrilex® - registration at ANVISA No. 10171110070 (not part of this product and not the object of this register), and allowing the use of another type of biomaterial, replacing cement Orthopedic Acrilex® since It can be injected into the femur and is regularly registered with ANVISA for this purpose.

**COUNTER-INDICATIONS**

By being a Kit made only by support instruments to the vertebroplasty or osteosynthesis procedure, there are no counter-indications referring to the Kit, but there are for the vertebroplasty or osteosynthesis procedure, which includes the utilization of the biomaterial that will be injected. So all counter-indications are related to the biomaterial used and to the vertebroplasty procedure, and not to the Vertebroplasty/Proximal Femur Kit object of this register.

In this case, all counter-indications shown in the instructions of use of the biomaterial to be used and the

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**RG ANVISA: 80128580088**  
B020 - REV15 - 18122023

# BIOMECANICA VERTEBROPLASTY KIT

counter-indications related to the vertebroplasty procedure shall be verified.

The vertebroplasty procedure shall be avoided in presence of infections, non-coagulated disturbances, hipper-sensibility to the biomaterial´s components, the vertebral body´s posterior wall destruction must be evaluated once the overflowing potential of biomaterial injected into the spinal channel can generate compression of the nervous structures.

The Vertebroplasty/Proximal Femut Kit must be used only for the vertebroplasty or osteosynthesis procedures, not being indicated to other purposes.

#### ADVERSE EFFECTS

The Adverse effects are related to the implanted biomaterial through the Vertebroplasty/Proximal Femur Kit, that is, the instructions of use of the biomaterial used shall be verified. By being a Kit made only of support instruments to the vertebroplasty or osteosynthesis procedure, there are no adverse effects related to the Kit object of this register, but to the total vertebroplasty or osteosynthesis procedure.

In spite of the low level of complications on the vertebroplasty or osteosynthesis procedure, attention shall be paid to the following risks:

A. Bleeding to the punch location, generally solved by compression to its local,
B. Temporary worsening of the pain and fever as reaction to the temperature produced by the biomaterial's polymerization;
C. Infection;
D. Roots or marrow injuries by unsuitable overflowing of contrast, with potential risk of paralysis and/or radiculopathy;
E. Liquorical Fistula;
F. Passage of the biomaterial to the venous system that can cause pulmonary embolism;
G.Passage of the biomaterial to the nervous system that can cause paraplegia.

#### WARNINGS

A. The vertebroplasty procedure neither correct the bone loss nor restores the original anatomy In case of tumorous metastasises the result is more restrict;
B. The Vertebroplasty/Proximal Femur Kit is a one-way use product otherwise it can occur biological contamination or transmission of viral diseases.
C. The Vertebroplasty/Proximal Femur Kit is disposable. A Kit must be used for each level of vertebra. Its reprocessing is expressly forbidden.

D. The vertebroplasty procedure must be monitored through X-rays.
E. The patient must be informed of its surgical procedure.
F. The use of other accessories to integrate this kit is prohibited.
G. The bone marrow cavity preparation results in medullar contents entering into the blood current. Enlarging a hole in the bone marrow cavity, effects can be similar to the arterial pressure with the introduction of the biomaterial.
H. A strange body inserted into the tissues increases the normal infection risk after surgical operations.
I. The Kit components shall not be used if the pack is violated, with holes; sterilization term is expired, with stains or failures. The product must be kept in its original package until the moment of use, according to the medical hospital asepasis procedures. If the package is violated or damaged, all the kit components must be discarded.
J. All storage, and product transportation conditions must be respected.
K. The surgeon doctor must be familiarized and must have sufficient knowledge of the implants fixation surgical technique and its limitation, including the pre and post-operative, adopted surgical procedure and potential precautions and risks.

#### THE PATIENT MUST BE INFORMED ABOUT:

A. Indications, counter-indications, adverse effects and warnings related to the vertebroplasty procedure and to the biomaterial to be injected;
B. The patient must be informed about the importance of the post-operative attendance. The lack of attendance obstructs the detection of post-surgical problems;
C. The patient must be informed about the biomaterial to be used (brand, code, name, share, etc.).

#### JOINT OPERATION WITH OTHER MEDICAL PRODUCTS

The Vertebroplasty/Proximal Femur Kit, is marketed as a complete and ready to use Kit, not being necessary other components to integrate this one.

In case the surgeon doctor, due to any eventuality in the surgical act, needs to use any other type of accessory, it is recommended that it is from the same company to avoid incompatibility on the material and engineering of the product.

For more information about compatibility of the Vertebroplasty/Proximal Femur Kit with other components, the producer must be consulted.

#### SPECIAL CONDITION OF STORAGE, CONSERVATION AND/OR MANIPULATION OF THE PRODUCT

 • Do not use if the package is damaged.

 • Fragile, handle it with care.

- Open the Package aseptically.
- The product must be stored and transported in dry and fresh place, with ambient temperature (+35°C), relative humidity around 30% to 70% and hidden from direct sun light.
- They cannot be stored directly on the floor (minimum height = 20cm). They cannot be in places too high, next to bulbs, which can cause package dryness or damage to the label.
- The transporter must be informed about the content and delivery date.
- They cannot be stored in places which contemned substances are used, such as, cleaning products, insecticide, pesticide, etc.

#### PRECAUTIONS BEFORE USING THE PRODUCT

Before using the Kit, the responsible, as designated by the hospital unit, shall verify the package physical conditions and later the product.

The responsible must verify if the cardboard of the package has a security seal and identified with the information referring to the material content in the pack.

After verifying this, the responsible shall remove the security seal in the external package and open the cardboard package.

The product packed in Double blister Pack shall be removed from the external package with care, in a surgical ambient The first blister shall be open by pulling the Tyvek at the edge of the blister. The second blister shall also be open by pulling the Tyvek at the edge of the blister.

With the product in hands, it shall be verified if the same does not show any notches, scratches or stains. In case any adverse effect occur affecting the user's safety, such as not working product, the responsible surgeon must communicate this adverse event to the competent sanitary organ and to Biomecanica through: e-mail: sac@biomecanica.com.br phone: +55 0xx(14) 2104-7900

In case of doubt, the responsible surgeon or health professional can communicate the adverse event through the Sanitarium Notification Vigilance System at the ANVISA website:

<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

#### PRODUCT DISCARD

All components present in the Vertebroplasty /Proximal Femur Kit should be unserviceable and discarded after use.

We recommend that the pieces are honed, bent or cut to their destruction, should in then be disposed according to the rules of the local health authority as they are considered medical waste (potentially contaminant product).destruction.

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