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of health care products - Requirements for validation and routine control - Industrial moist heat sterilization).

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ENGLISH EUA

Detailed description of the medical product, including the fundamentals of its functioning and its action, its content or composition, when applicable, as well as the list of accessories intended to integrate the product.

The Special Plates for Bone Synthesis are produced in stainless steel, according to ASTM F139 standards. The stainless steel used in the factoring of Special Plates for Bone Synthesis is analyzed in an outsourced laboratory to verify the conformity of the material.

'T' 4,5 mm Plate

Plates in 'T' 4,5 shape are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being used to support metaphyseal fractures of the proximal humerus, of the tibia medial plateau and distal of shin-bone. Its shape has a double fold allowing an adaptation on the indicated bone region. It has an oblong hole that allows a temporary fixation of the plate to the bone, also allowing some adjusting up and down.

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“T” 3,5 mm Plate

Plates in “T” 3,5 shape are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being indicated for setting in the distal radio region. Its shape has a double fold allowing an adaptation of the distal radio. It has an oblong hole that allows a temporary fixation of the plate to the bone, also allowing some adjusting the sagittal plan.

Left or Right 'L' Plates

'L' shaped plates are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being indicated for setting in the tibia proximal region. Its shape that has a double fold allows the adaptation to the tibia lateral plateau. It has an oblong hole that allows a temporary fixation of the plate to the bone, also allowing some adjusting the sagittal plan.

Malleable Plate

Malleable Plates are indicated for bone where your anatomy is complex, such as the pelvis around the acetabulum and distal humerus. They have holes that will receive the Cortical Screws for bone fixation.

Shermann Plate

Shermann Plates are indicated for intercondylar fractures of the humerus, pseudarthrosis of the humeral shaft and mandible. They have holes that will receive the Cortical Screws for bone fixation.

1/3 Tubular Plate

Tubular 1/3 format plates are used with 3.5 or 4.0 clamping. Its thickness of 1.5mm has limitations in its ability to provide stability. But it is useful in areas with minimal coverage of soft tissues, such as the lateral malleolus, olecranon and distal end of the ulna, in addition to fixation on the fibula, radius and metacarpal. The holes in the 1/3 tubular plates allow fixation with a Ø3.5 cortical screw or a Ø4.0 cancellous screw.

Composition: Special Plates for Bone Synthesis are manufactured in stainless steel, according to the ASTM F139 standard specifications.

Presentation Form: Implants are conditioned in plastic envelopes with a non-sterile tag and external label. The external label is made of a product tracking tag, in which contains the following information: a) Producer's Name; b) Commercial Name of the Product; c) Technical Name; d) Lot Number; e) ANVISA registration number; f) Product Code; g) Quantity; h) Description of the product that is included in the package (model); i) Production Date; j) Expiration Date: Undetermined Expiration Date; k) Description of the raw material used in the production of the product (Composition); l) Producer's Address; m) Name of the Responsible Technician; n) The Saying: "Non-sterile Product"; o) The Saying: "PROHIBID TO REPROCESS"; p) "Single Use Product, do not reutilize"; ; q) Before using, read the usage instructions; r) Warnings/Precautions/Special Cares/Symbology: Read the usage instructions.

Table 1 - Symbology, according to NBRISO 15223

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	Medical Product for unique use. Do not re-use.
	“Read Instructions”
	Fragile. Handle it carefully
	Keep away from sunlight
	Keep dry

Description of the packaging material

The non-sterile product is identified with a green tag of 20mm diameter with the phrase
“Non-Sterile”

Ø 20,0 mm

Figure 1 - Non-sterile product identification label

The Special Plates will be distributed individually according to the table below:

Table 2 - Distribution of Special Plates for Bone Synthesis

Model	Code	Qty.	Description	Raw material
'T' 4.5 Plate 	3653-03-000	01	'T' Plate – 4.5 of 3 holes by 68.5mm length	Stainless Steel F139
	3653-04-000	01	'T' Plate - 4.5 of 4 holes by 84mm length	Stainless Steel F139
	3653-05-000	01	'T' Plate - 4,5 of 5 holes by 100mm length	Stainless Steel F139
	3653-06-000	01	'T' Plate - 4,5 of 6 holes by 116mm length	Stainless Steel F139
	3653-07-000	01	'T' Plate - 4,5 of 7 holes by 132mm length	Stainless Steel F139
	3653-08-000	01	'T' Plate - 4,5 of 8 holes by 148mm length	Stainless Steel F139
	3653-09-000	01	'T' Plate - 4,5 of 9 holes by 164mm length	Stainless Steel F139
	3653-10-000	01	'T' Plate - 4,5 of 10 holes by 180mm length	Stainless Steel F139
	3653-11-000	01	'T' Plate - 4,5 of 11 holes by 196,5mm length	Stainless Steel F139
	3653-12-000	01	'T' Plate - 4,5 of 12 holes by 212.5mm length	Stainless Steel F139
Model	Code	Qty.	Description	Raw material
Semi Tubular Plate 	3520-03	01	Semi tubular plate 03 holes and 55mm length	Stainless Steel F139
	3520-04	01	Semi tubular plate 04 furos e 71mm length	Stainless Steel F139

	3520-05	01	Semi tubular plate 05 holes and 87mm length	Stainless Steel F139
	3520-06	01	Semi tubular plate 06 holes and 103mm length	Stainless Steel F139
	3520-07	01	Semi tubular plate 07 holes and 119mm length	Stainless Steel F139
	3520-08	01	Semi tubular plate 08 holes and 135mm length	Stainless Steel F139
	3520-09	01	Semi tubular plate 09 holes and 151mm length	Stainless Steel F139
	3520-10	01	Semi tubular plate 10 holes and 167mm length	Stainless Steel F139
Model	Code	Qty.	Description	Raw material
'T' 3.5 Plate 	3545-35-033	01	"T" plate 3.5 of 3 holes by 3 holes of 48mm length	Stainless Steel F139
	3545-35-034	01	T" plate 3.5 of 3 holes by 4 furos de 58mm length	Stainless Steel F139
	3545-35-035	01	T" plate 3.5 of 3 holes by 5 holes of 68mm length	Stainless Steel F139
	3545-35-043	01	T" plate 3.5 of 4 holes by 3 holes of 45,5mm length	Stainless Steel F139
	3545-35-044	01	T" plate 3.5 of 4 holes by 4 holes of 55mm length	Stainless Steel F139
	3545-35-045	01	T" plate 3.5 of 4 holes by 5 holes of 65mm length	Stainless Steel F139
	3545-35-046	01	T" plate 3.5 of 4 holes by 6 holes of 75mm length	Stainless Steel F139
Model	Code	Qty.	Description	Raw material
Left or Right 'L' Plate 	3652-03-E00		Left 'L' plate with 03 holes by 69mm in length	Stainless Steel F139
	3652-04-E00	01	Left 'L' plate with 04 holes by 85mm in length	Stainless Steel F139
	3652-05-E00	01	Left 'L' plate with 05 holes by 101mm in length	Stainless Steel F139
	3652-06-E00	01	Left 'L' plate with 06 holes by 117mm in length	Stainless Steel F139
	3652-07-E00	01	Left 'L' plate with 07 holes by 133mm in length	Stainless Steel F139
	3652-08-E00	01	Left 'L' plate with 08 holes by 149mm in length	Stainless Steel F139
	3652-03-D00	01	Left 'L' plate with 03 holes by 69mm in length	Stainless Steel F139
	3652-04-D00	01	Left 'L' plate with 04 holes by 85mm in length	Stainless Steel F139
	3652-05-D00	01	Left 'L' plate with 05 holes by 101mm in length	Stainless Steel F139
	3652-06-D00	01	Left 'L' plate with 06 holes by 117mm in length	Stainless Steel F139
	3652-07-D00	01	Left 'L' plate with 07 holes by 133mm in length	Stainless Steel F139
	3652-08-D00	01	Left 'L' plate with 08 holes by 149mm in length	Stainless Steel F139
Model	Code	Qty.	Description	Raw material
Malleable Plate 	3535-03-000	01	Malleable Plate - 03 holes per 55mm length	Stainless Steel F139
	3535-04-000	01	Malleable Plate - 04 holes per 71mm length	Stainless Steel F139
	3535-05-000	01	Malleable Plate - 05 holes per 87mm length	Stainless Steel F139
	3535-06-000	01	Malleable Plate - 06 holes per 103mm length	Stainless Steel F139
	3535-07-000	01	Malleable Plate - 07 holes per 119mm length	Stainless Steel F139

	3535-08-000	01	Malleable Plate - 08 holes per 135mm length	Stainless Steel F139
	3535-09-000	01	Malleable Plate - 09 holes per 151mm length	Stainless Steel F139
	3535-10-000	01	Malleable Plate - 10 holes per 167mm length	Stainless Steel F139
Model	Code	Qty.	Description	Raw material
Shermann Plate 	3512-00-003	01	Shermann plate with 03 holes for 44mm length	Stainless Steel F139
	3512-00-004	01	Shermann plate with 04 holes for 61mm length	Stainless Steel F139
	3512-00-005	01	Shermann plate with 05 holes for 78mm length	Stainless Steel F139
	3512-00-006	01	Shermann plate with 06 holes for 95mm length	Stainless Steel F139
	3512-00-007	01	Shermann plate with 07 holes for 122mm length	Stainless Steel F139
	3512-00-008	01	Shermann plate with 08 holes for 129mm length	Stainless Steel F139
	3512-00-009	01	Shermann plate with 09 holes for 146mm length	Stainless Steel F139
	3512-00-010	01	Shermann plate with 10 holes for 163mm length	Stainless Steel F139
	3512-00-011	01	Shermann plate with 11 holes for 180mm length	Stainless Steel F139
	3512-00-012	01	Shermann plate with 12 holes for 197mm length	Stainless Steel F139
Model	Code	Qty.	Description	Raw material
1/3 Tubular Plate 	3530-02-035	01	1/3 Tubular Plate with 02 holes for 26mm length	Stainless Steel F139
	3530-03-035	01	1/3 Tubular Plate with 03 holes for 38mm length	Stainless Steel F139
	3530-04-035	01	1/3 Tubular Plate with 04 holes for 50mm length	Stainless Steel F139
	3530-05-035	01	1/3 Tubular Plate with 05 holes for 62mm length	Stainless Steel F139
	3530-06-035	01	1/3 Tubular Plate with 06 holes for 74mm length	Stainless Steel F139
	3530-07-035	01	1/3 Tubular Plate with 07 holes for 86mm length	Stainless Steel F139
	3530-08-035	01	1/3 Tubular Plate with 08 holes for 98mm length	Stainless Steel F139
	3530-09-035	01	1/3 Tubular Plate with 09 holes for 110mm length	Stainless Steel F139
	3530-10-035	01	1/3 Tubular Plate with 10 holes for 122mm length	Stainless Steel F139
	3530-11-035	01	1/3 Tubular Plate with 11 holes for 134mm length	Stainless Steel F139
	3530-12-035	01	1/3 Tubular Plate with 12 holes for 146mm length	Stainless Steel F139

Attention

List of Instruments for the use of Special Plates for Bone Synthesis

In order to use the Special Plates for Bone Synthesis, it is necessary to use the instruments specified below.

The Instrumental Kit for Placing Plates and Screws for Bone Synthesis is registered with ANVISA and is not the object of this registration and is not part of this product.

The tables below describe the list of instruments that make up the trays of boxes 4.5 and 3.5.

Table 3 - List of Instruments that make up the Instrument tray for Placing Plates and Screws for Bone Synthesis 4.5

Tray of instruments for placing the Plates and Screws of the Bone Synthesis 4,5		
Qtd	Code	PRODUCT
01	6019-0G-000	Hohmann Retractor Large
01	6011-00-000	Drilling Guide
01	6003-32-150	Drill Ø3.2 x 150.0mm
01	6003-45-150	Drill Ø4.5 x 150.0mm
01	6015-00-000	Hexagonal Wrench
01	6007-00-000	Compressor Wrench
01	6005-00-000	Clamp Compressor
01	6016-00-000	Countersink
01	6013-00-000	Handle Double Guide Of Perforation N/ex/std
01	6013-01-000	Drill Guide 4.5mm - STD
01	6013-02-000	Drill Guide 4.5mm - N
01	6013-03-000	Drill Guide 4.5mm - EX
01	6014-00-000	Tap Protector Guide 4,5
01	6002-00-000	Tap For Cortical Screw Ø 4.5mm
01	6001-00-000	Tap For Cancellous Screw 6.5mm
01	6004-00-000	Cortical Measure
02	6043-00-000	Self Centering Tweezer
02	6042-00-000	Reduction Tweezer
01	6017-45-000	Bending Iron Plate (pair)

Table 4 - Relation of the instruments that compose the tray of the instrument for placing the Plates and Screws to the Bone Synthesis 3,5

Tray of instruments for placing the Plates and Screws of the Bone Synthesis 3,5		
QTD	REF	PRODUCT
01	6019-0M-000	Hohmann Retractor Medium (pair)
01	6019-0P-000	Hohmann Retractor Small (pair)
03	6003-20-150	Drill Ø2.0 / 2.5 / 3.5 x 150.0mm
01	6036-00-000	Hexagonal Wrench
01	6033-00-000	Small Countersink
01	6013-00-000	Double Drill Guide 4.5mm N/EX/STD
01	6030	Drill Guide STD/ N/ EX Ø3,5
01	6039-00-000	Drill Guide Ø2.5
01	6029-00-000	Drill Guide Ø3.5
01	6021-00-000	Tap For Cortical Screw
01	6024-00-000	Cortical Measure
01	6017-00-000	Bending Iron Plate (pair)
01	6006-00-000	Tweezer For Screw

Care for Handling and Transportation of the Medical Product:

- Do not use the product if it is damaged.
- The implants must be handled with care avoiding damage that can spoil the quality of the material and the safety of the patient.
- Preserve it in a cool, dry place, away from light and away from the action of weather.

- While transporting it, shocks and inadequate stacks must be avoided.

Note: Implant components from different manufacturers must not be used; therefore we recommend that the products have the same origin. Date of manufacture, expiration date and product batch: SEE LABEL.

Guarantee: the guarantee of the Special Plates for Bone Synthesis is related to the observance of this instruction for use. Misuse, non-compliance with the indications for use, information on use, contraindications, care and precautions, warnings, possible adverse effects, packaging, sterilization, cleaning and decontamination recommended in this instruction for use

Indication, purpose or use for which the medical product is destined, as indicated by the manufacturer.

Bone fixation device with holes for screw fixation to the bone. The special plates are implants that are fixed together to the bone with the purpose of giving stability to the fractures, osteosynthesis in general, member corrections, reconstructions, neutralization or compression, pseudoarthrosis and other procedures described in scientific literature. The special plates were developed for specific places. They have an anatomical shape, corresponding to the place they will be applied. The special plates for Bone Synthesis are available in several shapes, lengths and number of holes and are not reusable. The special plate's project is based on the AO/ASIF group special plates.

“T” 4,5 Plates: it has the purpose of acting as a thin cortex support or to prevent a collapse of the spongy bone defect, in metaphyseal fractures of the proximal humerus, of the tibia medial plateau and distal of shin-bone.

- Semi Tubular Plates: Semi Tubular Plates are easily deformable. The main indications to these plates are radio and proximal ulna fractures, specially comminuted fractures of the olecranon and top boot distal tibia fractures.

- “T” 3,5 Plates: Plates in “T” 3,5 shape are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being indicated for setting in the distal radio region.

- Left or Right 'L' Plates: They were designed based on the AO/ASIF group 'L' plates. The plates in L shape are plates destined to support a thin cortex or to prevent a collapse of the spongy bone defect and are relatively thick being indicated for a fixation on the proximal tibia region. Its shape has a double fold allowing an adaptation on the lateral tibia plateau.

- Malleable Plate: are indicated for bone in which its anatomy is complex, such as the pelvis around the acetabulum, distal humerus.

- Shermann plate: indicated for intercondylar fractures of the humerus, pseudarthrosis of the humeral shaft and mandible.

- 1/3 Tubular Plate: indicated in areas with minimal soft tissue coverage, such as the lateral malleolus, olecranon and distal end of the ulna, in addition to fixation on the fibula, radius and metacarpal.

Orthopedic Implants Recommendations

To reach better results in trauma is recommended:

1 - The use of orthopedic implants must be done only by enabled surgeons and trained to such procedure.

2 – There are surgical instruments and tests available to help in the surgical plate implantation. It is important that the surgical instruments used be the specific projected to the determined product. Variation on the project and surgical instrument dimensions and similar tests can compromise the critical required measures for a precise implantation.

3- In order to protect the doctor and its patient, in future problems, the hospital must be held responsible for the annotation on the patient's protocol, of the code and the lot numbers of the used implants. These are fundamental data to allow the tracking of the former.

4 – Before starting the surgery make sure that the implant collection and the respective instruments are whole and complete.

5 – Never mix implants from different manufacturers in the same patient. There may be differences in dimensions and tolerances in the slots of these implants, resulting in poor joints, precocious looseness or failure of the former.

6 – The orthopedic implants are classified as being “Single use products”, which means, they cannot be reused.

7 – The surgeon must also avoid producing creases or scratches on the orthopedic implants; because such damage can produce internal “stress” which can become focus of an eventual breach.

Responsibility Exemption: Biomechanica, as a manufacturer of this device, neither practice medicine nor recommend this or any other surgical technique for use in a particular patient. Biomechanica is not responsible for the selection of appropriate surgical technique to be used in each patient. The surgical techniques depend on the technique adopted by the surgeon, and to engage its choice of method, type and size of products to be placed.

Ancillary Components

The Ancillary components related below are not objects of this register, having to be bought separately, but are associated to the implantation process, without which the product cannot be implanted.

'T' 4,5 Plate specific ancillary components

Cortical screw Ø 4,5mm “ST” tip; total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);

Conventional cortical screw Ø 4,5mm, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);
Spongy screw Ø 6,5mm, total or partial spiral, manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product).

Semi Tubular specific ancillary components

Cortical screw Ø 4,5mm "ST" tip, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);
Conventional cortical screw Ø 4,5mm, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard(not object of this registration and not part of this product).

'T' 3,5 Plate specific ancillary components

Cortical screw Ø 3,5mm "ST" tip, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);
Conventional cortical screw Ø 3,5mm, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);

Left or Right 'L' Plate specific ancillary components

Cortical screw Ø 4,5mm "ST" tip, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);
Conventional cortical screw Ø 4,5mm, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);
Spongy screw Ø 6,5mm, total or partial spiral, manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product);

Specific ancillary components Malleable Plate

Cortical Screw Ø 4.5mm conventional "ST" tip, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)
Conventional Ø4.5 Cortical Screw, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Specific ancillary components Shermann plate

Cortical Screw Ø 4.5mm conventional "ST" tip, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)
Conventional Ø4.5 Cortical Screw, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Specific ancillary components 1/3 Tubular Plate

Cortical Screw Ø3.5mm conventional "ST" tip, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Conventional Cortical Screw Ø3.5mm without “ST” tip, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Ø4.0mm conventional sponge screw with “ST” tip, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Ø4.0mm conventional sponge screw without “ST” tip, full thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Ø4.0mm conventional sponge screw with “ST” tip, partial thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Ø4.0mm conventional sponge screw without “ST” tip, partial thread manufactured in austenitic stainless steel according to ASTM F138 standard (not object of this registration and not part of this product)

Accessories: Special Plates for Bone Synthesis do not have any accessory with the purpose of integrating the medical product.

Allowable Combinations with other Materials: The chart “Allowable Alloy of Contact” was based on the “NBRISO21534 – Implants for Non-active Surgery – Implants for Articulation Substitution – Particular Requirements” standard, but independent of the alloy mentioned on the “Of Contact” item, it is necessary that the combinations are from implants of the same brand, because they were projected for such combinations observing the superficial finish and superficial treatment and other demanded factors in projects that can interfere in such combinations. Because of that, metallic implants from different manufacturers are not recommended for chemical, physical, biological and functional incompatibility reasons.

Table 5 – Allowable Alloy of contact

Allowable Alloy of contact	
Implant	*De contato (quando não terá articulação)
Plate (Stainless Steel ASTM F138 - 18 Cromo-14 Nickel - 2,5 Molybdenum)	Stainless steel adapted 18 Cromo-14 Nickel -2,5 Molybdenum (NBRISO 5832-1 and ASTM F138) Stainless steel adapted of high nitrogen (NBRISO5832-9) Adapted league of titanium 6-aluminium 4-vanadium (NBRISO 5832-3 and ASTM F136) Adapted league of titanium 5-aluminium 2,5-iron (NBRISO5832-10) Adapted league of titanium 6-aluminium 7-niobium (NBRISO5832-11)

*For application where a metal or alloy is in contact with others and is not intended joint, provided that adequate attention is paid to the design, surface finish, surface treatment and metallurgical conditions.

The Special Plates for Bone Synthesis from Biomecanica brand must be used only with Biomecanica instruments.

Precautions, restrictions e warnings

PRECAUTIONS: SINGLE USE. DESTROY IT AFTER EXPLANTED. DO NOT REUSE THE PRODUCT. REPROCESS IS FORBIDDEN.

An orthopedic implant can only be used in one patient, one single time. Although it may appear undamaged, previous stress may create imperfections that may reduce the implantation success. Improper implant selection may cause uncommon stress on the implant and can result in subsequent fractures of the implant.

Due to its shape of presentation and to its characteristics, this product does not show side effects, however, it is suggested to monitor the patient's recovery through x-ray according to the International Safety Rules. Do not use magnetic resonance imaging, as the product is made of stainless steel causes interference in the examination harming the diagnosis.

The use of Special Plates for Bone Synthesis is conceived to be an assist element to fixation and boné stabilization on the superior and inferior members fracture treatment and not to substitute normal human skeleton structures.

The surgeon must be sure of the possibility of the patient's physical development, being able to have a necessity of reviewing the implant.

The surgeon must be familiarized and have sufficient knowledge of osteosynthesis and its limitations, including the pre and post-operative, adopted surgery technique, precautions and potential risks.

The surgeon must have attention to the bone preparation for the introduction of the special plates, because the stability and fixation success may be compromised.

During handling of the implant, there must be no scratches, nicks, or anything that may damage or mark the implant, because these defects are stress concentrators and may be sites of crack nucleation and decrease the corrosion resistance and may result in fractures on the implant. In certain situations it is necessary to shape the implant to fit it perfectly to the bone anatomy to be implanted. Should it prove necessary to modeling the implant to the bone anatomy is the limit of modeling.

To guarantee an adequate implantation, should be used in instruments manufactured by Biomechanica, because they are designed and manufactured specifically for the use of these implants. Using instruments from different manufacturers may compromise the surgery as well as not being in Harmony with the registration of the product in Agencia Nacional de Vigilância Sanitária; do not use the product if it is damaged;

NON-STERILE PRODUCT – STERILIZE BEFORE USE

SINGLE USE PRODUCT – IT IS HIGHLY FORBIDDEN TO REUSE THIS PRODUCT;

Special storage conditions: Store in a cool, dry place, away from light and away from the action of the weather.

Do not use the product in case the package is damaged;

Observation: Implants' components from different producers shall not be used; so we recommend that products have the same origin.

Manufacture Date, expiration date and product lot: SEE LABEL.

Contraindication: The use of this implant is contraindicated under the following situations:

- Sensitivity to metal;
- Fever or sign of local inflammation;
- Osteoporosis;
- Obesity;
- Inadequate tissue cover in the region to be operated;
- Pregnancy;
- Bone compromising from illness and/or infections;
- Elevated sedimentation rate unexplained by other diseases, raising the rate of WBC or suspected striking in the differential count
- Other infection.

Adverse Effects: mechanical slackening can be a result of defective fixation, unstable reconstruction of the fractured bone or hidden infection. Reactions to metal sensitivity in patients were rarely informed.

Implantation of strange material in tissues results in histological reactions involving macrophages and fibroblasts of several sizes. The chemical importance of this effect is uncertain, but similar changes can occur as a precursor, or during the healing process.

Information to be given the patient: The patient must be informed about the importance of the postsurgical attendance. The lack of attendance obstructs the detection of post-surgical problems such as components looseness or osteolysis occurrence. Not performing the review surgery when there is component looseness or osteolysis can result in progressive loss of periprosthetic bone stock. The patient must be informed as to the type of implant being used, brand, code, name, lot, surgeon's name and date of surgery.

Evaluation of the implanted Product: after the implantation, intra-operative professional must perform a radiological control to verify the correct positioning of the product. The responsible professional must perform, and it is his responsibility, clinical and radiological evaluations done after the surgical procedure at a frequency stipulated by him to check the state of the implant and its boné consolidation evolution. In case the product is found out of the correct positioning, or presents any nonagreement; it is the surgeon's responsibility to take the most appropriate corrective measure.

Useful Information to Avoid Recurrent Implantation Risks: to decrease the risks due to implantation this should be strictly followed: contraindications, instructions for use and all the information contained in the product's "Use Instructions".

Implantation Intrinsic Risks: the Special Plates for Bone Synthesis of materials are recognized biomedical use, which follow the requirements of the ASTM F139 Standard, Standard

Specification for Wrought 18Chromium-14Nickel-2.5 Molybdenum Stainless Steel Sheet and Strip for Surgical Implants.

Contamination Risks: there are risks of biological contamination and viral diseases conveyance like hepatitis and HIV, because the metallic components of the Special Plates for Bone Synthesis get in contact with body tissues and fluids. Explanted products must be treated as highly contaminating.

Sterilization: this product supplied non-sterile. It must be sterilized before use. We recommend the steam sterilization in Autoclave at a hospital (ISO 11134: 1994 Sterilization of health care products – Requirements for validation and routine control – Industrial moist heat sterilization).

Product Discard: No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the Special Plates for Bone Synthesis explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Tracking: together with the implantable component's package are 2 tracking tags containing the data of the implant used. This tag is the same tag contained in the external label of the product. The information that is in this tracking tag is related in the 'Presentation Form' item of this technical report. This tag must be attached to the patient's handbook and the other for hospital control. The hospital must register this information in its files and give this information to the patient. We inform that the patient's handbook must have this information register which allows the implanted product to be tracked. Attaching this tag at the handbook allows that tracking of the implanted product used. Among the main information we emphasize as essential the name of the manufacturer of the implant used, name of the implant used, code of the implant, lot number and register number of the product at ANVISA. This information is written on the tracking tags that goes with the product and also on the external label. Other information also must be considered essential, as the date of the surgery, name of the patient who received the implant, name of the surgeon, patient's weight, patient's age and Other information requested on the patient's handbook are also to be filled.

Marking: The Special Plates for Bone Synthesis contains the following information laser recorded to allow the tracking and patient's attendance at the post-operative: -Biomechanica's Logo; -

Manufacturing Lot Number; -Material Acronym (Stainless Steel) used to manufacture the product (S1); -Code; -Dimension.

Post Sell (Client's Complain): in case there is a necessity of doing any complaints about the Special Plates for Bone Synthesis related to any adverse effect that affects the user's safety, such as non working product, damage of the implantable metallic component, serious problems or death related to these components the responsible surgeon must communicate this adverse event to the competent sanitary agency and to Biomecanica through the e-mail sac@biomecanica.com.br or through the phone 00xx14 2104 7900. In case of doubt the responsible surgeon or health professional can communicate the adverse event through the Sanitary Vigilance Notification System at the ANVISA web-site: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Send of Material for the Manufacturer to Analyze: In case of sending the implant for the manufacturer to analyze, this must be cleaned at the hospital using a bactericide solution and large specter antifungal. Then it must be disinfected or steam sterilized in autoclave or in ethylene oxide. Must be send to Biomecanica inside full package, identified with the clean method, sterilization and data of the product.

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