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ENGLISH EUA

Quantity and Description of Labels and Traceability Cards that Accompany the Product:

- 6 Component identification and traceability adhesive labels;

	XXXXXX	✓ Product Lote Number
	XXXX-XX	✓ Manufacturing Date
	XXXX-XX	✓ Use by
	XXXX-XX	✓ Code

➤ **Expiration Date:** XX/XXXX

➤ **Sterile Product Identification:**

Steril

Sterilization in ETO Gas

➤ **Single use product indication:**

“SINGLE USE PRODUCT - PROHIBITED TO REPROCESS”

➤ **Special Conditions for Storage, Conservation and/ or Handling of the Product:**

Open the Packaging Aseptically	
 35°C °C Temperature Limits	Transport and store it in dry and cool places and at room temperature (maximum 35°C) and Relative humidity around 30% and 70%.
Do not store it in direct contact to the floor (minimum height = 20cm) and neither in very high places and close to lamps because that may lead to the dryness and damage to the label.	
Do not store it in places where contaminants substances are stored such as cleaning materials, insecticide, pesticide, etc..	

	Attention Instructions, warnings and precautions: Vide instruction for use
	Do not use it if the packaging is damaged
	Fragile, carefully handle it
	Keep it away from the Sun
	Keep it dry
	"VIDE INSTRUCTION FOR USE"

2.2. Product Presentation:

The product is packaged in double blister-type packaging, which is made of non-toxic PET (Polyethylene Terephthalate), allowing contact with food and hospital products.

The two blister packs are sealed with dye-free and puncture-free Tyvek using the heat seal process. The double sealed blister pack is also enclosed in an outer cardboard carton. The implants are sterilized in ETO with a validity of 5 years.

2.3. Traceability:

Along with the product, there are 6 labels containing the necessary data for product traceability. A label must be obligatorily affixed to the patient's clinical record, another to the report delivered to the patient, another to the tax documentation that generates the charge, on the AIH, in the case of a patient treated by the SUS, or on the sales invoice, in the case of patient assisted by the complementary health system, another available for the control of the supplier (distribution history record - RHD, another available for the control of the responsible surgeon (main) and the last one for medical insurance, if any.

The following information is essential in the medical record: name of the implant used; date of surgery; product code; lot number; name of the patient who received the implant; surgeon's name. All models of Femoral Component with Restriction - BM receive marking, through laser, following the sequence below within the limits allowed by the size of the model.

Biomechanica Logo - Manufacturing Batch Number - Material Symbol used in the manufacture of the product – Code – Size

3. Product Description:


The Femoral Component with Restriction - BM is a non-absorbable, non-active implant intended for both primary and revision surgery, for replacement or restitution of the femoro-tibial and femoro-patellar joint (knee) with sacrifice of the posterior cruciate ligament, since it has a restriction system.

Restricted Femoral Component Models - BM


3.5. Restricted Femoral Component Models - BM

3.5.1. Models Type BKS III

All Type BKS III models are manufactured from a cobalt-chromium-molybdenum cast alloy in accordance with ASTM F75. All models have external polished tracks for articulation with the Patellar Component - BM, BKS III models, not object of this registration, Biomecanica brand. The internal surface was designed for fixation in the distal portion of the femur by means of acrylic bone cement (not object of this record), Biomecanica brand. In the intercondylar region, there is a blocking system (restriction) that requires the sacrifice of the posterior cruciate ligament.

The only difference between the BKS III models intended for primary surgery and the BKS III models intended for revision surgery is that the models intended for revision surgery have a cylinder for connection with the Femur Revision Nail - BM, BKS III models, not object of this registration, brand Biomecânica. This rod is compressed inside the cylinder, which in turn has a threaded lock system, which keeps the rod connected to the Femoral Component with Restriction - BM, model BKS III. This thread lock is manufactured from titanium 6-aluminum 4-vanadium alloy and is shipped with the revision Femoral Component as part of this product.

When fitting the Revision Femur Stem - BM, models BKS III, in the Femoral Component with Restriction - BM, Models BKS III, the locking system must be tightened by means of a 3.5mm hexagon spanner. The Revision Femur Nail - BM, models BKS III, as well as the Femoral Component with Restriction - BM, Models BKS III, will be fixed to the femur by means of bone cement, therefore the locking system only serves to maintain the Revision Femur Nail - BM , models BKS III, fixed to the Femoral Component with Restriction - BM, models BKS III, until their cementation in the femur. After cementation, the locking system is no longer important since the bone cement will cement the two prostheses together in the bone. The Femoral Component with Restriction - BM, models BKS III, are commercialized in models for primary surgery and in models for revision surgery, in sizes Extra Large, Large, Medium, Small and Extra Small, and for each size it is subdivided into Right and Left, as shown in Table 1. This variation in dimensions allows the surgeon to choose the appropriate model, according to the patient's bone characteristics and the surgical technique adopted.

Table 1: List of BKS III type models

Code	Description	Surgery	Side
2340-11-000	BKS III - Extra Small Femoral Component - Right	Primary	Right
2340-01-000	BKS III - Small Femoral Component - Right		
2340-03-000	BKS III - Middle Femoral Component - Right		
2340-05-000	BKS III - Large Femoral Component - Right		
2340-07-000	BKS III - Extra Large Femoral Component - Right		
2348-11-000	BKS III - Femoral Component - Extra Small Right (Revision) with locking system	Revision	
2348-01-000	BKS III - Femoral Component - Small Right (Revision) with locking system		

2348-03-000	BKS III - Femoral Component - Middle Right (Revision) with locking system		
2348-05-000	BKS III - Femoral Component - Large Right (Revision)		
2348-07-000	BKS III - Femoral Component - Extra Large Right (Revision) with locking system		
2340-12-000	BKS III - Extra Small Femoral Component - Left	Primary	Left
2340-02-000	BKS III - Small Femoral Component - Left		
2340-04-000	BKS III - Middle Femoral Component - Left		
2340-06-000	BKS III - Large Femoral Component - Left		
2340-08-000	BKS III - Extra Large Femoral Component - Left		
2348-12-000	BKS III - Femoral Component - Extra Small Left (Revision) with locking system	Revision	
2348-02-000	BKS III - Femoral Component - Small Left (Revision) with locking system		
2348-04-000	BKS III - Femoral Component - Middle Left (Revision) with locking system		
2348-06-000	BKS III - Femoral Component - Large Left (Revision) with locking system		
2348-08-000	BKS III - Femoral Component - Extra Large Left (Revision) with locking system		

3.6. Ancillary Components

The Femoral Component with Restriction - BM, models type BKS III, have as ancillary components the following products according to Table 2 (Not objects of this registration and not part of this product).

Table 2: List of Ancillary Components of type BKS III Models

Description	Raw Material
Fixed Articular Insert with Stabilizer - BM (BKS III models)	Ultra High Molecular Weight Polyethylene (UHMWPE)
Tibial Component - BM (BKS III models)	Cast Cobalt-Chromium-Molybdenum Alloy (ASTM F75)
Patellar Component - BM (BKS III models)	Ultra High Molecular Weight Polyethylene (UHMWPE)
*Acrilex® Orthopedic Cement	Bone Cement (NBR ISO 5833)

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3.7. Accessories

The Femoral Component with Restriction - BM, models 2348-11-000, 2348-01-000, 2348-12-000, 2348-02-000, 2348-03-000, 2348-04-000, 2348-05-000, 2348-06-000, 2348-07-000 and 2348-08-000, intended for revision surgery, have the following products as accessories as per Table 4 (Not objects of this registration and not part of this product).

Table 4: List of Accessories for type BKS III Models for Revision Surgery

Description	Raw Material
Femur Review Nail - BM (BKS III models)	Cast Cobalt-Chromium-Molybdenum Alloy (ASTM F75)

3.8. Support Instruments

There are instruments (not objects of this record and not part of this product) that serve as support for the total knee arthroplasty procedure.

5221-00-000 - Flat Milling Cutter (Rough), 5321-00-000 - Alignment Bar, 5320-00-000 - Alignment Bar with Coupler, 5336-00-000 - Tibial Drill Guide, 5355-02-000 - Patellar Forceps 3 Handles, 5335-00-000 - Handle for Tibial Tray, 5352-00-000 - Tweezers to remove from the tibial plateal, 5334-01-000 - Tibial Tray - Small, 5334-02-000 - Tibial Tray - Medium, 5334-03-000 - Tibial Tray - Large, 5334-04-000 - Tibial Tray - Ex. Large, 5341-00-000 - Test Tibial Impactor, 5340-00-000 - Test Tibial Extractor, 5322-00-000 - Femoral A/P Gauge, 5337-00-000 - Small / Medium Tibial Reamer, 5338-00-000 - Large / Extra Large Tibial Reamer, 5319-00-000 - Angular Support, 5326-00-000 - Extender for Tibial Cutting Guide, 5313-00-000 - Femoral Alignment Guide, 5328-00-000 - Depth gauge for tibial cutting, 5339-00-000 - Handle for Tibial reamer, 5354-00-000 - Pin Extractor Pliers, 5318-00-000 - Distal Femoral Cutting Guide, 5311-00-000 - Drill Ø14, 0 mm, 5310-00-000 - Stepped Drill, 5324-00-000 - G/EG Intercondylar Guide, 5364-00-000 - Punch, 5325-00-000 - Tibial Guide Rod, 5344-01-000 - Test Tibia Small, 5344-02-000 - Medium Test Tibia, 5344-03-000 - Large Test Tibia, 5344-04-000 - Extra Large Test Tibia, 5351-00-000 - Tibial Cutting Guide, 5350-00-000 - Tibial Ang. Cutting Guide, 5348-00-000 - Femoral Impactor, 5349-00-000 - Tibial Impactor, 5345-08-000 - Plateau Teste Small Medium -08 mm, 5345-10-000 - Small Medium Test Plateau -10mm, 5345-12-000 -Small Medium Test Plateau -12 mm, 5345-15-000 - Small Medium Test Plateau -15 mm, 5346-08-000 - Large / Extra Large Test Plateau - 08 mm, 5346-10-000 - Large / Extra Large Test Plateau -10mm, 5346-12-000 - Large / Extra Large Test Plate -12 mm, 5346-15-000 - Large / Extra Large Test Plateau - 15 mm, 5341-01-000 - Small Right Test Femur, 5341-02-000 - Small Left Test Femur, 5341-03-000 - Medium Right Test Femur, 5341-04-000 -Femur Medium Left Test, 5341-05-000 - Large Right Test Femur, 5341-06-000 - Large Left Test Femur, 5341-07-000 - Extra Large Right Test Femur, 5341-08-000 - Extra Large Left Test Femur, 5316-00-000 - Anterior Femoral Support Guide, 5323-01-000 - Cutting Guide Small Femoral, 5323-02-000 - Medium Femoral Cutting Guide, 5323-03-000 - Large Femoral Cutting Guide, 5323-04-000 - Extra Large Femoral Cutting Guide, 5342-08-000 - Spacer Block 8, 0 mm, 5342-10-000 - Spacer Block 10.0 mm, 5342-12-000 - Spacer Block 12.0 mm, 5342-15-000 - Spacer Block 15.0 mm, 5315-00-000 - Cutting Guide for Anterior Femoral, 5327-00-000 - Cutting Guide Initial Tibial, 5312-00-000 - Rotational Femoral Alignment Guide, 5363-00-000 - Tibial Retractor, 5362-00-000 - Patellar Tendon Retractor, 5317-00-000 - Distal Femoral Cutting Guide, 5314-00-000 - Impactor, 5359-00-000 - Universal extractor, 5361-00-000 - Hohmann type retractor (24.5x252 mm), 5360-00-000 - Hohmann type retractor (24.5x252 mm), 5367-00-000 - Impactor extractor wrench Screw, 5228-00-000 - Osteotome 20 mm, 5307-00-000 -Quick Coupling, 5372-00-000 - Screw Tip, 5333-00-000 - Tibial Recession Level Guide, 5211-00-000 - Osteotome 22 mm, 6003-02-032 - Drill Ø3,2x150 mm with quick coupling, 5381-00-000 - Patella guide 3 pegs, 5383-00-000 - Test Femur Extractor, 5382-00-000 - Caliper, 5381 - 01-000 - Patellar Drill Ø7.0 mm, 5384-00-000 - Patellar Clamp, 5385-00-000 - Extramedullary Tibial Cutting Guide, 5386-0 0-000 - Plateau Introducer Forceps, 5313-12-000 - Bar for femoral alignment guide -120 mm, 5313-15-000 - Bar for femoral alignment guide -150 mm, 5313-20-000

- Bar for guide of femoral alignment –200 mm, 5324-10-000 - S/M intercondylar guide, 5353-28-000 - Test patella Ø28 mm (3 pegs), 5353-30-000 - Test patella Ø30 mm (3 pegs) s), 5353-32-000 - Test patella Ø32 mm (3 pegs), 5353-34-000 - Test patella Ø34 mm (3 pegs), 5353-36-000 - Test patella Ø36 mm (3 pegs) ´s), 5211-16-000 - Osteotome 16mm, 5211-25-000 - Osteotome 25mm.

4. Indications, Contraindications, Adverse Effects and Warnings

4.1. Indications

The Restricted Femoral Component - BM aims to replace or restore the femoro-tibial and femoro-patellar joints of the knee, together with their ancillary components, when these joints are no longer able to maintain their functional capacity due to deformity, pain and movement blockage. This procedure is known as “total knee arthroplasty”, and seeks to create an artificial joint that has the same functionality as the natural knee. and knee stabilization.

The Femoral Component with Restriction - BM has a restriction system (blocking), that is, it requires the sacrifice of the posterior cruciate ligament during surgery.

Pain caused by chondrocalcinosis and pseudogout in an elderly patient is an occasional indication for arthroplasty in the absence of complete loss of cartilaginous space. In cases in which the flexion contracture advances beyond 20°, gait becomes very impaired, and arthroplasty is also indicated. In cases of laxity, knee arthroplasty is also indicated, but the use of prostheses with restriction in these cases is only feasible if the laxity is severe, otherwise unrestricted prostheses are recommended. For elderly patients with advanced osteoarthritis, the best treatment is total knee replacement. In younger patients, total knee replacement is indicated only in cases of tricompartmental disease (medial and lateral tibiofemoral and patellofemoral), requiring future revision surgery. In cases of rheumatoid arthrosis, total knee replacement is indicated for all ages, due to the complications presented by the disease.

The Restricted Femoral Component - BM is indicated in total knee arthroplasty procedures, without preservation of the posterior cruciate ligament, due to:

- Osteoarthritis: direct friction between the bones caused by the wear and tear of the joint cartilage that covers them, generating an extensive picture of pain, movement limitation and deformity;
- Rheumatoid Arthritis: a rheumatic disease that causes joint pain and deformity in general;
- Consequences of comminuted fractures: many fragments of the knee where it is practically impossible to achieve good coaptation.

4.2. Contraindications

The contraindications below do not limit the use of the Femoral Component with Restriction - BM models, but must be analyzed by the surgeon:

- C. Gross weakness of the quadriceps;
- D. Systemic infections or diseases close to the joint;
- AA. Knee rearthrodesis after arthroplasty;
- BB. Prior patellectomy;
- CC. Presence of recurvatum deformity;
- DD. Neuropathic arthritis;
- EE. Severe discontinuity of the extensor mechanism;
- FF. Recurvatum deformity due to muscle weakness and the presence of a well-functioning painless knee arthrodesis;
- GG. In younger patients, total knee replacement is indicated Only in cases of tricompartmental disease (medial and lateral tibiofemoral and patellofemoral) requiring future revision surgery;
- HH. Morbid obesity;
- II. Fever or leukocytosis;
- JJ. Inability of the patient to receive anesthesia;
- KK. Atherosclerotic leg disease ;
- LL. Skin problems such as psoriasis;
- MM. Neuropathic arthropathy;
- NN. Recurrent urinary tract infections;
- OO. History of osteomyelitis near the knee;
- PP. Severe neuromuscular, muscular or vascular deficiencies in the affected region;
- QQ. Severe injury to the bone structures generating an unfeasible risk for the application of the prosthesis;
- RR. Improper state of the ligamentous system component model;
- SS. Bone tumors close to the implantation site of the component;
- TT. Muscle atrophy;
- UU. Bone immaturity;
- VV. Insufficiency of the knee extensor mechanism;
- WW. Structural soft.tissue failure of the knee;
- XX. Knee ankylosis;
- YY. Insufficient bone structure in the tibia, fêmur and patella;
- ZZ. Patients unwilling or unable to follow post-operative care and instructions.

4.3. Adverses Effects

- Q. Joint line migration;
- R. Change in gait pattern;
- S. Patellar tendo rupture;
- T. Decrease or increase in patelar height;
- U. Hematomas;

- V. Suture dehiscence;
- W. Localized and asymptomatic progressive bone resorption (osteolysis) around the component;
- X. Loosening migration or fracture of the medical device;
- Y. Superficial or deep infection;
- Z. Vascular disorders, including thrombosis and pulmonary embolism.
- AA. Limb shortening due to bone resorption;
- BB. Allergic reaction to the foreign body, which may result in histological reactions involving various types of macrophages and fibroblasts;
- CC. Pain, discomfort or abnormal sensations due to the presence or migration of the medical product, and also due to the surgical procedure;
- DD. Neural or neurological damage due to surgical trauma (including paralysis and soft tissue injuries);
- EE. Inability to perform normal daily physical activities;
- FF. Periprosthetic fractures.

4.4. Warnings

The Femoral Component with Restriction - BM brand Biomecanica must be used only with instruments brand Biomecanica. Implant components from different manufacturers should not be used, so we recommend that the products have the same origin.

- M. The Restricted Femoral Component - BM is intended to replace the damaged joint and not normal structures of the human body;
- N. The Femoral Component with Restriction - BM has a restriction system (blocking), which requires the sacrifice of the posterior cruciate ligament during surgery;
- O. Surgical implants should never be reused, and explanted implants should never be implanted again. Small defects and internal pressure patterns can occur, even if the implant appears intact, causing implant failure and rupture;
- P. The correct selection of the Femoral Component with Restriction - BM model is extremely important in accordance with the patient's bone characteristics, analyzed by the surgeon, and with the surgical technique adopted by the surgeon;
- Q. Wear of the ultra high molecular weight polyethylene (UHMWPE) ancillary components may occur due to poor positioning of the Restricted Femoral Component - BM, inadequate ligament balance or aseptic loosening of the prosthesis. For this reason, it is important that the surgeon has sufficient knowledge about the surgical technique adopted, as well as the implants and other instruments used;
- R. The infections often generated are linked to comorbidity;
- S. Aseptic failure of total knee arthroplasty can be caused by component loosening, wear of the polyethylene manufactured ancillary component, ligament laxity, periprosthetic fracture, arthrofibrosis, and patellofemoral complications;

- T. The results presented in revision surgeries are inferior compared to the results presented in primary surgeries;
- U. Looseness, both in flexion and extension, is generally caused by insufficient size of the prosthesis implanted too proximally to the femur;
- V. Patients who smoke are more likely to have longer operative times and higher costs;
- W. The results are more reserved in cases of absence of joint cartilage before surgery;
- X. Knee replacement has a finite expected survival, adversely affected by activity level.

4. The patient must be informed of:

- G. The need for periodic medical follow-up, in order to observe possible alterations in the status of the implant and the adjacent bone. Only monitoring can detect possible loosening of components or the occurrence of osteolysis;
- H. Failure to perform revision surgery when loosening components or osteolysis may result in progressive loss of periprosthetic bone stock;
- I. The patient must be fully understood and warned that the product does not replace and does not have the same performance as normal bone;
- J. Information related to topics: indications, contraindications, adverse effects and warnings;
- K. The patient must be instructed to inform that he has a metallic implant in case he undergoes the Magnetic Resonance exam, as the implant manufactured from metal causes interference in the exam, jeopardizing the diagnosis, and the exam may cause displacement of the implant, causing problems for the patient;
- L. The patient must be advised of the need to attend regular post-operative follow-up exams during the entire period in which the product remains implanted.

5. Joint operation with other medical products

According to the NBR ISO 21534 standard - Implants for Non-Active Surgery - Implants for Joint Replacement - Particular requirements, the molten cobalt-chromium-molybdenum alloy, raw material of the Femoral Component with Restriction - BM is acceptable for the manufacture of implants and can be articulated with ultra high molecular weight polyethylene (UHMWPE) implants, raw material for the Fixed Articular Insert with Stabilizer - BM. Still according to the standard, the titanium alloy 6-aluminum 4-vanadium, raw material for the threaded lock, is also acceptable for the manufacture of implants and can establish contact with cobalt-based alloys, raw material for the Femoral Component with Constraint - BM. A list of acceptable combinations was formulated with the Femoral Component without Restriction - BM (Table 4) and with the threaded lock (Table 5), models type BKS III (Table 3). Regardless of this relationship, the combinations must also be of implants of the same brand, as they were designed for such combinations, observing surface finish and surface treatment and other factors required in projects that may interfere with the combinations.

Table 4: Permissible Alloys for Contact with the Femoral Component with Restriction - BM, models type BKS III

Admissible contact alloys	
*Articulation	**Contact (when there will be no articulation)
- UHMWPE Ultra High Molecular Weight Polyethylene (NBR ISO 5834-1, NBR ISO 5834-2, ASTM F648); - Conformed Cobalt-Nickel-Chromium-Molybdenum Alloy (NBR ISO 5832-12).	- Titanium based alloy (NBR ISO 5832-3, NBR ISO 5832-11, ASTM F136); - Cobalt-based alloy (NBR ISO 5832-4, NBR ISO 5832-5, NBR ISO 5832-6, NBR ISO 5832-7, NBR ISO 5832-8, NBR ISO 5832-12, ASTM F75) - Stainless Steel (NBR ISO 5832-9)

**For articulation surfaces of implants for articulation replacement, they were considered acceptable in specific applications, provided that adequate attention was paid to the design, surface finish and surface treatment.*

***For applications where one metal or alloy is in contact with another and no articulation is intended, provided proper attention is paid to design, surface finish, surface treatment and metallurgical conditions.*

Table 5: Admissible Contact Alloys with the Locking System present in the Femoral Component with Restriction - BM models type BKS III

Admissible contact alloys	
*Articulation	**Contact (when there will be no articulation)
- UHMWPE Ultra High Molecular Weight Polyethylene (NBR ISO 5834-1, NBR ISO 5834-2, ASTM F648);	- Cobalt-based alloy (NBR ISO 5832-4, NBR ISO 5832-5, NBR ISO 5832-6, NBR ISO 5832-7, NBR ISO 5832-8, NBR ISO 5832-12, ASTM F75) - Stainless Steel (NBR ISO 5832-1 and NBR ISO 5832-9)

**For articulation surfaces of implants for articulation replacement, they were considered acceptable in specific applications, provided that adequate attention was paid to the design, surface finish and surface treatment.*

***For applications where one metal or alloy is in contact with another and no articulation is intended, provided proper attention is paid to design, surface finish, surface treatment and metallurgical condition.*

8. Precautions before using the product

Before using the product, the person in charge, as designated by the hospital unit, must check the physical conditions of the packaging and subsequently of the product.

The person in charge must verify that the outer cardboard packaging is sealed with a security seal and identified with information regarding the material present in the packaging. After this verification, the person in charge must remove the security seal and open the cardboard packaging. The product packaged in a double blister should be carefully removed from the outer packaging, in a surgical environment. The first blister should be opened by pulling the surgical paper on its edge at the edge of the blister. The second blister should also be opened by pulling the remaining surgical paper on the edge of the blister. With the product in hand, check that it does not have dents, scratches or stains.

In the event of any adverse effect that affects user safety, such as product not working, damage to implantable components, serious problems or death related to these components, the surgeon in charge must report this adverse event to the competent health agency and Biomechanics through:

e-mail: sac@biomecanica.com.br

phone: (14) 2104-7900

In cases of doubt, the surgeon in charge or the health professional may communicate the adverse event through the Health Surveillance Notification System on the ANVISA website:

<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

9. Product disposal

No surgical implant should be reused. Any metallic implant, after use, must be discarded. Even if it appears to be undamaged, it may have small defects that can lead to fatigue fractures in addition to contamination. Thus, the Femoral Component with Restriction - BM explanted from patients must be properly discarded by the hospital institution. The hospital institution is responsible for the complete mischaracterization of the implant, preventing its reuse. The hospital institution is responsible for the method used to mischaracterize the implant. Biomecânica recommends that the explanted implants be mechanically deformed with the aid of a hammer or impact press, and must then be identified with the saying "Improper for Use". As described in Resolution No. 2605, of 08/11/2006, implantable devices of any nature classified as single use are prohibited from being reprocessed.

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