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## ENGLISH EUA







**Detailed description of the medical device, including the fundamentals of its operation and its action, its content or composition, when applicable, and the list of accessories that integrate the product.**

BM Non-Cannulated Screws for Small and Large Fragments are implants that are attached to the bone in order to provide fixation, stabilization and correction of anatomic form of the bone (integrating the set: screw, plate and bone). They are made of stainless steel.

The chemical, mechanical and metallographic composition of these materials are defined by the standard ASTM F138. It has a threaded body and a cylindrical head with a hexagon underside to allow the introduction into the bone with the use of the wrench. Comprehensively, orthopedic screws can be classified as cortical screws, when the screw is fixed onto the outer and highest density region of the bone (cortical bone) or cancellous screws when the setting is predominantly in the inner and lower density region (spongy bone).

Cancellous screws have more deeply cut and more widely spaced threads compared to cortical screw. They are available in various sizes. BM Non-cannulated Screws for Large and Small Fragments are not reusable. They are divided into the following models:

**Table 1 - Screw Models**

Models				
	Full thread	Full thread	Partial thread	Partial thread
Fitting Types				

**Table 2 – BM Non-Cannulated Screws for Small and Large Fragment**

Produto	Produto
6.5mm Metallic Cancellous Screw (stainless steel) 	6.5mm Metallic Fully Threaded Cancellous Screw (stainless steel) 
Produto	Produto
4.5mm Metallic Malleolar Screw (stainless screw)	4.5mm (P-ST) Metallic Cortical Screw (stainless screw)

	
<b>Produto</b>	<b>Produto</b>
4.5mm Metallic Cortical Screw (stainless screw)	3.5mm (P-ST) Metallic Cortical Screw (stainless screw)
	
<b>Produto</b>	<b>Produto</b>
3.5mm Metallic Cortical Screw (stainless screw)	4.0mm Metallic Cancellous Screw (stainless steel)
	

### Composition


BM Non-Cannulated Screws for Small and Large Fragments are manufactured in stainless steel according to the standard ASTM F138. The stainless steel used in their manufacture is analyzed in outsourced laboratories to verify if the material complies with ASTM F138 for stainless steel.

### Product Presentation

The implants are packaged in plastic envelopes containing two labels: one of them informing it is a non-sterile product and the other containing the product information: a) Manufacturer b) Product name, c) Technical Product Name d) Batch number e) Registration Number at the Ministry of Health, f) Product Code g) Quantity; h) Product Description i) Manufacture date j) Expiration Date: Expiration Undetermined l) Description of the raw material used to manufacture the product, m) Manufacturer's Address (described in the plastic envelope); n) Name of the Responsible Technician o) and the statements: "Non-sterile"; "No Reprocessing"; "Product of single use, do not reuse it"; "Before use see instructions"

### Symbolism of the Traceability Label according to NBRISO 15223

Table 3 – Symbology according to NBR ISO 15223

Símbolo	Descrição
	Manufacturer date
	Do not use if package is damaged
	Expiration Date
	Medical Product for unique use. Do not re-use
	Read Instructions
	Fragile - Handle it carefully

	Keep away from sunlight
	Keep Dry

**Description of the packing material**  
**The non-sterile product is identified with a green label with the words “Non-sterile”**

Ø 20,0 mm

**Figure - Non-sterile product identification label**

Large Screws and small fragments, non-cannulated - BM will be distributed in a unitary form according to the table below:

**Table 4 - Distribution of large screws and small fragments, non-cannulated - BM**

Code	Qty	Description
3012-16-025	1	Cancellous screw Ø6.5mm x 25mm in length with 16mm thread in length
3012-16-030	1	Cancellous screw Ø6.5mm x 30mm in length with 16mm thread in length
3012-16-035	1	Cancellous screw Ø6.5mm x 35mm in length with 16mm thread in length
3012-16-040	1	Cancellous screw Ø6.5mm x 40mm in length with 16mm thread in length
3012-16-045	1	Cancellous screw Ø6.5mm x 45mm in length with 16mm thread in length
3012-16-050	1	Cancellous screw Ø6.5mm x 50mm in length with 16mm thread in length
3012-16-055	1	Cancellous screw Ø6.5mm x 55mm in length with 16mm thread in length
3012-16-060	1	Cancellous screw Ø6.5mm x 60mm in length with 16mm thread in length
3012-16-065	1	Cancellous screw Ø6.5mm x 65mm in length with 16mm thread in length
3012-16-070	1	Cancellous screw Ø6.5mm x 70mm in length with 16mm thread in length
3012-16-075	1	Cancellous screw Ø6.5mm x 75mm in length with 16mm thread in length
3012-16-080	1	Cancellous screw Ø6.5mm x 80mm in length with 16mm thread in length
3012-16-085	1	Cancellous screw Ø6.5mm x 85mm in length with 16mm thread in length
3012-16-090	1	Cancellous screw Ø6.5mm x 90mm in length with 16mm thread in length
3012-16-095	1	Cancellous screw Ø6.5mm x 95mm in length with 16mm thread in length
3012-16-100	1	Cancellous screw Ø6.5mm x 100mm in length with 16mm thread in length
3012-16-105	1	Cancellous screw Ø6.5mm x 105mm in length with 16mm thread in length
3012-16-110	1	Cancellous screw Ø6.5mm x 110mm in length with 16mm thread in length

3013-22-030	1	Cancellous screw Ø6.5mm x 30mm in length with 22.5mm thread in length
3013-27-035	1	Cancellous screw Ø6.5mm x 35mm in length with 27.5mm thread in length
3013-32-040	1	Cancellous screw Ø6.5mm x 40mm in length with 32mm thread in length
3013-32-045	1	Cancellous screw Ø6.5mm x 45mm in length with 32mm thread in length
3013-32-050	1	Cancellous screw Ø6.5mm x 50mm in length with 32mm thread in length
3013-32-055	1	Cancellous screw Ø6.5mm x 55mm in length with 32mm thread in length
3013-32-060	1	Cancellous screw Ø6.5mm x 60mm in length with 32mm thread in length
3013-32-065	1	Cancellous screw Ø6.5mm x 65mm in length with 32mm thread in length
3013-32-070	1	Cancellous screw Ø6.5mm x 70mm in length with 32mm thread in length
3013-32-075	1	Cancellous screw Ø6.5mm x 75mm in length with 32mm thread in length
3013-32-080	1	Cancellous screw Ø6.5mm x 80mm in length with 32mm thread in length
3013-32-085	1	Cancellous screw Ø6.5mm x 85mm in length with 32mm thread in length
3013-32-090	1	Cancellous screw Ø6.5mm x 90mm in length with 32mm thread in length
3013-32-095	1	Cancellous screw Ø6.5mm x 95mm in length with 32mm thread in length
3013-32-100	1	Cancellous screw Ø6.5mm x 100mm in length with 32mm thread in length
3013-32-105	1	Cancellous screw Ø6.5mm x 105mm in length with 32mm thread in length
3013-32-110	1	Cancellous screw Ø6.5mm x 110mm in length with 32mm thread in length
<hr/>		
3014-65-025	1	Cancellous screw Ø6.5mm x 25mm length with full thread
3014-65-030	1	Cancellous screw Ø6.5mm x 30mm length with full thread
3014-65-035	1	Cancellous screw Ø6.5mm x 35mm length with full thread
3014-65-040	1	Cancellous screw Ø6.5mm x 40mm length with full thread
3014-65-045	1	Cancellous screw Ø6.5mm x 45mm length with full thread
3014-65-050	1	Cancellous screw Ø6.5mm x 50mm length with full thread
3014-65-055	1	Cancellous screw Ø6.5mm x 55mm length with full thread
3014-65-060	1	Cancellous screw Ø6.5mm x 60mm length with full thread
3014-65-065	1	Cancellous screw Ø6.5mm x 65mm length with full thread
3014-65-070	1	Cancellous screw Ø6.5mm x 70mm length with full thread
3014-65-075	1	Cancellous screw Ø6.5mm x 75mm length with full thread
3014-65-080	1	Cancellous screw Ø6.5mm x 80mm length with full thread
<hr/>		
3010-20-000	1	Maleolar screw Ø4.5mm x 20mm length
3010-25-000	1	Maleolar screw Ø4.5mm x 25mm length
3010-30-000	1	Maleolar screw Ø4.5mm x 30mm length
3010-35-000	1	Maleolar screw Ø4.5mm x 35mm length
3010-40-000	1	Maleolar screw Ø4.5mm x 40mm length
3010-45-000	1	Maleolar screw Ø4.5mm x 45mm length
3010-50-000	1	Maleolar screw Ø4.5mm x 50mm length
3010-55-000	1	Maleolar screw Ø4.5mm x 55mm length

3010-60-000	1	Maleolar screw Ø4.5mm x 60mm length
3010-65-000	1	Maleolar screw Ø4.5mm x 65mm length
3010-70-000	1	Maleolar screw Ø4.5mm x 70mm length
<hr/>		
3007-45-014	1	Cortical screw Ø4.5mm with self-tapping tip x 14 mm length
3007-45-016	1	Cortical screw Ø4.5mm with self-tapping tip x 16 mm length
3007-45-018	1	Cortical screw Ø4.5mm with self-tapping tip x 18 mm length
3007-45-020	1	Cortical screw Ø4.5mm with self-tapping tip x 20 mm length
3007-45-022	1	Cortical screw Ø4.5mm with self-tapping tip x 22 mm length
3007-45-024	1	Cortical screw Ø4.5mm with self-tapping tip x 24 mm length
3007-45-026	1	Cortical screw Ø4.5mm with self-tapping tip x 26 mm length
3007-45-028	1	Cortical screw Ø4.5mm with self-tapping tip x 28 mm length
3007-45-030	1	Cortical screw Ø4.5mm with self-tapping tip x 30 mm length
3007-45-032	1	Cortical screw Ø4.5mm with self-tapping tip x 32 mm length
3007-45-034	1	Cortical screw Ø4.5mm with self-tapping tip x 34 mm length
3007-45-036	1	Cortical screw Ø4.5mm with self-tapping tip x 36 mm length
3007-45-038	1	Cortical screw Ø4.5mm with self-tapping tip x 38 mm length
3007-45-040	1	Cortical screw Ø4.5mm with self-tapping tip x 40 mm length
3007-45-042	1	Cortical screw Ø4.5mm with self-tapping tip x 42 mm length
3007-45-044	1	Cortical screw Ø4.5mm with self-tapping tip x 44 mm length
3007-45-046	1	Cortical screw Ø4.5mm with self-tapping tip x 46 mm length
3007-45-048	1	Cortical screw Ø4.5mm with self-tapping tip x 48 mm length
3007-45-050	1	Cortical screw Ø4.5mm with self-tapping tip x 50 mm length
3007-45-052	1	Cortical screw Ø4.5mm with self-tapping tip x 52 mm length
3007-45-054	1	Cortical screw Ø4.5mm with self-tapping tip x 54 mm length
3007-45-056	1	Cortical screw Ø4.5mm with self-tapping tip x 56 mm length
3007-45-058	1	Cortical screw Ø4.5mm with self-tapping tip x 58 mm length
3007-45-060	1	Cortical screw Ø4.5mm with self-tapping tip x 60 mm length
3007-45-062	1	Cortical screw Ø4.5mm with self-tapping tip x 62 mm length
3007-45-064	1	Cortical screw Ø4.5mm with self-tapping tip x 64 mm length
3007-45-066	1	Cortical screw Ø4.5mm with self-tapping tip x 66 mm length
3007-45-068	1	Cortical screw Ø4.5mm with self-tapping tip x 68 mm length
3007-45-070	1	Cortical screw Ø4.5mm with self-tapping tip x 70 mm length
<hr/>		
3000-45-014	1	Cortical screw Ø4.5mm x 14mm length
3000-45-016	1	Cortical screw Ø4.5mm x 16mm length
3000-45-018	1	Cortical screw Ø4.5mm x 18mm length
3000-45-020	1	Cortical screw Ø4.5mm x 20mm length
3000-45-022	1	Cortical screw Ø4.5mm x 22mm length
3000-45-024	1	Cortical screw Ø4.5mm x 24mm length
3000-45-026	1	Cortical screw Ø4.5mm x 26mm length
3000-45-028	1	Cortical screw Ø4.5mm x 28mm length
3000-45-030	1	Cortical screw Ø4.5mm x 30mm length
3000-45-032	1	Cortical screw Ø4.5mm x 32mm length
3000-45-034	1	Cortical screw Ø4.5mm x 34mm length
3000-45-036	1	Cortical screw Ø4.5mm x 36mm length
3000-45-038	1	Cortical screw Ø4.5mm x 38mm length
3000-45-040	1	Cortical screw Ø4.5mm x 40mm length
3000-45-042	1	Cortical screw Ø4.5mm x 42mm length

3000-45-044	1	Cortical screw Ø4.5mm x 44mm length
3000-45-046	1	Cortical screw Ø4.5mm x 46mm length
3000-45-048	1	Cortical screw Ø4.5mm x 48mm length
3000-45-050	1	Cortical screw Ø4.5mm x 50mm length
3000-45-052	1	Cortical screw Ø4.5mm x 52mm length
3000-45-054	1	Cortical screw Ø4.5mm x 54mm length
3000-45-056	1	Cortical screw Ø4.5mm x 56mm length
3000-45-058	1	Cortical screw Ø4.5mm x 58mm length
3000-45-060	1	Cortical screw Ø4.5mm x 60mm length
3000-45-062	1	Cortical screw Ø4.5mm x 62mm length
3000-45-064	1	Cortical screw Ø4.5mm x 64mm length
3000-45-066	1	Cortical screw Ø4.5mm x 66mm length
3000-45-068	1	Cortical screw Ø4.5mm x 68mm length
3000-45-070	1	Cortical screw Ø4.5mm x 70mm length
<hr/>		
3008-35-010	1	Cortical screw with self-tapping tip Ø3.5mm x 10 mm length
3008-35-012	1	Cortical screw with self-tapping tip Ø3.5mm x 12 mm length
3008-35-014	1	Cortical screw with self-tapping tip Ø3.5mm x 14 mm length
3008-35-016	1	Cortical screw with self-tapping tip Ø3.5mm x 16 mm length
3008-35-018	1	Cortical screw with self-tapping tip Ø3.5mm x 18 mm length
3008-35-020	1	Cortical screw with self-tapping tip Ø3.5mm x 20 mm length
3008-35-022	1	Cortical screw with self-tapping tip Ø3.5mm x 22 mm length
3008-35-024	1	Cortical screw with self-tapping tip Ø3.5mm x 24 mm length
3008-35-026	1	Cortical screw with self-tapping tip Ø3.5mm x 26 mm length
3008-35-028	1	Cortical screw with self-tapping tip Ø3.5mm x 28 mm length
3008-35-030	1	Cortical screw with self-tapping tip Ø3.5mm x 30 mm length
3008-35-032	1	Cortical screw with self-tapping tip Ø3.5mm x 32 mm length
3008-35-034	1	Cortical screw with self-tapping tip Ø3.5mm x 34 mm length
3008-35-036	1	Cortical screw with self-tapping tip Ø3.5mm x 36 mm length
3008-35-038	1	Cortical screw with self-tapping tip Ø3.5mm x 38 mm length
3008-35-040	1	Cortical screw with self-tapping tip Ø3.5mm x 40 mm length
3008-35-042	1	Cortical screw with self-tapping tip Ø3.5mm x 42 mm length
3008-35-044	1	Cortical screw with self-tapping tip Ø3.5mm x 44 mm length
3008-35-046	1	Cortical screw with self-tapping tip Ø3.5mm x 46 mm length
3008-35-048	1	Cortical screw with self-tapping tip Ø3.5mm x 48 mm length
3008-35-050	1	Cortical screw with self-tapping tip Ø3.5mm x 50 mm length
<hr/>		
3006-35-010	1	Cortical screw Ø3.5mm x 10mm length
3006-35-012	1	Cortical screw Ø3.5mm x 12mm length
3006-35-014	1	Cortical screw Ø3.5mm x 14mm length
3006-35-016	1	Cortical screw Ø3.5mm x 16mm length
3006-35-018	1	Cortical screw Ø3.5mm x 18mm length
3006-35-020	1	Cortical screw Ø3.5mm x 20mm length
3006-35-022	1	Cortical screw Ø3.5mm x 22mm length
3006-35-024	1	Cortical screw Ø3.5mm x 24mm length
3006-35-026	1	Cortical screw Ø3.5mm x 26mm length
3006-35-028	1	Cortical screw Ø3.5mm x 28mm length
3006-35-030	1	Cortical screw Ø3.5mm x 30mm length
3006-35-032	1	Cortical screw Ø3.5mm x 32mm length

3006-35-034	1	Cortical screw Ø3.5mm x 34mm length
3006-35-036	1	Cortical screw Ø3.5mm x 36mm length
3006-35-038	1	Cortical screw Ø3.5mm x 38mm length
3006-35-040	1	Cortical screw Ø3.5mm x 40mm length
3006-35-042	1	Cortical screw Ø3.5mm x 42mm length
3006-35-044	1	Cortical screw Ø3.5mm x 44mm length
3006-35-046	1	Cortical screw Ø3.5mm x 46mm length
3006-35-048	1	Cortical screw Ø3.5mm x 48mm length
3006-35-050	1	Cortical screw Ø3.5mm x 50mm length
3015-04-010	1	Cancellous screw Ø 4.0mm x 10mm with 5mm length thread
3015-04-012	1	Cancellous screw Ø 4.0mm x 12 mm with 6mm length thread
3015-04-014	1	Cancellous screw Ø 4.0mm x 14 mm with 7mm length thread
3015-04-016	1	Cancellous screw Ø 4.0mm x 16 mm with 7mm length thread
3015-04-018	1	Cancellous screw Ø 4.0mm x 18mm with 8mm length thread
3015-04-020	1	Cancellous screw Ø 4.0mm x 20mm with 9mm length thread
3015-04-022	1	Cancellous screw Ø 4.0mm x 22 mm with 10mm length thread
3015-04-024	1	Cancellous screw Ø 4.0mm x 24mm with 11mm length thread
3015-04-026	1	Cancellous screw Ø 4.0mm x 26mm with 12mm length thread
3015-04-028	1	Cancellous screw Ø 4.0mm x 28mm with 13mm length thread
3015-04-030	1	Cancellous screw Ø 4.0mm x 30mm with 14mm length thread
3015-04-035	1	Cancellous screw Ø 4.0mm x 35mm with 15mm length thread
3015-04-040	1	Cancellous screw Ø 4.0mm x 40mm with 15mm length thread
3015-04-045	1	Cancellous screw Ø 4.0mm x 45 mm with 16mm length thread
3015-04-050	1	Cancellous screw Ø 4.0mm x 50 mm with 16mm length thread
3015-04-055	1	Cancellous screw Ø 4.0mm x 55 mm with 16mm length thread
3015-04-060	1	Cancellous screw Ø 4.0mm x 60 mm with 16mm length thread

### Attention

For the use of large screws and small fragments, not cannulated - BM, it is necessary to use the instruments specified below:

The Instrumental Kit for Placing Large Screws and Small Fragments, Non-cannulated - BM is registered with ANVISA (not object of this registration and not part of this product).

**Table 5 - List of Instruments that make up the Instrument tray for Placing Plates and Screws for Bone Synthesis 4.5**

Tray of Instruments for Placement of Plates and Screws for Bone Synthesis 4.5		
Qty	Code	PRODUCT

01	6019-0G-000	Hohmann Retractor, large, 2 each
01	6011-00-000	Point Drill Guide
03	6170-32-160 6170-35-160 6170-45-160	Drill Bit Ø3.2 x 160mm Drill Bit Ø3.5 x 160mm Drill Bit Ø4.5 x 160mm
01	6157-00-000	T- Handle with quick coupling
01	6013-00-000	Double Drill Guide 4.5mm N/EX/STD
01	6150-00-000	T-Handle Hexagonal Screwdriver
01	6015-00-000	4.5mm Hexagonal Screwdriver
02	6183-00-000	Plate/bone Spacer
01	6182-45-000	Axis Positioner
01	6157-05-000	Quick Coupling - with broach connection
01	6152-00-000	Countersink
03	6013	Drill Guide 4.5mm - EX / M / STD
02	6181	Drill Guide 3,2 / 3,5
01	6014-00-000	Tap Protection Sleeve, 4.5
02	6180-45-000	Threaded Drill Guide, 4.5mm
01	6153-00-000	Cortical Tapper Ø4,5
01	6156-00-000	Cancellous Tapper Ø6,0
01	6151-00-000	Depth Gauge (Large Fragments)
02	6043-00-000	Self Centering Tweezer 4.5
02	6042-00-000	Reduction Tweezer 4.5
02	6006-00-000	Screw Holding Forceps
02	6017-45-000	Bending Plates

**Table 6 - List of Instruments that make up the Instrument Tray for Placing Plates and Screws for Bone Synthesis 3.5**

Tray of Instruments for Placement of Plates and Screws for Bone Synthesis 3.5		
QTD	REF	PRODUCT
01	6019-0M-000	Hohmann Retractor, medium, 2 each
01	6019-0P-000	Hohmann Retractor, small, 2 each
01	6170-25-120	Drill Bit - Ø2,5mm
01	6170-28-120	Drill Bit - Ø2,8mm
01	6170-35-110	Drill Bit - Ø3,5mm
01	6157-00-000	T- Handle with quick coupling
01	6175-00-000	T-Handle Hexagonal Screwdriver 3.5mm
01	8565-00-000	Screwdriver with Holding Sleeve
02	6183-00-000	Plate/bone Spacer
01	6182-35-000	Axis Positioner
01	6157-05-000	Quick Coupling - with broach connection
01	6174-00-000	Countersink, quick-coupling
01	6030-00-000	Double Drill Guide 3.5mm N/EX/STD
03	6030	Drill Sleeve 2.5mm - STD / M / EX
04	6181-25-053	Drill Sleeve Ø2,5 / Ø2,8
01	6172-00-000	Tap for cortical - Ø3.5mm
01	6173-00-000	Tap for Cancellous - Ø4.0mm
01	6179-00-000	Cortical Meter for SPS Plate
02	6041-00-000	Self Centering

02	6040-00-000	Reduction Forceps
02	6177-00-000	Reduction Forceps 130mm
02	6178-00-000	Reduction Forceps 140mm
01	8564-00-000	Plate Bending Plier
02	6017-05-000	Bending Plates
01	6039-00-000	Drill Sleeve 2.5
01	6029-00-000	Drill Sleeve 3.5

### **Care with Medical Product Handling and Transport**

- Do not use the product if it is damaged.
- The implants should be handled with care to avoid damage that could impair the quality of the material and patient safety.
- Store in a cool, dry place, away from light and away from the action of weather.
- When transporting it, crashes and inadequate stacks must be avoided.

### **Indication, purpose or use for which the medical product is intended, as indicated by the manufacturer.**

#### **Instructions for Use**

Warning: This implant is supplied non-sterile. It must be sterilized according to the instructions recommended in the instruction for use. Handle BM Non-cannulated Screws only in appropriate environments and with care. Only trained professionals should handle and implant these screws. Surgical techniques vary according to the surgeon's choice. The surgeon is responsible for choosing the method, type and size of the products to be used, including the criteria for evaluating the results of surgery. BM Non-Cannulated Screws for Small and Large Fragments should be applied and adapted in accordance with the requirements and appropriate surgical techniques. The mechanical contribution of the implant should be examined by observing the exact adjustment and compression, since a stable reconstruction of the fractured bone minimizes the load to be supported by the implant.

The torque to be applied during insertion will depend on the state and characteristic of the bone and the doctor must decide how much torque will be applied. Based on tests performed, the doctor should consider the levels of deployment. Weight, degree of activity and patient conditions can impact the performance of the implant. It should be considered the fact that no implant is as strong as natural bone and thus it has limitations in biomechanical applications.

BM Non-Cannulated Screws for Small and Large Fragments are manufactured to assist the orthopedic surgeon in internal fixation procedures. Do not expect these implants the same performance of a healthy bone or to support intense and continuous mechanical stress. The performance of the implant is related to factors such as body weight, physical activity, obedience to the instructions given by the doctor after surgery, especially in consolidation absence or delayed bone consolidation. Bone consolidation is entirely related to the mechanical contribution and biological reaction of the implant. The biological reactions will be determined, in most cases, by the fracture stability. The implant mechanical contribution should be examined by observing

the exact adjustment and compression, because a stable reconstruction of the fractured bone reduces the load to be supported by the implant. Another important factor is related to the implant placement. When the fixation of bone fragments is stable, it makes the plates and screws to meet their indication of bone fixation, since it reduces the load placed on the implants. The quality of screw fixation and torque will depend on bone density, the contact area of the bone screw, and bone quality: such as changes in the bone diameter, cortical thickness and characteristics of the trabecular sites during insertion and tightening of the screw.

### **Indication for Use**

The screws are used to fix plates or similar devices onto the bone, or like compression screws, to pull the pieces of the bone together. BM Non-Cannulated Screws for Small and Large Fragments are fixed to the bone in order to provide fixation, stabilization and correction of anatomic bone form. It can be used separately or in conjunction with nonrigid fixation plates. These screws are not reusable.

#### **6.5mm Cancellous Screw**

It has a head and a threaded body. The screw head is cylindrical with a hexagon underside to allow the introduction into the bone with the use of the wrench. It is partially threaded, with 16mm or 32mm threads. The 30mm and 35mm screws have different size threads.

The 30mm screw has a 22,5mm thread and the 35mm screw has a 27.5mm thread. The thread is specific to cancellous bone. The threads for cancellous bone are characterized by a relatively thin core and a wide and deep thread. This design provides a larger surface area for the spongy trabecular bone (characterized by the metaphyseal and epiphyseal areas of bone) to grip. Cancellous screws are indicated for use in diaphyseal bone fractures and for fixing plates.

#### **6.5mm Cancellous screw (Fully thread)**

It has a head and threaded body. The screw head is cylindrical with a hexagon underside to allow the introduction into the bone with the use of the wrench. It is fully threaded. The thread is specific to cancellous bone. The threads for cancellous bone are characterized by a relatively thin core and a wide and deep thread. This design provides a larger surface area for the spongy trabecular bone (characterized by the metaphyseal and epiphyseal areas of bone) to grip. Cancellous screws are indicated for use in diaphyseal bone fractures and for fixing plates.

#### **4.5mm Malleolar Screw**

It has a head and threaded body. The screw head is cylindrical with a hexagon underside to allow the introduction into the bone with the use of the wrench. Malleolar screws are indicated for the fixation of the medial malleolus. They have the same thread pitch of cortical screws and a trephine tip that allows them to cut their own path in cancellous bone. Because of their design, they are suitable for use in the fixation of bone fragments.

#### **4.5mm Cortical Screw**

It has a head and a fully threaded body, specifically for cortical bones. The screw head is cylindrical with a hexagon underside to allow the introduction into the bone with the use of the wrench. It can be used alone or in conjunction with Biomechanica Osteosynthesis Plates. Cortical screws are indicated for epiphyseal, metaphyseal and shaft fractures in order to provide stabilization and correction of the limb or anatomical shape of the bone (integrating the set: screw, plate and bone).

#### **3.5mm Cortical Screw**

It has a head and fully threaded body, specifically for cortical bones. The screw head is cylindrical with a hexagon underside to allow the introduction into the bone with the use of the wrench. It can be used alone or in conjunction with Biomechanica Osteosynthesis Plates. Cortical screws are indicated for epiphyseal, metaphyseal and shaft fractures in order to provide stabilization and correction of the limb or anatomical shape of the bone (integrating the set: screw, plate and bone).

#### **4.0mm Cancellous Screw**

It has a head and threaded body. The screw head is cylindrical with a hexagon underside to allow the introduction into the bone with the use of the wrench. It is partially threaded and has variety of thread lengths. The thread is specific to cancellous bone. The threads for cancellous bone are characterized by a relatively thin core and a wide and deep thread. This design provides a larger surface area for the spongy trabecular bone (characterized by the metaphyseal and epiphyseal areas of bone) to grip. Cancellous screws are indicated for use in diaphyseal bone fractures and for fixing plates.

#### **Disclaimer**

Biomechanica, as a manufacturer of this device, does not practice medicine or recommend this or any other surgical technique for use in a particular patient. The surgeon performing any implant procedure is responsible for defining and using the proper techniques for deploying the medical device in each patient. Biomechanica is not responsible for the selection of the appropriate surgical technique to be used for each patient.

#### **Ancillary components used in the implantation of BM Non-Cannulated Screws for Small and Large Fragments**

Non-Cannulated Screws for small and large fragments are used alone or together with non-rigid fixation plates. Below are the ancillary components that can be used together with the screws, they are not object of this record, and must be purchased separately, however they are associated with the implantation process, without which the product cannot be implanted. Screws not mentioned are used separately.

#### **Specific Ancillary Components for Ø6.5mm Screw Cancellous**

95° Autodynamic Plate for non-rigid fixation, made of austenitic stainless steel, according to ASTM F138 (not object of this registration and not part of this product).

135° Autodynamic Plate for non-rigid fixation, made of austenitic stainless steel, according to ASTM F138 (not object of this registration and not part of this product).

95° Autodynamic Plate for nonrigid fixation, made of austenitic stainless steel, according to ASTM F138 (not object of this registration and not part of this product).

#### **Specific Ancillary Components for Ø6.5mm Screw Cancellous (Fully Threaded)**

T Plate for non-rigid fixation, made of austenitic stainless steel according to ASTM F139 (not object of this registration and not part of this product).

#### **Specific Ancillary Components for 4.5mm Cortical Screw**

4.5mm Straight Plate - Broad for non-rigid fixation, made of austenitic stainless steel according to ASTM F138 (not object of this registration and not part of this product).

4.5mm Straight Plate - Narrow for non-rigid fixation, made of austenitic stainless steel according to ASTM F138 (not object of this registration and not part of this product).

L Plate - Right / Left made of austenitic stainless steel according to ASTM F139 (not object of this registration and not part of this product).

95° Autodynamic Plate for non-rigid fixation, made of austenitic stainless steel according to ASTM F138 (not object of this registration and not part of this product).

135° Autodynamic Plate for non-rigid fixation, made of austenitic stainless steel according to ASTM F138 (not object of this registration and not part of this product).

### Specific Ancillary Components for 3.5mm Cortical Screw

3.5mm Straight Plate (Small fragments) for non-rigid fixation, made of austenitic stainless steel according to ASTM F138 (not object of this registration and not part of this product).

3.5mm Oblique Plate - Left (not object of this registration and not part of this product) / Right for non-rigid fixation, made of austenitic stainless steel according to ASTM F139

1/3 Tubular Plate for non-rigid fixation, made of austenitic stainless steel according to ASTM F139 (not object of this registration and not part of this product).

3.5mm Reconstruction Plate for non-rigid fixation, made of austenitic stainless steel according to ASTM F138 (not object of this registration and not part of this product).

3.5mm T Plate for non-rigid fixation, made of austenitic stainless steel according to ASTM F139 (not object of this registration and not part of this product).

### Accessories

BM Non-Cannulated Screws for Small and Large Fragments have the following accessories:

#### Washer

Code	Diameter	Used with	Photo
3005-07	Ø7.0mm (Stainless Steel F-138)	Cortical Screw Ø3.5mm and Cancellous Screw Ø4.0mm	
3005-13	Ø13.0mm (Stainless Steel F-138)	Cortical Screw Ø4.5 and Cancellous Screw Ø6.5mm	

#### Allowable combinations with other materials

The table below "Allowable Contact Leagues" was based on the standard "NBRISO21534 – Non-active Surgical Implants – Joint Replacement Implants – Specific requirements," but despite the mentioned alloys in the item "In Contact" the combinations must be made with implants of the same brand because they were designed for such combinations observing surface finishing and surface treatment and other factors required in their projects that can interfere with the combinations. Therefore, metal implants from different manufacturers are not recommended because of chemical, physical, biological and functional incompatibility.

Tabela 7 – Allowable contact alloys

Allowable contact alloys	
Implant	* In contact (when there is no articulation)
Conformed Stainless steel Screws Wrought 18 Chromium – 14 Nickel - 2,5 Molybdenum	Stainless steel - NBRISO 5832-1 and ASTM F138) Wrought High nitrogen stainless steel (NBRISO5832-9) Wrought Titanium 6-aluminum 4-vanadium Alloy (NBRISO 5832-3 and ASTM F136) Wrought Titanium 6-aluminum 7-niobium alloy (NBRISO5832-11)

\* For applications when a metal or alloy is in contact with others and articulation is not intended, paying adequate attention to the design, surface finishing, surface treatment and metallurgical conditions.

BM Non-Cannulated Screws for Small and Large Fragments must be only used with instruments manufactured by Biomechanica.

### **Contraindications**

It is contraindicated the use of this implant in the following situations: -sensitivity to metal;-fever or signs of local inflammation; -osteoporosis -obesity; -inadequate tissue coverage in the region to be operated; -pregnancy; -bone endangerment caused by diseases and / or infections; - unexplained high sedimentation rate through other diseases, increase of the WBC rate or suspected striking changes in the count; - visible infection.

**Precautions, restrictions, warnings, special care and clarification on the use of medical products, as well as their storage and transport.**

### **Warnings and Precautions**

**SINGLE-USE PRODUCT. DESTROY AFTER REMOVAL. DO NOT REUSE THE PRODUCT.  
PROHIBITED TO REPROCESS IT.**

An orthopedic implant can be used in a single patient only once. Although it may appear undamaged, previous stress may create imperfections that may reduce the success of the implant. Improper selection of implant can cause unusual stress on it and may result in subsequent fracture of the implant.

Due to its form and characteristics, this product does not have side effect, however, it is suggested to monitor the patient's recovery through X-ray according to the International Security Standards. Do not use magnetic resonance imaging, because the product manufactured in stainless steel causes interference in the examination harming the diagnosis, and may also cause displacement of the implant causing trouble for the patient.

The use of BM Non-Cannulated Screws for Small and Large Fragments are designed to be an aid to the fixation and stabilization in the treatment of bone fractures in the lower and upper limbs, not to replace the normal structures of the human skeleton. The surgeon should be aware of the physical development of the patient, a revision of the implant may be needed. The surgeon should also be familiar with and hold enough knowledge of osteosynthesis and its limitations, including pre and post operative surgical technique used, precautions and potential risks. The surgeon should be attentive during the bone preparation for the introduction of the BM Non-Cannulated Screws for Small and Large Fragments, because the stability and success of fixation may be compromised;

During handling of the implant scratches, grooves, or anything that may damage or mark the implant must be avoided, because these defects are stress concentrators and may be sites of crack nucleation and decrease the corrosion resistance and may result in fracture of the implant.

To ensure proper implantation, Biomechanica instruments must be used, because they are designed and manufactured specifically for the placement of these implants. Using instruments from different manufacturers may compromise the surgery as well as not respect product registration in the National Health Surveillance.

**Do not use the product if it is damaged.**

**NON STERILE PRODUCT - STERILIZE PRIOR TO USE.**

**PRODUCT OF SINGE USE - IT IS STRICTLY PROHIBITED TO RE-USE THIS PRODUCT.**

**Special storage conditions:** Store it in a cool, dry place, protected from light and weather.

Do not use the product if the packaging is damaged.

**Note:** do not use implant components from different manufacturers. We, therefore, recommend that the products have the same provenance.

Manufacture Date, expiration date and product batch number: SEE LABEL.

#### **Product Conformation Limit**

While handling the implant, scratches, notches, or anything that could damage or mark the implant should be avoided, as these defects are stress concentrators and can be sites of crack nucleation and decrease corrosion resistance, which may result in fracture of the implant.

#### **Adverse Effects**

Mechanic loosening may be the result of a defective fixation, unstable reconstruction of the fractured bone or hidden infection. Metal sensitivity reactions in patients were rarely informed. Implantation of a strange material in tissues results in histologic reactions involving various sizes of macrophages and fibroblasts. The importance of this chemical effect is uncertain, but similar changes can occur as a precursor, or during the healing process.

#### **Information to be provided to the patient**

The following information must be provided to the patient: Indications, Contraindications, Information for use, Possible Adverse Effects, Precautions and Warnings. The patient must be informed about the importance of the postoperative follow-up. The lack of monitoring hampers the detection of problems after surgery. The failure to perform revision surgery when there is a component release or osteolysis may result in progressive loss of periprosthetic bone stock. The patient should be informed about the type of implant used, brand, code, name, batch number, surgeon's name and surgery date.

#### **Importante**

The surgeon must orient the patient about the care to be taken and appropriate physical support to be used during the time the implant is being used. When the fracture is perfectly stable, load can be applied. If the fracture is not stabilized, no load can be applied to the implanted system.

### **Implanted Product Evaluations**

After implantation, in the intraoperative, the professional in charge must perform a radiological control to verify the correct positioning of the product. It is this professional's responsibility to make the clinical and radiological evaluation after the surgical procedure in the frequency stipulated by him/her to check the implant condition and the bone development. If the product is in the wrong position or presents any non-compliance, it is the surgeon's responsibility to take the appropriate corrective action.

### **Useful information to avoid risks derived from the implantation**

To lower the implantation risks the information contained in the contraindications, instructions for use and any information contained in "Product Instructions for Use" must be followed strictly.

### **Implantation Intrinsic Risk**

BM Non-Cannulated Screws for Small and Large Fragments are made of biomedical use materials, following the requirements of ASTM F138 -00 - Standard Specification for Wrought 18 Chromium-14 Nickel-2, 5 Molybdenum Stainless Steel Bar and Wire for Surgical Implants (UNS S31673).

### **Contamination Risks**

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the metallic components of BM Non-Cannulated Screws for Small and Large Fragments come in contact with body tissues and fluids. Implanted products should be treated as highly polluting.

### **Decontamination, Cleaning and Sterilization**

To ensure the reduction of infection risk, the implants (which were not used in surgery, but were placed on the operating table or those that are not sterilized anymore) and surgical instruments must be cleaned after surgery and before sterilization.

Cleaning can be developed through manual or mechanical methods.

Cleaning Equipment:

- Ultrasonic Washer
- Sterilizing Washer
- Washer – disinfectant

Manual Cleaning: Is the procedure which the dirt is removed by means of physical action with the aid of detergent, water and artifacts such as sponge and brush. Manual cleaning is recommended for causing less harm to implants and surgical instruments. In manual cleaning water should be used at room temperature, using nylon brushes, never steel or abrasive brushes, because they can damage the protective layer of the material. In the decontamination process, antibacterial and

broad-spectrum antifungal solution should be used; aggressive cleaning agents, brushes, metal must be avoided not to damage the metallic implants.

In case of cleaning machines, implants and surgical instruments should be organized in drawers not to hit each other and avoid material damage. People responsible for cleaning the implants and surgical instruments should be aware of the type of cleaning product used and respect the concentration recommended by the manufacturer. Implants and surgical instruments must be thoroughly rinsed until all residues are removed and dried immediately after cleaning.

#### Types of cleaners

Enzyme Cleaners: consist mainly of enzymes, surfactants and solubilizers. A balanced combination of these elements means that the product can remove organic material in a short period of time.

Enzyme solutions: have excellent cleaning action, but do not have bactericidal and bacteriostatic activity.

Enzymes: are substances produced by living cells and which cause chemical reactions of the process. Once produced by the cells, an enzyme can be isolated and will maintain its catalytic properties, if certain conditions are maintained in its manufacture. Enzymes are classified into three major functional groups depending on the type of substrate that they will affect: proteases, amylases and lipases that act on proteinaceous substrates, fats and carbohydrates, which tend to solubilize and come off the products.

#### Desinfection:

It is a process where pathogenic or not microorganisms are killed, except for bacterial spores, by physical or chemical methods.

#### Levels of disinfection:

High level: destroys all microorganisms, with the exception of high numbers of bacterial spores => Glutaraldehyde 2% - 20 - 30 minutes.

Indication: preferably in hospitals

Intermediate level: eliminate vegetative bacteria, most viruses, fungi and microbacteria => Sodium hypochlorite 1% - 30 minutes.

Indication: day care, rest homes, nursing homes.

Low-level: eliminates most bacteria, some viruses and fungi, but does not eliminate microbacteria => Sodium hypochlorite 0.025%

Indication: nutrition.

### **Sterilization**

This product is supplied non-sterile. Before its use it should be sterilized. We recommend steam sterilization in autoclave in the hospital (ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, calibration and routine control of a sterilization process for medical devices)

## AUTOCLAVE

It is a sterilizing equipment that uses saturated steam at high pressures.

### OPERATING INSTRUCTIONS

To start the equipment, check:

- whether the circuit breaker is on;
- if the water mains valve is open;
- if the discharge register is closed.

1 - open the equipment door;

2 - properly accommodate the material to be sterilized;

3 - close the equipment door;

4 - select the desired cycle according to the material to be sterilized;

5 - turn on the main switch;

6 - the cycle will run automatically, in sequence;

7 - when turning on the "end of cycle" lamp, partially open the door for ten (10) minutes, approximately, for the material to cool.

It is recommended that the following physical sterilization parameters be applied in autoclaves (saturated steam):

**Table 8 - Parameters of physical sterilization in autoclaves**

Cycle	Temperatura	Exposure Time (minimum)	Drying Time (maximum)
Conventional (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Note: The time must be set when the heat in the sterilization chamber reaches the desired temperature.

### FOR CLARIFICATION CONSULT THE AUTOCLAVE MANUAL.

**Other sterilization method:** In addition to autoclave, the following sterilization method may be used as defined by the hospital institution:

- **Sterilization by ethylene oxide (E.T.O.)** - parameters and procedures established in the validation protocol and in EN550:1994 Sterilization of Medical Devices - Validation and Routine Control of Ethylene Oxide Sterilization.

### NOTE

The method of sterilization, equipment, controls, and instructions for sterilization used is under the hospital responsibility.

Care with the Sterilized Products: Storage Conditions of sterilized products:

The environment: it must be clean, airy, dry; it should be a restricted area.

The product: After the sterilization process, do not place it on a cold surface (stone or stainless steel), use baskets or containers until they cool down; mantle (woven cotton, non-woven fabric, paper surgical grade, corrugated cardboard, paper, film, TYVEK or perforated metal boxes) must remain intact and be handled slightly to prevent tearing packages or release the seal; be stored in closed cabinets with shelves, racks identified in order to facilitate the removal of material; material should be stored according to the maturity date of sterilization to facilitate the distribution and not be overcome in rolling stock; store separately from non-sterile to reduce the level of external contaminants.

### **Product Discard**

No surgical implant should be reused. Any metallic implant, after use, must be discarded. Even if it appears to be undamaged, it may have small defects that can lead to fatigue fractures in addition to contamination. In this way, large Screws and small fragments, non-cannulated - BM explanted from patients must be properly discarded by the hospital institution. The hospital institution is responsible for the complete mischaracterization of the implant, preventing its reuse. The hospital institution is responsible for the method used to mischaracterize the implant. Biomecânica recommends that the explanted implants be mechanically deformed with the aid of a hammer or impact press, and must then be identified with the saying "Improper for Use".

As described in Resolution No. 2605, of 08/11/2006, implantable devices of any nature classified as single use are prohibited from being reprocessed.

### **TRACEABILITY**

On the packaging of the implantable component there are 2 traceability labels containing the data of the implant used. These labels are the same label contained in the external product labeling. The information contained on this traceability label is related in the "Product Presentation" of this technical report. This label must be affixed to the patient's records and another is for the hospital control.

The hospital must have this information in its files and give them to the patient. We inform that it is mandatory to have the information in the patient's records that permits the traceability of the implanted product.

The attachment of this label in the patient's records enables the traceability of the product used. Among the main information we highlight as essential the name of the manufacturer of the implant used, name of the implant used, the implant code, batch number and registration number of the product in ANVISA. This information is described on the traceability labels that come with the product and is available on the external labeling. Other information must also be considered essential, such as, surgery date, name of patient who received the implant, surgeon's name, patient's weight, patient's age and other information requested in the patient's record must also be filled.

## **MARKING**

BM Non-Cannulated Screws for Small and Large Fragments contain the following information recorded with laser to trace and monitor the patient after surgery: - Biomecanica Logo; - Manufacturing Batch Number; - Raw Material Symbol; - Size

### **Post-Sales (Customer Complaint)**

If there is a need to make a complaint about BM Non-Cannulated Screws for Small and Large Fragments related to any adverse effect that affects the security of the user, such as problems with the performance of the product, damage to the implantable metal component, serious problems or deaths related to these components, the surgeon in charge must report this adverse event to the national health authority and Biomecanica by e-mail or by phone sac@biomecanica.com.br 0xx14 2104 7900. In cases of doubt, the surgeon in charge or health professional can make the communication of the adverse event through the Public Health Notification System at ANVISA website:

<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

### **Sending Material for the Manufacturer to Analyze**

If implants are sent to the manufacturer for analysis, they must be cleaned in the hospital using a broad-spectrum bactericidal and antifungal solution. Then it must be disinfected or sterilized by steam in an autoclave or ethylene oxide. They must be sent to Biomecanica in intact packaging, identified with the method of cleaning, sterilization and product data.

### **Performance provided for in the General Requirements of ANVISA's regulations that provide for the Essential Safety, Effectiveness and Safety Requirements applicable to products, described in Resolution - RDC 56, of 06/04/2001:**

Every implant is subject to breakage, wear or loosening and its sizing is determined by a number of requests over a certain period of time. The time foreseen in traumatology surgeries is until the fracture consolidates. These requests are related to the stress caused by normal activities such as walking, lifting or climbing stairs. So that we can minimize the exposed risks, the large screws and small fragments, not cannulated - BM are manufactured with high stability materials, using reputable raw materials such as ASTM F138, in addition to NBRISO14630 and NBRISO21534 standards. Large screws and small fragments, non-cannulated - BM were tested to find their resistance limits, thus having torque and pullout values.

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