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ENGLISH EUA

Medical product detailed description, including operation and action foundations, content or composition, if applicable, and listing of accessories intended to be part of the product.

Muller Acetabulum Model

Muller Acetabulum is made of polyethylene UHMWP 1900 (NBR ISO 5834-2). There are external grooves to assist fastening with acrylic cement in acetabular bone cavity and internal cavity to support metal femoral head. It also has a Ø1.5mm metal ring made of stainless Steel ASTM F138 to enable acetabulum position viewing in X-Ray.


Acetabulum Model CP3

PMMA centralizer/spacer which function is centralizing the acetabulum in acetabular bone cavity and assist cement homogeneous distribution

Ø1.5mm Metal ring to enable X-Ray acetabulum positioning viewing. Muller and CP3 acetabulum have metal ring made of Stainless Steel ASTM F138.

CP3 Acetabulum is made of polyethylene UHMWP 1900 (NBR ISO 5834-2). In external surface, apart from grooves, there are 4 PMMA centralizers/spacers which function is improving X-Ray acetabulum positioning viewing and centralizing acetabulum in acetabular bone cavity and assist in Acrilex Orthopedic Cement homogeneous distribution, which is not object of this registration and is not integral part of this product, Biomecânica band, ANVISA registration no. 10171110070.

It also has a Ø1.5mm metal ring made of Stainless Steel ASTM F138, which function is improving X-Ray acetabulum viewing. It has an internal cavity to support Biomecânica brand metal femoral head.

Composition: Muller Acetabulum model is made of polyethylene UHMWP 1900 (NBR ISO 5834-2). Acetabulum model CP3 is made of polyethylene UHMWP 1900 (NBR ISO 5834-2) and has centralizers/spacers made of PMMA (Polymethylmethacrylate). Muller and CP3 acetabulum have a metal ring made of Stainless Steel ASTM F138, which purpose is enabling X-Ray acetabulum positioning viewing.

Presentation Form: Components are stored in blister double packaging sealed with surgical paper, plus cardboard box external packaging and ready to be used labeling, already sterilized in ethylene oxide.


Product traceability labels include the following information:

- s) Manufacturer’s Name
- t) Product Trade Name
- u) Technical Name
- v) Batch Number
- w) ANVISA Registration Number
- x) Product Code
- y) Quantity
- z) Product description that contains the packaging (model)
- aa) Manufacturing Date
- bb) Valid Date Undefined Valid Date
- cc) Description of raw material used to manufacture the product
- dd) Manufacturer’s Address
- ee) Technical Responsible’s Name
- ff) The statements: **‘Sterile Product’**.
- gg) The statements: **‘REPROCESSING FORBIDDEN’**.
- hh) ‘Single use product, do not reuse it’
- ii) Before using, read label instructions
- jj) Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use

The product packaging provides Safety Symbols printed in cardboard unit boxes, according to the following model:

Table 1 - Symbols, as per NBR ISO 15223

Symbol	Description
	Instructions, warnings and precautions: see Instructions for Use
	E.T.O. Gas Sterilized

	Single-use hospital medical good
	Do not use if the packaging is damaged
	Fragile, handle with care
	Keep away from sun light
	Keep it dry

Polymer Cemented One-Component Acetabulum o BM will be distributed in unit form according to the following table:

Table 2 - Polymer Cemented One-Component Acetabulum - BM

Code	Description
2421-40-000	External Ø 40mm and internal Ø 28mm Muller Acetabulum
2421-42-000	External Ø 42 mm and internal Ø 28mm Muller Acetabulum
2421-44-000	External Ø 44 mm and internal Ø 28mm Muller Acetabulum
2421-46-000	External Ø 46 mm and internal Ø 28mm Muller Acetabulum
2421-48-000	External Ø 48 mm and internal Ø 28mm Muller Acetabulum
2421-50-000	External Ø 50 mm and internal Ø 28mm Muller Acetabulum
2421-52-000	External Ø 52 mm and internal Ø 28mm Muller Acetabulum
2421-54-000	External Ø 54 mm and internal Ø 28mm Muller Acetabulum
2421-56-000	External Ø 56 mm and internal Ø 28mm Muller Acetabulum
2421-58-000	External Ø 58 mm and internal Ø 28mm Muller Acetabulum
2424-28-040	Ø 28 x 40 mm – CP3 Acetabulum
2424-28-042	Ø 28 x 42 mm - CP3 Acetabulum
2424-28-044	Ø 28 x 44 mm - CP3 Acetabulum
2424-28-046	Ø 28 x 46 mm - CP3 Acetabulum
2424-28-048	Ø 28 x 48 mm - CP3 Acetabulum
2424-28-050	Ø 28 x 50 mm - CP3 Acetabulum
2424-28-052	Ø 28 x 52 mm - CP3 Acetabulum
2424-28-054	Ø 28 x 54 mm - CP3 Acetabulum
2424-28-056	Ø 28 x 56 mm - CP3 Acetabulum
2424-28-058	Ø 28 x 58 mm - CP3 Acetabulum
2424-28-060	Ø 28 x 60 mm - CP3 Acetabulum

Attention

Listing of instruments to assist Polymer Cemented One-Component Acetabulum - BM implantation

To use Polymer Cemented One-Component Acetabulum - BM, it is required to use instruments 'Instrument Kit to Place Prosthesis for Hip Arthroplasty', it is registered in ANVISA and is not object of this registration and is not integral part of this product. Instruments and relevant trays are provided as follows:

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Medical Product Handling and Transport Special Cautions: The medical product storage place must be clean, dry and illuminated, in order to keep optimal storage and transport conditions, as well as its physical and chemical integrity. The product must be stored transported in dry and ventilated place, with room temperature (max. 35°C) and relative humidity around 30%, and away from direct sun light.

Vibration, bump, corrosion, defective settling effects during motion and transport, inadequate stacking during storage, must be prevented.

They cannot be stored directly on the floor (minimum height = 20cm). They cannot be in too high places, close to lamps, as the packaging might get dried or the label might be damaged. They have to be stored in places where contaminant substances are stored, including, for instance, cleaning materials, insecticides, pesticides, etc.

The carrier must be notified on delivery contents and term.

Transport carefully, to prevent dropping and attrition, in order not to damage packaging and device, assuring thus material sterility.

Blister sealed packaging opening can be performed manually, in adequate location (surgical environment), pulling surgical grade paper through spares around the blister. This opening must be performed only in surgical act, to prevent product contamination.

The product cannot be used if the packaging is breached.

Remark: Different manufacturer implant components cannot be used; thus we recommend that the products have the same origin.

Manufacturing date, valid date and product batch: SEE LABEL.

Sterilized Material Handling: As the sterilized material is handled with aseptic technique, some standards have to be complied with to keep it sterile:

- it is essential to wash hands with water and soap before handling the sterilized material;
- use material with intact, dry, without stains packaging, with identification (material type and sterilization type);
- work facing the material;
- handle the material at waist level up;
- avoid coughing, sneezing, talking over the exposed material;
- do not move over the sterilized area;
- check the packaging valid date and adequacy;

- work at a clean, quite, dry environment, and without draft;
- keep a certain distance between the body and the material do be handled;
- follow other asepsis principles.

The nursing technique required to sterilized material handling is:

- open it, starting by the end opposite to the handler;
- protect material exposed with sterilized fields that involves it;
- touch with the hands only in packet external part;
- do not store as sterilized material a previously open package.

Warranty: Polymer Cemented One-Component Acetabulum warranty is linked to following these Instructions for Use.

Medical Product Indication, Purpose or Intended Use

It is indicated in cemented hip arthroplasties.

A full hip arthroplastry performance must be indicated to treat the following diseases:

- Arthritis
 - Rheumatoid
 - Juvenile Rheumatoid (Still´s Disease)
 - Ankylosing Spondylitis
 - Degenerative Joint Disease (osteoarthritis, hypertrophic)
 - Primary
 - Secondary
 - Femoral head epiphysis sliding
 - Hip congenital luxation/dysplasia
 - Flat thigh (Legg-Perthers Disease)
 - Paget´s Disease
 - Traumatic luxation
 - Acetabulum fracture
 - Hemophilia
- Vascular necrosis
 - Post-fracture or idiopathic luxation
 - Femoral head epiphysis sliding
 - Haemoglobinopathies (sickle-cell disease)
- Renal disease
 - Induced by corticosteroids
 - Alcoholism
 - Decompression sickness
 - Lupus
 - Gaucher´s Disease

Pseudarthrosis, femoral neck and trochanteric fractures with compromising of the head

- Pyogenic arthrodesis or osteomyelitis
Hematogenous
Postoperative
- Tuberculosis
- Sub-luxation or congenital luxation
- Hip arthrodesis and pseudarthrosis
- Failed reconstruction
Osteotomy
Dome arthroplasty
Femoral head prosthesis
Girdlestone Procedure
Total hip arthroplasty
Resurfacing arthroplasty
- Bone tumor involving proximal femur or acetabulum
- Hereditary diseases (e.g.: achondroplasia)

Polymer Cemented One-Component Acetabulum is available in the following measures:

Muller Acetabulum Model

Code	external Ø	internal Ø (mm)
2421-40-000	40.0	28.0
2421-42-000	42.0	28.0
2421-44-000	44.0	28.0
2421-46-000	46.0	28.0
2421-48-000	48.0	28.0
2421-50-000	50.0	28.0
2421-52-000	52.0	28.0
2421-54-000	54.0	28.0
2421-56-000	56.0	28.0
2421-58-000	58.0	28.0

Acetabulum Model CP3

Code	external Ø	internal Ø (mm)
2424-28-040	40	28.0
2424-28-042	42	28.0
2424-28-044	44	28.0
2424-28-046	46	28.0
2424-28-048	48	28.0
2424-28-050	50	28.0
2424-28-052	52	28.0
2424-28-054	54	28.0
2424-28-056	56	28.0
2424-28-058	58	28.0

2424-28-060	60	28.0
-------------	----	------

Femoral Cementing

Option 1) Using ACRILEX ORTHOPEDIC CEMENT (not object of this registration): the cement must be mixed in the container (supplied with the product) for 1 minute and 15 seconds approximately, and then poured in the cement drum that must be vertical position for 30 seconds. As the indicated time is waited, the nozzle is threaded in the drum. The cement is introduced in retrograde way (from distal to proximal region).

Option 2) Digital Application: Using ACRILEX ORTHOPEDIC CEMENT (not object of this registration), mix the components with the spatula supplied with the packaging for approximately 1 minute and 15 seconds. It is very important to minimize air input while it is mixed, as air presence damages cement efficacy and can prevent its mechanical performance. Handling process should continue from 15 seconds to 1 minute while the cement is ready to be applied to the bone. The application must last 1 minute and a half, after that time, the cement becomes too viscous to be inserted in the narrow cavity. The prosthesis is inserted when viscosity starts to increase. The goal is delaying rod insertion as much as possible, as during rod insertion canal interface pressures are directly linked to cement viscosity. Orthopedic Cement hardening final time (at 20°C) is usually 7 to 8 minutes.

Orthopedic Prosthesis Recommendations: To achieve the best hip arthroplasty results, the following is recommended:

- 1 - Using orthopedic prosthesis must be performed only by surgeons qualified and trained in such procedure. It is essential that a careful preoperative planning is developed, including by using product transparency or template.
- 2 - All care must be taken to prepare the receiving bone, to achieve optimal implant fit, preventing radiolucency and unwanted micro-motion appearance.
- 3 - Surgical instruments and tests are available to assist in orthopedic prosthesis surgical implantation. It is important that surgical instruments and tests are the ones specifically designed to the specific product. Surgical instrument and similar test project and dimension variations may compromise critical measures required to accurate implantation.
- 4 - In order to protect the physician and his/her patient, in occasional future problems, the hospital may be responsible for annotation in patient protocol of used implant batch code and number. Such data is critically relevant to enable traceability.
- 5 - Before starting the surgery, check if implant collection and relevant instruments are intact and complete.
- 6 - Never mix different manufacturer implants in the same patient. There might be fit dimension and tolerance differences of such implants, causing deficient junctions, early loosening or their failure.
- 7 - Orthopedic prostheses are classified as 'Single Use Products', i.e., they cannot be reused.

8 - In every case, postoperative established surgical practices must be followed. The patient must be warned on replacement or restitution limitations of a specific joint and such practice use recent history.

9 - The surgeon must also prevent creases or scratches in orthopedic prostheses, as such damage and/or faults can produce internal stress, that may make it a focus of eventual break.

10 - A joint surface should never be removed and reinserted during surgical procedure. Despite being apparently intact, it might have undergone dimensional modifications that impact its fastening and survival.

11 - Postoperative cautions, as well as patient capacity and willingness to follow instructions, are two of the most relevant aspects for arthroplasty success. The patient must be warned that non-compliance with postoperative instructions may lead to orthopedic prosthesis break or migration, requiring new surgery for revision or removal. 12 - All effort should be made to use compatible biomedical materials (metals and alloys), while using orthopedic prostheses, as mixing different materials, in the same patient, and component micro-motions can cause metalosis.

13 - A specific joint replacement or restitution success potential is incremented by adequate selection of size, form and project of orthopedic prosthesis. While adequate selection can enable minimizing risks, size and form, as well as human bone and adjacent soft tissue quality offer certain limitations concerning such implant mechanical strength.

14 - The orthopedic prosthesis is designed to replace a specific damaged joint and not human skeleton normal structures.

15 - Orthopedic prostheses cannot withstand activity levels and/or loads equal to the ones withstood in normal and healthy bones. This detail must be reported by the surgeon and understood by the patient.

16 - An obese or very heavy patient can produce overload in implant used, causing metal fatigue and leading to component loosening, break or failure.

17 - Excessive physical activity, and traumatismas that impact joint replaced have implied in orthopedic prosthesis early failure, whether by implant loss, fracture or wear. The patient must be warned to keep his/her activities according to the situation, protecting the joint replaced against excessive stress.

18 - In general, orthopedic prostheses are supplied in sterile condition, in double blister, in an external rigid cardboard box, properly identified by adhesive labels, with all product relevant legal information, that enable their full identification and traceability. The packaging must be intact at the time of receiving (do not use the product if the packaging is tampered). Check sterilization valid date (do not use the product if the sterilization term is expired).

Disclaimer: Biomecânica, as this device manufacturer, does not practice medicine and does not recommend any specific surgical technical to be used in a specific patient. The surgeon that performs any implant procedure is responsible for setting and using appropriate techniques to implant prosthesis in each patient. Biomecânica is not responsible for appropriate surgical

technique choice to be used in each patient. Surgical Techniques depend on the technique adopted by the surgeon, and the latter is liable to choose the method, type and dimension of the products to be placed.

Ancillary Components Ancillary components below must be purchased separately, as they are not integral part of this product.

- Femoral Metal Head is registered in ANVISA (it is not object of this registration and not integral part of this product).
- Cemented Modular Rod is registered in ANVISA (it is not object of this registration and not integral part of this product).
- Acrilex Orthopedic Cement is registered in ANVISA (it is not object of this registration and not integral part of this product).

Accessories: Polymer Cemented One-Component Acetabulum does not have any accessory with the purpose of integrating the medical product.

Admissible Combinations with other Materials: The following table 'Contact Admissible Alloys' was based on standard 'NBR ISO 21534 - Non-Active Surgery Implants - Implants to Replace Joints - Particular Requirements', but regardless of alloys provided in item of table '**Joint Contact Admissible Alloys***', it is required that the combinations are of implants of the same brand, as they have been designed to such combinations observing superficial finishing and superficial treatment, and other factors required in projects that can interfere in combinations. Thus, metal implants from different manufacturers are not recommended due to chemical, physical, biological and function incompatibility.

Implant	Joint contact admissible alloys
Polymer Cemented One-Component Acetabulum	Conformed Stainless Steel (ASTM F138, NBR ISO 5832-1) High Nitrogen Conformed Stainless Steel (NBR ISO 5832-9) Cobalt-chromium-molybdenum cast alloy (ASTM F75, NBR ISO 5832-4) Cobalt-chromium-tungsten-nickel conformed alloy (NBR ISO 5832-5) Cobalt-chromium-nickel-molybdenum-tungsten-iron forged and cold conformed alloy (NBR ISO 5832-7) Cobalt-nickel-chromium-molybdenum-tungsten-iron conformed alloy (NBR ISO 5832-8) titanium 6-aluminum 4-vanadium conformed alloy (ASTM F136, NBRISO5832-3) titanium 6-aluminum 7-niobium conformed alloy (NBR ISO 5832 -11) Alumina base ceramic materials (NBR ISO 6474) Zirconia base ceramic materials (NBR ISO 13356)

*For implant joint surfaces to replace the joint, superficial finishing and surface treatment have been considered acceptable in specific application, as long as adequate project attention is provided.

Biomecânica brand Polymer Cemented One-Component Acetabulum must be used only with Biomecânica brand instruments.

Precautions, Restrictions and Warnings

Precautions: SINGLE USE PRODUCT. DESTROY AFTER EXPLANTED. DO NOT REUSE THE PRODUCT. REPROCESSING FORBIDDEN

An implant cannot be used again. Even though damages are not apparent, previous tensions may have created imperfections that could reduce implant life.

Check if the product packaging is intact. Check the sterilization date. Open the packaging aseptically, upon checking if the dimension is the one chosen. Sterilization valid date: 5 years (indicated in product internal and external packagings).

The surgeon must be aware of patient physical development possibility, and implant revision may be required.

Contraindications: According to Harkess, total hip arthroplasty is a large size procedure and is linked to a significant complication number and death rate from 1% to 2%. Thus, as total hip arthroplasty is indicated, patients must be assessed carefully, especially concerning systemic diseases. Cardiac, pulmonary, hepatic, genitourinary or metabolic diseases, hypertension and malignancy must be treated and controlled before total hip arthroplasty. Other situations are also not advised in total hip arthroplasty surgeries, including:

Upper femur osteomyelitis, hip joint pyrogenic infection or evident infections, are absolute contraindications.

All effort must be made to rule preoperative asepsys possibility in a patient that has one or more of the following abnormalities:

- signs of fever, local inflammation;
- quick joint destruction, apparent bone absorption in roentgenogram;
- not explained sedimentation rate elevation through other diseases, WBC sum elevation, or a higher change in differential sum.

According to exceptional circumstances, low grade subacute infections, followed by joint total replacement, prosthesis removal can be treated by professionals, complete debridement, and new device reimplantation. However, the infection repetition has arthroplasty revision higher probability than in a recent case. Under no circumstances, it is recommended that the insertion is performed where there was bone extensive involvement; or involved soft tissue cannot be removed; or when organism is gram-negative, and not sensitive to *Staphylococcus aureus*, or where highly virulent organisms or mixed are involved.

Additionally, in infection distant focus, including genitourinary, pulmonary, skin, or other locations, are relatively not advised, due to hematogenesis transmission to implant location that may occur. Infection focus must be treated before, during and after implantation.

This implant use is not advised where there is musculature loss, or neuromuscular infection, compromising the affected member, what would make the procedure unjustifiable.

Warnings: Improper selection of implant, insertion, placement or excessive muscular activity that can cause unusual tensions in implant and can result in subsequent prosthesis fracture.

Whenever there is insufficiency or loss of bone acetabular reserve, grafts or other reinforcement procedures are advisable to provide acetabular support.

Postoperative caution acceptable techniques must be followed.

The patient must be warned on total joint reconstruction limitations and recent use history. Excessive physical activity and trauma that impact joint replaced were implied in reconstruction early failure by looseness, fractures, in implant use. The patient must be warned on total joint reconstruction limitations and recent use history. Excessive physical activity and trauma that impact joint replaced, and imply in reconstruction early failure by looseness, fractures, in implant use.

Joint excessive demands replaced through superactivity or non-control of body weight. It is suggested monitoring patient through X-Ray. Acetabulum has a metal ring in external part to enable X-Ray acetabulum positioning viewing. We recommend not using magnetic resonance, as this stainless steel ring causes interference in the exam, damaging the diagnosis, and it can also cause implant dislocation, causing problem to the patient.

Assays (Load Restrictions/Mechanical Strength) Roughness and sphericity requirements were observed, in order to eliminate lack of congruence between Polymer Cemented One-Component Acetabulum - BM and femoral metal head (not object of this registration).

Adverse Effects: Early loosening of total hip replacement components can occur. Early mechanical loosening can be the result of defective fastening or hidden infection; biological complication or mechanical problem early loosening (high tensions located). Peripheral neuropathies and heterotopic bone formation were informed in total hip arthroplasty. Subclinical nerve damage happens more frequently, possibly due to surgical trauma result.

Patient metal sensibility reactions, followed by total joint replacement, have been rarely reported. Sensibilization meaning and effects still required future evidence and clinical evaluation.

Foreign material implantation in tissues can result in histological reaction involving several macrophage and fibroblast sizes. Similar changes can occur, as a precursor, or during healing process.

Dislocations and sub-luxation have been reported, due to implant component improper positioning. Fibrous tissue muscle and slackness can also enable such conditions.

Informations to Be Supplied by Patient: Talk to your orthopaedist about any medical condition that can impact surgery. Total prostheses are successful surgeries in more than 90% of patients. Whenever there are complications, most of them can be treated successfully. Possible complications include:

Infection: the infection can occur in operative injury or can be deep (around the prosthesis). It can appear while the patient is at hospital or at home. It can even appear years later. Smaller infections in operative injury are usually treated with antibiotics. Greater or deep infections can require new surgery (for deep cleaning) or even prosthesis removal. Any infection in the body (bladder, throat, teeth, ears, etc.) can carry germs through the blood to prosthesis and cause infection.

Thrombosis: they are blood clots resulting from several factors, including mobility reduction, and make that the blood circulate more slowly in the veins, enabling thrombosis appearance. Blood clots can be suspect if there is pain and edema in your thigh or calf. If that happens, your orthopaedist must make tests and exams to assess your leg veins, and evaluation and follow-up by vascular surgeon can be required. Several measures can be used to reduce blood clot possibility:

- Measurement to 'thin' the blood (anticoagulants)
- Elastic stockings
- Exercise to increase leg circulation.
- Plastic boots that inflate and compress feet sole and calf, increasing venous return.

Even through use of such preventive measures, blood clots still must happen. If you notice an edema (swelling), redness or calf pain after hospital discharge, contact your orthopaedist.

Slackness: prosthesis slackness inside your bone can occur upon surgery. That can cause pain if the slackness is significant. A revision surgery (prosthesis change) may be required. New materials and new fastening methods can minimize such problem.

Luxation: eventually, upon a hip prosthesis, the prosthesis head can leave the acetabulum (cap). This is what we call luxation. In most cases, the hip can be repositioned in place, not requiring a new surgery. To prevent luxation, it is important that muscles are strong (do the exercises recommended by your orthopaedist) and that hip is not flexed above 90 degrees in the first months.

Wear: some kind of wear can be detected in any kind of prosthesis. Excessive wear can enable slackness and require a new surgery (revision).

Prosthesis Break: As they are implanted, biomaterials used in prostheses necessarily contact body fluids. Such fluids, even though they seem harmless, can, through time, degrade significantly most materials of considerable 'chemical inertia'. Additionally, a large share of orthopedic implants is submitted to static and/or cyclic mechanical stresses, that can cause material failure, i.e., prosthesis break.

With the materials currently used in prosthesis, implant fracture is very hard. However, if that happens, a new surgery will be required to change the fractured prosthesis.

Nervous lesions: your prostheses neighboring nerves can be injured during a surgery, even though that is very rare. This lesion can occur easily whenever the surgery involves major deformity corrections in pelvis or too short leg elongation due to wear. Through time, such nervous lesions generally improve and can fully recover. Eventually, involved nerve surgical exploration can be chosen.

Implanted Product Evaluations: Upon implantation, within intraoperative, the responsible professional must perform the radiological control to check the product correct positioning. The responsible professional must perform, and is responsible for, clinical and radiological evaluations submitted upon surgical procedure within the frequency set forth by him/her, to check the implant status and bone consolidation evolution. If the product is not in correct positioning, or has any non-compliance, the surgeon is responsible for taking the most appropriate corrective action.

Useful Information to Prevent Implantation-Related Risks: To decrease implantation-related risks, following strictly: contraindications, instructions for use and all information provided in product "Instructions For Use".

Implantation Intrinsic Risks: Polymer Cemented One-Component Acetabulum is manufactured with materials with acknowledged biomedical use. As per NBR ISO 5834-2 - Ultra High Molecular Weight Polyethylene Surgery Implants - Part 2: Molded Form Products, ASTM F138 - Standard Specification for Wrought 18 Chromium-14 Nickel-2,5 Molybdenum Stainless Steel Bar and Wire for Surgical Implants (Uns S31673) and NBR ISO 5833 - Surgery Implants - Acrylic Resin Cements.

Contamination Risks: There are biological contamination risks and viral disease transmission, including HIV and hepatitis, as Polymer Cemented One-Component Acetabulum - BM components contact tissue and body fluids. Explanted products must be treated as highly contaminant.

Sterilization: This product is sterilized by Ethylene Oxide Sterilization (E.T.O.), with the following valid date: 5 years. In case sterilization valid date is expire or the packaging is damage, the product must be returned to the responsible dealer or directly to Biomecânica.

Product Disposal: Implants that comprise Polymer Cemented One-Component Acetabulum explants from patients must be properly disposed by the hospital institution. The hospital institution is responsible for implant complete identity and usability destruction, preventing its reuse. The hospital institution is responsible for the method used to implant identity and usability destruction. Biomecânica recommends that explanted implants are deformed mechanically assisted by a hammer or impact press, and they must be identified afterwards with the statement "Unfit for Use".

As provided in Resolution no. 2605, of 08/11/2006, implanting devices of any nature rated as for single use, are forbidden to be reprocessed.

Support materials: Polymer Cemented One-Component Acetabulum - BM Instructions for Use.

TRACEABILITY: Along with the product, there are 6 labels providing data required to product traceability. A label must be necessarily stuck in patient clinical record, another one in the report delivered to the patient, another in fiscal documentation that generates collection, at AIH, in case of patient serviced by SUS, or in sale invoice, in case of patient serviced by complementary health system, another one provided to supplier control (distribution record history - RHD), another one provided to responsible surgeon control (main), and the last one to health insurance, if applicable. Information provided in these labels are listed in item 'Presentation Form' of these Instructions for Use.

Marking: Polymer Cemented One-Component Acetabulum contains the following laser-etched information, to enable traceability and patient follow-up in postoperative:

- Manufacturing Batch Number
- Dimension

After Sales (Customer Complaint): In case it is required to perform any complaint related to Polymer Cemented One-Component Acetabulum - BM related to any adverse effect that impacts the user safety, including non-operating product, implanting metal component damage, severe problems or death, related to these components, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7926. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website:

<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Submitting Material for Manufacturer Analysis: In case of submitting implants for manufacturer analysis, it must be sanitized using a wide spectrum bactericide and fungicide solution. Next, it must be disinfected or sterilized with steam in autoclave or ethylene oxide. They must be forwarded to Biomecânica in intact packagings, identified with cleaning method, sterilization and data of product.

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