discomfort or abnormal sensations due to the presence or migration of the medical device or due to the surgical procedure and migration, limb edema, pseudarthrosis; movement limitations; failure of fracture healing. Due to the incompatibility of the materials, design and dimension the use of components from different manufacturer is forbidden. The use of different manufacturers may cause problems to the product and bring risks. The company disclaims all responsibility of any adverse event that may occur. Do not use implants from different manufacturers. Therefore it is recommended that the products have the same provenance.

**Sterilization of the external fixators:** The external fixators are distributed not sterile. Before using the external fixators they should be sterilized.

It is recommended the steam sterilization in Autoclave at hospital (ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices)

**AUTOCLAVE:** It is a device for sterilization through saturated steam under pressure.

**INSTRUCCIONES:** Before starting the device check: if the switch is on, if the water valve is open and the flush valve is closed.

1 – open the door of the equipment; 2 – properly settle the material to be sterilized; 3 – close the door of the equipment; 4 – select the desired cycle according to the material to be sterilized; 5 – Turn on the switch; 6 – the cycle will run automatically in the sequence; 7 – when the lamp "cycle end" turns on partially open the door for about ten (10) minutes to cool down the material.

It is recommended to follow the parameters of physical sterilization in autoclave (saturated steam) showed below.

**Table – Parameters of physical sterilization in autoclave**

Ciclo	Temperature	Exposure Time (minimum)	Time to Drying (maximum)
Conventional (1 atm de pressure)	121°C (250°F)	30min.	-
High Vacuum	132°C (270°F)	6min.	15min.

**Note.:** Register the time when the heat of the sterilization chamber reaches the desired temperature.

**For further information check the instruction manual that comes with each autoclave.**

**Other method of sterilization:** Besides the autoclave the sterilization process may be in Ethylene Oxide (E.T.O.) – parameters and procedures established by the validation protocol and ISO 11135– Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

**NOTE:** The methods of sterilization, the equipments, the controls and the instructions on the applied sterilization under responsibility of the hospital. Storage conditions of the sterilized articles: - as the environment: it must be cleaned; ventilated and dried and it should be restricted to the team of the sector; - as the article: after the sterilization process do not settle them on cold surfaces (stone or stainless steel), use baskets or perforated recipients until they cool down; - wrapping (raw-cotton fabric, nonwoven fabric, surgical paper, crepe paper, film paper, tyvec or perforated metallic boxes) must remain intact and little handled in order to avoid the packages to tear or lose the seal; - store it in cabinets with shelves; - identify the shelves in order to make it easier for the removal of the material; the material should be stored according to the expiration date of the sterilization in order to make it easier for the distribution of the material and to prevent having expired material in the stock; store the sterile materials separately from the non-sterile ones in order to reduce the level of external contaminants.

**Storage and transport of the product:** The storage location should be clean, dry and illuminated in a way that the proper conditions of storage and transport are preserved as well as the physical and chemical integrity. The product must be stored and transported in dry and fresh environment at room temperature (Max. 35°C) and relative humidity around 30% to 70% and away from direct sunlight. The effects of vibration, shocks, corrosion, defective settlement during the handling and transport, improper piling in the storage must be avoided. The product must be carefully handled in order to avoid damages that may impair the quality of the material and the safety of the patient. The products can not be stored in direct contact with the floor (minimum height = 20cm), they can not remain in very high places close to lamps because it may respect the packaging or damage the label. They should not be stored with contaminant substances such as cleaning materials, insecticide, pesticide, etc. The transporter must be advised of the content and delivery deadline. Carefully transport the products in order to avoid fall and attrition in order not to damage the device.

**Identification:** The external fixators come with labels that present the following information: Data of the manufacturer; Technician in charge; nº of the register with ANVISA; Technical name of the product; Commercial name of the product; Commercial model of the product; Raw-material; Quantity of items per package; Lot number; Code of the product; Manufacturing date; Expiration date; Address of the manufacturer; The external fixators come with the following engraved information for: traceability; Biomechanica logo; Number of the manufacturing lot and code.

**Customer Service:** In case it is necessary to make a complaint about the use of the external fixators related to some adverse effect that may affect the safety of the user the responsible surgeon must inform the Health Authorities and Biomechanica through the email - e-mail sac@biomechanica.com.br – or by the phone - 0xx14 2104 7900. In case of doubts the responsible surgeon or the health Professional may inform about the adverse event through NOTIVISA on the webpage of ANVISA: <http://www.anvisa.gov.br/hotite/notivisa/index.htm>

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