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<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

ENGLISH EUA

Necessary information for the user to identify the product and its content as well as the description of the physical principles and bases of the technology of the product applied to its operation and action, available models and composition.

The bone plates with a device for Sliding fixation are metallic and anatomical implants used in large bones of the lower limbs of the human skeleton to provide alignment, reduction and temporary or permanent stabilization of different types of fractures. Part of the body of the Bone plates with a device for sliding fixation is straight with holes for the insertion of the screws. The other part of its body has an angled tube. The inner of this tube is oblong to receive the Sliding Screw II whose dimension is 7,9mm x 7,1mm.

The bone plates with a device for Sliding Fixation are used in conjunction with non-cannulated screws – BM - Ø4,5mm cortical and/or Ø6,5mm cancellous models. They are available in two options – with or without the locking ring. The bone plates with a device for sliding fixation are identified and packed in plastic packs.

Models

Figure 3 - Autodynamic Plate 95°

Figure 4 - Autodynamic Plate 135°

Composition: The bone plates with a device for sliding fixation are manufactured in stainless steel. The chemical, mechanical and metallographic composition of those materials is in accordance with ASTM F138 standards. They are available in two options, with or without the locking ring. The ring that comes with the Bone plates with a device for sliding fixation is in titanium in accordance with ASTM F136 standard. Those materials are proper combinations when articulation is not intended (NBRISO21534). The stainless steel and titanium used for the

manufacturing of the Bone plates with a device for sliding fixation are analyzed in outsourced laboratories for the verification of their compliance with ASTM F138 standard for stainless steel and ASTM F136 - for titanium alloy.

The Bone Plates with Sliding Fixation Device will be distributed in a unitary form according to the table below:

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PLACA ÓSSEA COM DISPOSITIVO DE FIXAÇÃO DESLIZANTE		
Code	Qty	Descrição
3970-02-000	1	Autodynamic plate 95° of 02 holes x 52mm length with Locking Ring
3970-03-000	1	Autodynamic plate 95° of 03 holes x 70mm length with Locking Ring
3970-04-000	1	Autodynamic plate 95° of 04 holes x 88mm length with Locking Ring
3970-05-000	1	Autodynamic plate 95° of 05 holes x 106mm length with Locking Ring
3970-06-000	1	Autodynamic plate 95° of 06 holes x 124mm length with Locking Ring
3970-07-000	1	Autodynamic plate 95° of 07 holes x 142mm length with Locking Ring
3970-08-000	1	Autodynamic plate 95° of 08 holes x 160mm length with Locking Ring
3970-09-000	1	Autodynamic plate 95° of 09 holes x 178mm length with Locking Ring
3970-10-000	1	Autodynamic plate 95° of 10 holes x 196mm length with Locking Ring
3970-02-000	1	Autodynamic plate 95° of 11 holes x 214mm length with Locking Ring
3970-12-000	1	Autodynamic plate 95° of 12 holes x 232mm length with Locking Ring
3970-02-000	1	Autodynamic plate 95° of 13 holes x 250mm length with Locking Ring
3970-14-000	1	Autodynamic plate 95° of 14 holes x 268mm length with Locking Ring
3970-15-000	1	Autodynamic plate 95° of 15 holes x 286mm length with Locking Ring
3970-16-000	1	Autodynamic plate 95° of 16 holes x 304mm length with Locking Ring
3970-17-000	1	Autodynamic plate 95° of 17 holes x 322mm length with Locking Ring
3970-18-000	1	Autodynamic plate 95° of 18 holes x 340mm length with Locking Ring
3970-19-000	1	Autodynamic plate 95° of 19 holes x 358mm length with Locking Ring
3970-20-000	1	Autodynamic plate 95° of 20 holes x 376mm length with Locking Ring
3970-21-000	1	Autodynamic plate 95° of 21 holes x 394mm length with Locking Ring
3970-22-000	1	Autodynamic plate 95° of 22 holes x 412mm length with Locking Ring
3730-02-000	1	Autodynamic plate self compression 95° of 02 holes x 51,5mm length without Locking Ring
3730-03-000	1	Autodynamic plate self compression 95° of 03 holes x 67,5mm length without Locking Ring
3730-04-000	1	Autodynamic plate self compression 95° of 04 holes x 83,5mm length without Locking Ring
3730-05-000	1	Autodynamic plate self compression 95° of 05 holes x 99,5mm length without Locking Ring
3730-06-000	1	Autodynamic plate self compression 95° of 06 holes x 114,5mm length without Locking Ring
3730-07-000	1	Autodynamic plate self compression 95° of 07 holes x 131,5mm length without Locking Ring
3730-08-000	1	Autodynamic plate self compression 95° of 08 holes x 147,5mm length without Locking Ring
3730-09-000	1	Autodynamic plate self compression 95° of 09 holes x 163,5mm length without Locking Ring
3730-10-000	1	Autodynamic plate self compression 95° of 10 holes x 179,5mm length without Locking Ring
3730-11-000	1	Autodynamic plate self compression 95° of 11 holes x 195,5mm length without Locking Ring
3730-12-000	1	Autodynamic plate self compression 95° of 12 holes x 211,5mm length without Locking Ring
3730-13-000	1	Autodynamic plate self compression 95° of 04 holes x 227,5mm length without Locking Ring
3730-14-000	1	Autodynamic plate self compression 95° of 14 holes x 243,5mm length without Locking Ring
3730-15-000	1	Autodynamic plate self compression 95° of 15 holes x 259,5mm length without Locking Ring
3730-16-000	1	Autodynamic plate self compression 95° of 16 holes x 275,5mm length without Locking Ring
3730-17-000	1	Autodynamic plate self compression 95° of 17 holes x 291,5mm length without Locking Ring

3730-18-000	1	Autodynamic plate self compression 95° of 18 holes x 307,5mm length without Locking Ring
3730-19-000	1	Autodynamic plate self compression 95° of 19 holes x 323,5mm length without Locking Ring
3730-20-000	1	Autodynamic plate self compression 95° of 20 holes x 339,5mm length without Locking Ring
3730-21-000	1	Autodynamic plate self compression 95° of 21 holes x 355,5mm length without Locking Ring
3730-22-000	1	Autodynamic plate self compression 95° of 22 holes x 371,5mm length without Locking Ring
3975-02-000	1	Autodynamic plate 135° of 02 holes x 62mm length with Locking Ring
3975-03-000	1	Autodynamic plate 135° of 03 holes x 80mm length with Locking Ring
3975-04-000	1	Autodynamic plate 135° of 04 holes x 98mm length with Locking Ring
3975-05-000	1	Autodynamic plate 135° of 05 holes x 116mm length with Locking Ring
3975-06-000	1	Autodynamic plate 135° of 06 holes x 134mm length with Locking Ring
3975-07-000	1	Autodynamic plate 135° of 07 holes x 152mm length with Locking Ring
3975-08-000	1	Autodynamic plate 135° of 08 holes x 170mm length with Locking Ring
3975-09-000	1	Autodynamic plate 135° of 09 holes x 188mm length with Locking Ring
3975-10-000	1	Autodynamic plate 135° of 10 holes x 220mm length with Locking Ring
3975-11-000	1	Autodynamic plate 135° of 11 holes x 224mm length with Locking Ring
3975-12-000	1	Autodynamic plate 135° of 12 holes x 242mm length with Locking Ring
3975-13-000	1	Autodynamic plate 135° of 13 holes x 260mm length with Locking Ring
3975-14-000	1	Autodynamic plate 135° of 14 holes x 278mm length with Locking Ring
3975-15-000	1	Autodynamic plate 135° of 15 holes x 296mm length with Locking Ring
3975-16-000	1	Autodynamic plate 135° of 16 holes x 314mm length with Locking Ring
3975-17-000	1	Autodynamic plate 135° of 17 holes x 332mm length with Locking Ring
3975-18-000	1	Autodynamic plate 135° of 18 holes x 350mm length with Locking Ring
3975-19-000	1	Autodynamic plate 135° of 19 holes x 368mm length with Locking Ring
3975-20-000	1	Autodynamic plate 135° of 20 holes x 386mm length with Locking Ring
3975-21-000	1	Autodynamic plate 135° of 21 holes x 404mm length with Locking Ring
3975-22-000	1	Autodynamic plate 135° of 22 holes x 422mm length with Locking Ring
3720-03-000	1	Autodynamic plate self compression 135° of 02 holes x 60mm length without Locking Ring
3720-03-000	1	Autodynamic plate self compression 135° of 03 holes x 76mm length without Locking Ring
3720-04-000	1	Autodynamic plate self compression 135° of 04 holes x 92mm length without Locking Ring
3720-05-000	1	Autodynamic plate self compression 135° of 05 holes x 108mm length without Locking Ring
3720-06-000	1	Autodynamic plate self compression 135° of 06 holes x 124mm length without Locking Ring
3720-07-000	1	Autodynamic plate self compression 135° of 07 holes x 140mm length without Locking Ring
3720-08-000	1	Autodynamic plate self compression 135° of 08 holes x 156mm length without Locking Ring
3720-09-000	1	Autodynamic plate self compression 135° of 09 holes x 172mm length without Locking Ring
3720-10-000	1	Autodynamic plate self compression 135° of 10 holes x 188mm length without Locking Ring
3720-11-000	1	Autodynamic plate self compression 135° of 11 holes x 204mm length without Locking Ring
3720-12-000	1	Autodynamic plate self compression 135° of 12 holes x 220mm length without Locking Ring
3720-13-000	1	Autodynamic plate self compression 135° of 13 holes x 236mm length without Locking Ring
3720-14-000	1	Autodynamic plate self compression 135° of 14 holes x 252mm length without Locking Ring
3720-15-000	1	Autodynamic plate self compression 135° of 15 holes x 268mm length without Locking Ring
3720-16-000	1	Autodynamic plate self compression 135° of 16 holes x 284mm length without Locking Ring
3720-17-000	1	Autodynamic plate self compression 135° of 17 holes x 300mm length without Locking Ring
3720-18-000	1	Autodynamic plate self compression 135° of 18 holes x 316mm length without Locking Ring
3720-19-000	1	Autodynamic plate self compression 135° of 19 holes x 332mm length without Locking Ring
3720-20-000	1	Autodynamic plate self compression 135° of 20 holes x 348mm length without Locking Ring
3720-21-000	1	Autodynamic plate self compression 135° of 21 holes x 364mm length without Locking Ring
3720-22-000	1	Autodynamic plate self compression 135° of 22 holes x 380mm length without Locking Ring

3755-00-050	1	Sliding Screw II 50 mm total length (22 mm thread length)
3755-00-055	1	Sliding Screw II 55 mm total length (22 mm thread length)
3755-00-060	1	Sliding Screw II 60 mm total length (22 mm thread length)
3755-00-065	1	Sliding Screw II 65 mm total length (22 mm thread length)
3755-00-070	1	Sliding Screw II 70 mm total length (22 mm thread length)
3755-00-075	1	Sliding Screw II 75 mm total length (22 mm thread length)
3755-00-080	1	Sliding Screw II 80 mm total length (22 mm thread length)
3755-00-085	1	Sliding Screw II 85 mm total length (22 mm thread length)
3755-00-090	1	Sliding Screw II 90 mm total length (22 mm thread length)
3755-00-095	1	Sliding Screw II 95 mm total length (22 mm thread length)
3755-00-100	1	Sliding Screw II 100 mm total length (22 mm thread length)
3755-00-105	1	Sliding Screw II 105 mm total length (22 mm thread length)
3755-00-110	1	Sliding Screw II 110 mm total length (22 mm thread length)
3755-00-115	1	Sliding Screw II 115 mm total length (22 mm thread length)
3755-00-120	1	Sliding Screw II 120 mm total length (22 mm thread length)
3755-00-125	1	Sliding Screw II 125 mm total length (22 mm thread length)
3755-00-130	1	Sliding Screw II 130 mm total length (22 mm thread length)
3700-0A	1	Compression screw with 27 mm total length

Attention


Biomechanica Instrument set is necessary for the implantation of Bone Plates with a device for sliding fixation.

Presentation

The Bone plates with a device for sliding fixation are composed of integrating or complementary parts such as II Sliding screws and Connection Screw. The implants are packed in plastic envelops which are labeled with a nonsterile tag and outer labeling. The sliding screws are distributed in separate packs to prevent them to collide with the Plates and also because of the variation in the dimension of the screws. The outer labeling consists of the traceability tags which provide the following information:

Name of the manufacturer; Commercial name of the product; Lot Number; Number of the register with ANVISA; Code of the product; Quantity; Description of the Product of the Pack (model); Manufacturing Date; Expiration Date; Undefined Expiration date; Description of the raw-material used for the manufacturing of the product (Composition); Manufacturer's address;
Name of the Technician in charge; And the sayings "Non-Sterile Product"; "DO NOT REPROCESS"; "Single Use Product. Do not re-use"; "Vide Instruction for Use before Use"; Warnings/Precautions/Cares/Symbology: Vide Instructions for Use.

The following graphical information comes on the packaging of the product: Safety symbols printed on the plastic envelops as showed in the model bellow:

Fragile, handle with care Keep dry Keep it Away from the sun
Figure 1 – Symbology according to NBR ISO 15223.

Special conditions for storage, preservation and/or handling of the product.

Cares with the handling and transport of the medical device:

Do not use the product in case it is damaged; the implants must be carefully handled in order to avoid damages that may impair the quality of the material or the safety of the patient; keep the product in ventilated and dry places away from the light and weather conditions. Shocks and improper pilling should be avoided during transportation.

Indication, purpose or use to which the medical device is intended according to the manufacturer's indication

Indication for Use

Device for Bone fixation with holes for the fixation of the screws to the bone. The plates are implants that are fixed to the bone in order to allow fixation. They are indicated for the large bones of the lower limbs of the human skeleton to provide alignment, reduction and temporary or permanente stabilization of the different types of fractures. They are available in several lengths and number of holes. The plates with a device for sliding fixation are not re-usable.

95° Autodynamic Plate

95° Autodynamic Plate is indicated for distal fixation of the femur, intercondylar, intertrochanteric and proximal subtrichanteric fractures. The plates are straight with a 95° angled tube in the distal region to be fixed to the fractured femoral condyle with the help of a II Sliding screw and a compression screw. They have holes that will receive the non-canalated screws – FB for the bone fixation. They are available in two options, with or without the locking ring (self-compression). The option with the locking ring allows the screw to be inserted into the bone in different positions (convergent, parallel and divergent at 360°) and the selfcompression holes for the insertion of self-compression screws. The locking ring in the holes of 95° Autodynamic Plates has an inner threaded cone for receiving the non-canalated screw – FB. The outside of this locking ring has a spherical surface that will connect with 95° Autodynamic Plate and which will also fix the set When tightened. The locking ring has a groove for its expansion when it receives a noncanalated screw – FB. The hole of 95° Autodynamic Plate is spherical similarly to the outside of the ring for receiving the ring and for the lock.

135° Autodyamic Plate

135° Autodynamic Plate is indicated for the proximal fixation of the femur in fractures of 31 A1, A2 and A3 types. It is straight with a 135° angled tube in the proximal region to be fixed to the fractured femoral neck with the help of a II Sliding screw and a compression screw. They have holes for receiving non-canalated screws – FB for the bone fixation. They are available in two

options, with or without the locking ring (self-compression). The option with the locking ring allows the screw to be inserted into the bone in different positions (convergent, parallel and divergent at 360°) and self-compression holes for the insertion of the self-compression screws. The locking ring in the holes of 135° Autodynamic Plates has an inner threaded cone to receive non-cannulated screw – FB. The outside of this ring has a spherical surface that will connect with 135° Autodynamic Plate and which will fix the set when tightened. The locking ring has a groove for its expansion when it receives non-cannulated screw – FB. The hole of 135° Autodynamic Plate is spherical similarly to the outside of the ring for receiving the set and for the lock.

Contraindication

The use of this implant is contraindicated in the following situations: metal sensitivity; fever or signs of local inflammation; osteoporosis; obesity; improper tissue coverage on the region to be operated; pregnancy; bone compromising due to disease and/or infection; unexplained increase in the sedimentation rate; increase in the WBC rate or strong suspicion in differential count; visible infection.

Warnings and Precautions

SINGLE USE PRODUCT. DESTROY AFTER EXPLANTED. DO NOT RE-USE THE PRODUCT. DO NOT REPROCESS

An orthopedic implant must be used for one single patient and for one single time. Although it seems to be intact, previous tensions may have created imperfections that may impair the success of the implant. An improper selection of the implant may promote uncommon stress on the implant which results in its subsequent failure. Due to its presentation and characteristics, this product does not present any side effect, therefore, it is recommended to monitor the patient's recovery through x-ray according to the International Safety Standards. Do not use magnetic resonance imaging because the product in stainless steel causes interference in the test and that impairs the diagnosis besides it may provoke the detachment of the implant which may bring problems to the patient. The Bone plates with a device for sliding fixation are designed to work as an auxiliary element for the bone fixation and for stabilization in the treatment of fractures of upper and lower limbs and not to replace healthy structures of the human skeleton. The surgeon must be aware of the chance of the patient's physical development which may ask for an implant revision. The surgeon must have enough knowledge and be familiar with osteosynthesis and its limitations besides the pre and postoperative, the adopted surgical techniques, precautions and potential risks; the surgeon must have attention to the bone preparation for the insertion of the bone plates with a device for bone fixation because the stability and success of the fixation may be compromised. In order to guarantee a proper implantation, Biomechanica instruments must be used because they were designed and manufactured specially for the use of those implants. The use of instruments from different manufacturers may compromise the surgery besides it is not in accordance with the register of the product in the National Health Surveillance Agency;

Do not use the product in case it is damaged;
NON-STERILE PRODUCT – STERILIZE BEFORE USE; SINGLE USE PRODUCT –
THE RE-USE OF THIS PRODUCT IS EXTREMELY FORBIDDEN;

Conditions of storage: Keep it in ventilated and dry places away from the light and weather conditions;

Do not use the product in case it is damaged;

Observation: Components of implants from different manufacturers should not be used together; therefore we recommend the products to have the same provenance.

Manufacturing Date, expiration date and product batch: VIDE LABEL.

Forming Limit of the Product

Scratches, notches or any other thing that might damage or mark the implant should be avoided. Those defects are stress concentration and may be sites of crack nucleation which would compromise the resistance to corrosion and lead to the fracture of the implant. In certain situations it is necessary to mold the implant in order to perfectly fit it to the treated bone. In case the modeling of the implant is necessary the bone anatomy is the limit for the it.

Warning

In terms of mechanical resistance of this design, there aren't specifications for the Plates for bone synthesis related to the load limit defined for each type of plate. The bone plates with a device for sliding fixation are manufactured to help the orthopedist in internal fixation procedures. It is not expected the implant to have the same performance of a healthy bone or to support intense or permanent mechanical stress. The performance of the implant is related to the patient's weight, level of physical activity, compliance with the instructions given by the surgeon for the postoperative and especially in cases of failure or retardation of healing. The bone healing is strictly linked to the mechanical contribution and biological reaction of the implant. The biological reaction will be determined in most of the cases through the stability of the fracture.

Side Effects

The mechanical looseness may be the result of a defective fixation, an unstable reconstruction of the fractured bone or occult infection. Metal sensitivity reactions in patients have been rarely reported.

The implantation of foreign material in the tissues results in histological reactions that involve several sizes of macrophage and fibroblasts. The chemical importance of this effect is uncertain as well as similar changes may come as a precursor or during the healing process.

Information to be given to the patient

The sayings of this instruction in the items: Indication – Contraindications – Information for Use – Potential Adverse Effect – Precaution and Warnings. The patient must be informed as the importance of a postoperative follow-up. The failure to perform a follow-up prevents the detection

of post-surgical problems such as the looseness of components or the incidence of osteolysis. The patient must be informed as the type of the used implant, the brand, code, name, lot, the surgeon's name and the date of the surgery.

Traceability

Two traceability tags with the data of the used implants come with the pack of the implantable component. Those tags are equal to the ones of the outer labeling of the product. The information present in those traceability tags is listed in item – “Presentation” in this instruction for use. One of those tags should be pasted on the patient's report and another one should be available for the control of the hospital. The hospital should register those data in its files and also forward them to the patient. It is mandatory that all information for the traceability of the implanted product is registered on the patient's report. The tags on the patient's report allow the traceability of the used product. Among the main information we highlight as indispensable – the name of the manufacturer of the used implant, the name of the used implant, the code of the implant, the lot number and number of the register of the product with ANVISA. That information is described in the traceability tags that come with the product and in the outer labeling. Other information such as the date of the surgery, the patient's name, the surgeon's name as well as the patient's weight and age should also be considered indispensable among other asked in the patient's report.

Performance foreseen in the General Requirements of ANVISA regulations that rules on the Essential requirements of safety and efficiency applicable to the products described in Resolution – RDC 56 from 04/06/2001:

The bone plates with a device for sliding fixation are manufactured in stainless steel according to ASTM F138 and the locking ring is manufactured according to ASTM F136. They are verified through laboratory reports of tests and quality certificates of the product. Those materials are acceptable for the manufacturing of the implants as defined in annex A of NBRISO21534 standard. According to NBRISO21534, in annex C.1, the metallic combination of non-articulating contact surfaces between the stainless steel NBRISO5832-1 (similar to ASTM F138) and titanium alloy NBR ISO 5832-3 (similar to ASTM F136) are considered proper combinations. The manufacturing and packaging of the Bone plates with a device for sliding fixation is based on specific internal procedures of cleaning and packaging. The bone plates with a device for sliding fixation are distributed with information that encompasses indications and enough information for the correct handling of the product as well as the safety of people involved with the handling.

Ancillary Components

The ancillary components for the implantation of the Bone plates with a device for sliding fixation are described below

Ancillary components specifically for 95° Autodynamic Plates without the locking ring

- Ø4,5mm cortical screws with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5mm Cancellous screws with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5 cancellous screw with 16mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5 cancellous screws with 32mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)

Ancillary components specifically for 95° Autodynamic plates with a locking ring

- Ø4,5mm cortical screws – FB with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5mm cancellous screws – FB with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5 cancellous screws – FB with 16mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5 cancellous screws – FB with 32mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)

Ancillary components specifically for 135° Autodynamic Plates without the locking ring

- Ø4,5mm cortical screws with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5mm cancellous screw with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5 mm cancellous screw with 16mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5 mm cancellous screws with 32mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)

Ancillary components specifically for 135° Autodynamic plates with the locking ring

- Ø4,5mm cortical screws – FB with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5mm cancellous screws – FB with full thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5mm cancellous screws – FB with 16mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)
- Ø6,5mm cancellous screws – FB with 32mm partial thread manufactured in austenitic stainless steel according to ASTM F138 standard. (not object of this registration and not part of this product)

The ancillary components described above are not part of this register therefore they have to be verified separately however they are associated with this process of implantation and without them the product can not be implanted.

Accessories

The bone plates with a device for sliding fixation do not present any accessory with the purpose of integrating with the medical device.

Evaluation of the implanted product

After the implantation and during the intraoperative the professional in charge must perform a radiological control in order to verify the correct positioning of the product. The Professional in charge should perform the clinical and radiological evaluations after the surgical procedure at a frequency stipulated by him to verify the conditions of the implant and the evolution of the bone healing. In case the product is out of its correct positioning or presents some non-conformity the surgeon must take the most appropriate corrective actions.

Useful information for avoiding risks due to implantation

In order to minimize the risks resulted from implantation the following should be strictly obeyed: contraindications, instructions for use and all information in the "Instruction for use" of the product.

Intrinsic risks of implantation

The bone plates with a sliding fixation are made of well known biomedical use materials that comply with the requirements of ASTM F138 -00- Standard Specification For Wrought 18 Chromium-14 Nickel-2,5 Molybdenum Stainless Steel Bar And Wire For Surgical Implants (Uns S31673) Standards in the manufacturing of the plates for the manufacturing of the plates and ASTM F136 - 02 - Standard Specification For Wrought Titanium-6aluminum-4vanadium Eli (Extra Low Interstitial) Alloy For Surgical Implant Applications (Uns R56401) for the manufacturing of the rings.

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the bone plates with a device for sliding fixation explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Additional information about the procedures before the use of the product; sterilization method

Risks of contamination

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the metallic component of the Bone plates with a device for sliding fixation gets in contact with body tissues and fluids. The explanted products must be treated as highly contaminant.

Decontamination, cleaning and Sterilization

The product comes in its original packaging previously cleaned and ready for sterilization. It should be carefully handled in order to avoid contamination. In case the product is submitted to conditions that compromise its cleaning state, it should be cleaned before sterilization. Bellow some suggestions of cleaning procedures:

The cleaning may be performed through manual or mechanical methods.

Devices for cleaning procedures:

- Ultrasound washer
- Washer – sterilization device
- Washer – disinfecting device

Manual cleaning: it is a manual procedure where the dirtiness is removed through physical actions with the help of detergent, water and tools such as sponge and brush. The manual cleaning is mostly recommended because it is less aggressive to the implants and surgical instruments. For the manual cleaning one should use water at room temperature, nylon brushes instead of steel or abrasive ones because they may harm the protective layer of the material.

For the decontamination one should use a bactericide and antifungal solution of broad spectrum and avoid aggressive cleaning agents such as metal brushes so that the quality of the metallic implants is not impaired.

In case of cleaning machines the implants and surgical instruments must remain positioned in drawers in a way that they do not collide with each other to avoid damage to the material.

People in charge of the cleaning of implants and surgical instruments should be aware of the type of cleaning product used while complying with the concentrations indicated by the manufacturer. The implants and surgical instruments should be carefully rinsed up to they are free of residue and immediately dried after the cleaning. Type of cleaners Enzymatic cleaners: they are basically composed of enzymes, surfactant and solubilizers. A balanced combination of those elements makes it possible to remove the organic matter from the material in a short period of time.

Enzymatic solutions: they present an excellent cleaning action however they do not have any bactericide and bacteriostatic.

Enzymes: they are substances from live cells and govern the chemical reactions of the process. Once the enzymes are produced by the live cells they may be isolated and will still keep their

catalytic properties if certain conditions are preserved in their production. The enzymes are classified in three major functional groups depending on the type of substratum that they will affect: protease, lipase and amylase which tend to solubilize and detach from the articles.

Disinfection:

It is a process that destroys microorganisms (pathogenic or not) of the articles except the bacterial spores through physical or chemical means.

Level of disinfection: High Level: it destroys all microorganisms except the ones with a high number of spores => 2% Glutaraldehyde – 20 – 30 minutes.

Indications: rather hospital

Medium level: it eliminates vegetative bacteria, most of the virus, fungi and microbacteria =>1% Sodium hypochlorite – 30 minutes.

Indication: for UBS, Day care, nursing home and shelter.

Low level: It eliminates most of the bacteria, some virus and fungi but it does not eliminate microbacteria => 0,025% of hypochlorite sodium. Indication: nutrition.

Sterilization

This product is provided non-sterile. It must be sterilized before use. We recommend steam sterilization in Autoclave at hospital (ISO 11134: 1994 Sterilization of health care products – Requirements for validation and routine control – Industrial moist heat sterilization).

AUTOCLAVE

It is a device for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS

Before starting the device, check:

if the switch is on;

if the water valve is open;

if the flush valve is closed.

1 – open the door of the device;

2 – properly set the material to be sterilized;

3 – close the door of the device;

4 – select the desired cycle according to the material to be sterilized;

5 – turn on the general key;

6 – the cycle will automatically run in the sequence;

7 – when the lamp “end of cycle” turns on partially open the door for about ten (10) minutes for the material to cool down.

The following parameters of physical sterilization in autoclave are recommended (saturated steam):

Table 2 - Parameters of physical sterilization in autoclave

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
Conventional (1atm of pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Obs.: The time must be recorded when the heat of the sterilization chamber reaches the desired temperature.

FOR A BETTER CLARIFICATION CHECK THE INSTRUCTION MANUAL THAT COMES WITH EACH AUTOCLAVE.

Other sterilization method: Besides the autoclave the following sterilization method can be adopted as defined by the hospital:

- **Sterilization by ethylene oxide (E.T.O)** - parameters and procedures established in the validation protocol and in EN550:1994 Sterilization of Medical Devices – Validation and Routine Control of Ethylene Oxide Sterilization.

NOTE: The methods of sterilization, the devices, the controls and the instructions for sterilization are under responsibility of the hospital.

Cares with the sterilized articles and storage conditions of the sterilized articles

- As the environment: it should be clean; ventilated and dry and limited to the staff of the sector.
- As the article: do not set the products on cold surfaces (Stone or stainless steel) after the sterilization process; keep them in perforated baskets or recipients up to they cool down.
- Wrapping (raw cotton fabric, surgical paper, crepe paper, film paper, tyvec or perforated metallic boxes) should remain intact and be little handled in order to avoid the packs to rip or lose their seal; They should be stored in closed cabinets with shelves;
- The shelves should be identified in order to make it easier for the removal of the material;
- The material must be stored according to the expiration date of the sterilization in order to make it easier for the distribution and prevent to have expired material in stock;
- Store the sterile products separately from the sterile ones in order to minimize the level of external contaminants.

Send the material to the manufacturer for analysis

In case the implant is sent to the manufacturer for analysis it must be sanitized at hospital with the help of a bactericide and antifungal solution of broad spectrum. Then it should be disinfected or stem sterilized in autoclave or ethylene oxide. They should be sent to Biomechanica in intact packs identified with the method of cleaning and sterilization and data of the product.

Post sale (Client complaint): In case a complaint about the Bone plates with a device for sliding fixation is necessary and such complaint is **related to some adverse event that affects the**

safety of the user such as if the product is not working, if the implantable metallic device is damaged, or in case of deaths or serious problems related to those components, the responsible surgeon must inform the competent health authorities and Biomecanica about this adverse event through the email sac@biomecanica.com.br or by the phone 0xx14 2104 7900.

In case of questions the responsible surgeon or the health Professional may report the adverse event through NOTIVISA on the website of ANVISA: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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