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ENGLISH EUA

Product Description:

Accessories for Femoral Components - BM are non-absorbable, non-active implants, accessories for the Femoral Component with Restriction - BM, models intended for revision surgery, Biomecanica brand.

Accessories Models for Femoral Components - BM

BKS III Femoral Revision Stem: They are manufactured from a cast cobalt-chromium-molybdenum alloy, in accordance with the ASTM F75 standard and were designed for fixation in the distal portion of the femur, inside the femoral canal, without flexion force, connected to the Restricted Femoral Component - BM intended for revision surgery, brand Biomecanica. These models are intended to provide greater stability for the Restricted Femoral Component - BM intended for revision surgery, Biomecanica brand, fixed to the distal portion of the femur by means of bone cement.

The surface of the component that will be in contact with the femoral canal is sandblasted with alumina oxide. The lower surface of the component was designed to be connected and locked in the cone present in the Femoral Component with Restriction - BM intended for revision surgery, brand Biomecanica. The locking system must be tightened using a 3.5mm hexagonal wrench.


Figure 1 - BKS III Femoral Revision Stem


Figure 2 - Models Femoral Revision Stem BKS III, connected to the Femoral Component with Restriction - BM, intended for revision surgery

The BKS III Femoral Revision Stem models are universal, that is, they can be used both with Femoral Component with Restriction – BM Right and Left. The BKS III Femoral Revision Stem models are available in different sizes, as shown in Table 1. This variation in dimensions allows the surgeon to choose the appropriate model, according to the patient's bone characteristics and the surgical technique adopted.

Table 1 - List of BKS III Femoral Revision Stem Models

Code	Description	
2326-10-060	BKS III - Femoral Revision Stem Ø10 x 60 mm	
2326-10-070	BKS III - Femoral Revision Stem Ø10 x 70 mm	
2326-10-080	BKS III - Femoral Revision Stem Ø10 x 80 mm	
2326-10-100	BKS III - Femoral Revision Stem Ø10 x 100 mm	
2326-12-060	BKS III - Femoral Revision Stem Ø12 x 60 mm	
2326-12-070	BKS III - Femoral Revision Stem Ø12 x 70 mm	
2326-12-080	BKS III - Femoral Revision Stem Ø12 x 80 mm	
2326-12-100	BKS III - Femoral Revision Stem Ø12 x 100 mm	
2326-14-060	BKS III - Femoral Revision Stem Ø14 x 60 mm	
2326-14-070	BKS III - Femoral Revision Stem Ø14 x 70 mm	
2326-14-080	BKS III - Femoral Revision Stem Ø14 x 80 mm	
2326-14-100	BKS III - Femoral Revision Stem Ø14 x 100 mm	

2327-10-060	BKS III - Decentered Femoral Revision Stem Ø10 x 60 mm	
2327-10-070	BKS III - Decentered Femoral Revision Stem Ø10 x 70 mm	
2327-10-080	BKS III - Decentered Femoral Revision Stem Ø10 x 80 mm	
2327-10-100	BKS III - Decentered Femoral Revision Stem Ø10 x 100 mm	
2327-12-060	BKS III - Decentered Femoral Revision Stem Ø12 x 60 mm	
2327-12-070	BKS III - Decentered Femoral Revision Stem Ø12 x 70 mm	
2327-12-080	BKS III - Decentered Femoral Revision Stem Ø12 x 80 mm	
2327-12-100	BKS III - Decentered Femoral Revision Stem Ø12 x 100 mm	
2327-14-060	BKS III - Decentered Femoral Revision Stem Ø14 x 60 mm	
2327-14-070	BKS III - Decentered Femoral Revision Stem Ø14 x 70 mm	
2327-14-080	BKS III - Decentered Femoral Revision Stem Ø14 x 80 mm	
2327-14-100	BKS III - Decentered Femoral Revision Stem Ø14 x 100 mm	

The BKS III - Decentered Femoral Revision Stem models present a differentiated anatomy for cases with irregular anatomy of the femoral canal.


Femoral Shims: Femoral Shims are used to correct asymmetric segmental defects caused by advanced arthrosis producing posteromedial or lateral bone loss.

Femoral Shims are manufactured from cobalt-chromium-molybdenum cast alloy in accordance with ASTM F75. The Femoral Shims were designed for fixation on the femur, fitted to the Femoral Component with Restriction - BM intended for revision surgery. These models are intended to provide greater stability for the Restricted Femoral Component - BM intended for revision surgery, in cases where the surgeon judges that there is inadequate bone stock for anatomical alignment and fixation with bone cement. The surface of the Femoral Shims is sandblasted in alumina oxide.

The Shims are presented in different sizes, formats and angulations. This variation in dimensions allows the surgeon to choose the appropriate model, according to the patient's bone characteristics and the surgical technique adopted.

Table 2 - List of Femoral Wedge Models

Code	Description	Illustration
2391-00-010	Distal Femoral Shim - Right Medial / Left Side - EP (10.0mm)	
2390-00-010	Distal Femoral Shim - Right Side / Left Medial - EP (10.0mm)	
2391-00-005	Distal Femoral Shim - Right Medial / Left Side - EP (5.0mm)	
2390-00-005	Distal Femoral Shim - Right Side / Left Medial - EP (5.0mm)	
2391-01-010	Distal Femoral Shim - Right Medial / Left Side - P (10.0mm)	
2390-01-010	Distal Femoral Shim - Right Side / Left Medial - P (10.0mm)	
2391-01-005	Distal Femoral Shim - Right Medial / Left Side - P (5.0mm)	
2390-01-005	Distal Femoral Shim - Right Side / Left Medial - P (5.0mm)	
2391-02-010	Distal Femoral Shim - Right Medial / Left Side - M (10.0mm)	
2390-02-010	Distal Femoral Shim - Right Side / Left Medial - M (10.0mm)	
2391-02-005	Distal Femoral Shim - Right Medial / Left Side - M (5.0mm)	
2390-02-005	Distal Femoral Shim - Right Side / Left Medial - M (5.0mm)	
2391-03-010	Distal Femoral Shim - Right Medial / Left Side - L (10.0mm)	

2390-03-010	Distal Femoral Shim - Right Side / Left Medial - L (10.0mm)	
2391-03-005	Distal Femoral Shim - Right Medial / Left Side - L (5.0mm)	
2390-03-005	Distal Femoral Shim - Right Side / Left Medial - G (5.0mm)	
2391-04-010	Distal Femoral Shim - Right Medial / Left Lateral - EG (10.0mm)	
2390-04-010	Distal Femoral Shim - Right Side / Left Medial - EG (10.0mm)	
2391-04-005	Distal Femoral Shim - Right Medial / Left Lateral - EG (5.0mm)	
2390-04-005	Distal Femoral Shim - Right Side / Left Medial - EG (5.0mm)	
2392-00-005	Posterior Femoral Shim - EP (5.0mm)	
2392-00-010	Posterior Femoral Shim - EP (10.0mm)	
2392-01-005	Posterior Femoral Shim - P (5.0mm)	
2392-01-010	Posterior Femoral Shim - P (10.0mm)	
2392-02-005	Posterior Femoral Shim - M (5.0mm)	
2392-02-010	Posterior Femoral Shim - M (10.0mm)	
2392-03-005	Posterior Femoral Shim - G (5.0mm)	
2392-03-010	Posterior Femoral Shim - G (10.0mm)	
2392-04-005	Posterior Femoral Shim - EG (5.0mm)	
2392-04-010	Posterior Femoral Shim - EG (10.0mm)	

The Femoral Shims are fitted to the Femoral Component with Restriction - BM intended for revision surgery as shown below:


Figure 3 - Femoral Shims connected to the Femoral Component with Restriction - BM, model BKS III intended for revision surgery

 XXXXXX	Product Batch Number
 XXXX-XX	Manufacturing date
 XXXX-XX	Expiration date

REF XXXX-XX | **Code**

➤ **Sterile Product Identification:**

Sterilized in ETO Gas

➤ **Single Use Product Indication:**

“SINGLE USE PRODUCT - REPROCESSING PROHIBITED”

Commercial Presentation of the Product:

The components are packaged in double blister-type packaging (primary and secondary) sealed with surgical grade paper, plus outer protective packaging such as a cardboard box and ready-to-use labeling, already sterilized in ethylene oxide. Inside the package are:

- 1 instruction for use
- 6 labels (traceability labels) to maintain product traceability.

The external labeling is the product traceability labels that contain the following information:

- q) Manufacturer’s name
- r) Commercial Name of the Product
- s) Technical name
- t) Batch number
- u) ANVISA registration number
- v) Product code
- w) Quantity
- x) Description of the product that contains the packaging (model)
- y) Date of manufacture
- z) Maturity date: Indeterminate maturity
- aa) Description of the raw material used to manufacture the product
- bb) Manufacturer’s address (described on the plastic envelope)
- cc) Name of the Technical Manager
- dd) The saying: “**Sterile product**”
- ee) The saying: “**PROHIBITED TO REPROCESS**”.
- ff) Safety symbols printed on the labeling labels:

XXXXXX Product Lot Number	XXXX-XX Date of Manufacture	Manufacturer Data
XXXX-XX Expiration Date	XXXX-XX Product Code	“SINGLE USE PRODUCT – REPROCESSING PROHIBITED”

Sterilized in ETO Gas

All Accessories for Femoral Components - BM are sold individually, ready for use, already sterilized in ethylene oxide with a validity of 5 years. Accessories for Femoral Components - BM will be distributed individually, each in its respective packaging, as shown in Figure 4 and Figure 5.

Figure 4 - Presentation of the distribution of the Revision Femoral Stems (box size: 300x140x35mm)

Figure 5 - Femoral Shims distribution presentation form (box measurement: 130x130x40mm)

Other Packaging Information

Accessories for Femoral Components - BM are packaged in double blister packs, which are made of non-toxic PET (Polyethylene Terephthalate). This raw material is indicated for the manufacture of blisters and thermoformed cradles (Vacum - Forming) in the medicine, food, toys, electronic components, cosmetics and automobile industries. The two blister packs are sealed with dye-free and puncture-free Tyvek using the heat seal process.

Accessories
for Femoral
Components -
BM

Figure 6 - Accessories for Femoral Components - BM Packaged in Double Blister Packing

The double sealed blister pack is also enclosed in an outer cardboard carton. In the case of the “Calços Femorais” models, the outer cardboard packaging measures 130x130x40mm. For the “BKS III Femoral Revision Stem” models, the outer cardboard package measures 300x140x35mm. There is no difference in the information contained in the two packages, the only difference between them being the size and arrangement of the information.


Figure 7 - Femoral Revision Stem in the outer cardboard packaging

Traceability:

Along with the product, there are 6 labels containing the necessary data for product traceability. A label must be obligatorily affixed to the patient's clinical record, another to the report delivered to the patient, another to the tax documentation that generates the charge, on the AIH, in the case of a patient treated by the SUS, or on the sales invoice, in the case of patient assisted by the complementary health system, another available for the control of the supplier (historical distribution record - HDR, another available for the control of the responsible surgeon (main) and the last one for medical insurance, if any.

The following information is essential in the medical record: name of the implant used; manufacturer's name, date of surgery; product code; lot number; name of the patient who received the implant; surgeon's name.

All models of Accessories for Femoral Components - BM receive marking, through laser, according to the sequence below, within the limits allowed by the size of the model.

Logo - Batch - Raw Material Acronym (C4) - Product Code – Dimension


Figure 8 - Laser marking on the Femoral Revision Stem - BKSIII Figure	Figure 9 - Laser marking on the Femoral Shims
--	---

Ancillary Components: Accessories for Femoral Components - BM have the following products as ancillary components, as shown in Table 3 (Not objects of this registration and not part of this product).

Table 3 - List of Ancillary Components of Accessories for Femoral Components - BM

Description	Raw Material
Fixed Articular Insert with Stabilizer - BM	Ultra High Molecular Weight Polyethylene (UHMWPE)
Tibial Component - BM	Fused Cobalt-Chromium-Molybdenum Alloy (ASTM F75)
Accessories for Tibial Components - BM	Fused Cobalt-Chromium-Molybdenum Alloy (ASTM F75)
Patellar Component - BM	Ultra High Molecular Weight Polyethylene (UHMWPE)
*Orthopedic Cement Acrilex®	Bone Cement (NBR ISO 5833)

*The use of another type of biomaterial is allowed as an ancillary component of Accessories for Femoral Components - BM, replacing Acrilex® Orthopedic Cement, provided it is suitable for fixing the prosthesis to the bone and is regularly registered with ANVISA for this purpose.

Support Instruments: The instruments for placing the knee prosthesis are used as support for the total knee arthroplasty procedure and must be purchased separately.

Table 4 - Support Instruments for Total Knee Arthroplasty

Código	Descrição	Matéria-Prima
5221-00-000	Plane reamer (Grosa)	Stainless steel 304
5321-00-000	Alignment Bar	Stainless steel 304
5320-00-000	Aligment Rod with coupler	Stainless steel 304
5336-00-000	Guide for Tibial Drill bit	Stainless steel 304
5355-02-000	Patellar fórceps 3 grips	Stainless steel 304
5335-00-000	handle for tibial tray guide	Stainless steel 304
5352-00-000	tibial Plateau extractor	Stainless steel 420
5334-01-000	Tibial tray – small	Stainless steel 304
5334-02-000	Tibial tray – medium	Stainless steel 304
5334-03-000	Tibial tray – large	Stainless steel 304
5334-04-000	Tibial tray – extra large	Stainless steel 304
5341-00-000	Tibial trail impactor	Stainless steel 304
5340-00-000	Tibial trail extractor	Stainless steel 304
5322-00-000	A/P femoral measuring device	Stainless steel 304
5337-00-000	Tibial reamer – small/médium	Stainless steel 304
5338-00-000	Tibial reamer – large/extra large	Stainless steel 304
5319-00-000	Angular support	Stainless steel 304
5326-00-000	Extensor for IM tibial cutting guide	Stainless steel 304
5313-00-000	Guide for femoral alignment	Stainless steel 304
5328-00-000	Tibial depth ressection gauge	Stainless steel 304
5339-00-000	Handle for tibial reamer	Stainless steel 304
5354-00-000	Extractor pin plier	Stainless steel 420
5318-00-000	Distal femoral trimming guide	Stainless steel 304
5311-00-000	Drill Ø14,0 mm	Stainless steel 420

5310-00-000	Staggered drill bit	Stainless steel 420
5324-00-000	Intercondylar guide L/EL	Stainless steel 304
5364-00-000	Puncture	Stainless steel 304
5325-00-000	Tibial stem guide	Stainless steel 304
5344-01-000	Tibia test- small	Stainless steel 304
5344-02-000	Tibia test- medium	Stainless steel 304
5344-03-000	Tibia test - large	Stainless steel 304
5344-04-000	Tibia test- extra large	Stainless steel 304
5351-00-000	Tibial trimming guide	Stainless steel 304
5350-00-000	Valgus tibial trimming guide	Stainless steel 304
5348-00-000	Femoral impactor	Politec
5349-00-000	Tibial impactor	Politec
5345-08-000	Plateau test small/medium – 08mm	Polyphenylene
5345-10-000	Plateau test small/medium – 10mm	Polyphenylene
5345-12-000	Plateau test small/medium – 12mm	Polyphenylene
5345-15-000	Plateau test small/medium – 15mm	Polyphenylene
5346-08-000	Plateau test large/extra large – 08mm	Polyphenylene
5346-10-000	Plateau test large/extra large – 10mm	Polyphenylene
5346-12-000	Plateau test large/extra large – 12mm	Polyphenylene
5346-15-000	Plateau test large/extra large – 15mm	Polyphenylene
5341-01-000	Femur test –small/right	Stainless steel 304
5341-02-000	Femur test -small/left	Stainless steel 304
5341-03-000	Femur test -medium/right	Stainless steel 304
5341-04-000	Femur test- medium/left	Stainless steel 304
5341-05-000	Femur test -large/right	Stainless steel 304
5341-06-000	Femur test – large/left	Stainless steel 304
5341-07-000	Femur test – extra large/left	Stainless steel 304
5341-08-000	Femur test – extra large/left	Stainless steel 304
5316-00-000	Anterior femoral cutting guide	Stainless steel 304
5323-01-000	Femoral cutting guide – small	Stainless steel 630
5323-02-000	Femoral cutting guide – médium	Stainless steel 630
5323-03-000	Femoral cutting guide – large	Stainless steel 630
5323-04-000	Femoral cutting guide – extra large	Stainless steel 630
5342-08-000	Spacer block 8,0 mm	Aluminum
5342-10-000	Spacer block 10,0 mm	Aluminum
5342-12-000	Spacer block 12,0 mm	Aluminum
5342-15-000	Spacer block 15,0 mm	Aluminum
5315-00-000	Anterior femoral cutting guide	Stainless steel 304
5327-00-000	Tibial cutting guide	Stainless steel 304
5312-00-000	Femoral rotational alignment guide	Stainless steel 304
5363-00-000	Tibial retractor	Stainless steel 304
5362-00-000	Patellar tendon retractor	Stainless steel 304
5317-00-000	Distal femoral cutting guide	Stainless steel 304
5314-00-000	Impactor	Stainless steel 304
5359-00-000	Universal extractor	Stainless steel 304
5361-00-000	Hohmann retractor (24.5x252 mm)	Stainless steel 304
5360-00-000	Hohmann retractor (24.5x252 mm)	Stainless steel 304

5367-00-000	Screw extractor	Stainless steel 420
5228-00-000	Osteotoma 20 mm	Stainless steel 420
5307-00-000	Quick coupling	Stainless steel 304
5372-00-000	Screw pointer	Stainless steel 304
5333-00-000	tibial resection guide	Stainless steel 420
5211-00-000	Osteotomo 22 mm	Stainless steel 420
6003-02-032	Drill bit Ø3,2x150 mm with quick coupling	Stainless steel 420
5381-00-000	Patellar guide 3 grips	Stainless steel 304
5383-00-000	Femur test extractor	Stainless steel 304
5382-00-000	Paquimeter	Stainless steel 304
5381-01-000	Patellar drill Ø7,0 mm	Stainless steel 420
5384-00-000	Patelar Clamp	Stainless steel 304
5385-00-000	Extramedullary tibial cutting guide	Stainless steel 304
5386-00-000	Plateau inserter	Stainless steel 420
5313-12-000	Bar for Femoral alignment guide 120mm	Stainless steel 304
5313-15-000	Bar for Femoral alignment guide 150mm	Stainless steel 304
5313-20-000	Bar for Femoral alignment guide 200mm	Stainless steel 304
5324-10-000	Intercondylar guide P/M	Stainless steel 304
5353-28-000	Patellar trial Ø28 mm (3grips)	Polyphenylene
5353-30-000	Patellar trial Ø30 mm (3grips)	Polyphenylene
5353-32-000	Patellar trial Ø32 mm (3 grips)	Polyphenylene
5353-34-000	Patellar trial Ø34 mm (3 grips)	Polyphenylene
5353-36-000	Patellar trial Ø36 mm (3 grips)	Polyphenylene
5211-16-000	Osteotomo16mm	Stainless steel 420
5211-25-000	Osteotomo25mm	Stainless steel 420

Indications, contraindications, adverse effects and warnings

Indications: The accessories for femoral components – BM have the purpose of correcting the asymmetric segment defects of the bone as well as giving more stability to the femoral component without preservation – BM designed for revision surgeries. They are used in revision surgeries for restitution or replacement of the femorotibial and patellofemoral joints (knee) in conjunction with their ancillary components. That procedure is known as "total knee arthroplasty" and tries to create an artificial articulation with the same functionality of the natural knee. The total knee arthroplasty has the purpose of eliminating the pain symptoms, correcting the defects, allowing a functional range of motion and stabilizing the knee.

Contraindications

The contraindications below do not limit the use of Accessories for Femoral Components - BM, but must be analyzed by the surgeon:

- C. Gross weakness of the quadriceps;
- D. Systemic infections or diseases close to the joint;
- W. Rearthrodesis of the knee after arthroplasty;
- X. Previous patellectomy;

- Y. Neuropathic arthritis;
- Z. Severe discontinuity or dysfunction of the extensor mechanism;
- AA. Recurvatum deformity due to muscle weakness and presence of a knee arthrodesis;
- BB. Morbid obesity;
- CC. Fever or leukocytosis;
- DD. Inability of the patient to receive anesthesia;
- EE. Atherosclerotic disease of the leg;
- FF. Skin problems such as psoriasis;
- GG. Recurrent urinary tract infections;
- HH. History of osteomyelitis near the knee;
- II. Severe neuromuscular, muscular or vascular deficiencies in the affected region;
- JJ. Severe injury to the bone structures generating an unfeasible risk for the application of the prosthesis;
- KK. Bone tumors close to the implantation site of the component;
- LL. Muscle atrophy;
- MM. Bone immaturity;
- NN. Insufficiency of the knee extensor mechanism;
- OO. Structural soft tissue failure of the knee;
- PP. Knee ankylosis;
- QQ. Insufficient bone structure in the tibia, femur and patella;
- RR. Patients unwilling or unable to follow post-operative care and instructions.

Adverse Effects

- P. Joint line migration;
- Q. Change in gait pattern;
- R. Decrease or increase in patellar height;
- S. Hematomas;
- T. Suture dehiscence;
- U. Localized and asymptomatic progressive bone resorption (osteolysis) around the component;
- V. Loosening, migration or fracture of the medical device;
- W. Superficial or deep infection;
- X. Vascular disorders, including thrombosis and pulmonary embolism.
- Y. Shortening of limb due to bone resorption;
- Z. Allergic reaction to the foreign body, which may result in histological reactions involving various types of macrophages and fibroblasts;
- AA. Pain, discomfort or abnormal sensations due to the presence or migration of the medical product, and also due to the surgical procedure;
- BB. Neural or neurological damage due to surgical trauma (including paralysis and soft tissue injuries);
- CC. Inability to perform normal daily physical activities;
- DD. Periprosthetic fractures.

Warnings

- I. Accessories for Femoral Components - BM help in stabilizing the Femoral Component with Restriction - BM, which in turn has the purpose of replacing the damaged joint and not normal structures of the human body;
- J. Surgical implants must never be reused, and explanted implants must never be implanted again. Small defects and internal pressure patterns can occur, even if the implant appears intact, causing implant failure and rupture;
- K. The correct selection of Femoral Components Accessories - BM models is extremely important in accordance with the patient's bone characteristics, analyzed by the surgeon, and with the surgical technique adopted by the surgeon;
- L. The infections often generated are linked to comorbidity;
- M. Aseptic failure of total knee arthroplasty can be caused by loosening of components, wear of the polyethylene manufactured ancillary component, ligament laxity, periprosthetic fracture, arthrofibrosis, and patellofemoral complications;
- N. The results presented in revision surgeries are inferior compared to the results presented in primary surgeries;
- O. Patients who smoke are more likely to have longer operative times and higher costs;
- P. Knee replacement has a finite expected survival, which is adversely affected by activity level.

The patient must be informed of:

- G. The need for periodic medical follow-up, in order to observe possible alterations in the status of the implant and the adjacent bone. Only monitoring can detect possible loosening of components or the occurrence of osteolysis;
- H. Failure to perform revision surgery when loosening components or osteolysis may result in progressive loss of periprosthetic bone stock;
- I. The patient must be fully understood and warned that the product does not replace and does not have the same performance as normal bone. Trials have shown that the use of this implant allows the patient to perform normal day-to-day activities, such as sitting, walking and climbing stairs. Activities such as sitting down to tie a sneaker or shoe should be done with caution by the patient.
- J. Information related to topics: indications, contraindications, adverse effects and warnings;
- K. The patient must be instructed to inform that he has a metallic implant in case he undergoes the Magnetic Resonance exam, as the implant manufactured from metal causes interference in the exam, jeopardizing the diagnosis, and the exam may cause displacement of the implant, causing problems for the patient;
- L. The patient must be advised of the need to attend regular post-operative follow-up exams during the entire period in which the product remains implanted.

Joint Operation with Other Medical Products: According to the norm "NBR ISO 21534 - Non-Active Surgical Implants - Joint Replacement Implants - Particular Requirements", the cast alloy of cobalt-


chromium-molybdenum, raw material of Accessories for Components Femoral - BM, is acceptable for the fabrication of implants.

A list of acceptable combinations with Accessories for Femoral Components - BM was formulated, according to the norm "NBR ISO 21534 - Implants for Non-Active Surgery - Implants for Joint Replacement - Particular Requirements" (Table 4). Regardless of this relationship, the combinations must also be of implants of the same brand, as they were designed for such combinations, observing surface finish and surface treatment and other factors required in projects that may interfere with the combinations. Therefore, metallic implants from different manufacturers are not recommended for reasons of chemical, physical, biological and functional incompatibility.

Table 4: Admissible Contact Alloys with Accessories for Femoral Components - BM

Admissible contact alloy (when there is no articulation)
- Titanium alloy (NBR ISO 5832-3, NBR ISO 5832-11, ASTM F136);
- Cobalt alloy (NBR ISO 5832-4, NBR ISO 5832-5, NBR ISO 5832-6, NBR ISO 5832-7, NBR ISO 5832-8, NBR ISO 5832-12, ASTM F75)
- Stainless steel (NBR ISO 5832-9)

Special conditions for storage, conservation and/or handling of the product (described in the external protective packaging of the product)

Open the Packaging Aseptically	
 35°C °C Temperature Limits	Transport and store it in dry and cool places and at room temperature (maximum 35°C) and Relative humidity around 30% and 70%.
Do not store it in direct contact to the floor (minimum height = 20cm) and neither in very high places and close to lamps because that may lead to the dryness and damage to the label.	
Do not store it in places where contaminants substances are stored such as cleaning materials, insecticide, pesticide, etc.	
	Attention Instructions, warnings and precautions: Vide instruction for use
	Do not use it if the packaging is damaged
	Fragile, carefully handle it
	Keep it away from the Sun
	Keep it dry

	Vide Instruction For Use
	Sterile Product
	

Precautions Before Using the Product: Before using the product, the person in charge, as designated by the hospital unit, must check the physical conditions of the packaging and subsequently of the product. The person in charge must verify that the outer cardboard packaging is sealed with a security seal and identified with the information regarding the material present in the packaging. After this check, remove the security seal and open the cardboard package.

The product packaged in a double blister should be carefully removed from the outer packaging, in a surgical environment. The first blister should be opened by pulling the surgical paper on its edge at the edge of the blister. The second blister should also be opened by pulling the remaining surgical paper on the edge of the blister. With the product in hand, check that it does not have dents, scratches or stains.

Load Bearable by the System Associated with the Weight of the Patient

Trials have shown that the Tibial Components - BM for revision can be used without restrictions in patients weighing up to 106 kg and in cases of patients weighing more than 80 kg, activities such as climbing stairs are not recommended.

After Sales Service (Customer Complaint): If there is a need to make a complaint about the Femoral Metal Heads - Biomechanics related to any adverse effect that affects the user's safety, such as product not working, damage to the implantable metal component, serious problems or death related to responsible for these components, the surgeon in charge must communicate this adverse event to the competent health agency and Biomecânica via email sac@biomecanica.com.br or by telephone 0xx14 2104 7926. In case of doubt, the surgeon in charge or the health professional can the communication of the adverse event through the Sanitary Surveillance Notification System on the ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Product Disposal: No surgical implant should be reused. Any metallic implant, after use, must be discarded. Even if it appears to be undamaged, it may have small defects that can lead to fatigue fractures in addition to contamination. Thus, the Accessories for Femoral Components - BM explanted from patients must be properly discarded by the hospital institution. The hospital institution is responsible for the complete mischaracterization of the implant, preventing its reuse. The hospital institution is responsible for the method used to mischaracterize the implant. Biomecânica

recommends that the explanted implants be mechanically deformed with the aid of a hammer or impact press, and must then be identified with the saying "Improper for Use". According to Resolution RE nº 2605, of 08/11/06, implantable devices of any nature classified as single use are prohibited from being reprocessed.

Traceability: Along with the product, there are 6 labels containing the necessary data for product traceability. A label must be obligatorily affixed to the patient's clinical record, another to the report delivered to the patient, another to the tax documentation that generates the charge, on the AIH, in the case of a patient treated by the SUS, or on the sales invoice, in the case of patient assisted by the complementary health system, another available for the control of the supplier (historical distribution record - RHD, another available for the control of the responsible surgeon (main) and the last one for medical insurance, if any.

Sending Material to the Manufacturer for Analysis: If implants are sent to the manufacturer for analysis, they must be cleaned in the hospital using a broad-spectrum bactericidal and antifungal solution. Then it must be disinfected or steam sterilized in an autoclave or ethylene oxide. They must be sent to Biomecanica in complete packaging, identified with the method of cleaning, sterilization and product data.

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