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ENGLISH EUA

Detailed description of the medical device including the basis of its functioning and action, content or composition, whenever applicable as well as the list of accessories designed to be used in conjunction with the product.

The Large and Small fragments screws, cannulated – BM are implants to be fixed to the bone with the purpose of providing fixation and stability. They are manufactured in stainless steel according to ASTM F138 standard. The cannulated screws are partially threaded and have an inner cannulated hole, cylindrical head with an inner hexagon that makes it easier for their insertion into the bone with the help of a hexagonal screwdriver and the guide wire. The cannulated screws are generally intended to cancellous bones, in other words, they are used in cases that the fixation happens predominantly in less dense inner regions (cancellous bone). The threads for cancellous screws are longer and deeper than the one of the cortical screws due to the characteristics of the cancellous bone. The advantage of the use of the cannulated screws is having a precise insertion of the screws into the cancellous bones in metaphyseal and epiphyseal areas through the guide wire. They are available in different sizes and are not reusable. The characteristics of the Large and small fragments screws, cannulated – BM are presented below:

Table 1 – Screw model



Model	Fitting Types
	
Partial Thread	Internal Hexagonal

Table 2 - Screws for large and small fragments, cannulated - BM

Product	Product
Ø4.5mm Cannulated Cancellous Screw – ST Tip 	
Produto Ø 6.5mm Cannulated Cancellous Screw - 16 mm Thread - ST Tip 	Produto Ø 6.5mm Cannulated Cancellous Screw –32 mm thread - ST Tip 
Produto Ø7.0mm Cannulated Cancellous Screw - 16 mm thread - ST Tip 	Produto Ø7.0mm Cannulated Cancellous Screw - 32 mm thread - ST Tip. 

Composition


The Large and Small Fragments screws, cannulated – BM are manufactured in stainless steel according to ASTM F138 standard. The stainless steel used for the manufacturing of the screws is analyzed in outsourced laboratories for the verification of its compliance with ASTM F138 standard.

Presentation of the Medical Product

The implants are packed in plastic envelopes that come with a tag of non-sterile product and outer labeling. The product comes with the following information: Manufacturer's name; Commercial Name of the product; Technical Name; Lot Number; Product registration Number with the Ministry of Health; Product Code; Quantity; Description of the product available in the packaging (model); Manufacturing date; Expiration Date; Undetermined expiration; Description of the raw-material used for the manufacturing of the product; Manufacturer's address (available on the plastic envelop); Name of the Responsible Technician; The sayings: “Non-sterile Product”, “Single Use Product, Do not reuse”; Read the instructions for use before use.

Symbology of the Traceability Tag (Label) according to NBRISO15223

Table 3 – Symbology of the Traceability Tag (label)

 XXXX	Date of Manufacture	 XXXX	Use by
 XXXXX	Batch Code		Manufacturer
	Non-Sterile	 XXXX-XX	Catalogue Number

The non-sterile product is identified with a green label that says “Non-Sterile”

NON STERILE

➤ **Symbology of the Outer Packaging:**

	Fragil, Handle with Care
	Keep away from the sunlight
	Keep dry

Screws for large and small fragments, cannulated - BM will be distributed in a unitary form according to the table below:

Table 4 - Distribution of Screws for large and small fragments, cannulated - BM

Código	Qtd.	Descrição
3019-45-020	1	Cannulated cancellous screw Ø 4.5 mm x 20 mm length with 7.0 mm of length
3019-45-022	1	Cannulated cancellous screw Ø 4.5 mm x 22 mm length with 7.0 mm of length
3019-45-024	1	Cannulated cancellous screw Ø 4.5 mm x 24 mm length with 8.0 mm of length
3019-45-026	1	Cannulated cancellous screw Ø 4.5 mm x 26 mm length with 9.0 mm of length
3019-45-028	1	Cannulated cancellous screw Ø 4.5 mm x 28 mm length with 9.0 mm of length
3019-45-030	1	Cannulated cancellous screw Ø 4.5 mm x 30 mm length with 10.0 mm of length
3019-45-032	1	Cannulated cancellous screw Ø 4.5 mm x 32 mm length with 11.0 mm of length
3019-45-034	1	Cannulated cancellous screw Ø 4.5 mm x 34 mm length with 11.0 mm of length
3019-45-036	1	Cannulated cancellous screw Ø 4.5 mm x 36 mm length with 12.0 mm of length
3019-45-038	1	Cannulated cancellous screw Ø 4.5 mm x 38 mm length with 13.0 mm of length
3019-45-040	1	Cannulated cancellous screw Ø 4.5 mm x 40 mm length with 13.0 mm of length
3019-45-042	1	Cannulated cancellous screw Ø 4.5 mm x 42 mm length with 14.0 mm of length
3019-45-044	1	Cannulated cancellous screw Ø 4.5 mm x 44 mm length with 15.0 mm of length
3019-45-046	1	Cannulated cancellous screw Ø 4.5 mm x 46 mm length with 15.0 mm of length
3019-45-048	1	Cannulated cancellous screw Ø 4.5 mm x 48 mm length with 16.0 mm of length
3019-45-050	1	Cannulated cancellous screw Ø 4.5 mm x 50 mm length with 16.0 mm of length
3019-45-052	1	Cannulated cancellous screw Ø 4.5 mm x 52 mm length with 17.0 mm of length
3019-45-054	1	Cannulated cancellous screw Ø 4.5 mm x 54 mm length with 18.0 mm of length
3019-45-056	1	Cannulated cancellous screw Ø 4.5 mm x 56 mm length with 19.0 mm of length
3019-45-058	1	Cannulated cancellous screw Ø 4.5 mm x 58 mm length with 20.0 mm of length
3019-45-060	1	Cannulated cancellous screw Ø 4.5 mm x 60 mm length with 21.0 mm of length
3019-45-064	1	Cannulated cancellous screw Ø 4.5 mm x 64 mm length with 21.0 mm of length
3019-45-068	1	Cannulated cancellous screw Ø 4.5 mm x 68 mm length with 23.0 mm of length
3019-45-070	1	Cannulated cancellous screw Ø 4.5 mm x 70 mm length with 24.0 mm of length
3023-16-030	1	Cannulated cancellous screw Ø 6.5 mm x 30 mm length with 16 mm of length
3023-16-035	1	Cannulated cancellous screw Ø 6.5 mm x 35 mm length with 16 mm of length
3023-16-040	1	Cannulated cancellous screw Ø 6.5 mm x 40 mm length with 16 mm of length
3023-16-045	1	Cannulated cancellous screw Ø 6.5 mm x 45 mm length with 16mm of length
3023-16-050	1	Cannulated cancellous screw Ø 6.5 mm x 50 mm length with 16 mm of length
3023-16-055	1	Cannulated cancellous screw Ø 6.5 mm x 55 mm length with 16 mm of length
3023-16-060	1	Cannulated cancellous screw Ø 6.5 mm x 60 mm length with 16 mm of length
3023-16-065	1	Cannulated cancellous screw Ø 6.5 mm x 65 mm length with 16 mm of length
3023-16-070	1	Cannulated cancellous screw Ø 6.5 mm x 70 mm length with 16 mm of length
3023-16-075	1	Cannulated cancellous screw Ø 6.5 mm x 75 mm length with 16 mm of length
3023-16-080	1	Cannulated cancellous screw Ø 6.5 mm x 80 mm length with 16 mm of length
3023-16-085	1	Cannulated cancellous screw Ø 6.5 mm x 85 mm length with 16 mm of length
3023-16-090	1	Cannulated cancellous screw Ø 6.5 mm x 90 mm length with 16 mm of length
3023-16-095	1	Cannulated cancellous screw Ø 6.5 mm x 95 mm length with 16 mm of length
3023-16-100	1	Cannulated cancellous screw Ø 6.5 mm x 100 mm length with 16 mm of length
3023-16-105	1	Cannulated cancellous screw Ø 6.5 mm x 105 mm length with 16 mm of length
3023-16-110	1	Cannulated cancellous screw Ø 6.5 mm x 110 mm length with 16 mm of length
3023-16-115	1	Cannulated cancellous screw Ø 6.5 mm x 115 mm length with 16 mm of length

3016-32-060	1	Cannulated cancellous screw Ø 7.0 mm x 60 mm length with 32 mm of length
3016-32-065	1	Cannulated cancellous screw Ø 7.0 mm x 65 mm length with 32 mm of length
3016-32-070	1	Cannulated cancellous screw Ø 7.0 mm x 70 mm length with 32 mm of length
3016-32-075	1	Cannulated cancellous screw Ø 7.0 mm x 75 mm length with 32 mm of length
3016-32-080	1	Cannulated cancellous screw Ø 7.0 mm x 80 mm length with 32 mm of length
3016-32-085	1	Cannulated cancellous screw Ø 7.0 mm x 85 mm length with 32 mm of length
3016-32-090	1	Cannulated cancellous screw Ø 7.0 mm x 90 mm length with 32 mm of length
3016-32-095	1	Cannulated cancellous screw Ø 7.0 mm x 95 mm length with 32 mm of length
3016-32-100	1	Cannulated cancellous screw Ø 7.0 mm x 100 mm length with 32 mm of length
3016-32-105	1	Cannulated cancellous screw Ø 7.0 mm x 105 mm length with 32 mm of length
3016-32-110	1	Cannulated cancellous screw Ø 7.0 mm x 110 mm length with 32 mm of length
3016-32-115	1	Cannulated cancellous screw Ø 7.0 mm x 115 mm length with 32 mm of length
3016-32-120	1	Cannulated cancellous screw Ø 7.0 mm x 120 mm length with 32 mm of length
3016-32-125	1	Cannulated cancellous screw Ø 7.0 mm x 125 mm length with 32 mm of length
3016-32-130	1	Cannulated cancellous screw Ø 7.0 mm x 130 mm length with 32 mm of length
3005-13-000	1	Washer Ø 13.0 mm stainless steel
3005-07-000	1	Washer Ø 7.0 mm stainless steel

Attention

For the implantation of the Large and Small fragments screws, cannulated – BM is necessary the use of the instruments presented in table below. Those instruments shall be acquired separately.

Table 5 - List of Instruments for placing Screws for large and small fragments, cannulated - BM

Qtd	Código	PRODUTO
02	6003-25-150	Drill Ø 2.0 x 150 mm
01	6086-00-000	Drill Guide
01	6059-01-000	Calibrated Wire Measurer
01	6082-00-000	Cannulated Reamer
01	6083-00-000	Cancellous Tap
01	6201-01-000	Cannulated Hexagonal Screwdriver
01	6016-00-000	Countersink
01	6087-00-000	Multiple Guide
01	6085-00-000	Soft Tissue Protector
02	6081-00-000	Guide Wire

Cares with the Storage, Conservation, Handling and Transport of the Medical Product

- Do not reuse the product in case it is damaged; - The implants shall be carefully handled in order to avoid damages that may impair the quality of the material and the patient's safety; - Keep it in a ventilated and dry place, away from the light and weather conditions; - During transportation it shall be avoided shocks and improper stacking; - Store and transport it in a dry and ventilated environment at room temperature (Max. 35°C) and around 30% of relative humidity; - Do not store it directly on the floor (minimum height = 20cm) nor in very high places close to lamps because that may resect the packaging or damage the label; - Do not store it with contaminants such as cleaning materials, insecticide, pesticide, etc.

Indication, purpose or use to which the medical product is intended based on the manufacturer's indication

Indication for Use

The large and small fragments screws, cannulated – BM are implants can be fixed with the help of a washer or not.

The use of cannulated screws is indicated in the following situations:

- When it requires a precise implantation of the screws into the cancellous bone in metaphyseal and epiphyseal areas. In that case a guide wire can be used for a higher precision implantation of the screw;
- For patients with dome high tibial valgus osteotomy;
- For femoral proximal epiphysiolisthesis treated with “in situ” fixation with only one cannulated screw because the results are the safest and the most effectiveness ones and that method provides a good fixation and stability;
- Fracture stabilization, osteosynthesis in general for the correction of limbs, fixation of ligaments, reconstruction, neutralization and treatment of pseudarthrosis;
- For surgeries with the need of bone fixation such as osteosynthesis fractures or osteotomies and for the treatment of ephysiolisis. Besides the screws are also indicated for the reduction and stabilization by compression and traction of the fractures;
- For the fixation of the slipped capital femoral epiphysis, for compression fixation in the metaphyseal areas of large bones such as the femoral neck, femoral condyles and tibial plateau and other areas similar to those ones;
- For high tibia osteotomy, stabilization of the epiphyseal slip and achievement of the epiphyseal plate closure;
- For upper limbs and dislocated fractures of the great tubercle of the humerus;
- For the ends of long bones of the lower and upper limbs such as the metaphysis of long bones, the distal radio, the distal humerus and the distal and proximal tibia.

Instructions for Use

This implant is provided non-sterile. It

shall be sterilized before use according to the instructions available in these instructions for use.

Exclusively handle the large and small fragments screws, cannulated – BM in proper environments and with the necessary care. Only qualified professionals shall handle and implant the screws. The surgical technique varies according to the choice of the surgeon who is also in charge of choosing the final method, type and dimension of the products to be implanted as well as the criteria for evaluation of the surgical results. The large and small fragments screws, cannulated – BM shall be implanted and adopted based on the surgical needs and techniques. The mechanical contribution of the implant shall be analyzed from the observance of a precise adaptation and compression, because a stable reconstruction of the fractured bone minimizes the load the implant will support, The torque to be applied during the insertion will depend on the

conditions and characteristics of the bone. The doctor shall decide the torque to be applied to. Based on the tests, the surgeon shall consider the implantation levels. The patient's weight, activity level and conditions may impact on the performance of the implant. It shall be considered that no implant is as strong as a natural bone and, therefore, it has limitations regarding the biomechanical stresses.

Disclaimer

Biomechanica, as the manufacturer of this device, does not practice or recommend this or other surgical technique to be used in the treatment of a certain patient. The surgeon that follows any implantable procedure is responsible for defining and using the appropriate techniques for implantation of a medical device in each patient. Biomechanica is not responsible for the choice of the appropriate surgical technique to be used in the treatment of each patient. The surgical techniques depend on the technique adopted by the surgeon who is also in charge of choosing the method, type and dimension of the products to be implanted.

Ancillary Components

The Large and small fragments screws, cannulated – BM do not use any ancillary components for their implantation.

Accessories (Integral Part) – Washer

The Large and small fragments screws, cannulated – BM have the following products as their accessories:

Product	Code	Diameter	Used with	Photo
Washer	3005-13-000	Ø13.0mm (Stainless Steel)	Cannulated Cancellous Screw Ø 6.5mm and Ø 7.0mm	
Washer	3005-07-000	Ø7.0mm (Stainless Steel)	Cannulated Cancellous Screw Ø 4.5mm	

Admissible combinations with other materials

The table “Admissible contacting alloys” is based on “NBRISO21534 standard – Implantes para Cirurgia não Ativos – Implantes para Substituição de Articulações – Requisitos Particulares”. However, independently of the alloys mentioned in item “Contact”, it is necessary that the combinations are between same brand implants because the implants were designed for such combinations based on the surface finishing and treatment and other factors required in the design that may interfere in the combinations. Therefore, the use of metallic implants from different manufacturers is not recommended due to chemical, physical, biological and functional incompatibility.

Table 6 – Admissible Contacting Alloysimplant

Admissible Contacting Alloysimplant
--

Implant	*Contact (whenever there is articulation)
Cannulated screw in Stainless Steel ASTM F138	Wrought 18Chromium-14Nickel-2.5 Molybdenum Stainless Steel (NBRISO 5832-1 and ASTM F138) Wrought high-nitrogen stainless steel (NBRISO5832-9) Wrought Titanium-6Aluminum-4Vanadium (NBRISO 5832-3 and ASTM F136) Wrought titanium 6-aluminum 7-Niobium alloy (NBRISO5832-11)

*For applications which the metal or alloy is in contact with other metal or alloy and articulation is not intended, on the condition that adequate attention is given to the design, surface finishing and treatment and metallurgical conditions.

Biomechanica Large and small fragments screws, cannulated – BM shall be only used with Biomechanica instruments.

Contraindications

The use of this implant in the following situations is contraindicated: pregnancy, osteoporosis; obesity; visible infection; poor tissue coverage; people with mental illness; patients that do not follow the postoperative instructions; foreign body sensitivity such as allergy; any health condition equivalent to a disease that limits the surgery; infection in the region to be operated or other infection histories; bone impairment due to illness and/or infections.

Precautions, restrictions, warnings, special cares and clarification on the use of the medical product as well as its storage and transport.

Warnings and Precautions

SINGLE USE PRODUCT. DESTROY IT AFTER EXPLANTED. DO NOT REUSE THE PRODUCT. DO NOT REPROCESS

An orthopedic implant may only be used in a single patient for a single time. Although it does not look damaged, previous stresses may have created imperfections that may impair the success of the implant. A poor selection of the implant may cause uncommon stresses to the implant and result in its fracture and subsequent breakage,

Due to its presentation and characteristics, this product does not have any side effects, however, it is suggested monitoring the patient's recovery through x-ray based on the International Safety Standards. Do not use magnetic resonance because the product in stainless steel causes interference in the exam and impairs the diagnosis, besides, it may also cause the looseness of the implant and bring problems to the patient. The use of Large and small fragments screws, cannulated – BM are designed to work as an auxiliary for the bone fixation and stabilization in the treatment of fractures of the upper and lower limbs and not to replace the healthy structures of the human skeleton.

The surgeon shall be aware of the chance of the patient's physical development, therefore it may be necessary the implant revision. The surgeon shall be familiar with and have enough knowledge of osteosynthesis and limitations, as well as the pre and postoperative, adopted surgical techniques, precautions and potential risks. The surgeon shall be careful with the bone

preparation for the insertion of the Large and small fragments screws, cannulated – BM, because the stability and success of the fixation may be impaired.

During the handling of the implant, scratches, grooves or any other thing that may damage or mark the implant shall be avoided because those defects are stress concentrations and may be sites for crack initiation and decrease the resistance to corrosion, what may lead to the breakage of the implant.

In order to guarantee a proper implantation, Biomechanica instruments shall be used because they were specifically designed and manufactured for the use of Biomechanica implants. The use of instruments from different manufacturers may impair the surgery, besides it is in disagreement with the product registration with the National Health Surveillance Agency (ANVISA).

Do not use the product in case it is damaged.

NON-STERILE PRODUCT – STERILIZE BEFORE USE.

SINGLE USE PRODUCT – THE REUSE OF THIS PRODUCT IS EXTREMELY PROHIBITED.

Special storage conditions: Keep it in ventilated and dry places, away from the light and weather conditions. Do not use the product in case the package is damaged.

Note: Implantable components from different manufacturers shall not be used. Therefore, it is recommended that the products have the same provenance. Manufacturing date, expiration date and product batch: VIDE LABEL.

Product Conformity

During the handling of the implant, it shall be avoided scratches, grooves or any other thing that may damage or mark the implant, because those defects are stress concentrations and may be sites for crack initiation and decrease the resistance to corrosion what may result in the breakage of the implant.

Adverse Events

The mechanical looseness may be the result of a poor fixation, unstable reconstruction of the fractured bone or occult infection. Sensitivity reaction to the metal in patients has rarely been reported. The implantation of a foreign material in tissues results in histological reactions that involve different sizes of macrophages and fibroblasts. The Chemical importance of that effect is uncertain. Similar changes may occur during or as a precursor to the healing process.

The poor implantation of the screws resulted from technical mistake may bring out many complications to the patient. The treatment also brings three risks: osteoarthritis, avascular necrosis and chondrolysis. The incorrect use of the crutches in the postoperative may cause complications with the implant. Neurological complications, bleeding, reactions to foreign bodies, bone fractures, postoperative pain and pain followed by abnormal sensations are also related to that.

Information to be provided to the Patient

The sayings of these instructions in items: Indications, Contraindications, Information for use, Potential Adverse Events, Precautions and Warnings. The patient shall be aware of the importance of a postoperative follow-up. The failure to have a follow-up prevents the detection of postoperative problems. The failure to undergo a revision surgery in case of the component loosening or osteolysis may result in a progressive loss of the periprosthetic bone stock. The patient shall know the type of used implant, the brand, code, name, batch, surgeon's name and the date of the surgery.

Important

The surgeon shall help the patient with the cares to be taken and the adequate physical support needed for the time of implant use. Before a perfectly stabilized fracture the load may be liberated. If the fracture is not stabilized, no load may be applied to the implanted system.

Evaluations of the Implanted Product

After implantation, in the intraoperative period, the responsible professional shall perform the radiological control to check the positioning of the implant. The responsible professional shall perform the clinical and radiological evaluations after the surgical procedures at a frequency previously stipulated in order to check the conditions of the implant and the evolution of the bone healing. In case the product is out of the correct positioning or presents any non-conformity, the surgeon is responsible for taking the most appropriate corrective actions.

Useful Information to Avoid Risks resulted from Implantation

In order to minimize the risks resulted from implantation the following shall be strictly followed: contraindications, instructions for use and all information available in the "Instructions for Use" of the product.

Intrinsic Implantation Risk

The Large and small fragments screws, cannulated – BM are manufactured with recognized materials for medical use that follow the requirements of ASTM F138 -00- Standard Specification For Wrought 18 Chromium-14 Nickel-2,5 Molybdenum Stainless Steel Bar And Wire For Surgical Implants (Uns S31673).

Contamination Risks

There are risks of biological contamination and transmission of viral diseases such as Hepatitis and HIV because the metallic components of the Large and small fragments screws, cannulated – BM get in contact with body tissues and fluids. The explanted products shall be treated as highly polluting products.

Decontamination, Cleaning and Sterilization

All models of the Large and small fragments screws, cannulated – BM are distributed non-sterile. Therefore, the institution that does the operation will be responsible for the sterilization and shall use regularly validated methods for that. The product is provided in its original packaging and it is previously cleaned. The packaging shall be checked for ruptures or damages before use. The product shall only be unpacked before sterilization in order to preserve intact the surface finishing and characteristics and it shall be handled the least as possible according to NBR ISO 8828 standard. In order to guarantee the minimization of the risks of infection, the implants (the ones that were not used in the surgery but were set on the surgical table or that have lost their sterility) and surgical instruments shall be cleaned before the sterilization.

The cleaning may be done through manual or mechanical methods.

Cleaning Equipments:

- Ultrasonic washer
- Washer – sterilizer
- Washer – Disinfectant

Manual Cleaning: it is a procedure that removes the dirtiness through physical action with the help of detergent, water and tools such as sponge and brush. The manual cleaning is mostly recommended due to it is less aggressive to the implants and surgical instruments. For the manual cleaning it shall be used water at room temperature, nylon brushes and never steel or abrasive ones because they may damage the protective layer of the material. For the decontamination it shall be used bactericidal and broad spectrum antifungal solutions and avoid aggressive cleaning agents and metal brushes so that the quality of the metallic implants is not impaired.

In case of cleaning washers, the surgical implants and instruments shall be set in drawers in order to avoid shocks and prevent damages to the material. People responsible for the cleaning of the implants and surgical instruments shall be aware of the type of cleaning product being used which shall obey the concentrations determined by the manufacturer. The implants and surgical instruments shall be carefully rinsed up to all residue is out and immediately dried after the cleaning.

✓ Type of Cleaners

Enzymatic Cleaners: they are basically composed of surfactant soluble enzymes. The balanced combination of those elements makes it possible for the product to remove the organic matter from the material in a short period of time.

Enzymatic Solutions: they present an excellent cleaning action; however, they do not present bactericidal and bacteriostatic activity.

Enzymes: they are substances produced by live cells that govern the chemical reaction of the process. Once an enzyme is produced by the cells, the enzyme may be isolated and it will still keep its catalytic properties if certain conditions are kept in its manufacturing. The enzymes are classified in three major functional groups depending on the type of substrate they will affect: proteases, lipases and amylases that act in protein substrates, fats and carbohydrates which tend

to solubilize and detach from the articles. Currently it is recommended the cleaning of complex characteristic articles for the cleaning guarantee.

✓ **Desinfection:**

It is a process that destroys microorganisms (pathogenic or not) of the articles except for the bactericidal spores by physical and chemical means.

✓ **Disinfection Levels:**

High Level: it destroys all microorganisms except for the high number of spores => Glutaraldehyde 2% - 20 – 30 minutes.

Indication: preferably hospital areas

Medium Level: it eliminated vegetative bacteria, most of the virus, fungi and microbacteria => Sodium hypochloride 1% - 30 minutes. Indication: UBS, Day cares, nursing home and almshouses.

Low Level: it eliminates most of the bacteria, some virus and fungi. However, it does not eliminate the microbacteria => Sodium Hypochloride 0,025%

Indication: Food

Sterilization

This product is provided non-sterile. It shall be sterilized before use. We recommend the moist heat sterilization in Autoclave at hospital (ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices).

AUTOCLAVE

It is an equipment for sterilization that uses saturated steam under pressure.

OPERATING INSTRUCTIONS

Before starting the equipment, check: if the switch is on; if the water valve is open; if the flush valve is closed.

- 1 – open the door of the equipment;
- 2 – properly set the material to be sterilized;
- 3 – close the door of the equipment;
- 4 – select the desired cycle according to the material to be sterilized;
- 5 – turn on the general key;
- 6 – the cycle will automatically run in the sequence;
- 7 – when the lamp “final of the cycle” turns on, partially open the door for about 10 (ten) minutes to cool down the material.

It is recommended the following physical parameters of sterilization in autoclave (saturated steam):

Table 7 - Physical parameters of sterilization in autoclave

Cyclo	Temperatura	Time of Exposure (minimum)	Time of Drying (maximum)
Conventional (1atm of pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Note: the time shall be recorded from the moment that the heat of the sterilization chamber reaches the desired temperature.

FOR FURTHER CLARIFICATION CHECK THE INSTRUCTIONS THAT FOLLOWS EACH AUTOCLAVE.

Other sterilization method: besides the autoclave the following method can be used as defined by the hospital institution:

- **Ethylene Oxide Sterilization (E.T.O.)** - parameters and procedures established in the validation protocol and by ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE

The sterilization method, the equipments, the controls and the instructions for sterilization used are under responsibility of the hospital institution.

Cares with the sterilized articles:

- Storage conditions for the sterilized articles:
- regarding the environment: it shall be clean; ventilated; dry; it shall be restricted to people of the sector.

Regarding the article: - after the sterilization process, do not set the device on cold surfaces (stainless steel or stone). Use baskets or perforated recipients up to the device cools down;

- The wrapping (raw cotton fabric, nonwoven fabric, surgical paper, crepe paper, paper with film, tyvek or perforated metallic boxes) shall remain intact and be little handled to avoid that the packs tear or lose their seal;
- It shall be stored in closed cabinets with shelves;
- The shelves shall be identified to facilitate the removal of the material;
- The material shall be stored considering the expiration date of its sterilization to facilitate the distribution and prevent having expired material in stock;
- Store the sterile materials separately from the non-sterile ones in order to minimize the level of external contaminants.

Product Discard

The Large and small fragments screws, cannulated – BM that are removed from the patients shall be properly discarded by the hospital institution. It is under responsibility of the hospital institution the complete mischaracterization of the implant in order to prevent the reuse and the method used for it. Biomechanica recommends that the explanted products are mechanically deformed with the help of a hammer or press and then it shall be identified it by the sayings “improper for

Use". According to resolution nº 2605, from 11/08/2006, implantable devices of any nature classified as single use products can not be reprocessed.

Traceability

Together with the implantable component there are 2 traceability tags with the data of the implant. Those tags are equal to the ones used for the outer labeling of the product. The information available in that traceability tag is available in the item "Presentation" of these instructions for use. One of those tags shall be pasted on the patient's report and the other one is for the control of the hospital. The hospital shall register and keep all information in its files and forward it to the patient. We inform that mandatorily the patient's chart shall have recorded all information that makes it possible for the traceability of the implanted product. The tag on the patient's chart allows the traceability of the used product. Among the main information, it is indispensable the name of the manufacturer of the implant, the name of the implant, the code of the implant, the batch number and the number of product registration with ANVISA. That information is available in the traceability tags that come with the product and in the outer label. Other information shall also be considered indispensable such as the date of the surgery, the name of the patient that had the implant implanted, the surgeon's name, the patient's weight and age and other information required by the patient's chart shall also be provided.

MARKING

For the large and small fragments screws, cannulated – BM, in case there is enough space, the following information is laser engraved in their body to allow the traceability and the patient's follow-up in the postoperative: Biomecanica logo; Manufacturing batch number; acronym of the raw-material used for the manufacturing of the product (S1); Code; Dimension.

If the space in the body of the screw is not enough, only the following information is laser engraved to allow the traceability and patient's follow-up in the postoperative: -Biomecanica Logo ; - Manufacturing Batch number; -Acronym of the material used for the manufacturing of the product (S1); - Dimension.

Post sales (Client's complaint)

In case it is necessary to make a complaint about the Large and small fragments screws, cannulated – BM related to some adverse event that affects the user's safety such as if the product does not work, in case of damage to the implantable metallic component, serious problems or death related to those adverse events, the responsible surgeon shall inform that adverse event to the authorized health body and Biomecanica by the e-mail sac@biomecanica.com.br or by the phone 0xx14 2104 7900. In case of questions the responsible surgeon or health Professional may inform the adverse event through the Health Surveillance Notification System on the website of ANVISA:

<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Send the Material to the Manufacturer for Analysis

In case the implant is sent to the manufacturer for analysis, it shall be sanitized at hospital with a bactericidal and broad spectrum antifungal solution. Then, it shall be disinfected or steam sterilized in autoclave or by ethylene oxide. It shall be sent to Biomechanica in intact packs and identified by the cleaning and sterilization method and the data of the product.

Performance provided in the General Requirements of ANVISA Regulation that rules on the Essential Requirements for Safety and Effectiveness applicable to the products according to Resolution - RDC 56, from 06/04/2001:

All implant is subjected to breakage, wear or looseness and its dimension is designed to be used for a certain time. The expected time for a traumatology surgery is up to the healing of the fracture. The use of the implant is related to the stress caused by regular activities such as walking, standing up or climbing the stairs. For us to minimize the exposed risks, the Large and small fragments screws, cannulated – BM are manufactured with high stability materials such as ASTM F138, besides following the Standards NBRISO14630 and NBRISO21534. The Large and small fragments screws, cannulated – BM were tested to find their resistance limits thereby providing torque and pullout values.

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