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
ENGLISH EUA

A detailed description of the medicinal product, including the basis for its operation and action, its content or composition, if any, and the list of accessories intended to integrate the device.

The family of non-cemented (metallic) acetabular capsules - BM is made of stainless steel according to ASTM F138. The inner part of the models that make up the Family of Non-Cemented Acetabular Capsules (metallic) - BM has a design to allow adjustment by means of an interference mechanism with the Unrestricted Polymeric Acetabular Core (not Cemented) - BM Registration in ANVISA no. 80128580099. Acetabular cores are available in the following models:

Acetabular Model Component Biomec III

The metal acetabular component made of stainless steel according to ASTM F138 standard with hemispherical shape and holes for insertion of the acetabular screw, are available with coating of:

- Titanium is bonded via Plasma Spray according to ASTM F 1580 specifications, which interacts with

adjacent bone, aiding in bone bonding. The porous surface of the plasma contains microscopic protrusion roughness, which penetrates the bone when the implant is placed. This roughness ensures the anchoring of the implant at the time of implantation, helping in the initial fixation of the implant. And the roughness of ssas (small bumps) during the healing period causes bone cells to continue to adhere and grow in these roughnesses, contributing to the stability and osseointegration of the implant.

The metal acetabular component is available in sizes Ø46, Ø48, Ø50, Ø52, Ø54, Ø56, Ø58, Ø60, Ø62, Ø64, Ø66 and Ø 68mm. This component internally has an appropriate space for pressure adjustment in the Acetabular Insert made of polyethylene not subject to this registration (registration in ANVISA no. 80128580099). Equatorial reinforcement increases this adjustment by snap adjustment, maximizing fixation. The acetabular component has 3 holes that allow the use of metal acetabular screws, manufactured according to ASTM F138 for additional fixing. These screws are accessories of the Family of Non-Cemented Acetabular Capsules (metallic) - BM (regulation nº 01 (03/02/2009 - item 1.7.1 of part 1 - hip joint)). These accessories have an asymmetrical thread, spherical head of ø8.0mm with low profile, hexagonal of 3.5mm to be used with the hexagonal key in its insertion or removal. It has a thread of Ø6.0mm, varying only in its length, with a trochanter tip that helps and facilitates the drilling and introduction of it.

The family of non-cemented (metallic) acetabular capsules - BM are available in the following models and sizes, as specified in the following tables:

Table - Model Acetabular Component Biomec III Gold Std

Gold Std	Code	Diameter	Raw material	Finishing	Sterilization
	2537-46-001	46 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-48-001	48 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-50-001	50 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-52-001	52 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-54-001	54 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-56-001	56 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-58-001	58 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-60-001	60 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-62-001	62 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-64-001	64 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-66-001	66 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
	2537-68-001	68 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
Gold Std	Code	Diameter	Raw material	Finishing	Sterilization
	2637-46-001	46 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-48-001	48 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-50-001	50 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-52-001	52 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-54-001	54 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-56-001	56 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-58-001	58 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-60-001	60 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-62-001	62 mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-64-001	64mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
	2637-66-001	66mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-68-001	68mm	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation	

Composition

Family of non-cemented acetabular capsules (metallic) - BM are manufactured in stainless steel according to ASTM F138 coated by titanium alloy by plasma spray according to ASTM F 1580 specifications.

List of instruments accompanying the family of non-cemented (metallic) acetabular capsules - BM

For the use of the Family of Non-Cemented Acetabular Capsules (metallic) - BM, it is necessary to use the Instrument Kit for the Placement of the Prosthesis for Hip Arthroplasty as specified in the table below. The Instrument Kits for the Placement of the Prosthesis for Hip Arthroplasty are not subject to this registration and are not part of this product and must be purchased separately.

Table - Instrumental Kit Relationship

Qty	Code	Product
01	5136-00-000	Acetabular screw clamp
01	2141-00-000	Cortical acetabular punch - uncemented shell
01	5609-42-000	Acetabular Scraper Component Ø42 - Press-Fit
01	5609-44-000	Acetabular Scraper Component Ø44 - Press-Fit
01	5609-46-000	Acetabular scraper component Ø46 - Press-Fit
01	5609-48-000	Acetabular scraper component Ø48 - Press-Fit
01	5609-50-000	Acetabular Scraper Component Ø50 - Press-Fit
01	5609-52-000	Acetabular scraper component Ø52 - Press-Fit
01	5609-54-000	Acetabular scraper component Ø54 - Press-Fit
01	5609-56-000	Acetabular scraper component Ø56 - Press-Fit
01	5609-58-000	Acetabular scraper component Ø58 - Press-Fit
01	5609-60-000	Acetabular scraper component Ø60 - Press-Fit
01	5609-62-000	Acetabular scraper component Ø62 - Press-Fit
01	5606-46-000	Acetabular component meter Ø46 - Press-Fit
01	5606-48-000	Acetabular component meter Ø48 - Press-Fit
01	5606-50-000	Ø50 Acetabular Component Meter - Press-Fit
01	5606-52-000	Acetabular component meter Ø52 - Press-Fit
01	5606-54-000	Ø54 Acetabular Component Meter - Press-Fit
01	5606-56-000	Ø56 Acetabular Component Meter - Press-Fit
01	5606-58-000	Acetabular component meter Ø58 - Press-Fit
01	5606-60-000	Ø60 Acetabular Component Meter - Press-Fit
01	5606-62-000	Ø62 Acetabular Component Meter - Press-Fit
01	5600-00-000	Flexible Drill Guide - Biomec III
01	5601-00-000	Depth gauge - Biomec III
01	5610-00-000	Flexible drill - Biomec III
01	7000-36-001	Stainless box for instrument acetabulum Biomec III
02	7000-36-002	Biomec III Acetabular Instrument Tray
01	5604-00-000	Extractor clamps insert
01	5602-10-000	Acetabulo Positioner / Impactor - Biomec III
01	5608-01-000	Shaving extender cable (1/4 chuck
01	5612-00-000	Spring-loaded gimbal wrench
01	5602-10-001	Fixed key 10/13
01	5607-10-000	Acetabular meter cable with stop
01	5121-10-000	Femoral Head Meter - Pair (38-62)
01	5605-10-000	Insertion Positioner / Impactor - Biomec III

Accessories

The family of non-cemented (metallic) acetabular capsules - BM according to regulation no. 01 (03/02/2009 - item 1.7.1 of part 1 - hip joint) has accessories that must be purchased separately in order to integrate the medical product as described in the following table:

Table - Accessories

Code	Description	Raw material	Type of sterilization	Illustration
2539-06-015	Acetabular screw Biomec III - Ø 6.0 X 15 mm	ASTM F138 STAINLESS STEEL	Ethylene oxide	
2539-06-020	Acetabular screw Biomec III - Ø 6.0 X 20 mm	ASTM F138 STAINLESS STEEL	Ethylene oxide	
2539-06-025	Acetabular screw Biomec III - Ø 6.0 X 25 mm	ASTM F138 STAINLESS STEEL	Ethylene oxide	
2539-06-030	Acetabular screw Biomec III - Ø 6.0 X 30 mm	ASTM F138 STAINLESS STEEL	Ethylene oxide	
2539-06-035	Acetabular screw Biomec III - Ø 6.0 X 35 mm	ASTM F138 STAINLESS STEEL	Ethylene oxide	
2539-06-040	Acetabular screw Biomec III - Ø 6.0 X 40 mm	ASTM F138 STAINLESS STEEL	Ethylene oxide	
Code	Description	Raw material	Type of sterilization	Illustration
2639-06-015	Acetabular screw Biomec III - Ø 6.0 X 15 mm	ASTM F138 STAINLESS STEEL	Gamma radiation	
2639-06-020	Acetabular screw Biomec III - Ø 6.0 X 20 mm	ASTM F138 STAINLESS STEEL	Gamma radiation	
2639-06-025	Acetabular screw Biomec III - Ø 6.0 X 25 mm	ASTM F138 STAINLESS STEEL	Gamma radiation	
2639-06-030	Acetabular screw Biomec III - Ø 6.0 X 30 mm	ASTM F138 STAINLESS STEEL	Gamma radiation	
2639-06-035	Acetabular screw Biomec III - Ø 6.0 X 35 mm	ASTM F138 STAINLESS STEEL	Gamma radiation	
2639-06-040	Acetabular screw Biomec III - Ø 6.0 X 40 mm	ASTM F138 STAINLESS STEEL	Gamma radiation	


Connecting components

Related hook components must be purchased separately, as they are not part of this product.

The following table illustrates the dimensional compatibility of implants:

- Unrestricted polymeric acetabular core (not cemented) - made of UHMWP 1900 polyethylene (NBRISO5834-2), registered with ANVISA (not subject to this registration and not members of this product).
- Femoral metal head - Biomecanica, made of stainless steel according to ASTM F138, registered with ANVISA (not subject to this registration and not included in this product).
- Femoral metal head - Biomecanica, made of stainless steel with high nitrogen content according to NBR ISO 5832-9, registered in ANVISA (not subject to this registration and not included in this product).

- Cemented modular rod - BM, made of stainless steel in forged state according to ASTM F621 standard, registered in ANVISA (not subject to this registration and not included in this product).
- Non-cemented modular rod - BM, manufactured in stainless steel according to ASTM F138 or manufactured in Cobalt-Chromium-Molybdenum according to ASTM F75 standard, registered in ANVISA (not subject to this registration and not included in this product).

Permitted combinations with other materials

The table "Permitted contact alloys" was based on the standard "NBRISO21534 - Implants for non-active surgery - Implants for joint replacement - Particular requirements", but regardless of the alloys mentioned in the item of the table "**Contact**" it is necessary that the combinations are of implants of the same brand, because they were designed for such combinations observing the surface finish and surface treatment and other factors required in projects that may interfere with Combinations. Therefore, metal implants from different manufacturers are not recommended for chemical, physical, biological and functional incompatibility.

Table - Permitted contact alloys

Permitted contact alloys	
Implant	*Contact (when there will be no articulation)
Family of non-cemented acetabular capsules (metallic) - BM in STAINLESS STEEL ASTM F138	Stainless steel formed 18 Chromium-14 Nickel -2.5 Molybdenum (NBRISO 5832-1 and ASTM F138) High nitrogen stainless steel (NBRISO5832-9) Titanium 4-vanadium aluminum 6-aluminum alloy (NBRISO 5832-3 and ASTM F136) Titanium 7-niobium 6-aluminum alloy (NBRISO5832-11) UHMWPE Polyethylene (ASTM F648, NBRISO5834-1, NBRISO5834-2)

*For applications where one metal or alloy is in contact with another and is not intended to be articulated, provided that appropriate attention is paid to design, surface finish, surface treatment and metallurgical conditions.

Commercial presentation of the product - Characteristics of the packaging

The family of non-cemented (metallic) acetabular capsules - BM is supplied in two types of sterilization, by ethylene oxide and gamma radiation. For each type of sterilization the containers are composed of:

Type of sterilization	Packaging description
Ethylene oxide	PVC or PET blister, crib, Tyvek, cardboard box and shrink PVC film
Gamma radiation	Non-toxic PET blister, crib, Tyvek, cardboard box and shrink PVC film

For sterilization in ethylene oxide gas, blister in PVC (polyvinyl chloride) or PET (polyethylene terephthalate) is used and PET (polyethylene terephthalate) containers are used for products sterilized by gamma radiation.

Both PET (Polyethylene Teeephthalate) and PVC (Polyvinyl Chloride) are raw materials indicated for the manufacture of ampoules and thermoformed cribs (Vacum - Forming) in the pharmaceutical, food, toy, electronic components, cosmetic and automotive industries. Both have as main characteristics: good impact resistance; optimal dimensional stability; non-toxic, allowing contact with food and hospital products; good moldability; Good printability and self-extinguishing. They are preferred for their mechanical properties, high transparency, virtually no defects, good hot fluidity and because they do not present any danger to human health.

The following figures illustrate a model of the family of non-cemented (metallic) acetabular capsules - BM and the acetabular screw packaged in the outer carton packaging. These boxes

are customized with the company logo and wrapped with shrink PVC film. Inside this cardboard outer packaging follows the double pack blister type sealed with Tyvek, implant traceability labels, instructions for use and patient card.

Figure - Model of the outer packaging of the family of non-cemented acetabular capsules (metallic) - BM

Figure - Model of the external packaging of the accessories (Biomec III acetabular screw)

Form of presentation of the product

External labeling is product traceability labels that contain the following information:

PRODUCTS STERILIZED IN ETHYLENE OXIDE GAS

Manufacturer details:

BIOMECANICA INDUSTRIA Y COMERCIO DE PRODUCTOS ORTOPÉDICOS LTDA.
 Street: Luiz Pengo, 145; 1st Industrial District; Postal Code: 17201-970; Jaú/SP Brazil
 Phone: +55 (14) 2104-7900 Telefax: +55 (14) 2104-7908
 CNPJ: 58.526.047/0001-73 / INSC EST: 401.042.207.113
www.biomecanica.com.br / biomecanica@biomecanica.com.br
 Aut. Func. ANVISA/MS nº X353Y38XY66Y
 BRAZILIAN INDUSTRY

- **Technical Manager**
- **ANVISA Registration Number**
- **Technical name of the product**
- **Trade name of the product**
- **Product Business Model**
- **Raw material**
- **Number of items per package**
- **Quantity and description of labels accompanying the product:**
 - 6 adhesive labels for identification and traceability of the product.
- **Say, "FORBIDDEN REPROCESSING"**

Symbology of the traceability label (Label) according to NBR ISO 15223

 xxx	 xxxx	 xxxx	 xxxx		
Product batch number	Date of manufacture	Use up to	Code	Identification of sterilized products in ETO Gas	"SINGLE-USE PRODUCT"

PRODUCTS STERILIZED IN GAMMA RADIATION


Manufacturer details:

BIOMECANICA INDUSTRIA Y COMERCIO DE PRODUCTOS ORTOPÉDICOS LTDA.
 Street: Luiz Pengo, 145; 1st Industrial District; Postal Code: 17201-970; Jaú/SP Brazil
 Phone: +55 (14) 2104-7900 Telefax: +55 (14) 2104-7908
 CNPJ: 58.526.047/0001-73 / INSC EST: 401.042.207.113
www.biomecanica.com.br /biomecanica@biomecanica.com.br
 Aut. Func. ANVISA/MS nº X353Y38XY66Y
 BRAZILIAN INDUSTRY

- **Technical Manager**
- **ANVISA Registration Number**
- **Technical name of the product**
- **Trade name of the product**
- **Product Business Model**
- **Raw material**
- **Number of items per package**
- **Quantity and description of labels accompanying the product:**
- 6 adhesive labels for identification and traceability of the product.
- **Say, "FORBIDDEN REPROCESSING"**

Symbology of the traceability label (Label) according to NBR ISO 15223

 xxxx	 xxxx	 xxxx	 xxx v		
Product batch number	Date of manufacture	Use up to	Code	Identification of sterile products in gamma radiation	"SINGLE-USE PRODUCT"

Information contained in the cardboard box (outer packaging)

Cardboard boxes (outer packaging) also have special conditions of Storage, Conservation and / or Handling of the Product with symbologies. The following table describes the information contained in this cardboard box.

Table - Special Conditions of Storage, Storage and/or Handling of the Product contained in the outer packaging of the product

Aseptically open the packaging

 <p>3°C °C</p> <p>Temperature limits</p>	<p>Store and transport in a dry and cool place, with ambient temperature (max. 35 °C) and relative humidity around 30% to 70%.</p>		
<p>Do not store directly on the floor (minimum height = 20cm) or in very high places near lamps, which could cause dryness of the packaging or damage to the label.</p>			
<p>Do not store in places where contaminants such as cleaning materials, insecticides, pesticides, etc. are stored.</p>			
	<p>Attention Instructions, warnings and precautions: See instructions for use</p>		<p>Fragile, handle with care</p>
	<p>Do not use if packaging is damaged</p>		<p>Keep out of the sun</p>
	<p>Keep dry</p>		<p>"SEE INSTRUCTIONS FOR USE"</p>
	<p>Identification of the sterile product</p>		

Commercial models of the non-cemented (metallic) acetabular capsule - BM family are individually distributed as described in the table below.

Table - Distribution of the family of non-cemented acetabular capsules (metallic) - BM

Code	Description	Raw material	Finishing	Sterilization
2537-46-001	Acetabular Component Biomec III Ø46mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-48-001	Acetabular Component Biomec III Ø48mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-50-001	Acetabular Component Biomec III Ø50mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-52-001	Acetabular Component Biomec III Ø52mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-54-001	Acetabular Component Biomec III Ø54mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide

2537-56-001	Acetabular Component Biomec III Ø56mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-58-001	Acetabular Component Biomec III Ø58mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-60-001	Acetabular Component Biomec III Ø60mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-62-001	Acetabular Component Biomec III Ø62mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-64-001	Acetabular Component Biomec III Ø64mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-66-001	Acetabular Component Biomec III Ø66mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2537-68-001	Acetabular Component Biomec III Ø68mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Ethylene oxide
2539-06-015	Acetabular screw Biomec III - Ø 6.0 X 15 mm	STAINLESS STEEL ASTM F 138	Polite	Ethylene oxide
2539-06-020	Acetabular screw Biomec III - Ø 6.0 X 20 mm	STAINLESS STEEL ASTM F 138	Polite	Ethylene oxide
2539-06-025	Acetabular screw Biomec III - Ø 6.0 X 25 mm	STAINLESS STEEL ASTM F 138	Polite	Ethylene oxide
2539-06-030	Acetabular screw Biomec III - Ø 6.0 X 30 mm	STAINLESS STEEL ASTM F 138	Polite	Ethylene oxide
2539-06-035	Acetabular screw Biomec III - Ø 6.0 X 35 mm	STAINLESS STEEL ASTM F 138	Polite	Ethylene oxide
2539-06-040	Acetabular screw Biomec III - Ø 6.0 X 40 mm	STAINLESS STEEL ASTM F 138	Polite	Ethylene oxide
Code	Description	Raw material	Finishing	Sterilization
2637-46-001	Acetabular Component Biomec III Ø46mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-48-001	Acetabular Component Biomec III Ø48mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-50-001	Acetabular Component Biomec III Ø50mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-52-001	Acetabular Component Biomec III Ø52mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-54-001	Acetabular Component Biomec III Ø54mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-56-001	Acetabular Component Biomec III Ø56mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-58-001	Acetabular Component Biomec III Ø58mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-60-001	Acetabular Component Biomec III Ø60mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-62-001	Acetabular Component Biomec III Ø62mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-64-001	Acetabular Component Biomec III Ø64mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-66-001	Acetabular Component Biomec III Ø66mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2637-68-001	Acetabular Component Biomec III Ø68mm - Gold Std	Stainless Steel - ASTM F138 with Titanium Alloy Coated by Plasma Spray - ASTM F 1580	Plasma spray	Gamma radiation
2639-06-015	Acetabular screw Biomec III - Ø 6.0 X 15 mm	STAINLESS STEEL ASTM F 138	Polite	Gamma radiation
2639-06-020	Acetabular screw Biomec III - Ø 6.0 X 20 mm	STAINLESS STEEL ASTM F 138	Polite	Gamma radiation
2639-06-025	Acetabular screw Biomec III - Ø 6.0 X 25 mm	STAINLESS STEEL ASTM F 138	Polite	Gamma radiation
2639-06-030	Acetabular screw Biomec III - Ø 6.0 X 30 mm	STAINLESS STEEL ASTM F 138	Polite	Gamma radiation
2639-06-035	Acetabular screw Biomec III - Ø 6.0 X 35 mm	STAINLESS STEEL ASTM F 138	Polite	Gamma radiation
2639-06-040	Acetabular screw Biomec III - Ø 6.0 X 40 mm	STAINLESS STEEL ASTM F 138	Polite	Gamma radiation

Indication, purpose or use for which the medical device is intended, as directed by the manufacturer.

Indication for use

The family of non-cemented (metallic) acetabular capsules - BM is indicated under the following morbid conditions*:

- IV. Primary osteoarthritis of the hip
- V. Osteoarthritis secondary to:
 - F. Rheumatic pathologies (rheumatoid arthritis, ankylosing spondylitis, acetabular protrusion of rheumatic origin)
 - G. Sequelae of childhood hip disease (developmental dysplasia of the hip, Legg-Calvé-Perthes disease, epiphosysis of the femoral head, non-rheumatic acetabular protrusion)
 - H. Osteonecrosis of the femoral head
 - I. Hip trauma (sequelae of femoral neck and head fractures and acetabular fracture)
 - J. Osteometabolic diseases (Paget's disease)
- VI. Hip fractures (neck, head and acetabulum)

* Due to the excellent risk and cost/effectiveness ratio of Hip Arthroplasty, at clinical discretion, other less frequent conditions that eventually evolve to joint destruction may be subject to the Total Hip Arthroplasty procedure, the contraindications and warnings of this document are observed.

Recommendations on orthopedic prostheses

To achieve better results in hip arthroplasties it is recommended:

- 1 - The use of orthopedic prostheses should be performed only by surgeons qualified and trained for such a procedure. It is essential to outline careful pre-operative planning, even with the help of transparency or product template.
- 2 - All care must be taken in the preparation of the recipient bone to obtain the perfect fit of the implant, avoiding radiolubricence and the undesirable appearance of micromovements.
- 3 - Instruments and surgical tests are available to assist in the surgical implantation of orthopedic prostheses. It is important that the surgical instruments and tests used are those designed specifically for the given product. Variation in the design and dimensions of surgical instruments and similar tests can compromise the critical measures required for precise implantation.
- 4 - To protect the doctor and his patient, in any future problems, the hospital must be responsible for the annotation in the patient's protocol, code and batch number of the implants used. This data is of fundamental importance to allow the traceability of this data.
- 5 - Before starting surgery, make sure the collection of implants and their instruments are intact and complete.

6 - Never mix implants from different manufacturers in the same patient. There may be differences in dimensions and tolerances in the fittings of these implants, causing poor joints, early release or implant failure.

7 - Orthopedic prostheses are classified as "Single-use products", that is, they cannot be reused.

8 - In all cases, the surgical practices consecrated in the postoperative period must be followed. The patient should be alerted to limitations of replacement or restitution of a given joint and recent history of use of this practice.

9 - The surgeon must also avoid producing folds or scratches on orthopedic prostheses, since these damages and / or malfunctions can produce internal "stress" that can become the focus of their eventual rupture.

10 - Postoperative care, as well as the patient's ability and willingness to follow instructions, are two of the most important aspects for the success of an arthroplasty. The patient should be advised that lack of attention to postoperative instructions may result in rupture or migration of the orthopedic prosthesis, requiring additional surgery for revision or removal.

11 - Every effort should be made to use compatible biomedical materials when using orthopedic prostheses, as mixing different materials in the same patient and micro movements of components can cause metallosis.

12 - The potential for success in replacing or restituting a given joint is increased by appropriate selection of the size, shape, and design of the prosthetic prosthesis. While proper selection can help minimize risks, the size and shape as well as the quality of human bones and adjacent soft tissues offer certain limitations regarding the size and mechanical strength of this implant.

13 - The orthopedic prosthesis is designed to replace a certain damaged joint and not the normal structures of the human skeleton.

14 - Orthopedic prostheses cannot withstand levels of activity and / or loads equal to those supported on normal and healthy bones. This detail must be communicated by the surgeon and understood by the patient.

15 - An obese or very heavy patient can produce overload in the implant used, causing fatigue in the metal and causing loosening, breakage or failure of the component.

16 - Excessive physical activity and trauma affecting the replaced joint have implied premature failure of orthopedic prostheses, either due to loss, fracture or wear of the implant. The patient should be warned to maintain his activities according to his situation, protecting the replaced joint against excessive "stress".

17 - In general, orthopedic prostheses are supplied in a sterile state, in double glazing, in an external box of rigid cardboard, duly identified by adhesive labels, with all the legal information relevant to the product, which ensure their complete identification and traceability. The packaging must be intact at the time of receipt (do not use the product if the packaging is tampered with). Check the validity of sterilization (do not use the product if the sterilization period is expired).

Resignation

Biomechanica, as the manufacturer of this device, does not practice medicine or recommend this or any other surgical technique for use in a given patient. The surgeon who performs any implant procedure is responsible for defining and using the proper techniques to implant the prosthesis in each patient. Biomechanica is not responsible for choosing the appropriate surgical technique to use in each patient. The surgical technique depends on the technique adopted by the surgeon, and it is up to him to choose the method, type and size of the products to be placed.

Contraindications

The only absolute contraindication to total hip arthroplasty is the presence of active infection at the surgical or systemic site. In this condition, implantation of a metallic material may perpetuate an existing infection or stimulate infection from other sites in the systemic infection. Patients with a remote history of infection in the hip joint may have a recurrence of the infectious process and therefore clinical, laboratory and imaging research is mandatory to avoid latent infections.

Precautions, restrictions, warnings, special care and clarification examinations on the use of the medical product, as well as its storage and transport.

Warnings

Although the results of total uncemented hip arthroplasty are well known and attest to its safety and efficacy, some conditions of use are not yet fully known or may still compromise the rescue of the patient's quality of life. Therefore, the following conditions should be considered:

6. **Use on immature hips**

7. **Inadequate bone support:**

Uncemented implants depend on their osseointegration into the adjacent bone, which occurs between 4-12 weeks after implantation. In this period, a primary mechanical stabilization of the implant is necessary, obtained by means of an adequate surgical technique and mainly of a bone support that allows a resistant fixation. The occurrence of proper stabilization can lead to a lack of integration and migration of the implant.

8. **Intraoperative fractures:** For immediate mechanical stability, contact pressure coupling is required, which means an implant of dimensions slightly larger than the cavity created intraoperatively, according to the surgical technique recommended by the manufacturer. Failure to observe this aspect can lead to the insertion of an oversized implant and thus cause fracture of the proximal femur.

9. **Early** loading:

Although with an adequate initial fixation (immediate) it is advisable that during the period of osseointegration, the load bearing is protected so as not to cause micromovement at the interface between implant and bone. The use of gait supports is recommended for 4-6 weeks after the operation.

10. **Implant cementation:**

Due to the mixture of plasma with titanium raw material, implants are indicated for use without bone cement. The cementation of these implants can predispose to corrosion in gaps in an environment without oxygen renewal (differential aeration cell) and thus cause catastrophic failures.

Precautions and warnings

SINGLE-USE PRODUCT. DESTROY AFTER EXPLANTING. DO NOT REUSE THE PRODUCT. PROHIBITED REPROCESSING

An orthopedic implant can only be used on a single patient once. Although it may appear undamaged, previous stresses can create imperfections that can reduce the success of the implant and lead to subsequent implant fracture.

The surgeon must be familiar with and have sufficient knowledge, including pre- and postoperative, surgical technique adopted, precautions and potential risks.

During implant handling, scratches, notches, or anything that could damage or mark the implant should be avoided, as these defects are stress concentrators and can be sites of crack nucleation and decrease corrosion resistance, which can result in implant fracture.

For implantation of the acetabular screw locked acetabulum, there is a clinically useful system for determining safe areas for screw placement. The system is based on two lines, one drawn from the anterosuperior iliac spine through the center of the acetabulum and the other drawn perpendicular to the first, thus creating four quadrants: anteroposterior, anteroinferior, posterior and posterior-inferior and posterior-inferior. The surgeon should be familiar with and have sufficient knowledge about this implantation technique. The rheterosuperior quadrant is the safest, and screws that exceed 25 mm can often be placed through the enlarged bone in this area. The anterosuperior quadrant should be avoided.

Make sure the packaging containing the product is intact. Check the sterilization date. Aseptically open the packaging after making sure the dimension is the one chosen. Sterilization validity: 5 years (indicated on the inner and outer packaging of the product). To ensure proper implantation, instruments manufactured by Biomechanica should be used, as these were designed and manufactured specifically for the use of these implants. The use of instruments from different manufacturers can compromise the surgery, in addition to being in disharmony with the registration of the product in the National Health Surveillance Agency.

The surgeon must be aware of the possibility of the patient's physical development, and there may be a need for revision of the implant.

Do not use the product if damaged.

SINGLE-USE PRODUCT: REUSE OF THIS PRODUCT IS EXPRESSLY PROHIBITED.

Do not use the product if the packaging is damaged.

Note: Implant components from different manufacturers should not be used. Therefore, we recommend that the products have the same origin.

Date of manufacture, shelf life and batch of the product: VIDE LABEL.

Adverse effects

Premature loosening of total hip replacement components may occur. Premature mechanical loosening may be the result of faulty fixation or hidden infection; Delayed release of biological complications or mechanical problems (high localized voltages). Peripheral neuropathies and heterotropic bone formation were reported in total hip arthroplasty. Subclinical nerve damage occurs most often, possibly due to the result of surgical trauma.

Metal sensitivity reactions in patients followed by total joint replacement were rarely reported. The importance and effects of sensitization await future clinical trials and evaluation.

The implantation of foreign material into tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. Similar changes can occur, such as a precursor, or during the healing process.

Displacements and subluxation can occur due to the result of incorrect positioning of the implant components. Muscle and sagging fibrous tissue can also contribute to these conditions.

Tell your orthopedist about any medical conditions that may affect surgery. Total prostheses are successful surgeries in more than 90% of patients. When complications arise, most can be treated successfully. Among the complications that may arise, we have:

Infection: Infection may occur in the surgical wound or may be deep (around the prosthesis). It can arise while the patient is in the hospital or at home. It may even appear years later. Minor infections in the surgical wound are usually treated with antibiotics. Major or deep infections may require additional surgery (for deep cleaning) or even removal of the prosthesis. Any infection in your body (bladder, throat, teeth, ears, etc.) can carry germs through the blood to your prosthesis and cause infection.

Thrombosis: These are blood clots resulting from several factors, including their reduced mobility, and cause blood to circulate more slowly in the veins, which can facilitate the appearance of thrombosis. Blood clots should be suspected if pain and edema appear in the thigh or calf. If this occurs, your orthopedist should do tests and tests to evaluate the veins in your leg, and it may be necessary to evaluate and follow up by a vascular surgeon. Several measures should be used to reduce the chance of blood clots, such as: - Medications to "thin" the blood (anticoagulants); - Elastic stockings; - Exercise to increase circulation in the legs and - Plastic boots that inflate and compress the soles of the feet and calf, increasing venous return.

Even with the use of these preventive measures, blood clots can still occur. If you notice edema (swelling), redness, or pain in your calf after discharge from the hospital, you should contact your orthopedist.

Loosening: Loosening of the prosthesis within the bone may occur after surgery. This can cause pain if loosening is significant. It is one of the indications for revision surgery (replacement of prostheses). New materials and new fastening methods should minimize this problem.

Wear: Some type of wear can be detected in any type of prosthesis. Excessive wear and tear can contribute to loosening and a revision of the surgery will be required.

Prosthesis breakage: Once implanted, the biomaterials used in the prostheses necessarily come into contact with body fluids. These fluids, although seemingly harmless, can, over time, significantly degrade most materials of considerable "chemical inertia." In addition, most orthopedic implants are subjected to static and/or cyclic mechanical stresses, which can cause material failure, i.e. rupture of the prosthesis.

Nerve injuries: Nerves in the vicinity of your prosthesis can be injured during surgery, although this is very rare. This injury is easier to occur when surgery involves corrections of large hip deformities or stretching of a very short leg due to wear and tear. Over time, these nerve injuries usually improve and can be fully recovered. Eventually, surgical exploration of the involved nerve may be chosen.

Information to be provided to the patient

To the statements of this instruction in the following elements: indications, contraindications, information of use, possible adverse effects, precautions and warnings.

The patient should be informed about the importance of post-surgical follow-up. The lack of follow-up prevents the detection of post-surgical problems and detects possible problems with the prosthesis.

The patient should be informed of the type of implant used, brand, code, name, batch, name of the surgeon and date of surgery.

Weight bearing, walking, muscle activity, among others, before the recommended time for bone consolidation can cause complications, and the patient should be aware.

In the postoperative care of an orthopedic surgical procedure, at will and the patient's ability to comply with instructions is extremely important. Children, elderly patients, patients with mental problems or chemical dependents, as they are prone to ignoring instructions and restrictions, can generate a high risk of failure of the procedure.

The patient should also be informed that the component does not replace and does not have the same normal bone performance. Therefore, it can deform, break or loosen due to excessive exertion or activities and early loading.

Evaluation of the implanted product

After implantation, intraoperatively, the responsible professional must perform a radiological control to verify the correct positioning of the product. The responsible professional must perform, and it is his responsibility, the clinical and radiological evaluations passed after the surgical procedure in the frequency stipulated by him to verify the status of the implant and the evolution of bone consolidation. If the product is out of the correct positioning, or has any non-compliance, it is the surgeon's responsibility to take the most appropriate corrective actions.

Useful information to avoid deployment risks

To reduce the risks arising from implementation, strictly follow: contraindications, instructions for use and all information contained in the "Instruction for use" of the product.

Intrinsic risks of implementation

The family of non-cemented (metallic) acetabular capsules - BM are manufactured with materials of recognized biomedical use, according to ASTM F138 -Standard Specification for Forged Stainless Steel Bars and Wires 18 Chromium-14 Nickel-2.5 Molybdenum for Surgical Implants (Uns S31673).

Contamination risks

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis, because components of the family of non-cemented (metallic) acetabular capsules - BM come into contact with tissue and body fluid. Explanted products should be treated as highly polluting products.

Sterilization

This product is supplied in sterile form. According to the customer's request, it can be supplied with the following sterilization processes:

Process	Validity
Ethylene oxide gas sterilization (ETO)	5 years
Sterilization by gamma radiation (gamma rays)	5 years

Note: The differentiation of the type of sterilization will be defined by the product code. ETO gas sterilized products have different gamma sterilized product codes.

If the product is expired due to sterilization and/or damaged packaging, it must be returned to the responsible distributor or directly to Biomecanica.

Medical Product Handling and Transportation Care

The storage place of the medical product must be clean, dry and illuminated to maintain ideal storage and transport conditions, as well as its physical and chemical integrity. The product should be stored and transported in a dry and cool place, with room temperature (max. 35 ° C), relative humidity about 30% to 70% and away from direct sunlight.

The effects of vibration, shock, corrosion, defective seats during handling and transport, incorrect stacking during storage should be avoided. Implants must be handled with care avoiding damage that may impair the quality of the material and patient safety.

They cannot be stored directly on the floor (minimum height = 20cm). They cannot be in very high places, near lamps, which could cause dryness of the container or damage to the label. They

should not be stored in places where contaminants such as cleaning materials, insecticides, pesticides, etc. are stored.

The carrier must be informed of the content and delivery time.

Transport carefully to avoid falls and friction so as not to damage the packaging and device, thus ensuring the sterility of the material.

The opening of the sealed blister can be done manually, in a suitable place (surgical environment), by pulling the Tyvek for the leftover. This opening should be carried out only in the surgical procedure to avoid contamination of the product.

Handling of sterilized material

When handling sterilized material with aseptic technique, some rules must be followed to keep it sterile:

- it is essential to wash hands with soap and water before handling sterile material;
- use material with intact, dry, unstained, identification packaging (type of material and date of sterilization);
- work in front of the material;
- handle the material at waist height up;
- avoid coughing, sneezing, talking about exposed material;
- do not make movements in the sterile area;
- ensure the validity and suitability of the packaging;
- work in a clean, quiet and dry environment without airflow;
- maintain a certain distance between the body and the material to be handled;
- Obey the other principles of asepsis.

The recommended nursing technique in the handling of sterilized material is:

- open it, starting at the opposite end of the manipulator;
- protect exposed material with the surrounding sterile field;
- touch with your hands only on the outside of the package;
- Do not store a previously opened package as sterile material.

Guarantee

The warranty of the Family of Non-Cemented (metallic) Acetabular Capsules - BM is related to the observance of this instruction for use. Where: indications and information for use, contraindications, care and precautions, warnings, possible adverse effects, packaging, sterilization, cleaning and decontamination recommended in this technical report and in the instructions for use.

Product disposal

The implants that make up the Family of Non-Cemented (metallic) Acetabular Capsules - BM explanted from patients must be properly disposed of by the hospital institution. It is the responsibility of the hospital institution to complete the decharacterization of the implant, avoiding its reuse and the method used for the decharacterization of the implant. Biomechanica recommends that the explanted implants be mechanically deformed with the help of a hammer or impact press and then identified as "Not suitable for use".

As described in Resolution No. 2605 of 08/11/2006, implantable devices of any nature framed as single-use are prohibited from being reprocessed.

Traceability

Next to the packaging of the implantable component follow 6 traceability labels containing the data of the implant used. This label is the same label on the external product labeling. The information on this traceability label is listed in the "Submission Form" item of this white paper. A label must be placed on the patient's medical records, another on the report delivered to the patient, another on the tax documentation that generates the collection, in the AIH, in the case of a patient attended by the SUS, or on the sales invoice, in the case of a patient attended by the complementary health system, another available for provider control (historical distribution record - RHD, another available for the control of the responsible surgeon (principal) and the last for medical insurance if any. We inform you that the patient's medical records must have information recorded to track the implanted product. The union of this label in the medical records allows the traceability of the product used. Among the main information we highlight as essential the name of the manufacturer of the implant used, name of the implant used, implant code, batch number and registration number of the product in ANVISA. This information is described in the traceability labels that follow with the product and in the external labeling. Other information should also be considered essential, such as the date of surgery, the name of the patient who received the implant, the name of the surgeon, the patient's weight, the patient's age, and other information requested in the patient's medical records should also be filled in.

Brand

The family of non-cemented (metallic) acetabular capsules - BM contains the following laser-etched information to allow traceability and follow-up of the patient in the postoperative period:

- Biomechanica logo
- Manufacturing batch number
- Acronym of the Material used in the manufacture of the product
- Code
- Dimension

The acetabular screws contain the following laser-etched information to allow traceability and follow-up of the patient in the postoperative period:

- Biomechanica logo

- Manufacturing batch number
- Acronym of the Material used in the manufacture of the product
- Dimension

After Sales (Customer complaint)

In case there is a need to reliably make any acetabular (metal) capsule - BM complaint related to any adverse effects affecting user safety, such as a non-functioning product, implantable metal component damage, serious problems or death related to these components, the responsible surgeon must report this adverse event to the competent health agency and Biomecanica via mail electronic sac@biomecanica.com.br or by Phone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health professional can communicate the adverse event through the Health Surveillance Notification System on the ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Sending material to the manufacturer for analysis

In case of sending implants to the manufacturer for analysis, it must be disinfected in the hospital using a broad-spectrum bactericidal and antifungal solution. It must then be disinfected or sterilized with autoclave steam or ethylene oxide. Biomecanica must be shipped in intact packaging, identified with the method of cleaning, sterilization and product data.

Description of the efficacy and safety of the medical product, in accordance with the anvisa regulations on the essential requirements of efficacy and safety of medical products.

The risks associated with the family of non-cemented (metallic) acetabular capsules - BM are in accordance with the essential requirements of safety and efficacy of the product according to RDC 56/01 for products falling within the scope of regulation 8, risk class 3 and risk management of products and processes are in accordance with NBRISO14971 as provided in procedure PRPD03 - Risk analysis management.

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