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Referência	Diâmetro (mm)	Comprimento (mm)	Indicação
6610-01-000	7,0	18,0	Hastes de Ø 1,5 a 2,5 mm
6610-02-000	8,0	25,0	Hastes de Ø 3,0 a 4,0 mm
6610-03-000	5,0	10,0	Hastes de Ø 1,5 a 2,0 mm
6610-04-000	6,5	12,0	Hastes de Ø 2,5 a 3,0 mm
6610-05-000	8,0	14,0	Hastes de Ø 3,5 a 4,0 mm
6610-06-000	7,0	18,4	Hastes de Ø 1,5 a 2,5 mm
6610-07-000	8,0	25,3	Hastes de Ø 3,0 a 4,0 mm
Construídas em aço inoxidável ASTM F138.			

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Referência	Diâmetro (mm)	Comprimento (mm)	Indicação
6611-01-000	7,0	18,0	Hastes de Ø 1,5 a 2,5 mm
6611-02-000	8,0	25,0	Hastes de Ø 3,0 a 4,0 mm
6611-03-000	5,0	10,0	Hastes de Ø 1,5 a 2,0 mm
6611-04-000	6,5	12,0	Hastes de Ø 2,5 a 3,0 mm
6611-05-000	8,0	14,0	Hastes de Ø 3,5 a 4,0 mm
6611-06-000	7,0	18,4	Hastes de Ø 1,5 a 2,5 mm
6611-07-000	8,0	25,3	Hastes de Ø 3,0 a 4,0 mm
Construídas em aço inoxidável ASTM F136.			

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## ENGLISH USA

### **Description of Available Models and Supplementary Technical Features.**

The Range of Flexible Intramedullary Nails – BM is comprised of metallic shafts whose purpose is providing stable fixation of bone fracture, such as tibia, fibula, femur, radius, ulna, humerus and clavicle, by means of intramedullary implantation.

This product's implantation procedure is less invasive, causing little or no further damage to the soft tissues affected by the surgery. The only accessory used in integration with the product is endcap, intended to prevent eventual sliding of the shaft after the surgery.

Figure 1 shows the Range of Flexible Intramedullary Nails - BM.

**Figure 1 – Range of Flexible Intramedullary Nails - BM.**

The implantable devices are single use, traded in units and non-sterile. The Range of Flexible Intramedullary Nails – BM is manufactured in two certifiably biocompatible materials – stainless steel according to the specification ASTM F138 and titanium alloy 6-aluminum 4-vanadium, according to the specification ASTM F136.

The range of products in this line comprises different diameters, ranging from 2 to 5 mm, with length ranging from 300 to 440 mm. Thus, the surgeon may define what the proper diameter is for the patient.

Details of the components of the Range of Flexible Intramedullary Nails – BM are provided in Tables 1 and 4.

**Tabela 2 – Presentation of the Flexible Intramedullary Nails F136 - Orion SP Flex**

Reference	Diameter (mm)	Length (mm)
3831-15-000	1,5	300
3831-20-000	2,0	440
3831-25-000	2,5	440
3831-30-000	3,0	440
3831-35-000	3,5	440
3831-40-000	4,0	440

Built of titanium alloy ASTM F136.

**Table 2 – Presentation of the Flexible Intramedullary Nail - F138 Orion SP Flex**

Reference	Diameter (mm)	Length (mm)
3832-15-000	1,5	300
3832-20-000	2,0	440
3832-25-000	2,5	440
3832-30-000	3,0	440
3832-35-000	3,5	440
3832-40-000	4,0	440
Built of stainless steel ASTM F138.		

**Table 3 – Presentation of Accessories (ENDCAPS) of the Flexible Intramedullary Nail - F138 Orion SP Flex**

Reference	Diameter (mm)	Length (mm)	Indication
6610-01-000	7,0	18,0	Nails with Ø 1,5 to 2,5 mm
6610-02-000	8,0	25,0	Nails with Ø 3,0 to 4,0 mm
6610-03-000	5,0	10,0	Nails with Ø 1,5 to 2,0 mm
6610-04-000	6,5	12,0	Nails with Ø 2,5 to 3,0 mm
6610-05-000	8,0	14,0	Nails with Ø 3,5 to 4,0 mm
6610-06-000	7,0	18,4	Nails with Ø 1,5 to 2,5 mm
6610-07-000	8,0	25,3	Nails with Ø 3,0 to 4,0 mm
Built of stainless steel ASTM F138.			

**Table 4 - Presentation of Accessories (ENDCAPS) of the Flexible Intramedullary Nail - F136 Orion SP Flex**

Reference	Diameter (mm)	Length (mm)	Indication
6611-01-000	7,0	18,0	Nails with Ø 1,5 to 2,5 mm
6611-02-000	8,0	25,0	Nails with Ø 3,0 to 4,0 mm
6611-03-000	5,0	10,0	Nails with Ø 1,5 to 2,0 mm
6611-04-000	6,5	12,0	Nails with Ø 2,5 to 3,0 mm
6611-05-000	8,0	14,0	Nails with Ø 3,5 to 4,0 mm
6611-06-000	7,0	18,4	Nails with Ø 1,5 to 2,5 mm
6611-07-000	8,0	25,3	Nails with Ø 3,0 to 4,0 mm
Built of stainless steel ASTM F136.			

### Composition

The shafts comprising the Range of Flexible Intramedullary Shafts – BM are manufactured in two metallic materials, to wit:

- 1) Stainless steel 18-Chrome 14-Nickel 2.5-Molybdenum – excellent biocompatibility material, successfully used for decades as orthopedic implantable devices, with chemical composition, microstructure and mechanical properties specified by ASTM F138.
- 2) Titanium alloy 6-Aluminum 4-Vanadium – titanium alloy alpha + beta, with chemical composition, microstructure and mechanical properties specified by ASTM F136. Also known in

use as orthopedic implant material, with lower density and rigidity as compared to stainless steel ASTM F138.

Table 6 contains a summary of the aforementioned information.

**Table 6 – Chemical rating of product raw material.**

Product	Material	Applied Standard
Flexible Intramedullary Nail F138 Orion SP Flex	Stainless steel 18-Chrome 14-Nickel 2.5- Molybdenum	ASTM F138*
Flexible Intramedullary Nail F136 Orion SP Flex	Titanium alloy Ti6Al4V	ASTM F136**
* Under standard ASTM F138 "Standard Specification for Wrought 18Chromium-14Nickel-2.5Molybdenum Stainless Steel Bar and Wire for Surgical Implants (UNS S31673)"		
** Under standard ASTM F136 "Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications"		

Attachment II – Testing and Analysis Report **shows the result of raw material characterization testing for manufacture of the Flexible Intramedullary Shafts in stainless steel 18-Chrome 14-Nickel 2.5-Molybdenum and titanium alloy Ti6Al4V. As reported in the inspection reports, the raw materials used meet the standards ASTM F138 and F136, respectively.**

**List of Tools to assist Implantation of the Flexible Intramedullary Shafts - BM**

For implantation of the Flexible Intramedullary Shafts – BM, using the tools specified in Table 7 below is required:

**Table 7 - List of Tools of the Range of Flexible Intramedullary Nails – BM**

Qtd	Código	DESCRIÇÃO	IMAGEM
01	4001-00-001	F Lever	
01	4001-00-002	Extraction Pliers	
01	4001-00-003	Flat Pliers	
01	4001-00-004	Drill Ø 2,7 x 125/ 100 mm	
01	4001-00-005	Drill Ø 3,2 x 195/ 170 mm	
01	4001-00-006	Drill Ø 4,5 x 195/ 170 mm	
01	4001-00-007	Hook Wrench	
01	4001-00-008	L Wrench	
01	4001-00-009	Cutter	

01	4001-00-010	Flexible Shaft Cutter - SPFlex	
01	4001-00-011	Flexible Shaft Insertion Guide - SPFlex	
01	4001-00-012	Double Drill Guide	
01	4001-00-013	Large Impactor for Flexible Shaft - SPFlex	
01	4001-00-014	Impactor for Flexible Shaft - SPFlex	
01	4001-00-015	Beveled Impactor	
01	4001-00-016	Insertion device for Flexible Shaft - SPFlex	
01	4001-00-017	Small-sized beveled impactor	
01	4001-00-018	Hammer	
01	4001-00-019	Combined hammer for Flexible Shaft - SPFlex	
01	4001-00-020	Sliding hammer	
01	4001-00-021	Wrench Piece	
01	4001-00-022	Puncture for Flexible Shaft - SPFlex	
01	4001-00-023	Angled Puncture	
01	4001-00-024	Wrench Pivot for insertion of buffer of Ø 1.5 to 2.5 mm	

The Flexible Intramedullary Shafts – BM, trademark Biomecanica, should be used only with Biomecanica tools. The tools comprise another record and should be acquired separately.

### Ancillary

The Range of Flexible Intramedullary Shafts - BM has no ancillary components.

### Admissible Combinations with other Materials

The shafts comprising the Range of Flexible Intramedullary Shafts – BM are used with no association with any other material. Thus, no assessment of admissible combinations with other materials is applied.

**Medical product's intended indication, purpose or use, as indicated by the manufacturer.**

**Indication of Use**


Products of the Range of Flexible Intramedullary Shafts - BM are indicated to reduce, align, stabilize and fix several types of fractures of the tibia, fibula, femur, radius, ulna, humerus and clavicle, by means of intramedullary implantation. These are devices which share the load and enable load support through the fracture site. According to the standard NBR ISO 14602 the non-active surgical implantable devices for osteosynthesis are used in corrective surgeries or for treatment of traumas, keeping reduction of fractured bones and stabilizing bone or adjacent structures, so as to provide bone merger or its correction, being removed or left *in situ* after compliance with the goals. The standard NBR ISO 15142-1 also reports use of intramedullary shafts in temporary stabilization of long bones with reduced strength due to fractures and/or diseases, also asserting that the shafts are usually removed when reaching their temporary stabilization purpose.

This range of shafts is aimed at children and teenagers. Therefore, it is different from the conventional ones due to their reduced dimensions and greater flexibility. Such features were defined so as to meet the specific requirements of such customers, whose bone structure is smaller and under growth. Thus, there is no locking in those shafts, because cross-fixation elements could cause pain and limit bone development.

**- INSTRUCTIONS FOR CHOOSING THE IMPLANT TO BE USED.**

The correct definition of the implant is essential for the successful surgery, the fracture fixation increases upon appropriate selection of the model. The implant dimensions are determined through X-Ray.

A shaft to be used is 30 to 40% of the diameter of the medullary cavity, use pre-surgery radiography to measure the isthmus and define the implant. Use of two Shafts with equal diameters is usually indicated to avoid varus or valgus in positioning. Pre-load of a shaft in the other should be in opposite directions as per Figure 2.

**Figure 2 - Opposite positioning of the Flexible Intramedullary Nails.**

**- LOAD RESTRICTIONS.**

Use of intramedullary shafts allows for the orthopedic surgeon to fix the bone fragments in the natural anatomic position. During the recovery of the patient, the surgeon controls the applied

load, gradually increasing as consolidation and the general psychological status of the patient occurs.

- FORECAST PERFORMANCE

- Total Fracture Fixation (Total Stability).
- Minimally Invasive Procedure.
- Avoid damage to soft tissues.

**Disclaimer**

Biomechanica, as the manufacturer of such device, neither practices medicine nor does it recommend this or any other surgical technique for use in certain patient. The surgeon under any implant procedure is responsible for defining and using the appropriate implant techniques in each patient. Biomechanica is not responsible for choosing the appropriate surgical technique to be used in each patient. The surgical techniques depend upon the technique adopted by the surgeon physician, and the latter shall choose the method, type and size of the products to be placed.

**Precautions, restrictions, warnings, special care and clarifications on medical device usage, as well as its storage and transportation.**

**Counter-indications**

The Flexible Intramedullary Shafts – BM should not be used in the following cases:

- Obliterated or blocked spinal cavity due to bone or implanted metal presence;
- Excessively narrow spinal cavity, not allowing for bone milling or natural accommodation of the implant;
- Any and all type of fracture which may not be properly fixed with the Lockable Intramedullary Shaft – Stainless steel;
- Severely curved femur or femur with severe distortion and/or deformity;
- Cases where the Lockable Intramedullary Shaft – Stainless steel may cross the open epiphyseal plate;
- Recent history of systemic or located infection;
- Local inflammation signs;
- Bone immaturity;
- Immunosuppressive pathologies;
- Mental disorder;
- Alcohol or drug abuse;
- Pregnancy;
- Severe Osteopenia;
- Severe instability due to lack of structural tissue integrity;

- Severe Osteoporosis;
- Morbid obesity;
- Fever or leukocytosis;
- Sensitivity or allergy to the metal in the medical product;
- All those known for use of bone fixation systems;
- Patients not willing or incapable of following the care and instructions after surgery;
- Any medical or surgical condition that could compromise successful surgical procedure;
- Every patient should be warned on the risks of walking without assistance of crutch or cane, before fracture consolidation, complete recovery and adaptation to the medical products used.

### **Warnings**

This product is single use and may not be reprocessed.

Don't mix components from different manufacturers. Use of components from different manufacturers inserts non-considered risks, removing any liability on the part of Biomecanica for adverse events that may occur.

Don't use components of materials different from those provided in the system.

### **Precautions and Warnings**

An orthopedic implant may only be used in a single patient, only once. Although it may seem to be undamaged, prior stress may create imperfections that may reduce successful implant. Improper implant selection may cause unusual stress to the implant and may result in subsequent implant fracture.

Due to its presentation and features, this product has no side effect, however, we advise monitoring the reestablishment of the patient by X-Ray as per International Safety Standards.

The Range of Flexible Intramedullary Shafts was designed to be a temporary stabilization element of fractures until bone consolidation and not to replace the regular structures of the human skeleton.

The surgeon should be familiarized and hold sufficient knowledge of osteosynthesis and its limitation, including pre- and post-surgery, adopted surgical technique, precautions and potential risks.

During implant handling, scratches, grooves, or anything that may damage or mark the implant should be avoided, because such faults are stress concentrators and may be crack centering sites and decrease resistance to corrosion, which may result in implant fracture.

In order to ensure an appropriate implantation, the tools manufactured by Biomecanica should be used, because they were designed and manufactured specifically for the use of such implantable devices. Using tools from different manufacturers may compromise the surgery in addition to being in disagreement to the product record at the National Sanitary Surveillance Agency.

**DON'T USE THE PRODUCT IF IT IS DAMAGED;**

**NON-STERILE PRODUCT – STERILIZE PRIOR TO USE;**

**SINGLE-USE PRODUCT – REUSING THIS PRODUCT IS EXPRESSLY PROHIBITED;**

**Special storage conditions:** Keep in ventilated, moist-free, sunlight-protected site away from weather inclemency.

Don't use the product if packaging is damaged;

**Note:** No implantable device components from different manufacturers should be used; therefore, we recommend the products to have the same origin.

Manufacturing date, due date and product batch: REFER TO LABEL.

### **Side effects**

Mechanical loosening may be the result of defective fixation, unstable reconstruction of the fractured bone or concealed infection. Reactions of metal sensitivity in patients were rarely reported.

Implantation of foreign matter in tissues results in histologic reactions involving several macrophage and fibroblast sizes. The chemical importance of such effect is uncertain, as well as similar changes may occur as a precedent, or during the healing process.

Other possible adverse events that could occur when implanting an intramedullary shaft are:

- surface or deep post-surgery infection;
- patients with exposed fracture have longer consolidation time as well as greater infection rate;
- vascular disorders, including thrombosis and lung embolism;
- angle diversions and limb or bone segment shortening operated as function of bone reabsorption;
- exposure of the synthesis material through skin necrosis, contributing towards installation of infectious processes, due to low quality of skin coverage;
- hemorrhage of blood vessels and/or hematomas;
- pain, discomfort or abnormal feelings due to the presence or migration of the medical product, and also due to the surgical procedure;
- neural or neurological damage as function of the surgical trauma (including paralysis and soft tissue injuries);
- premature or late failure of the medical product;
- loosening, migration, bending or fracture of the medical product;
- lack of consolidation or delayed bone consolidation which may lead to breaking of the medical product;
- vicious consolidation;
- loss of the correction degrees, height and/or reduction, obtained in the surgical procedure; e
- pseudo arthrosis (no joining).

### **Information to be Provided to the Patient or Responsible Person**

The statements of this instruction in the items:

Indications – Counter-Indications – Information of Use - Possible Side Effects – Precautions and Warnings.

The patient should be advised as to the importance of the post-surgical follow-up. Failure to monitor prevents detecting post-surgical issues, such as component loosening or osteolysis occurrence. Failure to carry out the revision surgery upon loosening of components or osteolysis may result in progressive loss of the periprosthetic bone stock.

The patient or the tutor should be reported as to the type of the implant used, trademark, code, name, batch, surgeon name and surgery date.

In the post-surgery care of an orthopedic surgical procedure, the patient's will and capacity to comply with the instructions is extremely important. Children, teenagers, or patients with mental disorders or drug addicts, since they tend to ignore the instructions and restrictions, may generate a high failure risk of the procedure.

The Patient or the tutor should also be informed that the component does not replace and does not have the same regular bone performance. Therefore, it may deform, break or loose due to excess effort or activities and premature load. The implant performance is related to factors, such as patient weight, physical activity level, compliance with the instructions provided by the physician post-surgery and mainly in case of delay or lack of consolidation.

### **Important**

The surgeon should instruct the patient on the care to be taken and the proper physical support to be used throughout use of the implant. In view of a perfectly stabilized fracture, the load may be released. If the fracture is not stabilized, no load may be given to the implanted system.

As recommended by the surgeon, remove the implant only after complete fracture consolidation. Removal of the implant does not require any specific tool, and the same tool is used for insertion, as well as for implant removal.

### **Product Shaping Limit**

When handling the implant, scratches, grooves, or anything that may damage or mark the implant should be avoided, because these flaws concentrate tensions and may be crack centering sites and decrease resistance, which may result in fracture or premature deformity of the implant.

### **Implanted Product Appraisals**

After implantation, during intra-surgery, the professional in charge should carry out the radiological control so as to check correct product positioning. The professional in charge should carry out, and it is under his liability, the clinical and radiological appraisals after the surgical procedure within the frequency set forth by him/her so as to check the implant status and bone

consolidation evolution. If the product is out of proper positioning, or has any non-compliance, the surgeon shall take the most appropriate corrective measure.

### **Useful Information to Avoid Risks Deriving from the Implantation**

In order to decrease the risks arising out of the implantation, the following shall be carefully followed: counter-indications, instructions of use and all pieces of information contained in the product's "Instructions of Use".

### **Contamination Risks**

There are biological contamination risks and viral disease transmission risks, such as HIV and hepatitis, because the metallic components of the Range of Flexible Intramedullary Shafts - BM contact body tissues and fluids. Explanted products should be treated as highly contaminant.

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All Models of the Flexible Intramedullary Shafts – BM are traded as non-sterile; therefore, the institution to carry out the surgery shall be liable for its sterilization shall use regularly validated methods.

The product is shipped in its previously clean original packaging. Prior to use, the packaging should be inspected aimed at checking that there are no cracks or damages. The product should only be taken out of its packaging prior to sterilization for the purpose of preserving the surface finishing and the set-up, handling it as little as possible, as recommended in the standard NBR ISO 8828.

In order to ensure reduction of the infection risks, the implantable devices (which were not used in the surgery, but were placed over the surgical table or which lost sterilization) and surgical devices should be cleaned after the surgery and prior to sterilization.

Cleaning may be developed through manual or mechanic methods.

Cleaning equipment:

- Ultrasonic washer
- Washer - sterilizer
- Washer – disinfectant

Manual cleaning: procedure in which dirtiness is removed through physical action upon assistance of detergent, water and articles such as sponge and brush. Manual cleaning is more recommended because it is less aggressive to the implantable and surgical devices. In manual cleaning, water at room temperature should be used, nylon brushes should be used, never steel or abrasive brushes, because it may impair the protective layer of the material.

A wide-range bactericidal and antifungal solution should be used in the decontamination, thus avoiding aggressive cleaning agents or metal brushes, so that quality of the implantable metallic devices is not impaired.

In case of cleaning machines, the implantable and surgical devices should be positioned inside drawers so that there is no shock between them, avoiding damages to the material.

Professionals in charge of cleaning the implantable and surgical devices should pay attention to the type of cleaning product used abiding by the concentrations recommended by the manufacturer. The implantable and surgical devices should be carefully rinsed until the entire waste is wiped-off and dried right after cleaning.

✓ Types of cleaners

Enzymatic cleaners: basically comprised of enzymes, surfactants and solubilizers. Balanced combination of such elements causes the product to remove the organic matter of the material in a short period of time.

Enzymatic solutions: excellent cleaning action, but no bactericidal and bacteriostatic activity.

Enzymes: substances produced by live cells and govern the process chemical reactions. Once produced by the cells, an enzyme may be isolated and shall keep its catalytic properties if certain conditions are kept in its manufacturing. Enzymes are rated in three larger functional groups, depending on the type of substrate they will affect: proteases, lipases and amylases which act in protein substrates, fats and carbohydrates, which tend to solubilize and loosen from the articles. Currently, cleaning of complex set-up articles is recommended so as to ensure cleaning.

✓ Disinfection:

Process which destroys pathogenic or non-pathogenic microorganisms of the articles, except for bacterial spores, by physical or chemical means.

✓ Disinfection levels:

High level: destroys all microorganisms, except for those at high number of spores => Glutaraldehyde 2% - 20 – 30 minutes.

Indication: preferably hospital area.

Medium level: remove vegetative bacteria, most virus, fungus and micro-bacteria => Sodium hypochlorite 1% - 30 minutes.

Indication: for UBS, nursery, rest homes, asylums.

Low level: removes most of bacteria, some virus and fungus, but does not remove micro-bacteria => Sodium hypochlorite 0,025%

Indication: nutrition.

**Sterilization:** This product is provided non-sterile. Prior to use, it should be sterilized. We recommend steam sterilization in Autoclave at the hospital (*ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

**AUTOCCLAVE:** Pressure saturated steam sterilization equipment.

**OPERATING INSTRUCTIONS:** For equipment start-up, check:

- that circuit breaker is ON;
- that water log is open;
- that the discharge log is closed.

- 1 – open equipment door;
- 2 – properly seat the material to be sterilized;
- 3 – close equipment door;
- 4 – select desired cycle according to the material to be sterilized;
- 5 – turn main switch ON;
- 6 – cycle automatically runs, in sequence;
- 7 – when “cycle end” lamp turns ON, partially open the door for approximately ten (10) minutes for material cooling.

The following autoclave physical sterilization parameters (saturated steam) is advisable to be applied:

**Table 8 - Parameters of physical sterilization in autoclaves**

Cycle	Temperature	Exposure Time (minium)	Drying Time (maximum)
Conventional (1atm of pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Note: The time shall be recorded when the heat of the sterilization chamber reaches the desired temperature.

**FOR FURTHER INFORMATION, PLEASE CHECK THE INSTRUCTION MANUAL FOLLOWING EACH AUTOCLAVE.**

**Another sterilization method:** to be used in addition to the autoclave as defined by the hospital institution there's the following:

- **Sterilization by ethylene oxide (ETO)** - parameters and procedures set in the validation protocol and in *ISO 11135-1 - Sterilization of health care products - Ethylene oxide -- Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.*

**NOTE**

The hospital institution is in charge of the sterilization method, equipment, the controls, and instructions of the sterilization used.

Care with the sterilized articles.

Storage conditions of the sterilized articles:

- as to the environment: should be clean; ventilated; dry; should be restricted to the sector team.
- as to the article: after the sterilization process, don't place it over cold surface (stone or stainless steel), use hollow baskets or containers until cooling;
- casing (raw cotton cloth, non-woven cloth, surgical degree paper, crepe paper, film paper, Tyvek or perforated metallic cartons) should remain complete and be little handled to avoid packs to tear or lose the seal;

- be stored in sealed wardrobes with shelves;
- shelves identified so as to make material removal easier;
- material should be stored according to sterilization due date to make distribution easier and not remaining due material in stock;
- store separately from non-sterile products towards reducing the outer contaminant level.

### **Precautions with Medical Device Storage, Conservation, Handling and Transportation**

- Don't use the product if it is damaged.
- The implantable devices should be carefully handled avoiding damages which may impair quality of the material and the patient's safety.
- Store in a vented, dry place, protected from sunlight and away from weather inclemency.
- When transporting, avoid shocks, and inadequate piling.
- Store and transport in dry and fresh site, with room temperature (Max. 35°C) and relative humidity ranging from 30% to 70%.
- Don't store directly on the floor (minimum height = 20cm) or in too high sites, near lamps, which could cause packaging desiccation or damage to the label.
- Don't store in sites where contaminant substances are stored such as, for example, cleaning materials, insecticides, pesticides, etc.

### **Product Disposal**

The implantable devices comprising the Range of Flexible Intramedullary Shafts - BM explanted from patients should be duly disposed of by the hospital institution. The hospital institution shall fully destroy the implant preventing its reuse and the method used for destruction of the implant. Biomecanica recommends that the implantable devices explanted are mechanically deformed upon assistance of a hammer or impact press and then they shall be identified with the sentence "Improper for Use".

As described in Resolution No. 2605, dated 08/11/2006, implantable devices of any nature characterized as single use, shall not be reprocessed.

### **Post-Sale (Customer Complaint)**

If you have any complaint as to the Range of Flexible Intramedullary Shafts - BM associated with any adverse event affecting user safety, such as product not operating, damage to the implantable metallic component, severe problems or death related to such components, the surgeon in charge is to communicate such adverse event to the applicable sanitary agency and Biomecanica by email [tecnovigilancia@biomecanica.com.br](mailto:tecnovigilancia@biomecanica.com.br) or phone 0xx14 2104 7900. When in doubt, the surgeon in charge or the healthcare professional may report the adverse event by Sanitary Surveillance Notice System at ANVISA website:

<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

## **Traceability**

6 traceability tags follow with the implantable component packaging, containing data of the implant used. Such tag is the same tag contained in the product's outer label. The information contained in the traceability tag is listed in item "**Presentation Form**" of this technical report. One tag should be placed to the patient's medical record, another one to the medical report delivered to the patient, another to the fiscal document generating the charge, at AIH, in case of a patient operated by SUS, or to the sale invoice, in case of a patient operated by the complementary healthcare system, another one available for supplier control (distribution historic record – RHD, another one available for main surgeon control (main) and the last one for medical covenant if any. The patient's medical record must contain information recorded which allows for tracking the implanted product. Placing of this tag to the medical record allows for tracking the product used. Among the main pieces of information, we highlight as essential the used implant manufacturer name, name of implant used, implant code, batch number and product record before ANVISA. This information is described in the traceability tags following the product and the outer label. Other information should also be deemed as essential, such as surgery date, name of patient who received the implant, surgeon name, patient weight, patient age and other information ordered in the patient's medical record should also be filled-out.

## **Warranty**

Warranty of the Range of Flexible Intramedullary Shafts - BM is related to compliance with the instructions of use. Being: indications and Information of Use, counter-indications, care and precautions, warnings, Possible Side Effects, packaging, sterilization, cleaning and decontamination recommended in these instructions of use.

## **Medical Device Presentation Method**

The implantable devices are stored in plastic envelopes with non-sterile tag and outer label distributed in unit form. The product contains the following information:

- a)** Manufacturer Name
- b)** Product Trade Name
- c)** Technical Name
- d)** Batch number
- e)** Record number at the Health Department
- f)** Product Code
- g)** Quantity
- h)** Product description containing the packaging (Model)
- i)** Manufacturing date
- j)** Due date: Indefinite due date
- k)** Description of raw material used to manufacture the product
- l)** Manufacturer Address (described in the plastic envelope)

- m) Name of expert in charge
- n) Statement: "Non-sterile product"
- o) Statement: "Reprocessing is prohibited"
- p) "Single-use product"
- q) Prior to use, read the instruction for use
- r) Packaging description and dimension: plastic envelope manufactured through low-density polyethylene, with 8 cm in width

**Symbols of the Traceability Tag (Label) as per NBRISO15223:**

	Manufacturing date		Due date
x) 	Product batch number		Manufacturer Data
	"Single-use product"	x) 	Code
	Non-sterile		

***The non-sterile product is identified with a green-colored tag containing "Non-sterile"***

NON-STERILE

➤ **Outer packaging symbol:**

				
Fragile, handle with care	Keep protected from sunlight	Keep dry	Caution	Don't use if packaging is damaged

Figure 4 shows the Tag Model used externally and internally in the packaging for identification and maintenance of the product traceability. In addition to the tag placed to the outer packaging, there are two other separate tags inside the packaging for traceability maintenance.

Figure 4 – Tag model

Figure 5 shows the plastic packaging model containing plate or screw (because they are traded in unit form and separated), instructions of use and label. The Flexible Intramedullary Nails - BM are to be traded in units.

Figure 5 – Product Packaging Illustration.

### Labeling

The Range of Flexible Intramedullary Nails - BM contains the following information recorded by laser to enable traceability and follow-up of the patient after the surgery, according to the space made available in each product:

- Biomecanica Logo
- Manufacturing Batch number
- Acronym of the raw material used in product manufacturing (T3 for Titanium Alloy).
- Code
- Dimension

### Material Shipment for Manufacturer's Analysis

If implantable devices are sent to manufacturer to enable analysis, it should be cleaned at the hospital using a wide-spectrum bactericidal and antifungal agent. Then, it should be disinfected or sterilized through steam in autoclave or ethylene oxide. It should be sent to Biomecanica in full packaging, identified with the cleaning, sterilization methods and product data.

**Description of efficacy and safety of the medical device, in compliance with ANVISA regulation providing for Core Medical Product Efficacy and Safety Requirements.**

Risks associated with the Range of Flexible Intramedullary Nails - BM comply with the Core Medical Product Efficacy and Safety Requirements as per RDC 56/01 for products listed in the standard 8, risk class 3 and Product and Process Risk Management is in accordance with NBR ISO14971 as per the procedure PRPD03 – Risk Analysis Management.

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