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Phone number: +55 (14) 2104-7900

CNPJ (Corporate Taxpayer Identification Number): 58.526.047/0001-73 / STATE REGIS-
TRATION: 401.042.207.113

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1.2. *Technical Manager:* Josiane da Silva Taietti - CRQ-SP
04272807

1.3. *ANVISA Registration no.:*

1.4. *Technical Product Name:* Bone Cement Applicator
and Mixer (Code: 1551630)

1.5. *Trade Product Name:* Cementation Kit - BM (Code:
4000-05-183)

1.6. Description of Kit Components, Quantity and Raw Material: Table 1 lists the components of the Cementation Kit - BM, quantities and raw materials. Table 2 shows the kit in its blister.

Table 1 - List of Components of Kit for Cementation - BM

Assembly Code and Description	Component Code	Item Description	Photo	Component Indication	Quantity	Raw Material
Cementation Kit - BM	084.0007A.1.1	Hole Nozzle		Connected to the hose, inject bone cement into the Cannulated Screw* (* not offered and not listed in this instruction)	3	Aluminum (ISO 209)
	084.0006H.1.8	Cover		Close the unused exits of three-way adapter	2	Polypropylene (PP)
	084.0006B.1.2	Hose		Carry the cement of the applicator to the surgical site	3	Policloreto de Vinila / Policarbonato (PVC / PC)
	140.00005.1.1	Three-way adapter		Used to adapt the applicator to up to 3 hoses.	1	Polyvinyl Chloride/ Policarbonato (PVC/ PC)

Cementation Kit - BM	140.00006.1.1	Funnel		Used to fill the applicator with bone cement	1	Polypropylene (PP)
	140.00004.1.1	Basis		Used to hold the applicator upright when handling bone cement or bone substitute	1	Polypropylene (PP)
	140.00008.1.1	Thread adapter		Used for adaptation	1	Polyamide (PA)
	140.00002.1.1	Applicator		Carry the mixture of orthopedic cement to the surgical site	1	Polyamide (PA) Polybutylene Terephthalate (PBT) Polypropylene (PP) Rubber Ethylene-Propylene-Diene (EPDM)
	140.00003.1.1	Applicator Thread		Used to move the plunger of applicator	1	Polyamide (PA)

Table 2 - Photo of the Cementation Kit - BM (Assembled on the blister)

Kit Code	Description	Illustration for Assembly in Blister
4000-05-183	Cementation Kit - BM	

1.7. Product ID:

Symbols of the Traceability Label (Label) according to NBR ISO 15223	
	Manufacturing date
	Expiration date
	Lot number
	Manufacturer data
	Code
	ETO Sterilization
	"Disposable product"
Reuse is forbidden	
Discard after Use	

1.8. Accessories:

No accessories.

1.9. Types of commercial presentation of the product:

Kit is doubled packaged in a blister manufactured in non-toxic PET (Polyethylene Terephthalate). Blister pack is heat sealed with puncture-free TYVEK paper and placed in a cardboard with a security seal.

2. Product Specifications

2.1 Operating Principle/Mechanism

Cementation Kit - BM is used to assist the surgeon during the preparation and cementation of the bone in surgeries. Kit contains a threaded applicator coupled with three-way adapter with hose connection. The hole nozzle is connected to each hose, enabling the injection of cement in the 3 ways simultaneously.

2.2 Indications

The Cementation Kit - BM is intended to aid in the fixation of bone fractures and bone reconstructions. When used for these indications, the Cementation Kit - BM can also be used to provide cement or bone substitute at the surgical site.

2.3 How to Use Product

It must be handled carefully in suitable environments. These devices must be handled only by skilled professionals. Surgical techniques may vary according to the surgeon's choice, and he will choose the method and the criteria for evaluating the results of the surgery. They must be applied and adapted according to the appropriate requirements and surgical techniques.

2.4 Expiration date

Expiration date: 5 years

When applicable, expiration date after opening: immediate use.

2.5 Sterile Product

Sterilization method: ETO.

2.6 Reprocessing

Forbidden.

2.7 Storage Conditions

The place must be clean, dry and clear to keep the ideal conditions of storage and transportation, as well as its physical and chemical integrity. The product must be stored and transported in a dry and cool place, in room temperature (Max. 35°C [95°F]), relative humidity around 30% to 70% and protected from direct sunlight. It cannot be stored directly on the floor (minimum height = 20cm). It cannot be placed in very high places, close to lamps, which could dry package or damage the label. It must not be stored in places where contaminating substances are stored, such as cleaning materials, insecticides, pesticides, etc.

2.8 Conditions for Transportation

transport it carefully to avoid falling and friction in order not to damage the packaging and the device. Avoid vibration, shock, corrosion, defective seating during handling and transportation, improper stacking during storage. The product must be handled carefully, avoiding damages which could impair the quality of the material and patient safety. Carrier must be informed about the content and delivery time.

2.9 Handling Conditions

The sealed blister pack can be opened manually in an appro-

ropriate place (surgical environment) by pulling Tyvek. Skilled personnel must open it only during surgery to avoid contamination of the product. During handling, avoid scratches, nicks, or anything which could damage or mark the components. Transport it carefully to avoid falling and friction in order not to damage the packaging and the component, thereby ensuring the sterility of the material.

2.10 Warning

- Cementation Kit - BM is a single-use product. Otherwise, biological contamination or transmission of viral diseases may occur.
- Cementation Kit - BM is disposable, and reprocessing is strictly forbidden.
- The use of additional accessories is forbidden.
- Kit components should not be used if the packaging is broken, with holes, expired sterilization period, stains or defects. The product must be kept in its original packaging until its use, according to the hospital medical asepsis procedures. If packaging is broken or damaged, all components of the Kit must be discarded.
- Conditions of storage and transportation of the product must be respected.

2.11 Caution

Surgeon must be familiar with and have enough knowledge on surgical fixation technique for spinal implants and its limitations, including pre and postoperative, used surgical technique, precautions and potential risks.

2.12 Contraindications

Not applicable.

As this Kit of instruments supports the procedure, there are no contraindications regarding the Kit. Contraindications are related to the used biomaterial and the surgical procedure. So, information on contraindications for bone cement used in surgery must be noticed.

2.13 Adverse Effects

Not applicable.

As this Kit of instruments supports the procedure, there are no adverse effects regarding the Kit. Adverse effects are related to the used biomaterial and the surgical procedure. So, information on how to use bone cement in surgery must be noticed.

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