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
Surgical Instrument in Nitinol – RG.: N° 80128580182

ENGLISH USA

Table 1: List of Surgical Instruments in Nitinol - BM

| CODE | DESCRIPTION | LINE | MATERIAL | INDICATION | DIMENSION | IMAGE |
|-------------|-----------------------------------|--------|---|--|-----------|---|
| 6472-00-008 | Lifter for Flexible Reamer 500 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 500 mm |  |
| 6472-00-009 | Lifter for Flexible Reamer 490 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 490 mm |  |
| 6472-00-010 | Lifter for Flexible Reamer 480 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 480 mm |  |

| | | | | | | |
|-------------|-----------------------------------|--------|---|--|--------|---|
| 6472-00-011 | Lifter for Flexible Reamer 470 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 470 mm |  |
| 6472-00-012 | Lifter for Flexible Reamer 460 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 460 mm |  |
| 6472-00-013 | Lifter for Flexible Reamer 450 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 450 mm |  |
| 6472-00-014 | Lifter for Flexible Reamer 440 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 440 mm |  |
| 6472-00-015 | Lifter for Flexible Reamer 430 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 430 mm |  |

| | | | | | | |
|-------------|-----------------------------------|--------|---|--|--------|---|
| 6472-00-016 | Lifter for Flexible Reamer 420 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 420 mm |  |
| 6472-00-017 | Lifter for Flexible Reamer 410 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 410 mm |  |
| 6472-00-018 | Lifter for Flexible Reamer 400 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 400 mm |  |
| 6472-00-019 | Lifter for Flexible Reamer 390 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 390 mm |  |
| 6472-00-020 | Lifter for Flexible Reamer 380 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 380 mm |  |

| | | | | | | |
|-------------|-----------------------------------|--------|---|--|--------|---|
| 6472-00-021 | Lifter for Flexible Reamer 370 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 370 mm |  |
| 6472-00-022 | Lifter for Flexible Reamer 360 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 360 mm |  |
| 6472-00-023 | Lifter for Flexible Reamer 350 mm | Trauma | Tube - Nitinol Alloy ASTM F2516 and Couplings - Stainless Steel 304 ASTM F899 | Connection instruments for reamers with coupling | 350 mm |  |

Traceability Tag Symbology (Label) according to NBR ISO 15223

| Manufacture Date | Expiration Date | Código do Lote | Manufacturer Data | Reference Number | Non-sterile |
|--|--|------------------|---|--------------------|---|
|  XXXX |  XXXX | LOT XXXXX |  | REF XXXX-XX |  |

TECHNICAL INFORMATION

Product Description: Surgical instruments can have different designs. Some have been developed specifically to meet the preferences of certain surgeons. Many bear the name of the surgeon for whom they were created. New instruments have been developed to improve various surgical procedures or even help to create new techniques. Although innovation means new designs and even new categories of instruments, in general instruments can be grouped by type or applications. Surgery techniques vary according to the knowledge and choice of the surgeon. It is up to the surgeon to make the final choice of the method, type and dimension of the instruments and implants to be used, as well as the criteria for evaluating post-surgical results. Biomechanica brand instruments must be used during the surgical procedure according to the purpose for which each instrument was developed. The use of these instruments is indicated for use by professionals only. Surgeons who manage the use of these products need to have a thorough understanding of the implant processes as well as the handling of implant instruments and components. Sharp instruments should only be used for their specific purposes. Misuse can cause blades and edges to become blunt, as well as affect their performance during the surgical procedure. Instruments that fall to the ground are invariably damaged. When this happens, the instrument must be set aside, carefully inspected and sent for repair. Counting the instruments is essential. Firstly, to protect the patient and, secondly, to avoid the possibility of being inadvertently sent to the laundry along with the surgical drapes, where they can be lost or damaged if they go unnoticed.

Manufacturing materials: Surgical instruments, on their structural part, are manufactured in Nitinol alloy, as specified in ASTM F2516; that establishes the chemical, physical, mechanical and metallurgical requirements for products composed of nickel-titanium alloys, used as a material in the manufacture of implants for surgery. Couplings are made of 304 stainless steel specified in ASTM F899.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Surgical instruments are indicated as an indispensable aid for the preparation of the bone segment that will receive the implant. Biomecânica instruments can only be used to help place Biomecânica implants, as they have been manufactured and designed to be used together. Using instruments from different manufacturers can compromise the surgery. With time of use, instruments tend to wear out naturally due to their regular use. Instruments should only be used for what they were specified. When reusing the Instruments in another patient, they should be washed, sanitized and then resterilized. This process must be performed whenever they are used in a new surgery. The reuse of instruments without going through a process of cleaning and re-sterilization is strictly prohibited. Worn and broken instruments can cause significant problems if they are not detected during inspection and testing. They may damage other instruments during the

cleaning/sterilization process. If they remain in the set of instruments, they can cause traumas in organic tissues to failure in critical moments during a surgical procedure. Identify problem instruments as quickly as possible. Separate them from the rest and put labels indicating that they must be replaced or repaired. The contraindication is linked to the implant to be used as described in the instructions for use accompanying the implants. Instruments should not be used for purposes other than those for which they were designed.

Cleaning: Proper care of surgical instruments begins with proper cleaning. This report mentions some methods of cleaning and sterilizing instruments. Proper care also means regular maintenance of instruments, including sharpening and adjustments. There is no standard schedule; maintenance will be determined by the frequency of use. In addition to the care of cleaning instruments, this instruction mentions several enemies of surgical instruments, such as: blood, tissue in general, surgical residues (which are the primary causes of pitting corrosion points), stains and discoloration of instruments. Water and moisture also have harmful effects, allowing these substances to dry or soak on your instruments causing undesirable stains. Other enemies used in washing instruments with improper solutions are: soap, bleach, disinfectants and other solutions not recommended. For the correct conservation of your instruments it is important to use recommended cleaning methods and to understand the causes of undesirable effects, such as stains. The spots appear orange or brown. The idea is to ensure the appropriate care to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, which lead to corrosion, appearing more often with an orange-brown color. If the instrument remains for an extended period of time (from one to four hours), in contact with these residues, marks and stains will appear on the instrument, especially if these residues dry together with the instrument. Instruments must be thoroughly cleaned and dried after use. Only sterilize a clean instrument. The high temperature of the autoclave will cause chemical reactions that can leave permanent stains on the instruments and/or their yellowing.

Sterilization of Instruments: Instruments are distributed non-sterile. Before use they must be sterilized.

We recommend steam sterilization in an Autoclave at the hospital (*ISO 17665-1: 1: Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*).

AUTOCLAVE: It is an equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: To start up the equipment, check: if the circuit breaker is turned on; whether the water supply register is open; if the discharge valve is closed.

1 - open the equipment door; 2 - properly accommodate the material to be sterilized; 3 - close the equipment door, 4 - select the desired cycle according to the material to be sterilized; 5 - turn on the main switch; 6 - the cycle will run automatically, in sequence;

7 - when turning on the “end of cycle” lamp, partially open the door for ten (10) minutes, approximately, to cool the material. It is recommended that the following parameters of physical sterilization be applied in autoclaves (saturated steam):

Table - Parameters of physical sterilization in autoclaves

| Cycle | Temperature | Exposure Time (minimum) | Dry Time (maximum) |
|------------------------------|---------------|-------------------------|--------------------|
| Conventional (1atm pressure) | 121°C (250°F) | 30 minutes | - |
| High Vacuum | 134°C (273°F) | 6 minutes | 15 minutes |

Note: The time must be marked when the heat from the sterilization chamber reaches the desired temperature.

For better clarification, consult the instruction manual that comes with each autoclave.

Another sterilization method that can be used in addition to the autoclave is **Ethylene Oxide Sterilization (E.T.O.)** - parameters and procedures established in the validation protocol and ISO 11135 - Sterilization of health care products - Ethylene oxide - Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: It is the responsibility of the hospital, the sterilization method, equipment, controls, and the sterilization instructions used. Care for sterile items.

Storage conditions for sterile items: - as for the environment: must be cleaned; airy; dry; must be restricted to the sector team; - as for the item: after the sterilization process, do not place it on a cold surface (stone or stainless steel), use empty baskets or containers until they cool; - wrapping (raw cotton fabric, non-woven fabric, surgical grade paper, creped paper, paper with film, tyvec or perforated metal boxes) must remain intact and be handled poorly to prevent packages from tearing or loosening the seal; - be stored in closed cabinets with shelves; - shelves identified in order to facilitate the removal of material; - material must be stored according to the expiration date of the sterilization to facilitate the distribution and not keep material expired in the stock; - store separately from non-sterile ones to reduce the level of external contaminants.

Storage: Non-sterile surgical instruments must be stored in such a way that the packaging is not affected. Packages should not be compressed, crushed, punctured or exposed to possible water damage. The storage area must not contain impurities, dust or pollutants of any kind and must not be exposed to extremes of temperature or humidity. Traffic in the area must be minimal. The instruments must be stored on trolleys or shelves, away from the floor, ceiling or walls. Store it in an airy, dry place, away from light and away from the action of weather. When transporting it, shocks and inadequate stackings should be avoided. Store and transport it in a dry and cool place, with room temperature (Max. 35°C) and relative humidity around 30% to 70%. Do not store it directly on the floor (minimum height = 20cm) or in very high places, near lamps, which could cause dryness of the packaging or damage to the label. Do not store it in places where contaminants such as cleaning materials, insecticides, pesticides, etc. are stored.


Identification: The Instruments follow with labels containing the information: Manufacturer's name; Product Commercial Name; Technical Name; Batch number; ANVISA registration number; Product code; Amount; Description of the product in the package (model); Manufacturing date; Expiration date; Description of the raw material used to manufacture the product (Composition); Manufacturer's address; Name of Technical Responsible; The words: "Non-sterile product"; Before using, see instructions for use; Warnings/Precautions/Special care/Symbology: See Instructions for Use.

Customer Service: If there is a need to make any complaints regarding the use of instrumentals related to any adverse effect affecting the safety of the user the surgeon responsible must communicate this adverse event to the competent health agency and Biomecânica through sac@biomecanica.com.br email or by phone 0xx14 2104 7900. In case of doubt, the surgeon responsible or the health professional may communicate the adverse event through the Health Surveillance Notification System at the ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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