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Polymer Blunt Non-Hinged Surgical Instruments - BM

RG.: 80128580123

**ENGLISH**  
**EUA**

**Table 1: Listing of Polymer Blunt Non-Hinged Surgical Instruments - BM**



Code	Description	Level (hip, knee, shoulder, etc.)	Raw Material
5112-28-C00	Test Interchangeable Head (Short) Ø28 mm Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5112-28-M00	Test Interchangeable Head (Medium) Ø28 mm Cone 12/14		
5112-28-L00	Test Interchangeable Head (Long) Ø28 mm Cone 12/14		
5390-08-000	Extra Small Test Plateau 8 mm	Knee	Polyphenylsulfone (PPSU)
5390-10-000	Extra Small Test Plateau 10 mm		
5390-12-000	Extra Small Test Plateau 12 mm		
5390-15-000	Extra Small Test Plateau 15 mm		
5374-08-000	Small / Medium Test Plateau 08 mm - PL	Knee	Polyphenylsulfone (PPSU)
5374-10-000	Small / Medium Test Plateau 10 mm - PL		
5374-12-000	Small / Medium Test Plateau 12 mm - PL		
5374-15-000	Small / Medium Test Plateau 15 mm - PL		
5375-08-000	Large / Extra Large Test Plateau 08 mm - PL		
5375-10-000	Large / Extra Large Test Plateau 10 mm - PL		
5375-12-000	Large / Extra Large Test Plateau 12 mm - PL		
5375-15-000	Large / Extra Large Test Plateau 15 mm - PL		
5396-08-000	Extra Small Test Plateau 08 mm - PL	Knee	Polyphenylsulfone (PPSU)
5396-10-000	Extra Small Test Plateau 10 mm - PL		
5396-12-000	Extra Small Test Plateau 12 mm - PL		
5396-15-000	Extra Small Test Plateau 15 mm - PL		
5345-15-000	Small / Medium Test Plateau 15 mm	Knee	Polyphenylsulfone (PPSU)
5345-12-000	Small / Medium Test Plateau 12 mm		
5345-10-000	Small / Medium Test Plateau 10 mm		
5345-08-000	Small / Medium Test Plateau 08 mm		
5346-15-000	Large / Extra Large Test Plateau 15 mm		
5346-12-000	Large / Extra Large Test Plateau 12 mm		
5346-10-000	Large / Extra Large Test Plateau 10 mm		
5346-08-000	Large / Extra Large Test Plateau 08 mm		
5353-26-000	Test Patella Ø 26 mm (3PEG'S)	Knee	Polyphenylsulfone (PPSU)
5353-28-000	Test Patella Ø 28 mm (3PEG'S)		
5353-30-000	Test Patella Ø 30 mm (3PEG'S)		
5353-32-000	Test Patella Ø 32 mm (3PEG'S)		
5353-34-000	Test Patella Ø 34 mm (3PEG'S)		
5353-36-000	Test Patella Ø 36 mm (3PEG'S)		
5345-17-000	BKS III - Small / Medium Test Plateau 17mm	Knee	Polyphenylsulfone (PPSU)
5346-17-000	BKS III - Large / Extra Large Test Plateau 17mm	Knee	Polyphenylsulfone (PPSU)

5374-17-000	BKS III - Small / Medium Test Plateau 17mm - PL	Knee	Polyphenylsulfone (PPSU)
5375-17-000	BKS III - Large / Extra Large Test Plateau 17mm - PL	Knee	Polyphenylsulfone (PPSU)
5390-17-000	BKS III - Extra Small Test Plateau 17mm	Knee	Polyphenylsulfone (PPSU)
5396-17-000	BKS III - Extra Small Test Plateau 17mm - PL	Knee	Polyphenylsulfone (PPSU)
5112-28-EL0	Test Interchangeable Head - EX. Long - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5139-28-C00	Test Interchangeable Head - Short - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5139-28-L00	Test Interchangeable Head - Long - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5139-28-M00	Test Interchangeable Head - Medium - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5139-32-C00	Test Interchangeable Head - Short - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5139-32-L00	Test Interchangeable Head - Long - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5139-32-M00	Test Interchangeable Head - Medium - Ø28mm - Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5606-00-000	Interchangeable Head Test Adapter - Biomec III	Hip	Polyphenylsulfone (PPSU)
5620-00-000	Interchangeable Head Test Adapter - CP3 Constricted	Hip	Polyphenylsulfone (PPSU)
5112 -28-ELP	Test interchangeable head (extralong plus)	Hip	Polyphenylsulfone (PPSU)
5409-10-000	Extra Small Test Plateau - 10	Knee	Polyphenylsulfone (PPSU)
5409-12-000	Extra Small Test Plateau - 12	Knee	Polyphenylsulfone (PPSU)
5409-14-000	Extra Small Test Plateau - 14	Knee	Polyphenylsulfone (PPSU)
5409-17-000	Extra Small Test Plateau - 17	Knee	Polyphenylsulfone (PPSU)
5409-20-000	Extra Small Test Plateau - 20	Knee	Polyphenylsulfone (PPSU)
5409-23-000	Extra Small Test Plateau - 23	Knee	Polyphenylsulfone (PPSU)
5410-10-000	Small / Medium Test Plateau - 10	Knee	Polyphenylsulfone (PPSU)
5410-12-000	Small / Medium Test Plateau - 12	Knee	Polyphenylsulfone (PPSU)
5410-14-000	Small / Medium Test Plateau - 14	Knee	Polyphenylsulfone (PPSU)
5410-17-000	Small / Medium Test Plateau - 17	Knee	Polyphenylsulfone (PPSU)
5410-20-000	Small / Medium Test Plateau - 20	Knee	Polyphenylsulfone (PPSU)
5410-23-000	Small / Medium Test Plateau - 23	Knee	Polyphenylsulfone (PPSU)
5411-10-000	Large / Extra Large Test Plateau - 10	Knee	Polyphenylsulfone (PPSU)
5411-12-000	Large / Extra Large Test Plateau - 12	Knee	Polyphenylsulfone (PPSU)
5411-14-000	Large / Extra Large Test Plateau - 14	Knee	Polyphenylsulfone (PPSU)
5411-17-000	Large / Extra Large Test Plateau - 17	Knee	Polyphenylsulfone (PPSU)
5411-20-000	Large / Extra Large Test Plateau - 20	Knee	Polyphenylsulfone (PPSU)
5411-23-000	Large / Extra Large Test Plateau - 23	Knee	Polyphenylsulfone (PPSU)
5475-00-000	Fracture Junction Tweezers Extension	Knee	Polyphenylsulfone (PPSU)
6413-00-000	Sliding Screw Ruler	Trauma	Polyphenylsulfone (PPSU)
6413-01-000	Proximal Wire Ruler	Trauma	Polyphenylsulfone (PPSU)
5112-28-C01	Test Interchangeable Head (Short) Ø28 mm Cone 12/14	Hip	Polyphenylsulfone (PPSU)
5112-28-M01	Test Interchangeable Head (Medium) Ø28 mm Cone 12/14		Polyphenylsulfone (PPSU)

5112-28-L01	Test Interchangeable Head (Long) Ø28 mm Cone 12/14		Polyphenylsulfone (PPSU)
5112-28-EL1	Test Interchangeable Head - EX. Long - Ø28mm - Cone 12/14		Polyphenylsulfone (PPSU)
5390-08-001	Extra Small Test Plateau 8 mm (Revision)	Knee	Polyphenylsulfone (PPSU)
5390-10-001	Extra Small Test Plateau 10 mm (Revision)		Polyphenylsulfone (PPSU)
5390-12-001	Extra Small Test Plateau 12 mm (Revision)		Polyphenylsulfone (PPSU)
5390-15-001	Extra Small Test Plateau 15 mm (Revision)		Polyphenylsulfone (PPSU)
5345-15-001	Small / Medium Test Plateau 15 mm (Revision)	Knee	Polyphenylsulfone (PPSU)
5345-12-001	Small / Medium Test Plateau 12 mm (Revision)		Polyphenylsulfone (PPSU)
5345-10-001	Small / Medium Test Plateau 10 mm (Revision)		Polyphenylsulfone (PPSU)
5345-08-001	Small / Medium Test Plateau 08 mm (Revision)		Polyphenylsulfone (PPSU)
5346-15-001	Large / Extra Large Test Plateau 15 mm (Revision)		Polyphenylsulfone (PPSU)
5346-12-001	Large / Extra Large Test Plateau 12 mm (Revision)		Polyphenylsulfone (PPSU)
5346-10-001	Large / Extra Large Test Plateau 10 mm (Revision)		Polyphenylsulfone (PPSU)
5346-08-001	Large / Extra Large Test Plateau 08 mm (Revision)		Polyphenylsulfone (PPSU)

### Traceability Label Symbols (Label) as per NBR ISO 15223

<b>Manufacturing Date</b>	<b>Valid Date</b>	<b>Batch No.</b>	<b>Manufacturer Data</b>	<b>Code</b>	<b>Non-Sterile</b>
					

### TECHNICAL INFORMATION

**Product Description:** Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Sharp components can only be used for their specific purposes. Misuse can cause blade and edge cutting loss, and impact their performance within surgical procedure. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

**Manufacturing materials:** Surgical instruments are made of polymer materials like Polyphenylsulfone (PPSU).

**Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications:** Polymer Blunt Non-Hinged Surgical Instruments - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Polymer Blunt Non-Hinged Surgical Instruments - BM are reused in another patient, it is required that it is cleaned, sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

**Cleaning:** Appropriate care with surgical instruments starts with appropriate cleaning. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

**Instrument Sterilization:** Instruments are supplied non-sterile. They have to be sterilized before use.

**We recommend autoclave steam sterilization at the hospital** (ISO 17665-1: 2006  
*Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

**AUTOCLAVE:** It is a piece of equipment for sterilization through saturated steam under pressure.

**OPERATING INSTRUCTIONS:** For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

**Table 2 - Autoclave physical sterilization parameters**

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

**For further clarification, refer to instruction manual supplied with autoclave.**

**Another sterilization method:** that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

**NOTE:** The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

**Storage:** Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls. Keep the in ventilated, dry place, away from light and severe weather action. As they are transported, inadequate bumps and stacking must be avoided. Store and transport in dry and ventilated place, with room temperature (max. 35°C) and relative humidity around 30% to 70%. Do not store directly on the floor (minimum height = 20 cm) or in very high places, close to lamps, as the packaging might get dried or the label might be damaged. Do not store in places where contaminant substances are stored, including, for instance, cleaning materials, insecticides, pesticides, etc.

**Identification:** Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product' ;'Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

**Consumer Service:** In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email [tecnovigilancia@biomecanica.com.br](mailto:tecnovigilancia@biomecanica.com.br) or telephone 0xx14 2104 7906. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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