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Blunt Non-Hinged Surgical Instruments (Aluminium)

RG.: 8012858127

**ENGLISH
EUA**

Table: Listing of Blunt Non-Hinged Surgical Instruments (Aluminum) - BM

Code	Description	Level (hip, knee, shoulder, etc.)	Raw Material	Indication for Use
6076-00-000	ONE-STEP Insertion Wrench Centering Guide	Trauma	Aluminum (NBR ISO 209:10)	Implant Introduction
6304-00-000	Guide Intramedullary Rod SP2 (Femur)	Trauma	Aluminum (NBR ISO 209:10)	Femur proximal screw drilling and placement guide
			Stainless Steel AISI 630 (NBR 13911)	
5620-38-000	Measuring Acetabular Comp. Ø38	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5620-39-000	Measuring Acetabular Comp. Ø39			
5620-60-000	Measuring Acetabular Comp. Ø60			
5620-62-000	Measuring Acetabular Comp. Ø62			
5620-64-000	Measuring Acetabular Comp. Ø64			
5620-66-000	Measuring Acetabular Comp. Ø66			
5620-40-000	Measuring Acetabular Component Ø 40	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5620-42-000	Measuring Acetabular Component Ø 42			
5620-44-000	Measuring Acetabular Component Ø 44			
5620-46-000	Measuring Acetabular Component Ø 46			
5620-48-000	Measuring Acetabular Component Ø 48			
5620-50-000	Measuring Acetabular Component Ø 50			
5620-52-000	Measuring Acetabular Component Ø 52			
5620-54-000	Measuring Acetabular Component Ø 54			
5620-56-000	Measuring Acetabular Component Ø 56			

5620-58-000	Measuring Acetabular Component Ø 58			
5606-66-000	Measuring Acetabular Comp. - Ø 66.0mm - Press-Fit	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5606-42-000	Measuring Acetabular Comp. Ø 42.0mm - Press-Fit	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5606-44-000	Measuring Acetabular Comp. Ø 44.0mm - Press-Fit	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5606-46-000	Measuring Acetabular Component Ø46 - Press-Fit	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5606-48-000	Measuring Acetabular Component Ø48 - Press-Fit			
5606-50-000	Measuring Acetabular Component Ø50 - Press-Fit			
5606-52-000	Measuring Acetabular Component Ø52 - Press-Fit			
5606-54-000	Measuring Acetabular Component Ø54 - Press-Fit			
5606-56-000	Measuring Acetabular Component Ø56 - Press-Fit			
5606-58-000	Measuring Acetabular Component Ø58 - Press-Fit			
5606-60-000	Measuring Acetabular Component Ø60 - Press-Fit			
5606-62-000	Measuring Acetabular Component Ø62 - Press-Fit			
5606-64-000	Measuring Acetabular Component Ø64 - Press-Fit			
5606-66-000	Measuring Acetabular Component Ø66 - Press-Fit			
5168-48-000	Measuring Acetabular Component 48 mm			
5168-46-000	Measuring Acetabular Component 46 mm			
5168-50-000	Measuring Acetabular Component 50 mm			
5168-52-000	Measuring Acetabular Component 52 mm			
5168-54-000	Measuring Acetabular Component 54 mm			
5168-56-000	Measuring Acetabular Component 56 mm			
5168-58-000	Measuring Acetabular Component 58 mm			
5168-60-000	Measuring Acetabular Component 60 mm			
5168-62-000	Measuring Acetabular Component 62 mm			
5170-54-000	Test Bipolar Head 54 mm	Hip	Aluminum (NBR ISO 209:10)	Check acetabulum size to select implant size
5170-50-000	Test Bipolar Head 50 mm			
5170-48-000	Test Bipolar Head 48 mm			
5170-44-000	Test Bipolar Head 44 mm			

5170-40-000	Test Bipolar Head 40 mm			
5170-42-000	Test Bipolar Head 42 mm			
5170-46-000	Test Bipolar Head 46 mm			
5170-52-000	Test Bipolar Head 52 mm			
5170-39-000	Test Bipolar Head 39 mm			
5170-41-000	Test Bipolar Head 41 mm			
5170-43-000	Test Bipolar Head 43 mm			
5170-45-000	Test Bipolar Head 45 mm			
5170-47-000	Test Bipolar Head 47 mm			
5170-49-000	Test Bipolar Head 49 mm			
5170-51-000	Test Bipolar Head 51 mm			
5170-53-000	Test Bipolar Head 53 mm			
5170-55-000	Test Bipolar Head - Ø 55mm	Hip	Aluminum (NBR ISO 209:10)	Check acetabulum size
5170-56-000	Test Bipolar Head - Ø 56mm	Hip	Aluminum (NBR ISO 209:10)	Check acetabulum size
8529-00-100	Moldable Bar - 100	Column	Aluminum (NBR ISO 209:10)	Check vertebra alignment
8529-00-200	Moldable Bar - 200			
8529-00-300	Moldable Bar - 300			
5342-08-000	Spacer Block 8.0 mm	Knee	Aluminum (NBR ISO 209:10)	Space Femur Tibia
5342-10-000	Spacer Block 10.0 mm			
5342-12-000	Spacer Block 12.0 mm			
5342-15-000	Spacer Block 15.0 mm			
5621-38-000	Acetabular Meter Ø 38 - (Moore/Thompson)	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5621-40-000	Acetabular Meter Ø 40 - (Moore/Thompson)			
5621-42-000	Acetabular Meter Ø 42 - (Moore/Thompson)			
5621-44-000	Acetabular Meter Ø 44 - (Moore/Thompson)			
5621-46-000	Acetabular Meter Ø 46 - (Moore/Thompson)			
5621-48-000	Acetabular Meter Ø 48 - (Moore/Thompson)			
5621-50-000	Acetabular Meter Ø 50 - (Moore/Thompson)			


5621-52-000	Acetabular Meter Ø 52 - (Moore/Thompson)			
5621-54-000	Acetabular Meter Ø 54 - (Moore/Thompson)			
5621-56-000	Acetabular Meter Ø 56 - (Moore/Thompson)			
5621-39-000	Acetabular Meter Ø 39 (Moore/Thompson)			
5621-41-000	Acetabular Meter Ø 41 (Moore/Thompson)			
5621-43-000	Acetabular Meter Ø 43 (Moore/Thompson)			
5621-47-000	Acetabular Meter Ø 47 (Moore/Thompson)			
5621-49-000	Acetabular Meter Ø 49 (Moore/Thompson)			
5621-53-000	Acetabular Meter Ø 53 (Moore/Thompson)			
5621-55-000	Acetabular Meter Ø 55 (Moore/Thompson)			
5621-45-000	Acetabular Meter - Ø 45 mm (Moore/Thompson)	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5621-51-000	Acetabular Meter - Ø 51 mm (Moore/Thompson)	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5191-40-000	Test Humeral Head - Ø 40 mm	Shoulder	Aluminum (NBR ISO 209:10)	Test implant instrumentation
5191-43-000	Test Humeral Head - Ø 43 mm			
5191-45-000	Test Humeral Head - Ø 45 mm			
2306-40-000	Test Offset Humeral Head Ø 40 mm	Shoulder	Aluminum (NBR ISO 209:10)	Test implant instrumentation
2306-43-000	Test Offset Humeral Head Ø 43 mm			
2306-45-000	Test Offset Humeral Head Ø 45 mm			
6500-01-024	Reconstruction plate template	Trauma	Aluminum (NBR ISO 209:10)	Instrument used as implant modeling reference
6500-02-021	Template	Trauma	Aluminum (NBR ISO 209:10)	Instrument used as implant modeling reference
6500-08-002	Pin meter	Trauma	Aluminum (NBR ISO 209:10)	Measurement instrument
6500-09-008	Pin meter	Trauma	Aluminum (NBR ISO 209:10)	Screw placement bone measurement instrument
6500-11-012	Diameter meter	Trauma	Aluminum (NBR ISO 209:10)	Measurement instrument
6500-18-007	Screw Meter	Trauma	Aluminum (NBR ISO 209:10)	Measure Screw
6500-05-017	Reconstruction plate template	Pelvis	Aluminum (NBR ISO 209:10)	Instrument used as implant modeling reference.
6059-00-000	Wire meter	Trauma	Aluminum (NBR ISO 209:10)	Measurement instrument
6059-01-000	Calibrated wire meter			

8529-10-000	Moldable Bar 10	Column	Aluminum (NBR ISO 209:10)	Check vertebra alignment.
8529-20-000	Moldable Bar 20			
8529-30-000	Moldable Bar 30			
5168-36-000	Measuring Acetabular Component 36mm	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5168-38-000	Measuring Acetabular Component 38mm			
5168-39-000	Measuring Acetabular Component 39mm			
5168-40-000	Measuring Acetabular Component 40mm			
5168-41-000	Measuring Acetabular Component 41mm			
5168-43-000	Measuring Acetabular Component 43mm			
5168-44-000	Measuring Acetabular Component 44mm			
5168-45-000	Measuring Acetabular Component 45mm			
5168-47-000	Measuring Acetabular Component 47mm			
5168-49-000	Measuring Acetabular Component 49mm			
5168-51-000	Measuring Acetabular Component 51mm			
5168-53-000	Measuring Acetabular Component 53mm			
5168-55-000	Measuring Acetabular Component 55mm			
5606-44-000	Measuring Acetabular Component 44.0mm - Press-Fit	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5620-41-000	Measuring Acetabular Component \varnothing 41	Hip	Aluminum (NBR ISO 209:10)	Measure acetabulum
5620-43-000	Measuring Acetabular Component \varnothing 43			
5620-45-000	Measuring Acetabular Component \varnothing 45			
5620-47-000	Measuring Acetabular Component \varnothing 47			
5620-49-000	Measuring Acetabular Component \varnothing 49			
5620-51-000	Measuring Acetabular Component \varnothing 51			
5620-53-000	Measuring Acetabular Component \varnothing 53			
5620-55-000	Measuring Acetabular Component \varnothing 55			
5620-57-000	Measuring Acetabular Component \varnothing 57			
5620-59-000	Measuring Acetabular Component \varnothing 59			
5620-61-000	Measuring Acetabular Component \varnothing 61			

4000-00-049	Conformation template	Maxilofacial Skull	Aluminum (NBR ISO 209:10)	Conform plates.
4000-00-050	Conformation template	Maxilofacial Skull	Aluminum (NBR ISO 209:10)	Conform plates.
4002-00-011	Humeral Head Proof - Small	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-012	Humeral Head Proof - Medium	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-013	Humeral Head Proof - Large	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-075	Glenoidal Component Proof - 19.0 mm	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-076	Glenoidal Component Proof - 22.0 mm	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-077	Glenoidal Component Proof - 24.0 mm	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-078	Glenoidal Component Proof - 26.0 mm	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
4002-00-079	Glenoidal Component Proof - 30.0 mm	Trauma	Aluminum (NBR ISO 209:10)	Motion size and amplitude test
6500-17-022	Moldable Bar (template)	Column	Aluminum (NBR ISO 209:10)	Distraction rod modeling test instrument.
6500-18-002	Drill Bit Handle	Trauma	Aluminum (NBR ISO 209:10)	Drill bit use assist.
5342-09-000	Spacer Block 9.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-11-000	Spacer Block 11.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-13-000	Spacer Block 13.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-14-000	Spacer Block 14.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-16-000	Spacer Block 16.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-17-000	Spacer Block 17.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-18-000	Spacer Block 18.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-19-000	Spacer Block 19.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-20-000	Spacer Block 20.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-21-000	Spacer Block 21.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-22-000	Spacer Block 22.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5342-23-000	Spacer Block 23.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment

5342-24-000	Spacer Block 24.0mm	Knee	Aluminum (NBR ISO 209:10)	Used to define tibial insert thickness and assess varus and valgus alignment
5474-00-001	Final Tightening Protection Tube	Knee	Aluminum (NBR ISO 209:10)	Used as base to tighten tibial component and stabilization rods (implant)
5123-00-003	Head Impactor	Hip	Aluminum (NBR ISO 209:10)/Polypropylene	Impact head in femoral rod

Traceability Label Symbols (Label) as per NBR ISO 15223

Manufacturing Date	Valid Date	Batch No.	Manufacturer Data	Code	Non-Sterile
					

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

Manufacturing materials: Surgical instruments are made of aluminum, as provided by NBR ISO 209:10; that characterize their chemical composition, mechanical properties and final behavior of aluminum, to produce surgical instruments. Some components are made of stainless steel and polypropylene.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Blunt Non-Hinged Surgical Instruments (Aluminum) - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Blunt Non-Hinged Surgical Instruments (Aluminum) - BM are reused in another patient, it is required that it is cleaned,

sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

Cleaning: Appropriate care with surgical instruments starts with appropriate cleaning. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

Instrument Sterilization: Instruments are supplied non-sterile. They have to be sterilized before use.

We recommend autoclave steam sterilization at the hospital (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a piece of equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general

switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

Table 1 - Autoclave physical sterilization parameters

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

For further clarification, refer to instruction manual supplied with autoclave.

Another sterilization method: that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

Storage: Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls. Keep the in ventilated, dry place, away from light and severe weather action. As they are transported, inadequate bumps and stacking must be avoided. Store and transport in dry and ventilated place, with room temperature (max. 35°C) and relative humidity around 30% to 70%. Do not store directly on the floor (minimum height = 20 cm) or in very high places, close to lamps, as the packaging might get dried or the label might be damaged. Do not store in places where contaminant substances are stored, including, for instance, cleaning materials, insecticides, pesticides, etc.

Identification: Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product' ;'Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

Consumer Service: In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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