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Blunt Non-Hinged Surgical Instruments – Generation 3 - BM

RG.: 80128580147

ENGLISH
EUA

Table: Listing of Blunt Non-Hinged Surgical Instruments - Generation 3 - BM

CODE	DESCRIPTION	Level	Size (mm)	Raw Material	Indication for Use
6501-00-000	ORION SP ROD EXTRACTOR	Trauma	370	STAINLESS STEEL AISI 420 (NBR 13911)	Extract rod
6502-35-000	HEXAGON WRENCH 3.5 MM - QUICK-COUPLING	Trauma	230	STAINLESS STEEL AISI 420 (NBR 13911)	Tighten lock screws
6502-50-000	HEXAGON WRENCH 3.5 MM - QUICK-COUPLING	Trauma	230	STAINLESS STEEL AISI 420 (NBR 13911)	Tighten lock screws
6503-00-000	EXTENSION WRENCH	Trauma	60	STAINLESS STEEL AISI 420 (NBR 13911)	Block implant rod in insertion guide
6510-00-000	DRILL BIT SLEEVE Ø6.5	Trauma	180	STAINLESS STEEL AISI 304 (NBR 13911)	Used as guide to drill bit sleeves and lock screw insertion
6504-00-000	LOCK SCREW SLEEVE	Trauma	190	STAINLESS STEEL AISI 304 (NBR 13911)	Used as guide to drill bit sleeves and lock screw insertion
6505-00-000	DRILL BIT SLEEVE Ø3.7 mm	Trauma	180	STAINLESS STEEL AISI 304 (NBR 13911)	Used to drill bits Ø3.7
6506-00-000	DRILL BIT SLEEVE Ø4.7 mm	Trauma	180	STAINLESS STEEL AISI 304 (NBR 13911)	Used to drill bits Ø4.7
6523-00-000	NAVIGATION ASSIST BUSHING	Trauma	30	STAINLESS STEEL AISI 304 (NBR 13911)	Used with distal drilling guides to identify implant distal drilling
6508-12-000	FEMORAL PROXIMAL DRILLING CUTTER GUIDE	Trauma	185	STAINLESS STEEL AISI 420 (NBR 13911) W/ RUBBER HANDLE	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
6508-11-000	TIBIAL PROXIMAL DRILLING CUTTER GUIDE	Trauma	185	STAINLESS STEEL AISI 420 (NBR 13911) W/ RUBBER HANDLE	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
6514-00-000	FEMORAL INITIAL DRILLING SLEEVE	Trauma	82	STAINLESS STEEL AISI 420 (NBR 13911)	Used with proximal milling cutter guide
6515-00-000	TIBIAL INITIAL DRILLING SLEEVE	Trauma	82	STAINLESS STEEL AISI 420 (NBR 13911)	Used with proximal milling cutter guide
6509-00-000	ROD IMPACTOR	Trauma	80	STAINLESS STEEL 630	Used to connect in implant insertion guide and receive impactions
6535-00-000	ORION SP ELITE GUIDE INTRAMEDULLARY ROD (FEMUR)	Trauma	225	CARBON FIBER, STAINLESS STEEL V630, STAINLESS STEEL; AISI 304 (NBR 13911)	Guides intramedullary rod and lock screws
6540-00-000	ORION SP ELITE GUIDE INTRAMEDULLARY ROD (TIBIA)	Trauma	255	CARBON FIBER, STAINLESS STEEL V630, STAINLESS STEEL; AISI 304 (NBR 13911)	Guides intramedullary rod and lock screws

6061-00-001	Angular Guide 135°	Trauma	93	STAINLESS STEEL AISI 304 (NBR 13911)	Bone Drilling Guide
6517-00-000	FEMORAL INTRAMEDULLARY ROD TEMPLATE	Trauma	500	STAINLESS STEEL AISI 304 (NBR 13911)	Used to indicate rod dimension to be implanted
6480-00-000	ROD EXTRACTOR ADAPTER	Trauma	88	STAINLESS STEEL AISI 420 (NBR 13911)	Used to extract the rod
6480-00-001	ROD EXTRACTOR ADAPTER	Trauma	88	STAINLESS STEEL AISI 304 (NBR 13911)	Used to extract the rod
6518-00-001	DISTAL SCREW METER	Trauma	160	STAINLESS STEEL AISI 304 (NBR 13911)	Used to size lock screw length
6516-00-000	ROD HEIGHT INDICATOR SLEEVE	Trauma	45	STAINLESS STEEL AISI 304 (NBR 13911)	Size lock screw length
6519-00-000	REDUCTION ASSIST WRENCH	Trauma	500	STAINLESS STEEL 301 STAINLESS STEEL 304 (NBR 13911)	Used to enable guide wire change
6533-00-000	FEMORAL GUIDE EXTENSION	Trauma	120	STAINLESS STEEL 630 (NBR 13911)	Used to fasten implant rods in drilling guide
6534-00-000	TIBIAL GUIDE EXTENSION	Trauma	100	STAINLESS STEEL 630 (NBR 13911)	Used to fasten implant rods in drilling guide
6541-00-000	ORION EXTENSION	Trauma	80	STAINLESS STEEL 630 (NBR 13911)	Used to fasten implant rods in drilling guide
6537-00-000	POINTER LOCK WRENCH HANDLE	Trauma	100	STAINLESS STEEL 304 (NBR 13911) POLYACETAL	Used to hold the pointers
6538-09-000	FEMORAL POINTER 9.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6538-10-000	FEMORAL POINTER 10.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6538-11-000	FEMORAL POINTER 11.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6538-12-000	FEMORAL POINTER 12.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6539-09-000	TIBIAL POINTER 9.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6539-10-000	TIBIAL POINTER 10.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6539-11-000	TIBIAL POINTER 11.0	Trauma	155	STAINLESS STEEL AISI 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6413-02-000	PROXIMAL SCREW RULER	TRAUMA	182 x 35	STAINLESS STEEL 304 (NBR 13911)	Used to size humerus proximal screw length
6540-00-001	ORION SP ELITE GUIDE INTRAMEDULLARY ROD (TIBIA)	TRAUMA	250x150	STAINLESS STEEL V630, STAINLESS STEEL AISI 304, ALUMINUM (NBR 13911)	Guides intramedullary rod and lock screws
6508-21-000	TIBIAL SOFT PART GUARD	TRAUMA	146 x 175	STAINLESS STEEL AISI 420 (NBR 13911) W/ RUBBER HANDLE	Used to protect sharp instrument soft parts
6508-22-000	FEMORAL SOFT PART GUARD	TRAUMA	136 x 175	STAINLESS STEEL AISI 420 (NBR 13911)	Used to protect sharp instrument soft parts

				W/ RUBBER HANDLE	
6542-00-000	FLEXIBLE WRENCH 3.5mm	TRAUMA	310 x 12	STAINLESS STEEL 420, STAINLESS STEEL 304, (NBR 13911)	Used to guide screws and lock
6542-10-000	FLEXIBLE WRENCH CANNULATED 3.5mm	TRAUMA	310 x 12	STAINLESS STEEL 420, STAINLESS STEEL 304 (NBR 13911)	Used to guide screws and lock
6502-40-000	HEXAGON WRENCH 4.0 MM - QUICK-COUPLING	Trauma	228 x 9.5	STAINLESS STEEL AISI 420 (NBR 13911)	Used to tighten lock screws
6538-09-100	FEMORAL POINTER 9.0 - CANNULATED	Trauma	200 x 9.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6538-10-100	FEMORAL POINTER 10.0 - CANNULATED	Trauma	200 x 10.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6538-11-100	FEMORAL POINTER 11.0 - CANNULATED	Trauma	200 x 11.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6538-12-100	FEMORAL POINTER 12.0 - CANNULATED	Trauma	200 x 12.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6539-09-100	TIBIAL POINTER 9.0 - CANNULATED	Trauma	170 x 9.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6539-10-100	TIBIAL POINTER 10.0 - CANNULATED	Trauma	170 x 10.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6539-11-100	TIBIAL POINTER 11.0 - CANNULATED	Trauma	170 x 11.0	STAINLESS STEEL 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
6039-15-000	DRILL BIT GUIDE Ø1.5	Trauma	90 x 80	STAINLESS STEEL 420, STAINLESS STEEL 304, (NBR 13911)	Used to protect soft parts and guide drill bit
6517-01-000	TIBIAL INTRAMEDULLARY ROD TEMPLATE	Trauma	500x25	STAINLESS STEEL AISI 304 (NBR 13911)	Used to indicate rod dimension to be implanted
6517-02-000	PROXIMAL INTRAMEDULLARY ROD TEMPLATE	Trauma	500x25	STAINLESS STEEL AISI 304 (NBR 13911)	Used to indicate rod dimension to be implanted
6517-03-000	HUMERAL INTRAMEDULLARY ROD TEMPLATE	Trauma	500x25	STAINLESS STEEL AISI 304 (NBR 13911)	Used to indicate rod dimension to be implanted
6519-01-000	HUMERAL REDUCTION ASSIST WRENCH	Trauma	500 x 120	STAINLESS STEEL 301 STAINLESS STEEL 304 (NBR 13911)	Used to enable guide wire change

Traceability Label Symbols (Label) as per NBR ISO 15223


Manufacturing Date

Valid Date


Batch No.

Manufacturer Data

Code

Non-Sterile

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be

grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Sharp components can only be used for their specific purposes. Misuse can cause blade and edge cutting loss, and impact their performance within surgical procedure. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

Manufacturing materials: Most surgical instruments, in their structural part, are made of stainless steels, as provided by NBR 13911:2013 or ASTM F899:2012; that characterize their chemical composition, mechanical properties and final behavior of stainless steel, to produce surgical instruments. The other parts are made of high strength polymers, that support sterilization.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Blunt Non-Hinged Surgical Instruments - Generation 3 - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Blunt Non-Hinged Surgical Instruments - Generation 3 - BM are reused in another patient, it is required that it is cleaned, sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

Cleaning: Appropriate care with surgical instruments starts with appropriate cleaning. This report mentions some instrument cleaning and sterilization methods. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies,

including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

Instrument Sterilization: Instruments are supplied non-sterile. They have to be sterilized before use.

We recommend autoclave steam sterilization at the hospital (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a piece of equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

Table - Autoclave physical sterilization parameters

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

For further clarification, refer to instruction manual supplied with autoclave.

Another sterilization method: that can be used apart from autoclave is **Ethylene Oxide**

Sterilization (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

Storage: Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls.

Identification: Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product' ;'Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

Consumer Service: In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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