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Blunt Non-Hinged Surgical Instruments (Carbon Fiber) – BM

RG.: 80128580159

ENGLISH
EUA

Table 1 - Blunt Non-Hinged Surgical Instrument (Carbon Fiber) Listing - BM

Code	Description	Level	Size (mm)	MP	Indication for Use
6521-00-000	Distal Guide Intramedullary Rod	Trauma	260	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in femur distal part.
6521-10-000	Distal Guide Intramedullary Rod (Retrograde)	Trauma	260	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in femur distal part.
6522-00-000	Femur Distal Guide Intramedullary Rod Decentralizing Device	Trauma	65	Carbon Fiber, Stainless Steel AISI 304	Used to hold femur distal guide intramedullary rod
6525-00-000	Anterior Support Intramedullary Guide Rod	Trauma	260	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in tibia distal part.
6525-10-000	Small Anterior Support Intramedullary Guide Rod	Trauma	200	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in tibia distal part.
6526-00-000	Anterior Support Intramedullary Guide Rod Decentralizing Device	Trauma	65	Carbon Fiber, Stainless Steel AISI 304	Used to holder anterior support intramedullary guide rod.
6529-00-000	Femur proximal distal guide	Trauma	400	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in femur distal part.
6530-00-000	Humerus Distal Guide	Trauma	260	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in tibia distal part.
6531-00-000	Femur proximal distal guide decentralizing device	Trauma	65	Carbon Fiber, Stainless Steel AISI 304	Used to hold femur distal guide intramedullary rod
6532-00-000	Humerus Distal Guide Decentralizing Device	Trauma	65	Carbon Fiber, Stainless Steel AISI 304	Used to hold femur distal guide intramedullary rod
6540-10-000	Tibial Oblique Drilling Proximal Guide	Trauma	255	Carbon Fiber, Stainless Steel V630, 304	It is a guide used to indicate screw position placement in oblique proximal part.

6540-20-000	Tibial Lateral Drilling Proximal Guide	Trauma	105	Carbon Fiber, Stainless Steel V630, 304	It is a guide used to indicate screw position placement in static and dynamic lateral proximal part.
6459-10-000	Humeral Oblique Drilling Proximal Guide	Trauma	140	Carbon Fiber, Stainless Steel V630, 304	It is a guide used to indicate screw position placement in oblique proximal part.
6540-20-001	Tibial Proximal Support	Trauma	100x30x30	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in oblique proximal and medial lateral part
6540-20-00D	Tibial Proximal Right Oblique Support	Trauma	100x30x30	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in oblique proximal and medial lateral part
6540-20-00E	Tibial Proximal Left Oblique Support	Trauma	100x30x30	Carbon Fiber, Stainless Steel AISI 304	It is a guide used to indicate screw position placement in oblique proximal and medial lateral part
6525-00-001	Tibial Distal Lateral Medial Support	Trauma	460 x 90	Carbon Fiber, Stainless Steel AISI 304	Used to support distal guide
6525-00-00D	Tibial Distal Lateral Medial Support - Right	Trauma	460 x 90	Carbon Fiber, Stainless Steel AISI 304	Used to support distal guide
6525-00-00E	Tibial Distal Lateral Medial Support - Left	Trauma	460 x 90	Carbon Fiber, Stainless Steel AISI 304	Used to support distal guide
6540-10-001	Tibial Distal Drilling Support	Trauma	306 x 164	Carbon Fiber, Stainless Steel AISI 304, Aluminum	It is a guide used to indicate screw position placement in lateral medial and distal antero-posterior part
6540-30-000	Tibial Distal Drilling Support	Trauma	167 x 189	Carbon Fiber, Stainless Steel AISI 304, Aluminum	It is a guide used to indicate screw position placement in lateral medial and distal antero-posterior part
6535-00-001	Orion SP Elite Guide Intramedullary Rod (Femur)	Trauma	225 x 158	Carbon Fiber, Stainless Steel V630, 304	Guides intramedullary rod and lock screws
6521-00-001	Orion SP Femur Lateral Support	Trauma	390 x 64	Carbon Fiber, Stainless Steel AISI 304	Used to support distal guide
6521-00-00D	Orion SP Femur Lateral Support - Right	Trauma	390 x 64	Carbon Fiber, Stainless Steel AISI 304	Used to support distal guide
6521-00-00E	Orion SP Femur Lateral Support - Left	Trauma	390 x 64	Carbon Fiber, Stainless Steel AISI 304	Used to support distal guide

6521-20-000	Femoral Distal Drilling Support	Trauma	197 x 189	Carbon Fiber, Stainless Steel AISI 304, Aluminum	It is a guide used to indicate screw position placement in lateral medial and distal antero-posterior part
6428-00-000	Proximal Screw Alignment Support	Trauma	310 x 220	Carbon Fiber, Stainless Steel AISI 304	Used to proximal screw alignment with femoral head
4005-00-000	Soft Part Guard - Bioclocking	Trauma	200 x 70	Carbon Fiber	Protect soft parts.
4005-00-001	Hohmann Small Retractor (Pair) 3.5	Trauma	260 x 30	Carbon Fiber	Retract soft parts.
4005-00-002	Hohmann Medium Retractor (Pair) 3.5	Trauma	280 x 30	Carbon Fiber	Retract soft parts.
4005-00-003	Hohmann Large Retractor (Pair) 4.5	Trauma	300 x 30	Carbon Fiber	Retract soft parts.
4005-00-004	BKS III - Tibial Retractor	Knee	320 x 70	Carbon Fiber	Retract soft parts.
4005-00-005	BKS III - Patellary Tendon Retractor	Knee	320 x 70	Carbon Fiber	Retract soft parts.
4005-00-006	Double and Wide Retractor	Knee	320 x 80	Carbon Fiber	Retract soft parts.
4005-00-007	Bent Retractor	Knee	320 x 60	Carbon Fiber	Retract soft parts.
4005-00-008	Single Retractor	Knee	320 x 40	Carbon Fiber	Retract soft parts.
4005-00-009	Double and Long Retractor	Knee	320 x 50	Carbon Fiber	Retract soft parts.
4005-00-010	5-point U Assist Retractor	Hip	500 x 500 x 60	Carbon Fiber	Retract soft parts.
4005-00-011	Bent Retractor with Guide - 67.5 mm	Hip	220 x 67.5	Carbon Fiber	Retract soft parts.
4005-00-012	Bent Retractor with Guide - 92.5 mm	Hip	240 x 92.5	Carbon Fiber	Retract soft parts.
4005-00-013	Langebeck Retractor	Hip	290 x 80	Carbon Fiber	Retract soft parts.
4005-00-014	Farabeufh Medium Retractor	Hip	220 x 60 x 60	Carbon Fiber	Retract soft parts.
4005-00-015	Hohmann Modified Retractor (Straight)	Hip	280 x 30	Carbon Fiber	Retract soft parts.
4005-00-016	Hohmann Modified Retractor (Bent)	Hip	280 x 70 x 30	Carbon Fiber	Retract soft parts.
4005-00-017	Hohmann Bent Retractor	Hip	280 x 60	Carbon Fiber	Retract soft parts.
4005-00-018	Hohmann Long Bent Retractor	Hip	280 x 40	Carbon Fiber	Retract soft parts.
4005-00-019	Soft Part Guard	Trauma	220 x 70	Carbon Fiber	Retract soft parts.
4005-00-020	Double and Narrow Retractor	Shoulder	280 x 70	Carbon Fiber	Retract soft parts.
4005-00-021	Posterior Femoral Retractor	Knee	260 x 50	Carbon Fiber	Retract soft parts.
4005-00-022	Retractor 2 - Fork Tip	Knee	280 x 70	Carbon Fiber	Retract soft parts.
4005-00-023	Retractor 1 - Single Tip	Knee	280 x 70	Carbon Fiber	Retract soft parts.
4005-00-024	Hofmann Retractor	Hip	280 x 80	Carbon Fiber	Soft Part Retractor

4005-00-025	Bent Retractor with Guide - 117.5mm	Hip	280 x 117.5	Carbon Fiber	Retract soft parts
4005-00-026	Bent Retractor with Guide - 142.5mm	Hip	280 x 142.5	Carbon Fiber	Retract soft parts
4005-00-027	Farabeufh Large Retractor	Hip	220 x 80 x 80	Carbon Fiber	Retract soft parts
4005-00-028	Retractor with Long Tooth - 75 mm	Hip	220 x 75	Carbon Fiber	Soft Part Retractor
4005-00-029	Retractor with Long Tooth - 125 mm	Hip	220 x 125	Carbon Fiber	Soft Part Retractor
4005-00-030	Glenoidal Component Retractor - Small	Trauma	280 x 50	Carbon Fiber	Soft tissue retraction to expose glenoid
4005-00-031	Glenoidal Component Retractor - Large	Trauma	280 x 70	Carbon Fiber	Soft tissue retraction to expose glenoid
4005-00-032	Humeral Retractor - Small	Trauma	280 x 30	Carbon Fiber/Silicone	Soft tissue retraction to expose humerus
4005-00-033	Humeral Retractor - Large	Trauma	280 x 40	Carbon Fiber/Silicone	Soft tissue retraction to expose humerus
4005-00-034	Small Deltoid Retractor	Trauma	280 x 30	Carbon Fiber	Soft tissue retraction to expose humerus
4005-00-035	Large Deltoid Retractor	Trauma	280 x 40	Carbon Fiber	Soft tissue retraction to expose humerus
4005-00-036	Fukuda Retractor	Trauma	260 x 50	Carbon Fiber	Retracts glenoid and humerus
4005-00-037	Metal Hook Retractor	Trauma	280 x 30	Carbon Fiber/Silicone	To prevent luxation
4005-00-038	Radiolucent sleeve	Trauma	170	Carbon Fiber	Guide the drill
4005-00-039	Radiolucent sleeve	Trauma	200	Carbon Fiber	Guide the drill

Traceability Label Symbols (Label) as per NBR ISO 15223


Manufacturing Date

Valid Date

Batch No.

Manufacturer Data

Code

Non-Sterile

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

Manufacturing materials: Most surgical instruments, in their structural part, and that contact the patient, are made of carbon fiber, as provided by ISO 14127; that characterize their compositions, to produce surgical instruments. Some instruments have stainless steel, aluminum and silicone.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Blunt Non-Hinged Surgical Instruments (Carbon Fiber) - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Blunt Non-Hinged Surgical Instruments (Carbon Fiber) - BM are reused in another patient, it is required that it is cleaned, sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure

within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

Cleaning: Appropriate care with surgical instruments starts with appropriate cleaning. This report mentions some instrument cleaning and sterilization methods. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

Instrument Sterilization: Instruments are supplied non-sterile. They have to be sterilized before use.

We recommend autoclave steam sterilization at the hospital (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a piece of equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

Table 2 - Autoclave physical sterilization parameters

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
-------	-------------	----------------------------	--------------------------

Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

For further clarification, refer to instruction manual supplied with autoclave.

Another sterilization method: that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135 - Sterilization of health care products - Ethylene oxide - Requirements for development, validation and routine control of sterilization process for medical devices.

NOTE: The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

Storage: Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls.

Identification: Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product' ;'Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

Consumer Service: In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7900. In case

of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website:
<http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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