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Blunt Non-Hinged Surgical Instruments – Generation 4 – BM

RG.: 80128580160

ENGLISH
EUA

Technical Name: Orthopedic Implant Instrument **Code:** 1551560

Trade Name: Blunt Non-Hinged Surgical Instruments - Generation 4- BM **NCM:** 90189099

Table 1 - Listing of Blunt Non-Hinged Surgical Instruments - Generation 4 - BM

Code	Description	Level	Size (mm)	Raw Material	Indication for Use
4004-00-000	Proximal SP2 Guide	Trauma	250 x 180	Stainless Steel 304; Stainless Steel 420; Stainless Steel 630 (NBR 13911); Aluminum*	Used to guide intramedullary rod and lock screws
4004-00-001	Proximal SP2 Guide	Trauma	250 x 180	Stainless Steel 304; Stainless Steel 420; Stainless Steel 630 (NBR 13911); Polymer*	Used to guide intramedullary rod and lock screws
4004-01-000	Proximal Right Distal Guide SP2	Trauma	400	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Polymer*	Used to guide screw drilling and placement in femur right distal part
4004-01-001	Proximal Right Distal Guide SP2	Trauma	400	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Aluminum*	Used to guide screw drilling and placement in femur right distal part
4004-02-000	Proximal Left Distal Guide SP2	Trauma	400	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Polymer*	Used to guide screw drilling and placement in femur left distal part
4004-02-001	Proximal Left Distal Guide SP2	Trauma	400	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Aluminum*	Used to guide screw drilling and placement in femur left distal part
4004-03-000	Proximal Distal Drilling Guide SP2	Trauma	280 x 180	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Aluminum*	Used to perform distal drilling

4004-03-001	Proximal Distal Drilling Guide SP2	Trauma	280 x 180	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Polymer*	Used to perform distal drilling
4004-04-009	Proximal Pointer SP2 - 9.0mm	Trauma	9 x 200	Stainless Steel 304 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
4004-04-091	Proximal Pointer SP2 - 9.0mm	Trauma	9 x 200	Stainless Steel 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
4004-05-010	Proximal Pointer SP2 - 10.0mm	Trauma	10 x 200	Stainless Steel 304 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
4004-05-101	Proximal Pointer SP2 - 10.0mm	Trauma	10 x 200	Stainless Steel 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
4004-06-011	Proximal Pointer SP2 - 11.0mm	Trauma	11 x 200	Stainless Steel 304 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
4004-06-111	Proximal Pointer SP2 - 11.0mm	Trauma	11 x 200	Stainless Steel 420 (NBR 13911)	Used to stabilize intramedullary rod distal drilling
4004-07-000	Proximal Extension SP2	Trauma	18 x 70	Stainless Steel 304 (NBR 13911)	Used to fasten implanting rods in drilling guide
4004-07-001	Proximal Extension SP2	Trauma	18 x 70	Stainless Steel 420 (NBR 13911)	Used to fasten implanting rods in drilling guide
4004-07-002	Proximal Extension SP2	Trauma	18 x 70	Stainless Steel 630 (NBR 13911)	Used to fasten implanting rods in drilling guide
4004-08-001	Proximal Screw Ruler SP2	Trauma	50 x 200	Stainless Steel 304 (NBR 13911)	Used to size proximal implant size
4004-09-001	Proximal Wire Ruler SP2	Trauma	50 x 200	Stainless Steel 304 (NBR 13911)	Used to size proximal implant size
4004-10-000	Proximal Fracture Reduction Wrench SP2	Trauma	380 x 70	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to femoral proximal part fracture reduction
4004-11-000	Proximal Tapered Milling Cutter Guide SP2	Trauma	200 x 90	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
4004-11-001	Proximal Tapered Milling Cutter Guide SP2	Trauma	200 x 90	Stainless Steel 304 (NBR 13911); w/ Rubber Handle	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
4004-11-002	Proximal Tapered Milling Cutter Guide SP2	Trauma	200 x 90	Stainless Steel 420 (NBR 13911); w/ Rubber Handle	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
4004-05-050	Proximal Tapered Milling Cutter Guide SP2	Trauma	200 x 90	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Aluminum*	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
4004-05-051	Proximal Tapered Milling Cutter Guide SP2	Trauma	200 x 90	Stainless Steel 304 (NBR 13911); Aluminum*	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth

4004-05-052	Proximal Tapered Milling Cutter Guide SP2	Trauma	200 x 90	Stainless Steel 420 (NBR 13911); Aluminum*	Used to protect sharp instrument soft parts and set cylindrical initial milling cutter depth
4004-12-000	Proximal Hexagon Wrench 8.0 SP2	Trauma	180 x 40	Stainless Steel 420 (NBR 13911)	Used to tighten lock screws
4004-14-000	Proximal Impactor Hammer SP2	Trauma	240 x 80	Stainless Steel 304; Stainless Steel 420 (NBR 13911)	Used to assist rod insertion and removal
4004-14-001	Proximal Impactor Hammer SP2	Trauma	240 x 80	Stainless Steel 304; Stainless Steel 630 (NBR 13911)	Used to assist rod insertion and removal
4004-15-000	Proximal Rod Extractor SP2	Trauma	340 x 25	Stainless Steel 304 (NBR 13911)	Used to extract the implant
4004-15-001	Proximal Rod Extractor SP2	Trauma	340 x 25	Stainless Steel 420 (NBR 13911)	Used to extract the implant
4004-17-000	Proximal Wire Guide SP2	Trauma	15 x 200	Stainless Steel 304 (NBR 13911)	Used to guide proximal guide wire and should be used connected in proximal guide
4004-17-001	Proximal Wire Guide SP2	Trauma	15 x 200	Stainless Steel 420 (NBR 13911)	Used to guide proximal guide wire and should be used connected in proximal guide
4004-18-000	Proximal Drill Bit Guide SP2 Ø 4.2 mm	Trauma	15 x 200	Stainless Steel 304 (NBR 13911)	Used to guide proximal drill bit and should be used connected in proximal guide
4004-18-001	Proximal Drill Bit Guide SP2 Ø 4.2 mm	Trauma	15 x 200	Stainless Steel 420 (NBR 13911)	Used to guide proximal drill bit and should be used connected in proximal guide
4004-19-000	Proximal Lock Screw Guide SP2	Trauma	15 x 200	Stainless Steel 304 (NBR 13911)	Proximal lock screw insertion guide
4004-19-001	Proximal Lock Screw Guide SP2	Trauma	15 x 200	Stainless Steel 420 (NBR 13911)	Proximal lock screw insertion guide
4004-20-000	Proximal Distal Drill Bit Guide SP2 Ø 4.2 mm	Trauma	15 x 200	Stainless Steel 304 (NBR 13911)	Used to guide drill bit and protect soft parts
4004-20-001	Proximal Distal Drill Bit Guide SP2 Ø 4.2 mm	Trauma	15 x 200	Stainless Steel 420 (NBR 13911)	Used to guide drill bit and protect soft parts
4004-21-000	Proximal Impactor SP2	Trauma	35 x 60	Stainless Steel 304 (NBR 13911)	Used to connect in implant guide and receive impactions
4004-21-001	Proximal Impactor SP2	Trauma	35 x 60	Stainless Steel 420 (NBR 13911)	Used to connect in implant guide and receive impactions

4004-21-002	Proximal Impactor SP2	Trauma	35 x 60	Stainless Steel 630 (NBR 13911)	Used to connect in implant guide and receive impactions
4004-23-000	Proximal Wire Guide SP2	Trauma	40 x 80	Stainless Steel 304; Stainless Steel 420 (NBR 13911)	Used to assist guide wire introduction.
4004-27-000	Proximal Lock Screw Extractor Wrench SP2	Trauma	8 x 250	Stainless Steel 304 (NBR 13911)	Used to extract proximal lock screw
4004-27-001	Proximal Lock Screw Extractor Wrench SP2	Trauma	8 x 250	Stainless Steel 420 (NBR 13911)	Used to extract proximal lock screw
4004-28-000	Proximal Fracture Reduction Wrench SP2	Trauma	80 x 340	Stainless Steel 304 (NBR 13911); Stainless Steel F138; w/ Rubber Handle	Used to femoral proximal part fracture reduction
4004-34-000	Proximal Reduction Assist Wrench SP2	Trauma	120 x 500	Stainless Steel 302; Stainless Steel 304 (NBR 13911)	Used to enable guide wire change
4004-35-000	Proximal Lock Screw Handle SP2	Trauma	14 x 120 x 80	Stainless Steel 304 (NBR 13911)	Used to hold and tighten the pointers
4004-36-000	Hexagon Wrench 4.0 mm	Trauma	300	Stainless Steel 420 (NBR 13911)	Used to tighten lock screws
4004-37-000	Proximal Soft Part Guard SP2	Trauma	220 x 140	Stainless Steel 304 (NBR 13911)	Used to protect and retract soft parts
4004-39-000	Proximal Femur Guide SP2	Trauma	350 x 200	Stainless Steel 304 (NBR 13911); Aluminum*	Used to align proximal screw alignment with femoral head
4004-40-000	Hexagon Wrench - Quick-Coupling	Trauma	320	Stainless Steel 420 (NBR 13911)	Used to tighten lock screws
4004-41-000	Hexagon Wrench	Trauma	320	Stainless Steel 420 (NBR 13911)	Used to tighten lock screws
6413-00-001	Sliding Screw Ruler	Trauma	50 x 250	Stainless Steel 304 (NBR 13911)	Used to size proximal implant size
6413-01-001	Proximal Wire Ruler	Trauma	50 x 250	Stainless Steel 304 (NBR 13911)	Used to size proximal implant size
4004-99-000	Head Extraction Lever	Arthroplasty	50 x 240	Stainless Steel 304 (NBR 13911)	Used to bipolar head extraction
4004-98-000	Test Screw with Head Insertion Punch	Arthroplasty	12 x 200	Stainless Steel 304 (NBR 13911)	Used to insert screw with head
6508-21-000	Tibial Soft Part Guard	Trauma	146 x 175	Stainless Steel 420 (NBR 13911); w/ Rubber Handle	Used to protect sharp instrument soft parts
6508-22-000	Femoral Soft Part Guard	Trauma	136 x 175	Stainless Steel 420 (NBR 13911); w/ Rubber Handle	Used to protect sharp instrument soft parts
4004-88-000	Soft Part Guard	Trauma	200 x 70	Stainless Steel 304 (NBR 13911)	Protect soft parts
4004-88-001	Soft Part Guard	Trauma	220 x 80	Stainless Steel 304 (NBR 13911)	Protect soft parts

4004-87-000	Milling Cutter Guide	Trauma	200 x 90	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to protect sharp instrument soft parts
4004-86-000	Milling Cutter Guide	Trauma	220 x 90	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to protect sharp instrument soft parts
4004-88-002	Soft Part Guard	Trauma	180 x 60	Stainless Steel 304 (NBR 13911)	Protect soft parts
4004-88-003	Soft Part Guard	Trauma	200 x 80	Stainless Steel 304 (NBR 13911)	Protect soft parts
4004-42-000	Humerus Reduction Wrench SP2	Trauma	400x120	Stainless Steel 302; Stainless Steel 304 (NBR 13911)	Used to enable guide wire change
4004-47-000	Hexagon Wrench - Quick-Coupling	Trauma	220x9.5	Stainless Steel 420 (NBR 13911)	Used to tighten lock screws
4004-48-000	Humerus Lock Screw Guide SP2	Trauma	Ø18x160	Stainless Steel 304 (NBR 13911)	Used to guide screw and protect soft parts
4004-48-001	Humerus Lock Screw Guide SP2	Trauma	Ø18x160	Stainless Steel 420 (NBR 13911)	Used to guide screw and protect soft parts
4004-49-000	Humerus Drill Bit Guide Ø4.2 mm SP2	Trauma	Ø18x160	Stainless Steel 420 (NBR 13911)	Used to guide drill bit and protect soft parts
4004-49-001	Humerus Drill Bit Guide Ø4.2 mm SP2	Trauma	Ø18x160	Stainless Steel 304 (NBR 13911)	Used to guide drill bit and protect soft parts
4004-50-000	Humerus Drill Bit Guide Ø4.2 mm SP2	Trauma	Ø18x160	Stainless Steel 420 (NBR 13911)	Used to guide drill bit and protect soft parts
4004-50-001	Humerus Drill Bit Guide Ø4.2 mm SP2	Trauma	Ø18x160	Stainless Steel 304 (NBR 13911)	Used to guide drill bit and protect soft parts
4004-52-000	Initial Drilling Wire	Trauma	Ø2x450	Stainless Steel 304 (NBR 13911)	Used as initial drilling guide
4004-54-000	Straight Handle - Humerus Quick-Coupling SP2	Trauma	164x30	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to tighten and loosen
4004-55-000	Screw Ruler	Trauma	182x35	Stainless Steel 304 (NBR 13911)	Used to size proximal implant size
4004-56-000	Humerus SP2 Extension	Trauma	Ø15x110	Stainless Steel 304 (NBR 13911)	Used to fasten implanting rods in drilling guide
4004-57-000	Wire Guide	Trauma	Ø33x75	Stainless Steel 304; Stainless Steel 420 (NBR 13911)	Used to assist guide wire introduction
4004-58-000	Humerus SP2 Impactor	Trauma	Ø25x72	Stainless Steel 304 (NBR 13911)	Used to connect in implant guide and receive impactions
4004-58-001	Humerus SP2 Impactor	Trauma	Ø25x72	Stainless Steel 420 (NBR 13911)	Used to connect in implant guide and receive impactions

4004-58-002	Humerus SP2 Impactor	Trauma	Ø25x72	Stainless Steel 630 (NBR 13911)	Used to connect in implant guide and receive impactions
4004-60-000	Humerus SP2 Impactor Extractor	Trauma	200x70	Stainless Steel 304; Stainless Steel 630 (NBR 13911)	Used to assist rod insertion and removal
4004-61-000	Humerus Rod Extractor Adapter SP2	Trauma	Ø12.7x88	Stainless Steel 304 (NBR 13911)	Used to extract the implant
4004-61-001	Humerus Rod Extractor Adapter SP2	Trauma	Ø12.7x88	Stainless Steel 420 (NBR 13911)	Used to extract the implant
4004-61-002	Humerus Rod Extractor Adapter SP2	Trauma	Ø12.7x88	Stainless Steel 630 (NBR 13911)	Used to extract the implant
4004-62-000	Milling Cutter Guide	Trauma	170x90	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to protect sharp instrument soft parts
4004-62-001	Milling Cutter Guide	Trauma	170x90	Stainless Steel 304 (NBR 13911); w/ Rubber Handle	Used to protect sharp instrument soft parts
4004-62-002	Milling Cutter Guide	Trauma	170x90	Stainless Steel 420 (NBR 13911); w/ Rubber Handle	Used to protect sharp instrument soft parts
4004-62-003	Milling Cutter Guide	Trauma	170x90	Stainless Steel 304; Stainless Steel 420 (NBR 13911); Aluminum*	Used to protect sharp instrument soft parts
4004-62-004	Milling Cutter Guide	Trauma	170x90	Stainless Steel 304 (NBR 13911); Aluminum*	Used to protect sharp instrument soft parts
4004-62-005	Milling Cutter Guide	Trauma	170x90	Stainless Steel 420 (NBR 13911); Aluminum*	Used to protect sharp instrument soft parts
4004-63-000	Quick-Coupling T Handle	Trauma	130x110	Stainless Steel 340; Stainless Steel 420 (NBR 13911) w/ Rubber Handle	Used to tighten and loosen
4004-64-000	Extension Wrench	Trauma	Ø15.8x60	Stainless Steel 304 (NBR 13911)	Used to connect in handle and extension
4004-64-001	Extension Wrench	Trauma	Ø15.8x60	Stainless Steel 420 (NBR 13911)	Used to connect in handle and extension
4004-64-002	Extension Wrench	Trauma	Ø15.8x60	Stainless Steel 630 (NBR 13911)	Used to connect in handle and extension
4004-65-000	Humerus Guide SP2	Trauma	200x160	Stainless Steel 304; Stainless Steel 420; Stainless Steel 630 (NBR 13911); Aluminum*	Used to guide intramedullary rod and lock screws
4004-65-001	Humerus Guide SP2	Trauma	200x160	Stainless Steel 304; Stainless Steel 420; Stainless Steel 630 (NBR 13911); Polymer*	Used to guide intramedullary rod and lock screws
4004-66-000	Rod Extractor	Trauma	Ø33x300	Stainless Steel 304 (NBR 13911)	Used to extract the implant
4004-66-001	Rod Extractor	Trauma	Ø33x300	Stainless Steel 420 (NBR 13911)	Used to extract the implant

2326-50-000	Anchor Screw Applicator 5.0	Arthroscopy	35x23x265 mm	Stainless Steel 304 (NBR 13911)/Polypropylene*	Used to insert the implant
2326-35-000	Anchor Screw Applicator 3.5	Arthroscopy	35x23x265 mm	Stainless Steel 304 (NBR 13911)/Polypropylene*	Used to insert the implant
3126-27-000	Anchor Screw Applicator 2.7	Arthroscopy	35x23x265 mm	Stainless Steel 304 (NBR 13911)/Polypropylene*	Used to insert the implant
3126-20-000	Anchor Screw Applicator 2.0	Arthroscopy	35x23x265 mm	Stainless Steel 304 (NBR 13911)/Polypropylene*	Used to insert the implant
*Every instrument that has Aluminum, Polymer and Polypropylene in its composition does not contact the patient, only parts made of Stainless Steel contact the patient.					

Traceability Label Symbols (Label) as per NBR ISO 15223


Manufacturing Date

Valid Date


Batch No.

Manufacturer Data

Code

Non-Sterile

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Sharp components can only be used for their specific purposes. Misuse can cause blade and edge cutting loss, and impact their performance within surgical procedure. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

Manufacturing materials: Most surgical instruments, in their structural part, are made of stainless steels, as provided by NBR 13911:2013 or ASTM F899:12; that characterize their chemical composition, mechanical properties and final behavior of stainless steel, to produce surgical instruments. The other parts are made of high strength polymers, that support sterilization.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Blunt Non-Hinged Surgical Instruments - Generation 4 - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Blunt Non-Hinged Surgical Instruments - Generation 4 - BM are reused in another patient, it is required that it is cleaned, sanitized and then re-sterilized. This process has to be

performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in package inserts supplied with implants. Instruments should not be used for different purposes than intended.

Cleaning: Appropriate care with surgical instruments starts with appropriate cleaning. This report mentions some instrument cleaning and sterilization methods. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

Instrument Sterilization: Instruments are supplied non-sterile. They have to be sterilized before use.

We recommend autoclave steam sterilization at the hospital (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a piece of equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

Table 2 - Autoclave physical sterilization parameters

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

For further clarification, refer to instruction manual supplied with autoclave.

Another sterilization method: that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

Storage: Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls.

Identification: Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product'

;Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

Consumer Service: In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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