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Surgical Instruments NANC Aluminium - BM

RG.: 80128580162

ENGLISH

EUA



Table 1 - List of Surgical Instruments NANC Aluminum - BM

CODE	DESCRIPTION	MP	SIZE	INDICATION
4000-03-051	Guide to proximal humerus plate	Aluminium	200 X 90 MM	Drilling guide for umeral plate
4000-03-052	Guide to proximal right humerus plate	Aluminium	200 X 90 MM	Drilling guide for umeral plate
4000-03-053	Guide to proximal plate of left humerus	Aluminium	200 X 90 MM	Drilling guide for umeral plate
4000-03-054	Adapter for one-eral connection	Aluminium	150 X 80 MM	One-plate adapter guide
4000-03-055	Adapter for right oneel connection	Aluminium	150 X 80 MM	One-plate adapter guide
4000-03-056	Adapter for left oneel connection	Aluminium	150 X 80 MM	One-plate adapter guide
4000-03-057	Umero3 Holes Proximal Plate Template	Aluminium	200 X 30 MM	Used to indicate implant dimension
4000-03-058	Umero4 Holes Proximal Plate Template	Aluminium	210 X 30 MM	Used to indicate implant dimension
4000-03-059	Umero5 Holes Proximal Plate Template	Aluminium	220 X 30 MM	Used to indicate implant dimension
4000-03-060	Umero6 Holes Proximal Plate Template	Aluminium	230 X 30 MM	Used to indicate implant dimension
4000-03-061	Umero7 Holes Proximal Plate Template	Aluminium	240 X 30 MM	Used to indicate implant dimension
4000-03-062	Umero8 Holes Proximal Plate Template	Aluminium	250 X 30 MM	Used to indicate implant dimension
4000-03-063	Guide to distal femur plate L.D.	Aluminium	460 X 60 MM	Drilling guide for femoral plate
4000-03-064	Guide to distal femur plate L.E.	Aluminium	460 X 60 MM	Drilling guide for femoral plate
4000-03-065	Adapter for femoral connection	Aluminium	60 X 90 MM	Femoral plate adapter guide
4000-03-066	Adapter for femoral right connection	Aluminium	60 X 90 MM	Femoral plate adapter guide
4000-03-067	Adapter for femoral left connection	Aluminium	60 X 90 MM	Femoral plate adapter guide
4000-03-068	Guide to proximal tibia plate L.D.	Aluminium	460 X 80 MM	Proximal tibia drilling guide.
4000-03-069	Guide to proximal tibia plate L.E.	Aluminium	460 X 80 MM	Proximal tibia drilling guide
4000-03-070	Adapter for tibial connection	Aluminium	60 X 90 MM	Tibial plate adapter guide
4000-03-071	Adapter for right tibial connection	Aluminium	60 X 90 MM	Tibial plate adapter guide
4000-03-072	Adapter for tibial left connection	Aluminium	60 X 90 MM	Tibial plate adapter guide
4000-03-073	Template Distal Plate of Right Femur 04 Holes	Aluminium	280 X 35 MM	Used to indicate implant dimension
4000-03-074	Template Distal Plate of Right Femur 05 Holes	Aluminium	300 X 35 MM	Used to indicate implant dimension
4000-03-075	Template Distal Plate of Right Femur 06 Holes	Aluminium	320 X 35 MM	Used to indicate implant dimension
4000-03-076	Template Distal Plate of Right Femur 07 Holes	Aluminium	340 X 35 MM	Used to indicate implant dimension
4000-03-077	Right Femur Distal Plate Template 08 Holes	Aluminium	360 X 35 MM	Used to indicate implant dimension
4000-03-078	Template Distal Plate of Right Femur 09 Holes	Aluminium	380 X 35 MM	Used to indicate implant dimension
4000-03-079	Right Femur Distal Plate Template 10 Holes	Aluminium	400 X 35 MM	Used to indicate implant dimension
4000-03-080	Right Femur Distal Plate Template 11 Holes	Aluminium	420 X 35 MM	Used to indicate implant dimension
4000-03-081	Right Femur Distal Plate Template 12 Holes	Aluminium	440 X 35 MM	Used to indicate implant dimension
4000-03-082	Right Femur Distal Plate Template 13 Holes	Aluminium	460 X 35 MM	Used to indicate implant dimension
4000-03-084	Left Femur Distal Plate Template 04 Holes	Aluminium	280 X 35 MM	Used to indicate implant dimension
4000-03-085	Left Femur Distal Plate Template 05 Holes	Aluminium	300 X 35 MM	Used to indicate implant dimension
4000-03-086	Left Femur Distal Plate Template 06 Holes	Aluminium	320 X 35 MM	Used to indicate implant dimension
4000-03-087	Template Distal Plate of Left Femur 07 Holes	Aluminium	340 X 35 MM	Used to indicate implant dimension

4000-03-143	Right Tibia Proximal Plate Template	Aluminium	460 X 35 MM	Used to indicate implant dimension
4000-03-144	Left Tibia Proximal Plate Template	Aluminium	460 X 35 MM	Used to indicate implant dimension
4000-03-145	Proximal Plate Template of Tibia	Aluminium	460 X 35 MM	Used to indicate implant dimension
4000-03-180	Markers for minimally invasive guide	Aluminium	Ø 12 X 30 MM	It serves as guidance for the surgeon to know which drilling of the minimally invasive guide was used
4000-03-232	DRILLING GUIDE FOR PROXIMAL PLATE HUMERUSWITH STOP	Aluminium	60 X 30 MM	Guide to proximal drilling.
4000-13-232	DRILLING GUIDE FOR PROXIMAL PLATE Humerus	Aluminium	60 X 30 MM	Guide to proximal drilling.
4000-03-233	DRILLING GUIDE FOR DISTAL PLATE OF FEMUR L.D.	Aluminium	60 X 30 MM	Guide to proximal drilling.
4000-03-234	DRILLING GUIDE FOR DISTAL PLATE OF FEMUR L.E.	Aluminium	60 X 30 MM	Guide to proximal drilling.
4000-03-235	DRILLING GUIDE FOR PROXIMAL PLATE OF TIBIA LARGA L.D.	Aluminium	60 X 30 MM	Guide to proximal drilling.
4000-03-236	DRILLING GUIDE FOR PROXIMAL PLATE OF TIBIA LARGA L.E.	Aluminium	60 X 30 MM	Guide to proximal drilling.
4000-03-237	Left distal femur osteotomy guide	Aluminium	60 X 30 MM	Cutting guide
4000-03-238	Guide to right distal femur osteotomy	Aluminium	60 X 30 MM	Cutting guide
4000-03-239	Guide to medial tibia osteotomy	Aluminium	60 X 30 MM	Cutting guide
4000-03-240	Left proximal tibia osteotomy guide	Aluminium	60 X 30 MM	Cutting guide
4000-03-241	Right proximal tibia osteotomy guide	Aluminium	60 X 30 MM	Cutting guide
4000-03-242	Template for straight reconstruction plate 20 Holes	Aluminium	210MM	Conform to the plates.
4000-03-243	Template for straight reconstruction plate 19 Holes	Aluminium	200 MM	Conform to the plates.
4000-03-244	Template for straight reconstruction plate 18 Holes	Aluminium	190 MM	Conform to the plates.
4000-03-245	Template for straight reconstruction plate 17 Holes	Aluminium	180 MM	Conform to the plates.
4000-03-246	Template for straight reconstruction plate 16 Holes	Aluminium	170 MM	Conform to the plates.
4000-03-247	Template for straight reconstruction plate 15 Holes	Aluminium	160 MM	Conform to the plates.
4000-03-248	Template for straight reconstruction plate 14 Holes	Aluminium	150 MM	Conform to the plates.
4000-03-249	Template for straight reconstruction plate 13 Holes	Aluminium	140 MM	Conform to the plates.
4000-03-250	Template for straight reconstruction plate 12 Holes	Aluminium	130 MM	Conform to the plates.
4000-03-251	Template for straight reconstruction plate 11 Holes	Aluminium	120 MM	Conform to the plates.
4000-03-252	Template for straight reconstruction plate 10 Holes	Aluminium	110 MM	Conform to the plates.
4000-03-253	Template for straight reconstruction plate 09 Holes	Aluminium	100 MM	Conform to the plates.
4000-03-254	Template for straight reconstruction plate 08 Holes	Aluminium	90 MM	Conform to the plates.
4000-03-255	Template for straight reconstruction plate 07 Holes	Aluminium	80 MM	Conform to the plates.
4000-03-256	Template for straight reconstruction plate 06 Holes	Aluminium	70 MM	Conform to the plates.
4000-03-257	Template for straight reconstruction plate 05 Holes	Aluminium	60 MM	Conform to the plates.
4000-03-258	Template for straight reconstruction plate 04 Holes	Aluminium	50 MM	Conform to the plates.
4000-03-259	Template for straight reconstruction plate 03 Holes	Aluminium	40 MM	Conform to the plates.
4000-03-260	Template for reconstruction plate R75 20 Holes	Aluminium	210 MM	Conform to the plates.
4000-03-261	Template for reconstruction plate R75 19 Holes	Aluminium	200 MM	Conform to the plates.
4000-03-262	Template for reconstruction plate R75 18 Holes	Aluminium	190 MM	Conform to the plates.
4000-03-263	Template for reconstruction plate R75 17 Holes	Aluminium	180 MM	Conform to the plates.
4000-03-264	Template for reconstruction plate R75 16 Holes	Aluminium	170 MM	Conform to the plates.
4000-03-265	Template for reconstruction plate R75 15 Holes	Aluminium	160 MM	Conform to the plates.

4000-03-313	Template for reconstruction plate R108 03 Holes	Aluminium	40 MM	Conform to the plates.
4000-03-314	Template for plate for calcaneus	Aluminium	100 X 40 MM	Conform to the plates.
6431-00-000	Shard Stabilization Guide	Aluminium / Stainless Steel 304* / Polyacetal*	230 x 120 MM	Used to stabilize fragments of the femoral head

Symbology of traceability label (Label) according to NBR ISO 15223

Manufacturing Date	Validity	Lot No.	Manufacturer Data	Code	Non-Sterile
					

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been specifically developed to meet the preferences of certain surgeons. Many bear the name of the surgeon they were raised by. New instruments were being developed to perfect various surgical procedures or even help create new techniques. While innovation means new designs and even new categories of instruments, instruments can generally be grouped by type or applications. Surgery techniques vary according to the knowledge and choice of the surgeon physician. It is up to the surgeon to choose the final method, type and size of the instruments and implants to be used, as well as the criteria for evaluating the post-surgical results. Biomechanica-branded Instruments should be used during the surgical procedure according to the purpose for which each instrument was developed. The use of these instruments is indicated for use only by professionals. Surgeons who manage the use of these products need to know perfectly the implant processes as well as the handling of instruments and components for implants. Instruments that fall to the ground are invariably damaged. When this happens, the instrument must be placed apart, carefully inspected and sent for repairs. The counting of the instruments is essential. First, to protect the patient himself and, secondly, to avoid the possibility of being inadvertently sent to the laundry room along with the surgical fields, where they can be lost or damaged if they go unnoticed.

Manufacturing materials: Most surgical instruments, in their structural part, are manufactured in aluminum, as specified in the NBR ISO 209:10 standard, which characterize their chemical composition, mechanical properties and the final steel behavior of surgical instruments.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Surgical Instruments are indicated as an indispensable aid for the preparation of the bone segment that will receive the implant. The instrumental brand Biomechanica can only be used to support the placement of biomechanica brand implants, as these were manufactured and designed to be used together. Using instruments from different manufacturers can compromise surgery. With the time of use, the instruments tend to wear naturally due to their regular use. Instruments should be used only for what was specified. The SURGICAL INSTRUMENTS NANC Aluminum - BM, when reusing it in another patient, it is necessary that it be washed,

sanitized and then resterilized. This process should be performed whenever it is used in a new surgery. Reuse of it without going through a process of hygiene and reesterilization is expressly prohibited. Worn and broken instruments can cause significant problems if not detected during inspection and testing. They may damage other instruments during cleaning/sterilization processes. If they remain in the instrument set they can cause from trauma to organic tissues to malfunction at critical times during a surgical procedure. Identify the problem instruments as soon as possible. Separate them from the or any other ones and place labels indicating that they should be replaced or repaired. The contraindication is linked to the implant to be used as described in the instructions for use accompanying the implants. Instruments should not be used for purposes other than those for which they were designed.

Cleaning : Proper care with surgical instruments begins with proper cleaning. This report cites some methods of cleaning and sterilization of the instruments. Proper care also means regular maintenance of instruments by predicting edging and adjustments. There is no standard programming; maintenance will be determined by the frequency of use. In addition to the care of cleaning the instruments, this instruction cites several enemies of surgical instruments, such as: blood, tissue in general, surgical residues (which are the primary causes of pitting corrosion points), stains and discoloration of instruments. Water and moisture also have harmful effects, allowing these substances to dwell or get drunk on their instruments causing unwanted stains. Other enemies used in the washing of instruments with improper solutions are: soap, bleaches, disinfectants and other solutions not recommended. For the correct preservation of your instruments it is important to use recommended cleaning methods and understand the causes of undesirable effects, such as stains. The spots appear with an orange or brown coloring. The idea is to ensure proper care to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, which lead to corrosion, appearing more often with an orange-brown color. If the instrument remains for an extended period of time (one to four hours), in contact with these residues, marks and stains will appear on the instrumental, especially if these residues dry together with the instrumental. The instruments should be thoroughly cleaned and dry after use. Just sterilize a clean instrumental. The elevated temperature of the autoclave will cause chemical reactions that can leave permanent stains on the instruments and/or yellowing of the instruments.

Instrument Sterilization: Instruments are distributed non-sterile. Before use they should be sterilized.

We recommend steam sterilization in Autoclave in the hospital (ISO 17665-1: 2006 *Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a equipment for sterilization through saturated steam under pressure.

OPERATIONAL INSTRUCTIONS: For starting the equipment, check: if the circuit breaker is on; whether the water network register is open; whether the discharge register is closed.

and not to facilitate distribution and not become overdue material in stock; - store separately from non-sterile ones to reduce the level of external contaminants.

Storage : Non-sterile surgical instruments should be stored in such a way that the packaging is not affected. Packages should not be compressed, crushed, punctured or exposed to possible water damage. The storage area must not contain impurities, dust or pollutants of any kind and should not be exposed to extremes of temperature or humidity. Traffic in the area should be minimal. Instruments should be stored in carts or shelves, far from the floor, ceiling or walls.

Identification: The Instruments follow with labels that contain the information: Manufacturer's Name; Commercial Name of the Product; Technical Name; Lot Number; ANVISA registration number; Product Code; Quantity; Description of the product containing the packaging (model); Date of manufacture; Expiration date; Description of the raw material used for the manufacture of the product (Composition); Manufacturer's address; Name of the Technical Officer; The sayings: "Non-sterile product"; Before use, see instructions for use; Warnings/Precautions/Special Care/Symbology: See Instructions for Use.

Customer Service: If there is a need to make any complaint regarding the use of instruments related to any adverse effect that affects the safety of the user, the surgeon responsible must report this adverse event to the competent health agency and Biomecanica by e-mail sac@biomecanica.com.br or by phone 0xx14 2104 7900. In case of doubts the responsible surgeon or the health professional can communicate the adverse event through the Health Surveillance Notification System on the ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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