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Radel Blunt Surgical Instruments - BM

RG.: 80128580170

ENGLISH
EUA

Table 1 - Listing of Radel Blunt Surgical Instruments - BM

CODE	DESCRIPTION	MP	SIZE	INDICATION
4000-05-051	HUMERUS PROXIMAL PLATE GUIDE	RADEL	220 x 180 mm	Humeral Plate Drilling Guide
4000-05-052	RIGHT HUMERUS PROXIMAL PLATE GUIDE	RADEL	220 x 180 mm	Humeral Plate Drilling Guide
4000-05-053	LEFT HUMERUS PROXIMAL PLATE GUIDE	RADEL	220 x 180 mm	Humeral Plate Drilling Guide
4000-05-063	RH FEMUR DISTAL PLATE GUIDE.	RADEL	480 x 150 mm	Femoral Plate Drilling Guide
4000-05-064	LH FEMUR DISTAL PLATE GUIDE.	RADEL	480 x 150 mm	Femoral Plate Drilling Guide
4000-05-068	RH TIBIA PROXIMAL PLATE GUIDE.	RADEL	480 x 150 mm	Tibial Proximal Drilling Guide.
4000-05-069	LH TIBIA PROXIMAL PLATE GUIDE.	RADEL	480 x 150 mm	Tibial Proximal Drilling Guide
4000-05-192	DRILL BRIT GRADUATED GUIDE Ø2.0	RADEL	150 mm	Drilling and Measurement Guide
4000-05-388	HUMERAL WIRE METER	RADEL	150 mm	Drilling and Measurement Guide
4000-05-389	TIBIAL WIRE METER	RADEL	150 mm	Drilling and Measurement Guide
4000-05-390	FEMORAL WIRE METER	RADEL	150 mm	Drilling and Measurement Guide
4000-05-391	WIRE METER 1.5	RADEL	150 mm	Drilling and Measurement Guide
4000-05-392	WIRE METER 2.0	RADEL	150 mm	Drilling and Measurement Guide
4000-05-393	WIRE METER 2.5	RADEL	150 mm	Drilling and Measurement Guide
4000-05-394	WIRE METER 3.0	RADEL	150 mm	Drilling and Measurement Guide
4000-05-395	WIRE METER 3.5	RADEL	150 mm	Drilling and Measurement Guide
4000-05-193	KIRSCHNER WIRE METER	RADEL	150 mm	Measure the depth the wire was inserted.

4000-05-232	HUMERUS PROXIMAL PLATE DRILLING GUIDE WITH STOP	RADEL	80 x 20 mm	Proximal Drilling Guide.
4000-15-232	HUMERUS PROXIMAL PLATE DRILLING GUIDE	RADEL	80 x 20 mm	Proximal Drilling Guide.
4000-05-233	LH FEMUR DISTAL PLATE DRILLING GUIDE.	RADEL	80 x 20 mm	Distal Drilling Guide
4000-05-234	RH FEMUR DISTAL PLATE DRILLING GUIDE.	RADEL	80 x 20 mm	Distal Drilling Guide.
4000-05-235	RH WIDE TIBIA PROXIMAL PLATE DRILLING GUIDE.	RADEL	60 x 20 mm	Proximal Drilling Guide.
4000-05-236	LH WIDE TIBIA PROXIMAL PLATE DRILLING GUIDE.	RADEL	60 x 20 mm	Proximal Drilling Guide.
4000-05-396	RIGHT RADIUS DISTAL PLATE DRILLING GUIDE	RADEL	60 x 20 mm	Distal Drilling Guide.
4000-05-397	LEFT RADIUS DISTAL PLATE DRILLING GUIDE	RADEL	60 x 20 mm	Distal Drilling Guide.
4000-05-237	LEFT DISTAL FEMUR OSTEOTOMY GUIDE	RADEL	60 x 20 mm	Cutting guide
4000-05-238	RIGHT DISTAL FEMUR OSTEOTOMY GUIDE	RADEL	60 x 20 mm	Cutting guide
4000-05-239	TIBIAL MEDIAL OSTEOTOMY GUIDE	RADEL	60 x 20 mm	Cutting guide
4000-05-240	LEFT PROXIMAL TIBIAL OSTEOTOMY GUIDE	RADEL	60 x 20 x 40 mm	Cutting guide
4000-05-241	RIGHT PROXIMAL TIBIAL OSTEOTOMY GUIDE	RADEL	60 x 20 x 40 mm	Cutting guide
4004-05-000	FÊNIX PROXIMAL GUIDE (RADEL)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	350 x 200 x 40 mm	Guides intramedullary rod and lock screws
4000-05-398	TIBIAL OBLIQUE DRILLING PROXIMAL GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	400 x 20 x 200 mm	Used to perform oblique proximal drilling
4000-05-399	TIBIAL LATERAL DRILLING PROXIMAL GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	200 x 80 mm	Used to perform static and dynamic lateral proximal drilling
4000-05-400	ORION SP HUMERUS GUIDE (ELITE)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	200 x 80 x 160 mm	Guides intramedullary rod and lock screws

4000-05-401	HUMERAL OBLIQUE DRILLING PROXIMAL GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	200 x 80 mm	Used to perform oblique proximal drilling
4000-05-402	ORION SP PROXIMAL GUIDE (ELITE)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	350 x 40 mm	Guides intramedullary rod and lock screws
4000-05-403	FRAGMENT STABILIZATION GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	350 x 40 mm	Used to stabilize femoral head fragments
4000-05-404	FEMORAL DISTAL DRILLING GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	80 x 60 x 40 mm	Used to identify implant distal drilling assisted or not by monitoring devices
4000-05-405	TIBIAL DISTAL DRILLING GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	80 x 60 x 40 mm	Used to identify implant distal drilling assisted or not by monitoring devices
4000-05-406	HUMERAL DISTAL DRILLING GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	80 x 60 x 40 mm	Used to identify implant distal drilling assisted or not by monitoring devices
4000-05-407	FEMORAL PROXIMAL DISTAL DRILLING GUIDE	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	80 x 60 x 40 mm	Used to identify implant distal drilling assisted or not by monitoring devices
4000-05-408	DISTAL GUIDE INTRAMEDULLARY ROD (FEMUR)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to guide screw drilling and placement in femur distal part

4000-05-409	DISTAL GUIDE INTRAMEDULLARY ROD (RETROGRADE)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to guide screw drilling and placement in femur distal part
4000-05-410	ANTERIOR SUPPORT INTRAMEDULLARY GUIDE ROD	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	380 x 40 mm	Used to guide screw drilling and placement in tibia distal part.
4000-05-411	SMALL ANTERIOR SUPPORT INTRAMEDULLARY GUIDE ROD	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	380 x 40 mm	Used to guide screw drilling and placement in tibia distal part.
4000-05-412	TIBIAL PROXIMAL SUPPORT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	120 x 40 mm	Used to perform oblique proximal and lateral medial drilling
4000-05-413	TIBIAL PROXIMAL RIGHT OBLIQUE SUPPORT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	120 x 40 mm	Used to perform oblique proximal and lateral medial drilling
4000-05-414	TIBIAL PROXIMAL LEFT OBLIQUE SUPPORT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	120 x 40 mm	Used to perform oblique proximal and lateral medial drilling
4000-05-415	TIBIAL DISTAL LATERAL MEDIAL SUPPORT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to support distal guide support
4000-05-416	TIBIAL DISTAL LATERAL MEDIAL SUPPORT – RIGHT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to support distal guide support

4000-05-417	TIBIAL DISTAL LATERAL MEDIAL SUPPORT - LEFT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to support distal guide support
4000-05-418	TIBIAL DISTAL DRILLING SUPPORT	RADEL ALUMINUM, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	400 x 20 x 200 mm	Used to perform lateral medial and distal antero-posterior drilling
4000-05-419	TIBIAL DISTAL DRILLING SUPPORT	RADEL ALUMINUM, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	220 x 60 x 120 mm	Used to perform lateral medial and distal antero-posterior drilling
4000-05-420	ORION SP ELITE GUIDE INTRAMEDULLARY ROD (FEMUR)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	350 x 40 mm	Guides intramedullary rod and lock screws
4000-05-421	ORION SP FEMUR LATERAL SUPPORT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to support distal guide support
4000-05-422	ORION SP FEMUR LATERAL SUPPORT – RIGHT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to support distal guide support
4000-05-423	ORION SP FEMUR LATERAL SUPPORT – LEFT	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 mm	Used to support distal guide support
4000-05-424	FEMORAL DISTAL DRILLING SUPPORT	RADEL ALUMINUM, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	220 x 60 x 120 mm	Used to perform lateral medial and distal antero-posterior drilling.

4000-05-425	PROXIMAL SCREW ALIGNMENT SUPPORT	RADEL ALUMINUM, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	480 x 40 x 300 mm	Used to proximal screw alignment with femoral head
4000-05-426	ORION SP ELITE GUIDE INTRAMEDULLARY ROD (TIBIA)	RADEL STAINLESS STEEL V630, STAINLESS STEEL AISI 304 (NBR 13911) (NBR 13911)*	400 x 20 x 200 mm	Guides intramedullary rod and lock screws
4000-05-427	FEMUR DISTAL GUIDE INTRAMEDULLARY ROD DECENTRALIZING DEVICE	RADEL STAINLESS STEEL AISI 304 (NBR 13911)*	200 x 200 mm	Used to hold femur distal guide intramedullary rod
4000-05-428	TIBIA DISTAL GUIDE INTRAMEDULLARY ROD DECENTRALIZING DEVICE	RADEL STAINLESS STEEL AISI 304 (NBR 13911)*	200 x 200 mm	Used to hold femur distal guide intramedullary rod
4000-05-429	FEMUR PROXIMAL DISTAL GUIDE DECENTRALIZING DEVICE	RADEL STAINLESS STEEL AISI 304 (NBR 13911)*	200 x 200 mm	Used to hold femur distal guide intramedullary rod
4000-05-430	HUMERUS DISTAL GUIDE DECENTRALIZING DEVICE	RADEL STAINLESS STEEL AISI 304 (NBR 13911)*	200 x 200 mm	Used to hold humeral distal guide intramedullary rod
4000-05-431	RADIOLUCENT DRILL BIT SLEEVE	RADEL	Ø 13 x 220 mm	Used with guide to drill bit sleeves
*Every instrument that has Stainless Steel and Aluminum in its composition does not contact the patient, only parts made of Radel contact the patient.				

Traceability Label Symbols (Label) as per NBR ISO 15223

Manufacturing Date	Valid Date	Batch No.	Manufacturer Data	Code	Non-Sterile
 XXXX	 XXXX	 XXXXX		 XXXX-XX	

TECHNICAL INFORMATION

Product Description: Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Sharp components can only be used for their specific purposes. Misuse can cause blade and edge cutting loss, and impact their performance within surgical procedure. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

Manufacturing materials: Most surgical instruments, in their structural part, are made of Radel, as provided in PPSU ; that characterize their chemical composition, mechanical properties and final behavior of Radel, to produce surgical instruments. Some instruments also have stainless steel and aluminum in their composition.

Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications: Surgical Instruments are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Radel Blunt Surgical Instruments - BM are reused in another patient, it is required that it is cleaned, sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues

to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

Cleaning: Appropriate care with surgical instruments starts with appropriate cleaning. This report mentions some instrument cleaning and sterilization methods. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

Instrument Sterilization: Instruments are supplied non-sterile. They have to be sterilized before use.

We recommend autoclave steam sterilization at the hospital (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

AUTOCLAVE: It is a piece of equipment for sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS: For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves

(saturated steam):

Table 2 - Autoclave physical sterilization parameters

Cycle	Temperature	Exposure Time (minimum)	Drying Time (maximum)
Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

For further clarification, refer to instruction manual supplied with autoclave.

Another sterilization method: that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

Storage: Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls.

Identification: Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product' ;'Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

Consumer Service: In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report

such adverse event to the competent health surveillance agency and Biomecânica through email sac@biomecanica.com.br or telephone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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