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## ENGLISH

**DESCRIPTION OF THE PRODUCT**

The System of Implantation for Setting of Column Plates BM was projected for use in the Column Cervical and Lumbar Toraco. It is a set of modular metallic components you implant, compatible between anatomical itself and to the bone, used for setting anterior (cervical), lateral or anterolateral (toraco-lumbar), promoting an immediate stability and protecting the colunal lace. Finishing of implantations is by mechanical polishing after that this product goes through a process of pickling, that consists of chemical attack of the polished surface of the metal that removes the residual impurities of the manufacturing process, such as oils, greases and other manufacturing materials then this product goes through an ionization process, that is, a superficial treatment in which the purpose is to fortify the superficial layer of the products manufactured in titanium.

**MANUFACTURING MATERIAL**

The products are made using the league 6Al titanium 4V ELI NBR ISO 5832-3 or ASTM F-136. It is a biocompatible league and it is resistant to the corrosion and fatigue.

**PURPOSE**

The System of Implantation for Setting of Column Plates BM is indicated to promote the fusing of the cervic-thoracic cervical column in situations as atlanto-axial fracture with instability, cervical luxation, fracture/luxation, cervical pain originated in the disc with discal degeneration, espondilolistese, estenose of the column, tumors. Its use can be posterior setting of the cervical column using the cervical plates, where it will promote an immediate stability. Treatment of instability of the thoracic, resulting from breaking of column, tumor, degeneration of the disc, escolioses, Lordoses, estenoses or imperfections of a previous surgery in the vertebral column being used for fractures from explosion with the affected canal and tumors in the part of the vertebral body.

**PRESENTATION FORMS**

All implantation of the System of Implantation for Setting of Column - Plate BM is packed in plastic envelopes, identified with non barren label, external labeling, containing 1 unit of plantation, 1 set of instructions (Dulla seal) and 1 internal label to be fixed in the handbook of the patient. The hospital institution is responsible in assuring the identification and tracking of the medical product, when used, through the data and the necessary information, contained in the labeling (label) of the medical product. The labeling contains the following information:

Data of the Manufacturer/Description of the content of the packing/Measures/Raw material/ Reference Unit/ Non Barren/Method of the Sterilization/Lot/Amount/Date of Manufacturing/Validity/hospice medical item for unique use/Instructions, Warnings and Precautions: To see use instructions/Responsible Technician/Registration of the ANVISA Nº/Brazilian Industry/and the saying: “Product of Only Use - Not to reuse”

**ANCILLARY COMPONENTS**

The System of implantation for Setting of Column Plates BM does not possess ancillary components.

**ACCESSORIES**

The System of Implantation for Setting of Column Plates BM does not possess any accessory with the intention to integrate the medical product, granting to this product a function or complementary technical characteristic.

**SUPPORT MATERIAL**

The support materials that go with the implantations of the System of Implantation for Setting of Column

Plates BM, are:

- instruction for use

**INDICATIONS OF USE FOR WHAT THE PRODUCT IS MADE**

The System of Implantation for Setting of Column Plates BM is indicated to promote the fusing of the cervic-thoracic cervical column in situations as atlanto-axial fracture with instability, cervical luxation, fracture/luxation, cervical pain originated in the disc with discal degeneration, espondilolistese, estenose of the column, tumors. Its use can be posterior setting of the cervical column using the cervical plates, where it will promote an immediate stability and protection of the lace of the column and in the previous setting through the toraco-lumbar plates, in the vertebral column being used for fractures from explosion with the affected canal and tumors in the part of the vertebral body.

NOTE: The implantations indicated here had been projected specifically for above described use, being under option of the surgeon for other applications.

**USE INFORMATION**

- For the correct use of the metallic components that you implant of the System of Implantation of Setting of Column Plates BM, it is necessary:

1.Deep anatomical knowledge of each vertebral segment to suffer surgery and knowledge of the pathology to be treated;
2.The use of the System of Implantation for Setting of Column - Plates BM must be done by surgeons qualified and trained for this procedure;
3.Analysis of the bony structure of the patient verifying if the same one stands the dimensions of the implantations of the System of Implantation for Setting of Column –Plates BM;
4.Close attention is necessary in the access: in the bony perforation and putting of the screws. We recommend for better security the use of imageintensifier;
5.The hospital must fill the handbook of the patient writing down the code of the product and number of the lot of the implantation. Biomecânica includes 1 extra adhesive label to be glued in the handbook of the patient,where all data of the implantation mentioned above is;
6.The use of surgical instruments of the same origin, that is, from the same manufacturer is important, because such instruments were projected specifically for use with this system of implantation;
7.The implantations of the System of Implantation for Setting of Column Plates BM are of unique use. The reuse is COMPLETELY FORBIDDEN.

**PRODUCT DISCART**

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the metallic components of the System of Implantation for Setting of Column Plates BM explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecânica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings “Improper for use”. According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

**CONTRAINDICATIONS**

This implantation is not indicated in the following situations:

- sensitivity to the metal;
- fever or sign of local inflammation;
- osteoporose;
- obesity (heavy load in the spinal system produced by obese patients can compromise the setting of the implantation for column);
- inadequate tessue covering in the area to be operated on;
- pregnanacy;
- bony compromise by illnesses and/or infections;
- unexplained rise of sedimentation tax through other illnesses, rise of the WBC tax, or impressing suspicion in distinguishing counting;
- visible infection;
- precocious deambulation;
- other situations that might prevent the success of spinal implantation surgery.
\* WBC: White cells

**CARE AND PRECAUTION**

Biomecânica fixation system for cervical spine should be implanted with the help of Biomecânica instruments because they were specially designed and manufactured for those implants and that guaranteees a proper implantation. - the surgeon must be aware of the possibility of the physical development of the patient, possibly having the necessity of verification of the implantation;- the surgeon must be familiar with and withhold enough knowledge of the surgical technique of spinal setting and its limitation, including pre-and post-operative, adopted surgical technique, precaution and potential risk;- the surgeon must be attentive in the bony preparation for the introduction of the screws and placing of the plates, because the stability and the success of the setting can be compromised;- Scratches, carvings and inverse packings in the implantations during the modeling are to be avoided, because these defects are concentrative of tension and can be small farms of nucleation of cracks and diminish the resistance to corrosion, being able to result in breaking of the implantation;- to guarantee an adjusted implantation, the instruments of column manufactured by the Biomecânica must be used, because they have been projected and manufactured specifically for these implantations. The use of instruments by different manufacturers can compromise the surgery;
- Do not use the product if damaged;
- PRODUCT MUST STERILE, STERILIZE BEFORE USE;
- PRODUCT OF ONLY ONE USE, REUSE IS FORBIDDEN;- Special conditions for storage: Keep the product in aired, dry place, away from the light and far from inclementer agencies;
- Do not use the product if the packing is damaged and/or date expired;
- Note: Component of implantations of different manufacturers MUST NOT be used; So we recommend that the products have the same origin.

**WARNINGS**

Product of Only Use. Should be destroyed after removed. Do not reuse the product. An orthopedic implantation can be used in a unique patient, only once. Although it may not seem to be damaged, previous tension can create imperfections that can reduce the success of the implantation. Improper choice of implantation can cause uncommon tensions in the implantation and can result in subsequence breaking of the implantation. Due to its form of presentation and its characteristics, this product does not present collateral effect, however, monitoring the recovery of the patient through X-ray, according to International Norms of Security, is suggested. Attention when using magnetic resonance in patients who already have implantations of column manufactured in stainless steel, because it causes interference in the examination disturbing the diagnosis. The use of the System of Implantation for Setting of Column - Plates BM is conceived to promote the fusing of the cervical column and the cervic-thoracic junction and not to substitute the normal structures of the human skeleton.

Obesity has effect on the load in the spinal system being able to compromise the setting of the implantation for column. Extreme physical activity and traumas directly related to the thoracic and cervical segments can compromise the success of the implantation prematurely. We recommend that the activities are dosed preventing extreme “stress” in the operated area. The system of implantation for Setting of Column - Plates BM cannot tolerate loads or levels of similar activities supported by normal and healthy bones. The mixture of different non compatible and/or similar biomedical materials can also increase the corrosion risk, taking into consideration that despite the implantation being manufactured by biocompatible material when implanted in the patient the same is subject to constant ambient changes promoting the concentration of salts, acids and alkalis that can promote the corrosion. The clinical results and durability of the implantation depend on the existence of an accurate surgical technique.

Manufacturing date, stated period of validity and lot of the product: SEE LABEL.

**POSSIBLE ADVERSE EFFECT**

a) - Reactions of sensitivity to the metal in patients have rarely been informed. Implantation of strange material in tessues results in histological reactions that involve various sizes of macrophages and fibroblasts. The chemical importance of this effect is uncertain, as similar changes may occur as a precursor, or during the curative process;
b) - Delay or non bony consolidation may lead to breaking of the implantation;
c) - Mechanical loosening, component torsion, dismemberment, breaking of components – may be the result of defective setting or occult infection;
d) Superficial and/or deep infection;
e) - Shortening of the segment and the vertebral column because of bony reabsorption;
f) - Abnormal pain, discomfort or sensations due to the setting of implantable metallic component or the migration of the implantable mechanical components, thrombosis;
g) Pseudo-artrose;
h) Fracture;
i) Paralysis;
j) - Hematomas and hemorrhage in the sanguineous vessels;
l) - Pulmonary Embolism;
m) – Neural or neurological problems caused by the surgical traumas like: gastrointestinal disturbance,
n) - Precocious or delayed failure of the implantable metallic components;
o) - Scabs

# SYSTEM OF IMPLANTATION FOR SETTING OF COLUMN PLATES BM

#### PACKING

The system of implantation for Setting of Column - Plates BM is packed in plastic envelopes, duly identified guaranteeing its tracking. In case of isolated packing, this product must not be used. The adhesive label that goes together with the implantable metallic component must be glued to the handbook of the patient. The hospital must register these information in its archives and repass these information to the patient. This label helps tracking of the implantable metallic component.

#### SPECIAL CARE IN THE STORAGE AND TRANSPORT OF THE MEDICAL PRODUCT

Do not use the product if it is damaged./ The implantations must be dealt with care preventing damages that can harm the quality of the material and the security of the patient./ Keep it in aired, dry place, away from the light and far from imclemancies./When shipping it, shocks, vibrations and high temperatures (above 45°) and inadequate pilings up must be prevented.

\* The caution to be taken of during the shipping is described in the collective packing.

#### STERILIZATION

This product is supplied non sterile. Before the use they must be esterilized. We recommend the vapor sterilization the in Sterilizer in hospital (ISO 11134: 1994 Sterilization of health care products. Requirements for validation and routine control industrial moist heat sterilization).

#### Sterilizer

It is an equipment for sterilization through vapor saturated under pressure.

#### Operational Instructions

To start the equipment, verify: if the breaker is on; if the register of the water net is opened; if the discharge register is closed.

1 - open the door of the equipment; 2 - place the material to be esterilized adequately; 3 - close the door of the equipment; 4 - select the cycle desired in accordance with the material to be esterilized; 5 – turn on general key; 6 - the cycle will run automatically,in the sequence; 7 - when light is on "end of cycle", partially open the door for about ten (10) minutes, for cooling of the material.

#### FOR BETTER INFORMATION CONSULT THE INSTRUCTION MANUAL THAT GOES WITH STERILIZER.

Another method of sterilization: But the following method of sterilization as defined by the hospital institution may be used besides the sterilizer.- Sterilization by ethylene oxide (E.T.O.) - parameters and procedures established in the protocol of validation and the EN550:1994 Sterilization of Medical Devices Validation and Routine Control of Ethylene Oxide Sterilization.

#### NOTE

The method of sterilization, the equipment, the controls, and the instructions of sterilization used is responsibility of the hospital institution.

#### RECOMMENDED CLEANNES AND DECONTAMINATION

General recommendations

To guarantee the reduction of infection risks, the surgical implantations and instruments should be cleaned after the surgery and before the sterilization. These products must be decontaminated using bactericidal and antifungal solutions of wide specter, avoiding agents of aggressive cleanness, metal brushes so that the quality of the metallic implantations is not affected. The manual cleanness is best recommended by attacking the surgical implantations and instruments less. In manual cleanness water in natural temperature must be used, use nylon brushes, never steel or abrasive ones, because it may damage the protective layer of the material. In case cleanness machines are used the implantations and surgical instruments must be placed in drawers so that one does not hit the other preventing damage in the material. Placing material of different composition together during the washing in cleanness machines is not recommended. The person responsible for the cleanness of the implantations and the surgical instruments must be aware of the type of the cleaning product used respecting the concentrations indicated by the manufacturer. The surgical implantations and instruments must be carefully rinsed until all residue has been wiped and that they are dried immediately after the cleanness.

#### RESTERILIZATION

In case of reesterilization, all the above care must be taken care of.

#### WARRANTY

The warranty of the System of Implantation for Setting of Column - Plates BM is related to the observance of this instruction of use. Bad use, disrespect of the indications of use or information of use, contraindications, care and precautions, warnings, possible adverse effect, recommended packing, sterilization, cleanness and decontamination in this instruction of use affects the warranty of this product.

#### COMPLAINT

In case there is any complaint of the System of Implantation for Setting of Column Plates BM it should be reported to the authorized distributor of Biomecánica. In case the complaint is related to some adverse effect that affects the security of the user, as product not functioning, serious damage to the implantable metallic component, serious problems or death related to these components Biomecánica must be informed.

#### SENDING OF MATERIAL FOR MANUFACTURER ANALYSIS

In case of sending implantations for manufacturer to analyse, it must be disinfected in the hospital using a bactericidal and antifungal solution of wide specter. Later it must be disinfected or sterilized by vapor in sterilizer or oxidate of ethylene. They must be sent to Biomecánica in intact packings, identified with the method of cleanness, sterilization and data of the product.

#### USE INSTRUCTION

The System of implantation for Setting of Column Plates BM has the purpose or the indication to promote the fusing of the cervical column and the cervico-toracica junction in situations as atlanto-axial breaking with instability, cervical luxation, fracture/luxation, cervical pain of discogenic origin with discal degeneration, espondiloliste, estenose of the column, tumors. Its use can be posterior setting of the cervical column using the cervical plates, where it will promote an immediate stability and protection of the colunal lace and in the previous setting through the toraco-lumbar plates, in the vertebral column being used for fracture by explosion with the affected canal and tumors in the part of the vertebral body.

The surgical techniques depend on the technique applied by the surgeon, in which the method; type and dimension of the product to be used is chosen. The mixture of different non compatible and/or similar biomedic material may increase the risk of corrosion. The clinical results and durability of implantation depend on the existence of a precise surgical technique.

This product must only be manipulated by qualified people. Components of different manufacturers shall not be used. So we recommend that products have the same origin. All implantations is ONLY USE. So, when removed from the patient, it must be cleaned and nos shall be reused.

#### ASSEMBLING AND DISSASSEMBLING

The setting of the System of Implantation for Setting of Column Plates BM is done with the help of biomecánica´s column instrumental kit.

#### THE PATIENT MUST BE INFORMED ABOUT:

The sayings in this instruction in the items:

Indications,contraindications,information of use – possible adverse effects,precautions and warnings. The patient must be informed of the importance of an after surgical follow up. The lack of follow up impedes the detection of post-surgical problems. The patient must be informed of the type of implantation used, brand, code, name,lot, surgeon's name and date of surgery.

#### PERFORMANCE FORESEEN IN THE REQUIREMENTS ESSENCIALS TO SECURITY AND EFFECTIVENESS OF THE PRODUCT

Limits of modeling or mechanical conformation of implantable metallic components.

For the success of the surgery is essential the correct election of the implantations. It is of extreme importance the correct manuscript of the implantations. During the manuscript of the implantations scratches must be avoided, notchs and inverse packings must be prevented during the a modeling, therefore these defects are concentrative of tension and can be small farms of nucleation of cracks and diminish the resistance of the corrosion, being able to result in breaking of the implantation.

#### CONTAMINATION RISKS

Exist the risk of biological contamination and transmission of illnesses you capsize such as HIV and hepatitis

contact the metallic components of the System of Implantation for Setting of Column Plates BM gets in contact with fabric and corporal fluids. Explanted products must be treated as highly agents that can contaminate.

#### TRACKING

Next to the packing of the implantable component it follows 1 label more l contend the data of the used implantation. This label must be glued in the handbook of the patient. The hospital must register these information in its archives and repass these information to the patient. In the handbook the following information are essential: name of the used implantation / dates of the surgery / code of the product / lot number / name of the patient who received the implantation / name from the surgeon.

#### MARKING

All the implantations mensioned in this use instruction will have the following information recorded with laser to allow tracking and follow up of the patient in the post-surgery: brand logo of the manufacturer / size / lot number / reference

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