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## ENGLISH

### DESCRIPTION OF THE PRODUCT

This device is constituted of a component modulate humeral proximal and one shell of the glenoid. The component of the humeral stem has three brims or flanges positioned directly below of the neck to resist the rotation. Each brim is perforated to allow the relocation of the soft tissues or fragment of the bone. The device is slotted distally to enhance the fixation of the cement. The devices are manufactured in the diameters of 7,5/10/12mm and the humeral heads in the diameters of 40 /43 /45 mm centered or eccentric, in agreement with the surgical technique adopted by surgeon.

### FORM OF PRESENTATION

Double packing blister stamped with surgical paper, plus external protection packing type cardboard box. The commercial name, the validity and sterilization process and the production lots are suitable in the internal and external packing of the product according to effective legislation.

### INDICATIONS

The replacement of the prosthesis with this device should be indicated to alleviate severe pains or significant inability degenerative or traumatic disease of the articulation gleno-humeral, disunite fractures of humeral head of long duration, vascular necrosis of the humeral head or other complicated clinical management's where the arthrodesis or arthroplasty resections are not acceptable.

### CONTRAINDICATIONS

Visible infection is an absolute contraindication. All effort should be made to govern the possibility of assepsis preoperative in a patient that has one or more of the following abnormalities: (1) fever or indication of local inflammation; (2) fast destruction of the articulations or apparent absorption of the bone in roentgenogram; and (3) inexplicable elevation of the sedimentation rate for other diseases, elevation in the WBC count, or suspect of transfer in differential count of WBC. Additionally, distant focus of infection as: genitourinary, lung, skin, and other places are relative contraindications, because can occur hematogenesis which is spilled for the place of the implants. Thus, infection focus should be treated before, during, and after the implantation. The use of this implants is contraindicated in the presence of significant damage to the superior brachial plexus, paralysis of the nerve axial, suspect of bone loss, non operation of the deltoid or external muscles rotator, or where disease neuro-muscular, commits the affected member what would make the procedure unjustifiable (for example, neuropathy in the articulation).

### WARNINGS

INAPPROPRIATE SELECTION OF IMPLANTS, INSERT, PLACEMENT, OR EXCESSIVE MUSCULAR ACTIVITY CAN CAUSE UNCOMMON TENSIONS IN IMPLANTS AND CAN RESULT IN SUBSEQUENT FRACTURE OF PROSTHESIS. In whole case, surgical practices should be followed in postoperative care. The patient should be alerted on the limitations of reconstruction of the total articulation and its recent history of use. Excessive physical activity and trauma affecting the substituted articulation was implied in premature failure in the reconstruction by looseness, fractures, and the quota or use of the implants

The patient should be disposed to proceed with the demanded rehabilitation program and be cautious to govern its activities even to certain necessary point, to protect the substituted articulation of unnecessary tensions.

It is suggested to monitoring the reestablishment of patient through X Ray. Do not utilize magnetic resonance due interferences.

The patient should be informed of the importance of a post-surgical follow-up. The failure of that follow-up prevents the detection of post-surgical problems. The failure of submitting the patient to a revision surgery in case of the component looseness or osteolysis may result in a progressive loss of the periprosthetic bone stock.

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## PRECAUTIONS

Biomechanical shoulder prostheses should only be implanted with the help of Biomechanical instruments because those instruments were specifically designed and manufactured for the use of those implants for the guarantee of a proper implantation. An orthopedic implant may only be used in a single patient for one single time. Although it does not look damaged, previous stresses may have created imperfections and that may impair the success of the implant. The improper selection of the implant may cause uncommon stresses on the implant and result in its subsequent fracture. Verify if the pack that contains the product is entire. Verify the sterilization date. Open aseptically the packing, after certify that is the dimension that was chosen.

After the opening of the packaging at the surgical Center even if the implant was not used it must not be re-used because it was exposed to the environment and that may cause its contamination. The implant should be cleaned and sent to the manufacturer to be sterilized again.

The surgeon must be aware of physical development of patient, having the possibility of revision on implant

## ADVERSE EFFECTS

Were informed components loosening, displacement, iatrogenic fractures of the glenoid, and traumatic fractures under the component humeral followed by to the use of similar devices.

Sensitivity reactions of the metal in patients, followed by the substitution of the total articulation, it was rarely informed. The meaning and sensibility effects, expect clinical evidence and future evaluation.

Implantation of strange material in tissues can result in histological reactions that involve several macrophages sizes and fibroblast. The current clinical meaning of this effect is uncertain because similar changes can happen as a precursor or during the normal healing process of the wound.

## SPECIAL CARES IN THE STORAGE AND TRANSPORT OF THE MEDICAL PRODUCT:

- Do Not use the product if it seems to be damaged
- The implants should be manipulated carefully ,avoiding damages to harm the quality of the material and the patient's safety
- Conserve at airy place, and dry, to the shelter of the light and far away from the action of bad weather
- When transporting him, it should be avoided shocks, vibrations and high temperatures (above 45?) and inadequate piling up
- The cares to be taken during the transport are described in the packing.

## STERILIZATION

Product sterilized by Oxide of Ethylene E.T.O. Validity 3 years.

In the eventuality of sterilization is expired or surgical paper damaged, this product must be returned to manufacturer.

## RECOMENDATIONS

Skilled people should only manipulate this product.

In case if this product is withdraw from patient, he must be cleaned and discarded to avoid contamination.

\*WBC: white blood cells

We recommended the completion of the product certificate of the quality, and the devolution to Biomecánica.

## PRODUCT DISCARD

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the shoulder prostheses explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

## CONTAMINATION RISKS

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

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