Português

CENTRALIZADOR

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Thais de Paula Busquim - CRQ nº 04263005

English

The centralizer is an implant made of PMMA (Polymethylmethacrylate). The PMMA used is one of the components of the orthopedic cement *Acrilex*, neither an object of this record nor part of this product, made by Biomecânica, ANVISA Registration No 10,171,110,070. The Orthopedic Cement *Acrilex* is a self-solidification cement (NBR ISO 5833), also the basis of PMMA, used for setting the metal prosthesis to the bone in arthroplasty

procedures, and according to the standard NBR ISO 21534 is considered an acceptable material for the manufacture of implants.

The centralizer is inserted along with Cemented Modular femoral Stem, made by Biomecânica, named CP3, which is neither an object of this record nor part of this product; and after cementation, the centralizer is embedded in the cement, becoming part of it.

The graphical information presented below is illustrative only, and refer to the centralizers.

The centralizers will come in a double blister packaging sealed with surgical paper, and an extra external packaging for protection made of cardboard box and labelling. The centralizers are sold in a three-unit set per package, which has a unique sale code. All the centralizers available in the package have the same inner diameter, and different outer diameter; therefore each unit of the set has a unique manufacturing code, as shown in "Table1" below. During the surgery, the surgeon will decide between one of the three centralizers of the set, which will depend on the reaming of the femoral canal, and the two remaining centralizers should be unused and discarded.

Name of the product set	Code of the set	Units available in the set	Description of the units
Centralizer 4,0	2571-04-001	03	This set contains One 4.0 Centralizer ext: 8.0 mm (code: 2571-04-008), One 4.0 Centralizer ext: 10.0 mm (code: 2571-04-010); One 4.0 Centralizer ext: 12.0 mm (code: 2571-04-012).
Centralizer 5,0	2571-05-001	03	This set contains one 5.0 Centralizer ext: 8.0 mm (code: 2571-05-008), One 5.0 Centralizer ext: 10.0 mm (code: 2571-05-010); one 5.0 Centralizer ext: 12.0 mm (code: 2571-05-012).
Centralizer 7,0	2571-07-001	03	This set contains one 7.0 Centralizer ext: 10.0 mm (code: 2571-07-010), one 7.0 Centralizer ext: 12.0 mm (code: 2571-07-012); one 7.0 Centralizer ext: 14.0 mm (code: 2571-07-014).

Centralizer 9,0	2571-09-001	03	This set contains one 9.0 Centralizer ext: 12.0 mm (code: 2571-09-012), one 9.0 Centralizer ext: 14.0 mm (code: 2571-09-014); one 9.0 Centralizer ext: 16.0 mm (code: 2571-09-016).
Centralizer 11,0	2571-11-001	03	Conjunto conteniendo 01 unidad de Centralizador 11,0 Ø ext: 12,0mm (código: 2571-11-012); 01 unidad de Centralizador 11,0 Ø ext: 14,0mm (código: 2571-11-014); 01 unidad de Centralizador 11,0 Ø ext: 16,0mm (código: 2571-11-016).
Centralizer 13,0	2571-13-001	03	This set contains one 13.0 Centralizer ext: 14.0 mm (code: 2571-13-014), one 13.0 Centralizer ext: 16.0 mm (code: 2571-13-016), one 13.0 Centralizer ext: 18.0 mm (code: 2571-13-018).

All centralizers are made of Polymethylmethacrylate (PMMA) and the only difference between them is the diameter (inner and outer diameter).

The difference of the centralizer presentation is related to the size of the Cemented Modular Femoral Stem, the CP3. The choice of which Centralizer to use will depend on the assessment of the surgeon in accordance with the Femoral Rod to be used.

The centralizers are made of PMMA (Polymethylmethacrylate), the same component of the orthopedic cement *Acrilex*, neither an object of this record nor part of this product, made by Biomecânica, ANVISA Registration No 10,171,110,070, according to the standard “NBR ISO 5833 – Surgical Implants - Acrylic Resin Cement” "evidenced by laboratory testing reports and certification of quality of the product. This material is acceptable for the manufacture of implants as defined in Annex A of the standard “NBR ISO 21,534 - Non-active surgical implants – joint replacement implants - Particular requirements." The manufacturing and packaging of the centralizers are based on specific procedures. The centralizers come in a double sterile blister packaging sealed with surgical paper, and an extra external protection made of cardboard box containing sufficient information and indications. It also has sufficient information for the correct handling of the product including the safety of those involved with its handling. The symbols present in the packaging and labels are in accordance with the standard “NBR ISO 15223 - Products for health - symbols to be used on labels, labeling and information to be supplied with products for health."

SPECIAL CONDITIONS FOR STORAGE, CONSERVATION AND / OR MANIPULATION OF THE PRODUCT

A. Carry with care to avoid its falling and friction in order not to damage the package and the device, ensuring the sterility of the material;

B. Keep it in a dry, aired place, and protected from light and away from the action of the weather;

C. shocks, vibration, high temperatures (over 45°) and inadequate stacks must be avoided.

D. Do not stored them directly on the floor (minimum height = 20cm). Do not keep them in very high places, near the lamps, which could cause package dryness or label damage.

E. Do not stored them in places where contaminants are used or kept, such as cleaning materials, insecticides, pesticides, etc.

F. Further care to be taken during transportation are described in the collective package

INDICATION OF USE

The Orthopedic Cement *Acrilex*, neither an object of this record and nor part of this product, made by Biomecânica, Anvisa registration No 10171110070, is used to set the prosthesis to the bone in the process of joint arthroplasty, and vertebroplasty surgeries, filling in the bone space of the vertebrae.

The centralizers are indicated in hip surgeries, in the placement of Cemented Modular Femoral Stems, neither an object of this record nor part of this product, made by Biomecânica, named CP3, as indicated:

Code of the product set	Name of the product set	Indication
2571-04-001	4,0 Centralizer	4.0 mm CP3 Femoral Stem - Cone 12/14 (2575-04-000); brand Biomecânica.
2571-05-001	5,0 Centralizer	5.0 mm CP3 Femoral Stem - Cone 12/14 (2575-05-000); brand Biomecânica.
2571-07-001	7,0 Centralizer	7.0 mm CP3 Femoral Stem - Cone 12/14 (2575-07-000); brand Biomecânica.
2571-09-001	9,0 Centralizer	9.0 mm CP3 Femoral Stem - Cone 12/14 (2575-09-000); brand Biomecânica .
2571-11-001	11,0 Centralizer	11.0 mm CP3 Femoral Stem - Cone 12/14 (2575-11-000); brand Biomecânica.
2571-13-001	13,0 Centralizer	13.0 mm CP3 Femoral Stem - Cone 12/14 (2575-13-000); brand Biomecânica.

When the bone cement is inserted for cementing the femoral stem to the bone the Centralizer incorporates the bone cement becoming part of it.

The purpose of the centralizer is to centralize the femoral stem distally at the time of its insertion through the femoral canal, thus avoiding the misplacement of the stem The Centralizer allows, therefore, an adequate positioning of the femoral stem in the femoral canal, and permits the distal migration of the stem within

Stems, made by Biomecanica, neither an object of this record nor part of this product, the CP3, as shown in "Table 2".

The alert against the use of this product is related to the use of bone cement and the femoral stem.

It is recommended to read the contra-indications contained in the leaflet of the bone cement before the use of centralizers.

The centralizers are contraindicated:

- A) In the presence of an infection;
- B) In not compensated clotting disorders.

WARNINGS AND PRECAUTIONS

For the proper use of centralizers the surgeon must have expertise in the use of bone cement and the placement of cemented femoral stems. The clinical findings clearly indicate the need of proper and strict surgical aseptic technique.

The centralizers must be used with Cemented Modular Femoral Stems only, made by Biomecanica, not an object of this record, model CP3.

Properly cemented prostheses are stable and durable, however, there is a possibility of loosening and fracture of the cement, fracture of the prosthesis, or both, caused by trauma, inadequate cementing technique, mechanical failure of the materials or infection. A long-term post-operative monitoring is important for all patients.

After the surgery the patient must be warned that if an intercurrent infection occurs, the doctor must be seen in order to reduce the risk of implant infection.

Read the instructions for use before handling the product. The centralizers are products for single use. The reuse is strictly prohibited. The patient must be informed about the surgical procedure.

RECOMMENDATIONS

- A. Single-use product – Do not reuse it;
- B. The centralizers provided are sterile - do not resterilize them;
- C. Use it immediately after opening the sterilization seal;
- D. In the surgery, the surgeon must choose one of three centralizers of the package, and the other two centralizers must be unused and discarded, according to the medical institution's procedures;
- E. Discard and do not use opened or damaged devices. Only use the devices that are in sealed and undamaged packages;
- F. Do not use it in case of loss of sterility of the device;
- G. The patient should be informed of all information regarding the implantation of the centralizer.

INFORMATION TO BE PROVIDED TO THE PATIENT

A. Indications, contraindications, usage information, possible adverse effects, precautions and warnings;
B. The patient must be informed about the importance of post-surgery monitoring. The lack of monitoring prevents the detection of post-surgical problems;
C. The patient must be informed about the brand, code, name of product, and batch of the femoral stem to be implanted, the bone cement used, and the Centralizer. The surgeon's name and date of surgery must be informed to the patient as well.

TRACEABILITY

The centralizers do not have a batch/code recorded on them because of the strength of the material. So, the batch/code should be described on the internal and external package label.

The surgeon and his team can make use of the 5 labels given, available in the product packaging by pasting them into the patient's medical records to maintain the traceability of the material. It is recommended that one of the 5 labels is given to the patient, so that the patient has information of the implant used.

The following information is essential in the patient's medical record:
-name of the implant used
-/date of surgery
-/product code / batch number
-name of the patient who received the implant
-/name of the surgeon.

PERMISSIBLE COMBINATIONS WITH OTHER MATERIALS

The centralizers should only be used with the Orthopedic Cement Acrilex (ANVISA Registration No. 10171110070) and the Cemented Modular Femoral Stems, neither one nor the other is part of this record, made by Biomecanica, and must be purchased separately. The centralizer should not be used with materials from other manufacturers for reasons of chemical, physical, biological and functional incompatibility.

ANCILLARY COMPONENTS (Neither an object of this record nor part of this product)

- :: Specific Ancillary components for the 4.0 Centralizer Set(2571-04-001);
- 4.0mm CP3 Femoral Stem - Cone 12/14 (2575-04-000), made by Biomecânica;
- Orthopedic Cement Acrilex (2005-00-000), made by Biomecânica;

- :: Specific Ancillary Components for the 5.0 Centralizer Set (2571-05-001)
- 5.0mm CP3 Femoral Stem - Cone 12/14 (2575-05-000), made by Biomecânica;
- Orthopedic Cement Acrilex (2005-00-000), made by Biomecânica;

- :: Specific Ancillary components for the 7.0 Centralizer (2571-07-001)
- 7.0mm CP3 Femoral Stem - Cone 12/14 (2575-07-000), made by Biomecânica;
- Orthopedic Cement Arilex (2005-00-000), made by Biomecânica;

- :: Specific Ancillary Components for the 9.0 Centralizer (2571-09-001)
- 9.0mm CP3 Femoral Stem - Cone 12/14 (2575-09-000), made by Biomecânica;
- Orthopedic Cement Acrilex (2005-00-000), made by Biomecânica;

- :: Specific Ancillary Components for the 11.0 Centralizer (2571-11-001)
- 11.0mm CP3 Femoral Stem - Cone 12/14 (2575-11-000), made by Biomecânica;
- Orthopedic Cement Acrilex (2005-00-000), made by Biomecânica;

- :: Specific Ancillary Components for the 13.0 Centralizer (2571-13-001)
- 13.0mm CP3 Femoral Stem - Cone 12/14 (2575-13-000), made by Biomecânica;
- Orthopedic Cement Acrilex(2005-00-000), made by Biomecânica;

ACCESSORIES

The Centralizers do not hold any accessory with the aim of integrating the medical product.

EVALUATION OF IMPLANTED PRODUCT

After the implantation, during the surgery, the responsible person should perform a radiological control to verify the correct positioning of the product. It is the responsible person's responsibility to make the clinical and radiological evaluations after the surgery, in the frequency prescribed by him/her, to check the status of the implant and follow the evolution of bone consolidation. If the product is not in the correct position, or shows any non-compliance, it is the surgeon's responsibility to take the most appropriate corrective action.

USEFUL INFORMATION TO AVOID RISK DUE TO THE IMPLANTATION

To reduce the risks of implantation, contraindications, instructions for use and all information contained in the "Instructions for Use" of the product must be strictly followed.

INTRINSIC RISK OF IMPLANTATION

The implantation of a strange material in organic tissues can result in inflammation, which can occur, for example, in the presence of implant debris (such as metal debris or polyethylene),

which can cause a histiocytic response, as a granuloma, and thus causing bone destruction associated or not with implant loosening.

The patient must also be informed that any circumstances that lead to chronic tissue damage may be oncogene. Cancerous tissue found around the implants may be related to factors that are not related directly to the implant as metastases from primary tumors of the lung, breast, digestive system and others, or because of the implantation of cancer cells that may occur during surgical or diagnostic procedures, such as the biopsy or arising from the progression of Paget's disease.

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ADDITIONAL INFORMATION ABOUT THE PROCEDURE BEFORE USING THE PRODUCT, STERILIZATION METHOD

The centralizers provided are sterile through ethylene oxide (ETO) sterilization - parameters and procedures established in the validation protocol, and EN550: 1994 Sterilization of Medical Devices - Validation and Routine Control of Ethylene Oxide Sterilization - valid for 5 years.

IMPLANTS SENT TO MANUFACTURER TO BE ANALYZED


If an implant is sent to the manufacturer in order to be analyzed, it must be cleaned at the hospital using a bactericide and antifungicide broad spectrum. Then, they must be disinfected or sterilized in a steam autoclave or ethylene oxide. They must be shipped to Biomecanica in integrated packaging, identified by the cleaning and sterilizing methods, and product data.

POST SALE (CUSTOMER CLAIM)

If there is a need of making a complaint about the centralizers, related to any adverse effect that affects the safety of the user, such as if the product is not working, damage of the implantable metal component, serious problems or death related to these components, the responsible surgeon must report this adverse event to the national health authority and Biomecanica through the e-mail: sac@biomecanica.com.br –or phone: (14) 21047900
In case of doubts, the responsible surgeon or the health professional can report the adverse event through the System of Notifications in Health Surveillance at ANVISA's homepage: http://www.anvisa.gov.br/hotsite/notivisa / index.htm

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