COMPONENTE PATELAR

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
PATELLAR COMPONENT

All BKS III - Patellar Components are manufactured from Polyethylene Ultra High Molecular Weight (UHMWPE).

Ancillary Components
The BKS III - Patellar Component has the following ancillary components (see below). But only Biomecnica BKS III - Femoral Component with or without Preservation will form an articular contact with the BKS III - Patellar Component. Biomecnica Orthopedic Cement, Acrlie®®, will fix the BKS III - Patellar Component to the patella, forming the prosthesis + bone cement + bone set. The other ancillary components will be part of total knee arthroplasty, but will not be in contact with the BKS III - Patellar Component.

List of ancillary components of BKS III Patellar Component
BKS III - Femoral Component - Cobalt-Chromium-Molybdenum Alloy Casting(ASTM F75); BKS III - Femoral Component with Preservation - Cobalt-Chromium-Molybdenum Alloy Casting(ASTM F75); BKS III - Tibial Component - Cobalt-Chromium-Molybdenum Alloy Casting(ASTM F75); BKS III - Tibial Plateau - Ultra High Molecular Weight Polyethylene (UHMWPE); BKS III - Tibial Plateau with Preservation - Ultra High Molecular Weight Polyethylene (UHMWPE); BKS III - Revision Femoral Stem - Cobalt-Chromium-Molybdenum Alloy Casting(ASTM F75); BKS III - Revision Tibial Stem - Cobalt-Chromium-Molybdenum Alloy Casting(ASTM F75); Assembly Screw - Wrought Titanium 6-aluminum 4-vanadium Alloy (ASTM F136) * Acrlie®® Orthopedic Cement - Acrylic Resin Cement (ISO 5833).

* It is allowed the use of another type of biomaterial component as an ancillary component of the BKS III - Patellar Component, replacing Biomecnica Acrlie®® Orthopedic Cement, if it is appropriate to fix the prosthesis to the bone and is regularly registered at ANVISA for this purpose.

Accessories

The BKS III - Patellar Component does not have any accessories in order to integrate the medical product, giving this product a valuable complementary role or feature.

Instruments for support

There are instruments, manufactured by Biomecnica, which serve as a support for the total knee arthroplasty procedure as shown below: 5211-00-000 - BKS III - Flat Reamer - 304 stainless steel; 5321-00-000 - BKS III - Alignment Rod - 304 Stainless Steel; 5320-00-000 - BKS III - Alignment Rod with Coupler - 304 Stainless steel; 5336-00-000 BKS III - Guide for Tibial Drill Bit - 304 stainless steel; 5355-02-000 - BKS III - 3-Peg Patellar Forcpes - 304 Stainless steel; 5335-00-000 - BKS III Handle for Tibial Tray Guide - 304 Stainless Steel; 5352-00-000 - BKS III Tibial Plateau Extractor - 304 stainless steel; 5334-01-000 - BKS III - Tibial Tray - Small - 304 Stainless Steel; 5334-02-000 - BKS III - Tibial Tray - Medium - 304 Stainless Steel; 5334-03-000 - BKS III - Tibial Tray - Large - 304 Stainless Steel; 5334-04-000 - BKS III - Tibial Tray Extra Large - 304 Stainless Steel 304; 5341-00-000 - BKS III - Tibial Tray Impactor - 304 stainless steel; 5340-00-000 - BKS III - Tibial Tray Extractor - 304 stainless steel; 5322-00-000 - BKS III - A/P Femoral Measuring Device - 304 stainless steel; 5337-00-000 - BKS III - Tibial Reamer Small/Medium - 304 stainless steel; 5338-00-000 - BKS III - Tibial Reamer Large/Extra Large - 304 stainless steel; 5319-00-000 - BKS III - Angular Support - 304 Stainless Steel; 5326-00-000 - BKS III - Extensor for IM Cutting Guide - 304 Stainless Steel; 5313-00-000 - BKS III - Femoral Alignment Guide - 304 Stainless Steel; 5328-00-000 - BKS III - Tibial Depth Resection Gauge - 304 stainless steel; 5339-00-000 - BKS III - Handle for Tibial Reamer - 304 stainless steel; 5354-00-000 - BKS III - Extractor Pin Pliers - 304 stainless steel; 5318-00-000 - BKS III - Distal Femoral Blemming Guide - 304 Stainless steel; 5311-00-000 - BKS III - Drill Bit 14mm - 320 Stainless steel; 5310-00-000 - BKS III - Staggered Drill Bit - 420 Stainless Steel; 5324-00-000 - BKS III - Intercondryl Guide G/E - 304 stainless steel; 5364-00-000 - BKS III - Puncture - 304 Stainless Steel; 5325-00-000 - BKS III - Tibial Stem Guide - 304 Stainless Steel 304; 5344-01-000 - BKS III - Tibia Trial - Small - 304 Stainless Steel; 5344-02-000 - BKS III - Medium Tibia Trial - 304 Stainless Steel; 5344-03-000; BKS III - Large Tibia Trial - 304 stainless steel; 5344-04-000 - BKS III - Extra Large Tibia Trial - 304 Stainless Steel; 535-00-000 - BKS III – Tibial Trimming Guide - 304 Stainless steel; 5350-00-000 - Ang. Tibial Trimming Guide - 304 Stainless Steel; BKS III Femoral Impactor - Politec; 5349-00-000 - BKS III - Tibial Impactor - Politec; 5345-08-000 - BKS III – Trial Plateau Small/Medium - 08 mm - 304 stainless steel; 5345-10-000 - BKS III – Trial Plateau Small/Medium - 10 mm - 304 stainless steel; 5345-12-000 - BKS III - Trial Plateau Small/Medium - 12 mm - 304 stainless steel; 5345-15-000 - BKS III - Trial Plateau Small/Medium - 15 mm - 304 stainless steel; 5346-08-000 - BKS III - Plateau Large/Extra Large - 08 mm - 304 Stainless Steel 304; 5346-10-000 - BKS III - Trial Plateau Large/Extra Large - 10 mm - 304 stainless steel; 5346-12-000 - BKS III - Trial Plateau Large/Extra Large - 12 mm - 304 stainless steel; 5346-15-000 - BKS III - Trial Plateau Large/Extra Large - 15 mm - 304 stainless steel; 5341-01-000 - BKS III – Small/Right Femoral Trial - 304 stainless steel; 5341-02-000 - BKS III – Femoral Trial Small/Left - 304 stainless steel; 5341-03-000 - BKS III - Femoral Trial Medium/Right - 304 Stainless Steel; 5341-04-000 - BKS III - Femoral Trial Medium/Left - 304 stainless steel; 5341-05-000 - BKS III - Large Right Femur Test - 304 stainless steel; 5341-06-000 - BKS III – Femoral Trial Large/Left - 304 stainless steel; 5341-07-000 - BKS III - Femoral Trial Extra Large/Right - 304 stainless steel 304; 5341-08-000 - BKS III – Femoral Trial Extra Large/Left - 304 stainless steel; 5316-00-000 - BKS III – Anterior Femoral Cutting Guide - 304 stainless steel; 5323-01-000 - BKS III – Small Femoral Cutting Guide - 630 Stainless Steel; 5323-02-000 - BKS III – Medium Femoral Cutting Guide - 630 Stainless Steel; 5323-03-000 - BKS III - Large Femoral Cutting Guide - 630 Stainless steel; 5323-04-000 - BKS

stabilizes the knee.

In the case of elderly patients with advanced osteoarthritis, the best treatment is total knee replacement. In younger patients, the total knee replacement is indicated only in cases of tricompartmental illness (medial and lateral tibiofemoral, and patellofemoral), with the need of a future revision surgery. In cases of rheumatoid arthritis, total knee replacement is recommended for all ages because of the complications presented by the disease.

The BKS III - Patellar component is indicated in the procedures for knee arthroplasty due to:

- Osteoarthritis:** Known as the “wear-and-tear” kind of arthritis, it is a chronic condition characterized by the breakdown of the joint's cartilage, causing the bones to rub against each other, causing stiffness, pain and loss of movement in the joint.
- Rheumatoid arthritis:** is an autoimmune disease that causes chronic inflammation of the joints. It can also cause inflammation of the tissue around the joints, as well as in other organs in the body.
- Comminuted fracture:** a fracture in which the bone involved in the fracture is broken into several pieces, making it difficult to achieve a good coaptation.

Contraindications

The contraindications below do not limit the use of BKS III - Patellar Components, but must be examined by the surgeon:

- Quadricep weakness, infections or systemic diseases near joints; previous patellectomy; neuropathic arthritis; in younger patients, a total knee replacement is indicated only in case of tricompartmental illness (medial and lateral tibiofemoral, and patellofemoral), and requires a future revision surgery, Morbid obesity, fever or leukocytosis; intense neuromuscular, muscular or vascular deficiencies in the affected region;

Severe risk of bone structures making it infeasible the application of the prosthesis; inappropriate ligament system state for the type of component chosen; bone tumors near the site where the component will be placed; muscular atrophy; bone immaturity; failure of the knee extensor mechanism; structural failure of the knee soft tissues; Ankylosis of the knee; insufficient bone structure in the tibia, femur and patella; sensitivity or allergy to metal medical devices; Patients who are unwilling or unable to follow the care and instructions in the postoperative period.

Adverse Effects

Vascular damage during knee replacement is extremely rare and potentially destructive, Asymptomatic and progressive local bone resorption (osteolysis) around the component may occur, Loosening, migration or fracture of the medical product; Superficial or deep infection, vascular disorders, including thrombosis and pulmonary embolism; Shortening of limb because of bone resorption; material sensitivity or allergic reaction to foreign body that may result in histological reactions involving various types of macrophages and fibroblasts; Pain, discomfort or abnormal sensations due to the presence or migration of the medical product and surgical procedure; nerve or neurological damage due to surgical trauma (including soft tissue paralysis and injuries); inability to perform normal daily physical activities; hemorrhage of blood vessels and / or bruising;

Warnings

Biomecánica patellar component should be implanted with the help of Biomecánica instruments. Components of implants from different manufacturers should not be used together, therefore it is recommended that the products have the same provenance.

The BKS III - Patellar Component aims to replace the damaged joint and not normal structures of the human body. The BKS III - Patellar Component will be sold already sterilized with ethylene oxide, and therefore it is of single-use, and its reprocessing is forbidden; Open the packaging aseptically, the models containing three pegs offer greater fixation since it presents a larger contact area between the implant, the patella and the bone cement. The patient's weight, activity level and attention to the instructions regarding the limitations of these conditions, take effect on load and number of cycles to which the medical device is exposed. The surgeon must be thoroughly acquainted with both medical and surgical aspects of the medical device. The correct selection of BKS III - Patellar Component model is extremely important; the BKS III - Patellar component can not support activity levels and / or loads equal to those incurred in normal and healthy bones; Replacement of patellofemoral joint, in cases of osteoarthritis, is not necessary.

The replacement of the patellofemoral joint can lead to possible complications, such as the creation of two more interfaces (bone-cement and cement-prosthesis), changes in bone vascularization by the decreasing patellar resistance and patellar fracture caused by osteonecrosis. However, the replacement of the patellofemoral joint in rheumatoid arthritis patients should be performed because of the disease evolution. Patellar luxation may occur even where there was no joint replacement due to rotational bad positioning of the tibia, loosening of the reconstruction of medial capsule to the patella, comorbidity, and severe osteoporosis. Infections, often generated, are linked to comorbidity.

The patient should be informed of:

There is a need of a periodic medical monitoring to observe possible changes in the state of the implant and adjacent bone. Only the monitoring can detect possible components loosening or the occurrence of osteolysis; Failure to perform revision surgery, related to osteolysis or components loosening, may result in progressive loss of prosthetic bone stock.

The patient should completely understand and alert that the product does not replace or have the same performance of the normal bone and, therefore, can break, deform or loosen due to excessive efforts or activities, over the permitted ones and if applying excessive load, etc. information regarding the topics: indications, contraindications, side effects and warnings; The patient must be advised of the need to attend regular follow-up examinations after surgery during the period in which the product remains deployed.

Allowable Combination with Other Medical Products

A list of acceptable combinations with BKS III - Patellar Component in accordance with standard NBR ISO 21534 - Implants for Surgery non-active surgical implants – Joint Replacement Implants - Particular requirements was formulated. However, independently of this combination, the implants must be of the same brand because they were designed for such combinations observing the surface finishing, surface treatment and other factors required on projects that may interfere with the combinations. Therefore, metal implants from different manufacturers are not recommended because of chemical, physical, biological and functional incompatibility.

Allowable Contact Leagues

Implant - BKS III - Patellar Component:

Articulation (joint)

- Stainless Steel Conformed (NBR ISO 5832-1 or similar (ASTMF 138);
- Wrought High Nitrogen Stainless Steel (ISO 5832-9);
- Cobalt-Chromium-Molybdenum Casting Alloy (ISO 5832-4) or similar (ASTM F75);
- Wrought cobalt-chromium-tungsten-nickel alloy (ISO 5832-5);
- Forgeable and cold-formed cobalt-chromium-nickel-molybdenum-iron alloy (ISO 5832-7);

COMPONENTE PATELAR

- Wrought cobalt-nickel-chromium-molybdenum-tungsten-iron (I S O 5832-8);
- Wrought titanium 6-aluminum 4-vanadium (ISO 58 323) or similar (ASTM F136);
- Wrought titanium6-4-aluminium 7-niobium alloy (ISO 5832-11);
- Ceramic Materials Based on Alumina (ISO 6474);
- Ceramics Based on Zirconia (ISO 13 356); Wrought cobalt-nickel-chromium-molybdenum alloy (ISO 5832-6);
- Cobalt-chromium-molybdenum casting alloy (ISO 5832-4)

* For implant articulating surfaces of joint replacement, they were considered acceptable in specific applications; however, adequate attention must be given to the project, surface finishing and surface treatment.

In accordance with standard NBR ISO 21534 - Non-active surgical implants -- Joint replacement implants - Particular requirements, the Ultra High Molecular Weight Polyethylene (UHMWPE) is acceptable for the manufacture of implants and can be combined with Cobalt-Chromium-Molybdenum Casting Alloy, raw material of the BKS III - Femoral Component with or without Preservation.

Open Packaging Aseptically

Transport and store it in a cool dry place, with room temperature (max 35 °C) and relative humidity around 30%. Do not store directly on the floor (minimum height = 20cm) and in very high places, close to the bulbs, which could cause dryness or damage to the packaging label. Do not store in places where contaminants, such as, cleaning materials, insecticides, pesticides, are stored etc.

	Do not use if package is damaged
	Keep out of the sun
	Keep dry
	Fragile, handle with care

Product Care Before Use

Before using the product, the responsible person designated by the hospital should check the physical condition of the packaging and product. It is important to verify if the external cardboard packaging is wrapped in transparent PVC film and identified with the product information contained in the package. After this verification the transparent PVC film can be removed and the cardboard package can be opened. The product packed in a double blister should be carefully removed from the external packaging in a surgical environment. The blisters should be opened by pulling the medical grade paper. With product in hand, it is necessary to check if it has any grooves, scratches or stains.

If it occurs any adverse effect that affects the user's security, such as: product is not working, implantable components are damaged, serious problems or death related to these components, the surgeon should report this adverse event to the national health authority and Biomecánica by e-mail: sac@biomecánica.com.br – or phone: (14) 2104-7900 In cases of doubt, the surgeon or health professional can make the communication through the adverse event reporting system on the website of the Health Surveillance ANVISA: http://www.anvisa.gov.br/hotsite/notivisa/index.htm

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the BKS III - Patellar Component explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecánica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings “Improper for use”. According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

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