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

## ENGLISH

**Product Presentation**

The product is packed in double blister packs, manufactured in nontoxic PET (Polyethylene Terephthalate), allowing contact with food and hospital products. The two blister packs are sealed with Tyvek, without dyes, free of holes, through the heat sealing process. The sealed double blister packs are also placed into a cardboard packaging. The implantes are sterilized in ETO with 5 years validity.

**Traceability**

On the packaging of the implantable component there are 6 labels containing the necessary data for implant traceability. All the models of Biomecânica Cruciate-Retaining Femoral Components contain the following information recorded with laser (when possible): Biomecânica Logo - Manufacturing Batch Number - Material Symbol used in the manufacture of the product - Code - Size Symbology, labels and traceability NBRISO 15223

	xxxxxx	Product Lote Number		xxxx-xx	Manufacturing Date
	xxxx-xx	Use by		REF XXXX-XX	Code
		Sterile		STERILE EO	Sterilization in ETO Gas
		Product of Single Use - Do not reprocess			

**Product Description:**

Biomecânica Cruciate-Retaining Femoral Component is a non-absorbable and non-ative implant, for both primary surgery as well as revision surgery, for the femoral-tibial and patelofemoral (knee) replacement or reconstruction with the preservation of the posterior cruciate ligament.

**Biomecânica Cruciate-Retaining Femoral Component**

**BKS III Models**

All Biomecânica BKS III Cruciate-Retaining Femoral Component have polished external surface for articulating with Biomecânica BKSIII Patellar Component. In the intercondylar region, there is not a locking system that requires the substitution of the posterior cruciate ligament during surgery. The inner surface of the component, which will be in contact with the bone cement is sprayed. The outer surface of the component is polished to improve the articulation with Biomecânica BKS III Patellar Component.

Biomecânica BKS III Cruciate-Retaining Femoral Components are sold in several dimensions and in the following sizes: Extra Large, Large, Medium, Small, Extra small. There is a specific model for the right knee and another for the left. This variation enables the surgeon to choose the appropriate model, according to the bone characteristics of each patient and surgical technique used.

**List of Biomecânica BKS III Cruciate-Retainig Femoral Components**

2355-11-000 BKS III Cruciate-Retaining Femoral Component, Extra Small, Right
2355-01-000 BKS III Cruciate-Retaining Femoral Component, Small, Right
2355-03-000 BKS III Cruciate-Retaining Femoral Component, Medium, Right
2355-05-000 BKS III Cruciate-Retaining Femoral Component, Large, Right
2355-07-000 BKS III Cruciate-Retaining Femoral Component, Extra Large, Right
2355-12-000 BKS III Cruciate-Retaining Femoral Component, Extra Small, Left
2355-02-000 BKS III Cruciate-Retaining Femoral Component, Small, Left
2355-04-000 BKS III Cruciate-Retaining Femoral Component, Medium, Left
2355-06-000 BKS III Cruciate-Retaining Femoral Component, Large, Left
2355-08-000 BKS III Cruciate-Retaining Femoral Component, Extra Large, Left

All BKS III models are manufactured from cobalt-chromium-molybdenum molten alloy, according to ASTM F75.

**Ancillary Components**

Biomecânica BKS III Cruciate-Retaining Femoral Components have, as ancillary components, the following items listed below.

Description	Raw Material
Biomecânica BKS III Cruciate-Retaining Articular Insert	Ultra-high-molecular-weight polyethylene (UHMWPE)
Biomecânica BKS III Patellar Component	Ultra-high-molecular-weight polyethylene (UHMWPE)
Biomecânica BKS III Tibial Component	Cobalt-Chromium-Molybdenum Alloy Castings (ASTM F75)
* Acrillex® Orthopedic Cement	Bone Cement (NBR ISSO 5833)

\* It is permitted the use of another type of biomaterial as an ancillary component for Biomecânica BKS III Cruciate-Retaining Femoral Component for primary surgery, replacing the Acrillex® Orthopedic Cement, only if it is appropriate to fix the prosthesis to the bone and is regularly registered at ANVISA for this purpose.

**Complementary Instruments**

presence of an knee arthrodiesis without pain; In younger patients the total knee replacement is indicated only in case of tricompartmentar disease (medial and lateral tibiofemoral, and patellofemoral), and it demands further revision surgery; Morbid obesity, fever or leukocytosis; patient inability to receive anesthesia; leg atheroelctic disease; Skin problems such as psoriasis; neuropathic arthropathy; recurrent urinary tract infections; prevailing cases of osteomyelitis around the knee; neuromuscular, muscular or vascular intense deficiencies in the affected area; bone structures seriously damaged creating severe risk for the application of the prosthesis; Ligamentous system unable to fit the selected component model; bone tumors near the area where the component will be implanted; muscular atrophy, bone immaturity; knee extensor mechanism failure; Structural failure of knee soft tissues; Knee ankylosis; tibia, femur and patella insufficient bone structure; Patients who are unwilling or unable to follow post surgery cares and instructions.

#### Adverse Effects

Joint line migration; Walking standard change; patellar tendon rupture; patellar height decrease or increase, hematomas, wound dehiscence, progressive and asymptomatic bone resorption (osteolysis) around the component; Loosening, migration or fracture of the medical product; Superficial or deep infection, vascular disorders, including thrombosis and pulmonary embolism; Limb shortening because of bone resorption; Allergic reaction to foreign body that may result in histological reactions involving various types of macrophages and fibroblasts; Pain, discomfort or abnormal sensations depending on the presence or migration of the medical product, and also according to the surgical procedure; nerve or neurological damage due to surgical trauma (including paralysis and soft tissue injuries); inability to perform normal daily physical activities; periprosthetic fractures.

#### Warnings

Biomechanica femoral component with preservation should be implanted with the help of Biomechanica instruments. Components of implants from different manufacturers should not be used together; therefore it is recommended that the products have the same provenance.

Biomechanica Cruciate-Retaining Femoral Component is intended to replace damaged joint and not the normal structures of the human body; it does not have a restrained system (locking) that requires the preservation of the posterior cruciate ligament during surgery. Surgical implants must never be reused, and the explanted implants should never be implanted again. Small defects of internal pressure patterns may occur, even if the implant appears to be intact, causing rupture and failure of the implant;

The correct selection of Biomechanica Cruciate- Retaining Femoral Component is extremely important and must be in accordance to the patient's bone characteristics examined by the surgeon. Wear of the ancillary ultra-high-molecular-weight (UHMWPE) polyethylene components may occur due to the poor positioning of Biomechanica Cruciate-Retaining Femoral Component, inadequate ligament balance or aseptic loosening of the prosthesis. For this reason, it is important that the surgeon have sufficient knowledge about the surgical technique used as well as implants and other instruments used; Infections often generated are linked to comorbidity; Aseptic failure of total knee arthroplasty can be caused by components loosening, wear of the ancillary component made of polyethylene, laxity, periprosthetic fracture, arthrofibrosis and patellofemoral complications; The results presented in revision surgeries are lower when compared with the results presented in primary surgeries; laxity, in both flexion and extension, is usually caused by small prosthesis, since they are placed in the very proximal extremity of the femur; Patients who smoke are most likely to go through a longer operative period and have higher costs; the results are not satisfactory in case of articular cartilage absence before surgery. Knee replacement has a finite expected survival that is adversely affected by activity level.

#### The patient should be informed of:

The need for periodic medical monitoring to observe possible changes in the condition of the implant and adjacent bone. Only by monitoring, possible loosening of components or the occurrence of osteolysis can be detected; Not performing revision surgery when the component releases or osteolysis occurs may result in progressive loss of periprosthetic bone stock. The patient must fully understand and be warned that the product does not replace and does not have the same performance of the natural bone. The patient must be advised of the need to attend regular examinations during postoperative follow-up throughout the period in which the product remains in place.

#### Allowable Combination with Other Medical Products

According to NBR ISO 21534 - Non-active surgical implants - Joint replacement implants - Particular requirements, the cast cobalt-chromium-molybdenum alloy, raw material of Biomechanica Cruciate-Retaining Femoral Component, is acceptable for the manufacturing of implants and can be used with implants made of ultra high molecular weight polyethylene (UHMWPE), the raw material of Biomechanica Cruciate Retaining Articular Insert.

#### Allowable Contact Alloys with Biomechanica BKS III Cruciate-Retaining Femoral Component

<b>*Joint</b>	<b>**In contact (when there is no articulation)</b>
- Ultra high molecular weight polyethylene-titanium-based alloy (NBR ISO 5832-3, NBR ISO 5832-11, ASTM F136);	- Cobalt-based Alloy (NBR ISO 5832-4, NBR ISO 5832-5, NBR ISO 5832-6, NBR ISO 5832-7, NBRISO 5832-8, NBR ISO 5832-12, ASTM F75)
- Wrought cobalt-chromium-molybdenum alloy (NBR ISO 5832-12)	- Stainless Steel (NBR ISO 5832-9)

\* For articulating surfaces of joint replacement, they were considered acceptable in specific applications; however, adequate attention must be given to the project, surface finishing and surface treatment.

\*\* For applications where a metal or alloy is in contact with others and the design is not intended, considering that adequate attention is given to the design, surface finishing, surface treatment and metallurgical conditions.

#### Product Care Before Use

Before using the product, the responsible person designated by the hospital should check the physical condition of the packaging and product. The person in charge must verify that the outer cardboard packaging is sealed with a security seal and identified with the information regarding the material present in the packaging. After this verification is sealed with a security seal can be removed and the cardboard package can be opened. The product packed in a double blister should be carefully removed from the external packaging in a surgical environment. The blisters should be opened by pulling the medical grade paper. With product in hand, it is necessary to check if it has any grooves, scratches or stains.

If any adverse effect that affects the user's security occurs, such as: problems with the performance of the product, damage to the implantable component, serious problems or deaths related to these components, the surgeon in charge must report this adverse event to the Health Authority and Biomechanica by e-mail [sac@biomecanica.com.br](mailto:sac@biomecanica.com.br) or by phone +55 014 2104-7926. In cases of doubt, the surgeon or health professional can make the communication of the adverse event through the Brazilian Health Notification System at ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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#### Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little

defects that may lead to contamination or its breakage due to fatigue. This way Biomechanica Cruciate-Retaining Femoral Component explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolución nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

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