CABEÇA BIPOLAR PAR / ÍMPAR

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ENGLISH

Detailed description of the medical product, including the fundamentals of its operation and its action, its content or composition, when applicable, and the list of accessories designed to integrate the product.

Biomecânica Bipolar Head is a unique product that consists of a spherical concave (inner) surface intended to articulate with a femoral component of a total or partial hip joint. It also has a convex (outer) spherical surface which is intended to articulate with the biological acetabulum. The Ø22mm inner diameter is suitable for the perfect fit of Biomecânica Metallic Femoral Head. Its outer diameters vary from 38mm to 55mm.

Biomecânica Bipolar Head surface is made of stainless steel, and it's available in two models: with locking ring, in order to prevent dislocation of Biomecânica Metallic Femoral Head and without the locking ring (interference). The locking ring is made of polyethylene UHMWP 1900 (ISO5834-2), the same composition of the inner part of Biomecânica bipolar head.

Composition: The outer surface of Biomecânica Bipolar Head is made of stainless steel according to ASTM F139. The inner surface is coated with polyethylene UHMWP 1900 (ISO5834-2).

PRODUCT PRESENTATION: The components are packed in double blister packs, sealed with surgical grade paper. The sealed double blister packs are also placed into a cardboard packaging. The products are ready to be used, since they are sterilized with ethylene oxide.

The product traceability labels contain the following information: a) Manufacturer b) Product Commercial Name c) Technical name, d) Batch Number, e) Registry Number at ANVISA, f) Product Code g) Quantity; h) Description of the product, i) Manufacture Date, j) Expiration Date; k) Description of the raw material used to manufacture the product, l) Manufacturer Address; m) Name of Responsible Technician, n) The statement: "Sterile product. **DO NOT REPROCESS!**" Product for Single Use Only; do not reuse"
q) Before using, read the label; r) Warnings / Precautions / Special Care / Symbols; See Statement of Use

The product packaging contains the Safety symbols according to the model below:

 Do not use if package is damaged

 Keep out of the sun

 Keep dry

 Fragile, handle with care

Biomecânica Bipolar Head is sold individually with the following models, according to the table below:

2541-38-000 – Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 38 mm – Stainless steel and polyethylene/ 2541-40-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 40 mm – Stainless steel and polyethylene/ 2541-42-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 42 mm – Stainless Steel and polyethylene/ 2541-44-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 44 mm - Stainless Steel and polyethylene / 2541-46-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 46 mm - Stainless Steel and polyethylene / 2541-48-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 48 mm - Stainless steel and polyethylene / 2541-50-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 50 mm - Stainless Steel and polyethylene/ 2541-52-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 52 mm - Stainless Steel and polyethylene/ 2541-54-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 54 mm - Stainless Steel and polyethylene / 2541-56-000 - Bipolar Head with locking ring inner Ø 28.0mm and outer Ø 56mm - Stainless Steel and polyethylene/ 2540-39-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 39 mm - Stainless Steel and polyethylene / 2540-41-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 41 mm - Stainless Steel and polyethylene / 2540-43-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 43 mm - Stainless Steel and polyethylene/ 2540-45-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 45 mm - Stainless Steel and polyethylene/ 2540-47-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 47 mm - Stainless Steel and polyethylene/ 2540-49-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 49 mm - Stainless Steel and polyethylene/ 2540-51-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 51 mm - Stainless Steel and polyethylene/ 2540-53-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 53 mm - Stainless Steel and polyethylene/ 2540-55-000 - Bipolar Head Interference inner Ø 28.0mm and outer Ø 55 mm - Stainless Steel and polyethylene.

BIPOLAR HEAD ODD/EVEN

Accessories: Biomecanica Bipolar Head does not have any accessory with the purpose of integrating the medical product.

Allowable Combination with Other Medical Products

The description below “allowable combinations with other materials” was based on the standard “NBR ISO 21534 - Non-active metallic implants - Joint replacement implants - Particular requirements”, but regardless of the alloys listed in “allowable alloys”, the implants must also be made by the same manufacturer because they were designed for such combinations observing surface finishing, surface treatment and other factors required on projects that may interfere with the combinations. Therefore, metallic implants from different manufacturers must not be used together because of chemical, physical, biological and functional incompatibility.

Allowable Contact Alloys with Biomecanica Bipolar Head

Wrought Stainless steel (ASTM F138, NBRI505832-1); Wrought high nitrogen stainless steel (NBRI505832-9); Cobalt-chromium-molybdenum casting alloy (ASTM F75, NBRI50 5832-4); Wrought cobalt-chromium-tungsten-nickel alloy (NBRI50 5832-5); Forgeable and cold-formed cobalt-chromium-nickel-molybdenum-rhenium alloy (NBRI505832-7); Wrought cobalt-nickel-chromium-molybdenum-tungsten-iron (NBRI50 5832-8); Wrought titanium 6-aluminum 4-vanadium (ASTM F136, NBRI505832-3); Wrought titanium 6-aluminum 7-niobium alloy (NBRI50 5832-11); Ceramic Materials Based on Alumina (NBRI50 6474); Ceramics Based on Zirconia (NBRI50 13356); Wrought cobalt-nickel-chromium-molybdenum alloy (NBRI505832-9); Cobalt-chromium-molybdenum casting alloy (NBRI505832-4)

* For articulating surfaces of joint replacement, they were considered acceptable in specific applications; however, adequate attention must be given to the project, surface finishing and surface treatment.

Warnings, Restrictions and Precautions
PRODUCT FOR SINGLE-USE. DESTROY AFTER REMOVAL. DO NOT REUSE THE PRODUCT. PROHIBITED TO REPROCESS IT.

An orthopedic implant may be used in a single patient only once. Even if they appear undamaged, they may have small defects and stress patterns which may lead to fatigue failure. Check the sterilization date and if the product packaging is intact. Open the packaging aseptically, after certifying that the dimension is the one chosen. Sterilization validity: 5 years (indicated in internal and external product packaging). The surgeon should be aware of the physical development of the patient, a revision of the implant may be needed.

Contraindications: Osteomyelitis of the upper femur, pyogenic infection of the hip joint, or evident infections are absolute contraindications. Every effort should be made to guarantee the asepis preoperatively in patients who have one or more of the following abnormalities: fever, local inflammation, rapid joint destruction, apparent bone absorption shown in roentgenogram, unexplained high sedimentation rate through other diseases, increase of the WBC rate, or suspected striking changes in the count. Additionally, in distal site infection, such as genitourinary, lung, skin or in other organs, it is relatively contraindicated due to transmission of the infection to the implant site. The focus of infection must be treated before, during and after implant implantation. The use of this implant is contraindicated when there is loss of musculature or neuromuscular infection that compromises the affected limb, which would make the procedure unjustifiable. The endanger of the acetabular cartilage in cases of coxarthrosis is not a contraindication for this procedure.

Warning: improper implant selection, insertion, placement, or excessive muscular activity may cause unusual strain on the implant and may result in subsequent fracture of the prosthesis.

Where there is acetabular failure or loss of acetabular bone stock it is advisable the use of grafts or other procedures to provide support to the acetabular component of the liner. For safe and effective use of this system, the surgeon must be thoroughly familiar with the implantation procedure. Acceptable techniques in postoperative care must be followed.

The patient should be advised about the limitations of partial joint reconstruction and its recent history of use. Excessive physical activity and trauma that affect the partially replaced joint implicate in premature failure of the reconstruction caused by loosening and fracture in the use of implants. Excessive demands on the implant may occur because of excess of activities or for not controlling body weight. Due to its form and characteristics, this product does not have side effect, however, it is suggested to monitor the patient's recovery through X-ray according to the International Society Standards. Do not use magnetic resonance imaging, because the product manufactured in stainless steel causes interference in the examination impairing the diagnosis, and may also cause displacement of the implant causing trouble for the patient.

Adverse Effects: early release of components of partial hip replacement may occur. Premature mechanical loosening may be the result of impaired fixation or hidden infection; late loosening resulting from biological complications or mechanical problems (localized high strain). Reactions caused by metal sensitivity in patients were rarely reported.

The meaning and effect of raising expected future evidence and clinical evaluation. Implantation of foreign material in tissues may result in histological reactions involving various sizes of macrophages and fibroblasts. Similar changes may occur as a precursor, or during the healing process. Dislocations and subluxation were informed, as the result of improper positioning of implant components. Muscle and fibrous tissue laxity can also contribute to these conditions.

Information to be provided to the Patient: Talk to your orthopedic surgeon about any medical condition that may affect the surgery. Among the complications of partial arthroplasty that may arise are:

Infection: Infection may occur in the wound or it can be deep (around the prosthesis). It may appear while patients are in hospital or at home. It may appear even years later. Small surgical wound infections are usually treated with antibiotics. Large or deep infections may require further surgery (for deep cleansing) or even the removal of the prosthesis. Any infection in your body (bladder, throat, teeth, ears, etc.) may take germs through the blood to your prosthesis and cause infection.

Thrombosis: is the formation of a blood clot inside a blood vessel, obstructing the flow of blood through the circulatory system. Blood clots are the result of several factors, including its decreased mobility, and cause blood to flow more slowly in the veins, which may facilitate the development of thrombosis. A common sign of blood clots is pain and edema in the thigh or calf. If this occurs, your orthopedist will do tests and exams to evaluate the veins of your leg and may require monitoring and evaluation of a vascular surgeon. Several measures should be used to reduce the possibility of blood clots - Medication to thin the blood (anticoagulants) - Elastic Stockings - Exercise to improve circulation in the legs - plastic boots that inflate and compress the soles of the feet and calf, increasing venous return. Even with the use of these preventive measures, blood clots may still occur. If edema (swelling), redness or pain in the calf after leaving hospital is noticed, you should contact your orthopedic surgeon.

Loosening: Prosthesis loosening may occur after surgery. This can

cause pain if loosening is significant. A revision surgery (implant exchange) may be required. New materials and new methods of attachment should minimize this problem.

Luxation: Eventually, after a hip prosthesis, the head of the prosthesis may come out (dislocate) of the acetabulum. Most of the times, the dislocated femoral ball may be pulled back in place without operation. To prevent dislocation, muscles must be strong (do exercises recommended by your orthopedic surgeon) and the leg must not be flexed above the hip (over 90 degrees) in the first month. **Wear:** Some type of wear may be detected in any type of prosthesis. Excessive wear may contribute to loosening and may require further surgery (Revision Surgery).

Prosthesis Breakage: Once implanted, the biomaterials used in prostheses necessarily come into contact with bodily fluids. These fluids, while seemingly harmless, can, over time, significantly deteriorate the majority of materials of considerable “chemical inertness.” Additionally, most orthopedic implants are subjected to static and / or cyclic mechanical strain, which may cause the failure of the material, that is, the breakage of the prosthesis. With the materials currently used in implants, fracture of the implant is rare. However, if this occurs, further surgery to replace the fractured prosthesis will be required.

Implanted Product Evaluations: After implantation, in the intraoperative, the professional in charge must perform a radiological control to verify the correct positioning of the product. It is this professional's responsibility to make the clinical and radiological evaluation after the surgical procedure in the frequency stipulated by him/her to check the condition of the implant and the bone development. If the product is in the wrong position or presents any non-compliance, it is the surgeon's responsibility to take the appropriate corrective action. Useful information to avoid risks derived from implantation: To reduce the implantation risks the instructions contained in the Contraindications, Instructions for use and “Product Instructions for Use” must be strictly followed.

Implantation Intrinsic Risk: Biomecanica Bipolar Heads are made of biomedical materials. They are: Stainless Steel, according to the requirements of ASTM 139 and UHMWP 1900 polyethylene, according to ISO5834-2.

Sterilization: This product is sterilized by Ethylene Oxide Gas ETO. Its sterilization expires in 5 years. If the validity of sterilization is overdue or packaging is damaged, the product should be returned to the distributor or directly to Biomecanica.

Product discard
No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way Biomecanica Bipolars Heads explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings “Improper for use”. According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classed as single use.

Contamination Risks

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

TRACEABILITY: On the packaging of the implantable component there are 6 labels containing the necessary data for implant traceability. The information contained on this traceability label is related in the “Product Presentation” of this technical report.

MARKING: Biomecanica Bipolar Heads contain the following information recorded with laser to trace and monitor the patient after surgery: Biomecanica Logo / Manufacturing Batch Number / Symbol of the raw material used in the manufacture of the product. Code / Size

Post-Sales (Customer Complaint): If there is a need to make a complaint about Biomecanica Bipolar Head related to any adverse effect that affects the security of the user, such as problems with the performance of the product, damage to the implantable metal component, serious problems or deaths related to these components, the surgeon in charge must report this adverse event to the national health authority and Biomecanica by e-mail sac@biomecanica.com.br or by phone 55 14 2104-7926. In cases of doubts, the surgeon in charge or health professional can make the communication of the adverse event through the Public Health Notification System at ANVISA website: <http://www.anvisa.gov.br/hotstite/notivisa/index.htm>

Sending the Product for the manufacturer's Analysis: when the implant is sent for the manufacturer's analysis, it should be cleaned at the hospital using a bactericidal and broad-spectrum antifungal solution. Then, it should be disinfected or sterilized in a steam autoclave or with ethylene oxide. It should be sent to Biomecanica in an intact packaging, identifying the method of cleaning, sterilization, and product data.

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