cardboard box. The boxes are properly identified by adhesive tags that contain all legal information that is pertinent to the product for the guarantee of a complete identification and traceability of the product. The package must be intact at its receipt (do not use the product in case the packaging is violated). Check the expiration date of the sterilization (do not use the product in case the expiration date of the sterilization is expired).

Disclaimer: Biomecanica, as the manufacturer of this device, does not practice medicine nor recommends it or any other surgical technique to be applied to a certain patient. The surgeon who performs any implant procedure is responsible for determining and applying the proper techniques for the implantation of the prosthesis in each patient. Biomecanica is not responsible for the choice of a proper surgical technique to be applied to each patient. The surgical techniques depend on the techniques adopted by the surgeon who is in charge of selecting the method of implantation, the type and dimensions of the medical devices.

Anclary Components: The anclary components listed below must be bought separately because they do not integrate with the product in question. The anclary components listed below must be bought separately because they do not integrate with the product in question.

Accessories: Unrestricted polymeric liners (cementless) – BM do not have any accessory with the purpose of integrating with the medical device.
Allowable combinations with other materials: The table Bellow "Allowable contact implants" was based on "NBRISO21534 – Non-active surgical implants – Joint replacement implants – Particular requirements", but independently of the alloys mentioned in table "Contact allowable alloys of joint", the combinations must be between implants of same brand because they have been designed for such combination based on the surface treatment and finishing and other factors or combinations required by the design that may interpose the combinations. Therefore, metallic implants from different manufacturers are not recommended to be used together due to chemical, physical, biological and functional incompatibility.

Contact allowable alloys of joint
Wrought stainless steel (ASTM F138, NBRISO5832-1)
Wrought High Nitrogen stainless steel (NBRISO5832-9)
Cobalt-Chromium-molybdenum casting alloy (ASTM F75, NBRISO5832-4)
Wrought cobalt-chromium-tungsten-nickel alloy (NBRISO5832-5)
Forgeable and cold-formed cobalt-chromium-nickel-molybdenum-iron alloy (NBRISO5832-7)
Wrought cobalt-nickel-chromium-molybdenum-tungsten-iron alloy (NBRISO5832-8)
Wrought titanium 6-aluminum 4-vanadium alloy (ASTM F136, NBRISO5832-3)
Wrought titanium 6-aluminium 7-niobium alloy (NBRISO5832-11)
Ceramic materials based on alumina (NBRISO6474)
Ceramic materials based on zirconia (NBRISO 13356)
Wrought cobalt-nickel-chromium-molybdenum alloy (NBRISO5832-6)
Cobalt-chromium-molybdenum casting alloy (ASTM F75, NBRISO5832-4)
Wrought cobalt-chromium-molybdenum alloy (NBRISO5832-12)

*For Joint surfaces of Implants and joint replacement they were considered acceptable to specific applications since full attention is given to the project, surface finishing and surface treatment.

Unrestricted Polymeric Liner (cementless) – BM should be only used with Biomecanica instruments.

Precautions, restrictions and warnings
Precautions: SINGLE USE PRODUCT. DESTROY AFTER EXPLATED. DO NOT RE-USE THE PRODUCT. DO NOT REPROCESS.
An implant should never be re-used. Although it seems not to be damaged, previous stress may have created imperfections that may reduce the useful life of the implant.
Check if the packaging of the product is intact. Verify the sterilization date. Aseptically open the packaging after certifying that the dimension is the one that was chosen. Sterilization Validity: 5 years (indicated in the inner and outer packaging of the product).
The surgeon must be aware of the possibility of physical deformation and that may require a revision of the implant.

Contraindications: The total hip arthroplasty is a powerful procedure and is linked to a large number of implications and to 1% to 2% of the mortality rate. This way, when a total hip arthroplasty is indicated, the patients must be carefully evaluated specially in cases of Systemic disorders. Heart and liver disease, pulmonary, genitourinary or metabolic disease, hypertension and malignancy must be treated and controlled before the total hip arthroplasty. Total hip arthroplasty is also contraindicated in other situations such as: Osteomyelitis of the upper femur, hip joint pyogenic infection or evident infection are absolutely contraindicated. All effort is necessary in order to govern the possibility of preoperative asepsis in patients with one or more than the following abnormalities: - signs of fever and local inflammation; - rapid joint destruction and apparent bone absorption in roentgenograms; - increase of the sedimentation rate - unexplained through other diseases; increase of the WBC count or an accentuated change in differential count. Uncommon circumstances, such acute or low-level infection followed by the replacement of total joint or removal of the prosthesis may be treated by professionals, complete debridement and reimplantation of a new device. However, there is a higher probability of recurrence of infection in arthroplasty revision than in recent cases. In no circumstances the insertion is recommended if: -there has been an extensive bone involvement; the involved soft tissues can not be removed; where the body is gram-negative and not sensitive to Staphylococcus aureus or where highly virulent or mixed bodies are involved.
Additionally, the implantation in distant focus of infection such as the genitourinary system, lung, skin or other places is relatively contraindicated due to the chance of hematogenous transmission in the local of implantation. The focus of infection may be treated before, within or after the implantation.
The implantation of this device in areas of muscle loss or neuromuscular infection compromising the affected limb is contraindicated because it would make the procedure unjustifiable.

Warnings: Improper implant selection, insertion, placement or excessive muscle activity may cause uncommon stress to the implant and may result in subsequent fracture of the prosthesis. In areas where there is insufficiency or loss of acetabular bone stock, graft or other reinforcement procedures are recommended to give support to the shell and retention of the liner.
Due to laboratory data for the evaluation of the rates of use unavailable for other materials, it is recommended that only acetabular components molded in UHMWP polyethylene are used with this device. For the security and effective use of this system, the surgeon must be familiar with the implant procedure.
Acceptable techniques must be followed in the postoperative cares. The patient must be aware of the limitations of the total joint reconstruction and recent history of use. Excessive physical activity and trauma affect the replaced joint and lead to the early failure of the reconstruction due to looseness, fracture and use of the implants.

Overload on the replaced joint through overactivity or no control of the body weight.

Adverse effects: Potential early looseness of the components in the total hip replacement; early mechanical looseness that may be the result of defective fixation or occult infection; Biological implications or mechanical problems (high loaded stress); Peripheral neuropathy and formation of heterotrophic bone in Total hip arthroplasty and damage to subclinal nerve which is more frequent and possibly due to a surgical trauma.

As has rarely been reported sensitivity reactions to metal in patients followed by total joint replacement. The meaning and effects of sensitivity wait for evidences and future clinical evaluations.

The implantation of foreign tissues may result in historical reactions involving several types of macrophage and fibroblasts. Similar changes may occur during or as a precursor to the healing process.
Displacement and subluxation have occurred due to poor positioning of the implants. Muscles and looseness of fibrous tissues may also contribute to those conditions.

Information to be provided to the Patient: Inform your orthopedist of any medical condition that may affect the surgery. Total prosthesis surgeries have been successful in over 90% of the patients. When implications occur, they are usually successfully treated. Among the potential implications, we have:
Infection: the infection may occur in the surgical wound or it may be deep (around the prosthesis). It may occur while the patient is in hospital or at home. It may also occur years later. Minor infections in the surgical wound are usually treated with antibiotics. Major or deep infections may require a new surgery (for a deep cleaning) or the prosthesis removal. An infection in the body (Urinary bladder, throat, teeth, inner ear, etc.) may allow the germs to get to the prosthesis through the blood and cause infection.
Thrombosis: Formation of blood clots resulted from several factors such as their mobility reduction that make the blood to slowly circulate in the vessels which facilitates the appearance of thrombosis. Blood clots should be suspected in case of pain and edema in the thigh or calf. If that occurs an orthopedist must make tests and inspections in order to evaluate the vessels of the legs. The evaluation and follow up by a vascular surgeon may be necessary. Several measurements must be taken to minimize the chance of blood clots: - Blood-thinning medicines (anticoagulant); - Elastic stockings; - Exercises to improve the circulation in the legs; - Plastic boots that inflate and compress the sole of the foot and calf improving the venous return.
Even with those preventive measures, blood clots may still occur. In case of edema (swelling), redness or pain in the calf after leaving hospital an orthopedist must be informed.

Looseness: The looseness of the prosthesis inside the bone may occur after the surgery and that may cause pain if the looseness is significant and a revision surgery (change the prosthesis) may be necessary. New materials and methods of fixation should minimize this problem.
Luxation: eventually, after a hip prosthesis, the head of the prosthesis may get out of the acetabular cup (hat). That is called Luxation. In most of the cases the hip may be replaced with no need of a revision surgery. In order to prevent luxation it is important that the muscles are strong (Practice the exercises recommended by the orthopedist) and that the hip is not flexed more than 90 degrees during the initial months.

Wear: wear may be detected in any type of prosthesis. An excessive wear may contribute to the looseness and that may require a revision surgery.
Prevention Breakage: Once the biomaterials used in the prosthesis are implanted, they necessarily get in contact with the body fluids. Even though those fluids seem

to be harmless they can significantly degrade most of considerable "chemical inertness" materials over the years. Additionally, a large number of the orthopedic materials is submitted to mechanical, static and/or cyclic stress that may lead to the failure of the material, in other words, the breakage of the prosthesis. The fracture of the implant with the material used nowadays is unlikely. However, in case it occurs, a new surgery to replace the fractured prosthesis will be necessary.

Nerve injury: nerves adjacent to the prosthesis may be injured during a surgery, however it is very rare. Such injury is more likely to occur when the surgery involves resections of large deformities in the pelvis or elongation of a shorter leg due to wear over the years those nerve injuries get better and are totally recovered. Eventually it is possible to choose for the surgery of the nerve involved.

Product discard
No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way Unrestricted Polymeric Liner (cementless) – BM explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Useful information to avoid Risks due to implantation: In order to minimize the risks due to implantation, the following must be strictly followed: contraindications, instructions for use and all information in the "instruction for Use" of the product.

Intrinsic Risks of Implantation: Unrestricted polymeric liners (cementless) – BM are made of materials of well known biomedical use. And NBR ISO 5834-2 – Implants for Surgery – Ultra-high-molecular-weight Polyethylene-Part2: Molded forms.

Risks of Contamination: There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because Unrestricted Polymeric Liners (cementless) – BM get in contact with the tissue and body fluids. The explanted products must be treated as highly contaminants.

Sterilization: This product is sterilized in Ethylene Gas - ETO and the sterilization is valid for 5 years. In case the sterilization date is expired or the packaging is damaged, the product should be returned to the distributor or directly to Biomecanica.

Product Discard: The Implants part of Unrestricted Polymeric Liner (cementless) – BM and which are explanted must be properly discarded by the hospital. The complete mischar-acterization of the implant is under responsibility of the hospital in order to prevent its re-use. The method used for the mischaracterization of the implant is under the responsibility of the hospital. Biomecanica recommends the implants to be mechanically deformed with the help of a hammer or press and then identified with the sayings "Improper for Use".
As described in Resolution n° 2605 of 11/08/2006, implantable devices of any origin classified as single use ones can not be reprocessed.

Traceability: There are 4 tags that come with the product. Those tags contain the necessary data for the traceability of the product. One tag is mandatorily pasted on the patient's report; one should be handed to the patient; one on the tax documentation that generated the charging – AIH- in case the patient is supported by SUS or on the invoice in case the patient is supported by a complementary health system; one should be available for the control of the supplier (distribution historical record – RHD); one should be available for the control of the surgeon in charge (the main one) and the last one for the medical insurance in case there is one.

Marking: Unrestricted acetabular Liner (cementless) – BM comes with the following engraving information by machining in order to allow the traceability and follow up of the patient in the postoperative: - BM acronym; Number of the manufacturing lot; -Dimension.

Post Sale (Client Complaint): In case it is necessary to make a complaint about Unrestricted Acetabular Liner (cementless) – BM related to some adverse event that affects the safety of the user such as if the product is out of order or in case the implantable component is damaged or the occurrence of serious problems or death related to those components the surgeon in charge must inform the health authorities and Biomecanica of the adverse event through the email sac@biomecanica.com.br or by the phone -0xx14 2104 7900. In case of doubts the surgeon in charge or the health professional may report the adverse event through the System of notification and Health surveillance on the website of ANVISA: http://www.anvisa.gov.br/hotsite/notifica/index.htm

Sending of Material to the Manufacturer for analysis: In case the implant is sent to the manufacturer for analysis it must be sanitized at the hospital with the help of bactericide and antifungal solution of broad spectrum. Then it must be disinfected or steam sterilized in autoclave or ethylene oxide. They must be taken to Biomecanica in intact packaging and identified with the method of cleaning and sterilization and the data of the product.

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