FIXADOR EXTERNO

EXTERNAL FIXATOR

Nome Técnico: Fixadores Externos

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ENGLISH

Technical Name : External Fixator

Commercial Name: Bio Tub External Fixator Line – BM – Non-sterile

Product Description

The external fixators have two basic functions: a) Fuse the fractured parts of diaphyseal or epiphyseal bone tissue. That will promote dynamic compression or distraction or static stabilization between the fragments. b) Promote the elongation of the bone tissue in case of shortening of it. The external fixator is a medical product which is not implantable and it is designed to help the bone fixation and intended to: minimizing, stabilizing and fixing several types of fractures in the lower and upper limbs as well as in the human skeleton according to the indication for the use. It also has the purpose of transferring the load applied to the fractured bone to the device which keeps the bone free of load and able to regenerate after the fracture. The external fixators present a structure designed exclusively for tubes and connectors which are connected by wires and pins. The external fixators are available in different designs and with proper methods of application besides general principles that direct their applications. The wires and pins are ancillary components for the external fixators, therefore they are validated in another process of registration and sold separately.

Manufacturing Materials

The external fixators are made of materials that are well known in the medical field and that meet ASTM Standards for the steels and NBR Standards for the aluminum. Two types of raw-materials are used – ASTM F899 for 304 typed stainless steel; ASTM F899 for 420 typed stainless steel and NBR 6834 for 6065-T5 aluminum alloy. The table below brings the codes, the components of the set and the raw-material of each component of the fixator.

Table - List of external fixators

Code	Components
50064-00-000	Bar-Pin connector set
50067-00-000	Bar-bar connector set
50068-00-000	Squeeze key
50069-27-000	Drill Ø2,7mm x 220mm
50070-10-000	Tube Ø9,5mm x 100mm
50070-11-000	Tube Ø9,5mm x 110mm
50070-12-000	Tube Ø9,5mm x 120mm
50070-13-000	Tube Ø9,5mm x 130mm
50070-14-000	Tube Ø9,5mm x 140mm
50070-15-000	Tube Ø9,5mm x 150mm
50070-16-000	Tube Ø9,5mm x 160mm
50070-17-000	Tube Ø9,5mm x 170mm
50070-18-000	Tube Ø9,5mm x 180mm
50070-19-000	Tube Ø9,5mm x 190mm
50070-20-000	Tube Ø9,5mm x 200mm
50070-21-000	Tube Ø9,5mm x 210mm
50070-22-000	Tube Ø9,5mm x 220mm
50070-23-000	Tube Ø9,5mm x 230mm
50070-24-000	Tube Ø9,5mm x 240mm
50070-25-000	Tube Ø9,5mm x 250mm
50070-26-000	Tube Ø9,5mm x 260mm
50070-27-000	Tube Ø9,5mm x 270mm
50070-28-000	Tube Ø9,5mm x 280mm
50070-29-000	Tube Ø9,5mm x 290mm
50070-30-000	Tube Ø9,5mm x 300mm
50070-31-000	Tube Ø9,5mm x 310mm
50070-32-000	Tube Ø9,5mm x 320mm
50070-33-000	Tube Ø9,5mm x 330mm
50070-34-000	Tube Ø9,5mm x 340mm
50070-40-000	Tube Ø9,5mm x 400mm
50071-00-000	Schanz pin guide
50071-01-000	Drill guide Ø3,3mm
50069-35-000	Drill Ø3,5mm x 220mm
50072-00-000	Wrench 10/13

Symbology of the Traceability tag (Label) according to NBRISO15223

	No Sterile		REF XXXXX Code		Validity
	LOT XXXXX Lot n°		M Manufacturing date		Manufacturing data

TECHNICAL INFORMATION

Indications, Precautions, Restrictions/Warnings, Adverse effects and Contraindications: The external fixator is indicated for bones of lower and upper limbs in corrective surgeries or in the treatment of traumas because the fixators keep the reduction of the fractured bone and stabilizes the bone or adjacent structures in order to promote or correct the bone fusion. The surgeon must be qualified and familiar with osteosynthesis and its limitations including the pre and postoperative, surgical techniques to be adopted, precautions and potential risks. Do not use the product if it is damaged. The pre-operative and surgical procedures as well as the knowledge about the surgical technique and proper selection and implantation of the implant are important considerations for a successful use of the provisional devices for external fixation. Do not use the product in case the packaging is damaged. Some potential adverse effects that may occur in the use of the external fixators are: postoperative superficial or deep infection; pain;

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B046B - REV03 - 28022022

Modular Fixator

discomfort or abnormal sensations due to the presence or migration of the medical device...

Sterilization of the external fixators: The external fixators are distributed not sterile. Before using the external fixators they should be sterilized.

AUTOCLAVE: It is a device for sterilization through saturated steam under pressure.

INSTRUCCIONES: Before starting the device check: if the switch is on, if the water valve is open and the flush valve is closed.

1 - open the door of the equipment; 2 - properly settle the material to be sterilized; 3 - close the door of the equipment; 4 - select the desired cycle according to the material to be sterilized...

Table - Parameters of physical sterilization in autoclave

Table with 4 columns: Ciclo, Temperature, Exposure Time (minimum), Time to Drying (maximum). Rows include Conventional (1 atm of pressure) and High Vacuum.

Note.: Register the time when the heat of the sterilization chamber reaches the desired temperature.

For further information check the instruction manual that comes with each autoclave.

Other method of sterilization: Besides the autoclave the sterilization process may be in Ethylene Oxide (E.T.O.) - parameters and procedures established by the validation protocol and ISO 11135-1.

NOTE: The methods of sterilization, the equipments, the controls and the instructions on the applied sterilization under responsibility of the hospital. Storage conditions of the sterilized articles: - as the environment: it must be cleaned; ventilated and dried and it should be restricted to the team of the sector...

Storage and transport of the product: The storage location should be clean, dry and illuminated in a way that the proper conditions of storage and transport are preserved as well as the physical and chemical integrity.

Identification: The external fixators come with labels that present the following information: Data of the manufacturer; Technician in charge; nº of the register with ANVISA; Technical name of the product; Commercial name of the product; Commercial model of the product...

Customer Service: In case it is necessary to make a complaint about the use of the external fixators related to some adverse effect that may affect the safety of the user the responsible surgeon must inform the Health Authorities and Biomecnica through the email - e-mail: sac@biomecnica.com.br...

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