COMPONENTE TIBIAL CIMENTADO - BM

PORTUGUÊS

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CEMENTED TIBIAL COMPONENT - BM

In one study, the reasons that led to a revision surgery were observed, and these:

- Pain and disability due to aseptic loosening of the tibial component, inappropiate ligament balancing and positioning of the component, polyethylene fracture, or dislocation of the patella;
- Dislocation of the patella due to wrong positioning of the tibial component or rupture of the reinserting the patella medial vastus;
- Fracture of the patella;
- Acute or chronic joint infection.

Use Instructions: Note: This implant is provided sterile. It should be handled in appropriate environments and with due care. Only trained professionals should handle and deploy these devices. Surgical techniques vary according to the choice of the surgeon, leaving him the final choice of method, type and size of products to be used, and the criteria for evaluating the results of surgery. Cemented Tibial Component - BM should be applied and adapted in accordance with the requirements and appropriate surgical techniques. The mechanical contribution of the implant should be examined by observing the exact adaptation and compression, since a stable reconstruction of the fractured bone minimizes the load to be borne by the implant. Weight, activity level and patient condition may impact the performance of the implant. It must be considered the fact that no implant is as strong as the natural bone and therefore has limitations as to the biomechanical demands.

Contraindications

Contraindications below do not limit the use of Cemented Tibial Component - BM models, but it is considered the use of reorganized biological use, which follows: Infections or systemic diseases or near the joint; Knee reartrosis after arthroplasty; prior patellectomy; Presence of recurrent deformity; Neuropathic arthritis; Discontinuity or severe dysfunction of the extensor mechanism; In younger patients the total knee replacement is indicated only in the case of rheumatoid arthritis tricompartmental and lateral tibiofemoral, and patellofemoral, requiring a future revision surgery; Morbid obesity, fever or leukocytosis; Inability of the patient receiving anesthesia; Atherosclerotic disease of the leg; Skin problems such as, psoriasis; Neuropathic arthropathy; Recurrent urinary tract infections; History of osteomyelitis near the knee; Neurovascular, muscle or vascular intense disabilities in the affected region; Severe risk of bone structure generating infeasible risk for application of the implant, when the ligament system is severely affected; Choice of component model; Bone tumors near the site of component implantation; Muscular atrophy, bone immaturity; Failure of the extensor mechanism of the knee; Crash of structural knee soft tissue; Knee ankylosis; Insufficient bone structure in the tibia, femur and patella; Patients unwilling or unable to follow postoperative care and instructions.

Precautions, restrictions, warnings, special care and clarification about the use of medical products, as well as its storage and transportation.

Precautions and Warnings

PRODUCT FOR ONE USE ONLY. DESTROY AFTER EXPLANTED. DO NOT REUSE THE PRODUCT. REPROCESSING PROHIBITED.

An orthopedic implant can only be used on a single patient only once. Although it may not seem to be damaged, prior tensions can create imperfections that may reduce the success of the implant. Improper selection of implant can cause unusual tension on the implant and may result in subsequent fracture of the implant. The surgeon must be familiar with and hold sufficient knowledge, including pre and postoperative surgical technique used, precautions and potential risks. During handling of the implant should be avoided scratches, nicks, or anything that may damage or mark the implant, because these defects are tension concentrators and can be sites of crack nucleation and decrease the corrosion resistance and can result in fracture of the implant. Check if the package containing the product is intact, Check the date of sterilization. Aseptically open the package and make sure that the size is one that is chosen. Validity of sterilization: 5 years (indicated in the inner and outer product packaging). To ensure proper implant, the instruments used must be manufactured by Biomechanics, because they were designed and manufactured specifically for the use of these implants. Using instruments from different manufacturers may compromise the surgery besides being in disharmony with product registration in the National Agency for Sanitary Vigilance.

The surgeon must be aware of the possibility of the patient's physical development. Do not use the product if it is damaged.

PRODUCT FOR ONE USE ONLY - IS EXPRESSLY PROHIBITED RE-USE OF THIS PRODUCT.

Special storage conditions: Store in a cool, dry place away from light and away from the action of weathering. Do not use the product if the package is damaged.

Note: Not to be used implant components from different manufacturers. We therefore recommend that the products have the same origin.

Manufacturing date, expiry date and product batch: SEE LABEL.

Adverse Effects: Adverse effects that may occur with the use of cemented tibial components are: Migration of the joint line; Alteration of gait pattern; Patellar tendon rupture; Decrease or increase in patellar height; Bruising; Suture dehiscence; Progressive bone reabsorption (osteolysis) localized and asymptomatic around the component; Loosening, migration or fracture of the medical product; Superficial or deep infection; Vascular thrombosis and pulmonary embolism; Allergic reactions. Shortening of limb due to bone reabsorption; Allergic reaction to foreign body and may result in histological reactions involving various types of macrophages and fibroblasts; Pain, discomfort, or abnormal sensations due to the presence or migration of medical product, and also due to the procedure surgical; Neurological or neural damage due to surgical trauma (including paralysis and soft tissue lesions); inability to perform normal daily physical activities; periprosthetic fractures.

Warnings: Family of Cemented Tibial Component - BM aims to replace the damaged joint and not normal structures of the human body. Surgical implants should never be reused, and explanted implants should never be implanted again. Small defects and internal pressure changes can occur, even if the implant appear intact, causing failure and implant disruption. The correct selection of models that comprise the family of Cemented Tibial Component - BM is extremely important according to bone characteristics of the patient, examined by a surgeon, and the surgical technique adopted by the surgeon. Infections often generated are linked to comorbidity. The aseptic failure of total knee arthroplasty can be caused by components relaxation, attrition of ancillary manufactured component of polyethylene, ligamentous laxity, periprosthetic fracture, arthrofibrosis and patellofemoral complications. The results presented in revision surgeries are lower compared to the results presented in primary surgeries. The loosens, both in flexion and in extension, is usually caused by insufficient size of the implanted prosthesis very close to the femur. Smokers is very likely to have longer operation times and higher costs. The results are more reserved in case of absence of articular cartilage prior to surgery. A knee replacement has a finite expected survival, adversely affected by the level of activity.

Patient should be informed about: The need for periodic medical consultation, to observe possible changes in the status of the implant and the adjacent bone. Only monitoring can detect possible loosening of components or occurrence of infection. Not performing revision surgery when the component loosening or osteolysis can result in progressive loss of periprosthetic bone stock. Should be made fully understand and warn the patient that the product does not replace and does not have the same performance as normal bone. Related information on the topics: indications, contraindications, side effects and warnings. The patient should be instructed to inform the surgeon when symptoms occur. The patient should be instructed to inform the surgeon when symptoms occur. Not performing revision surgery when the component loosening or osteolysis can result in progressive loss of periprosthetic bone stock. Should be made fully understand and warn the patient that the product does not replace and does not have the same performance as normal bone. Related information on the topics: indications, contraindications, side effects and warnings. The patient should be instructed to inform the surgeon when symptoms occur. The patient should be instructed to inform the surgeon when symptoms occur.

Disclaimer: Biomechanics, the manufacturer of this device, does not practice medicine or recommend this or any other surgical technique for use in a particular patient. The surgeon performing any implant procedure is responsible for setting and using the proper technique to implant the prosthesis in each patient. Biomechanics is not responsible for choosing the appropriate surgical technique to be used in each patient. Surgical techniques depend on the technique used by the surgeon, being this the choice of method, type and size of products to be placed.

Evaluations of Product Implemented: After implant, in intraoperatively, the responsible professional must perform radiological control to check the correct positioning of the product. The responsible professional shall make, and it is your responsibility, clinical and radiographic examinations passed after the surgical procedure in the frequency stipulated by him to check the status of the implant and the evolution of bone. If the product is out of the correct position, or provide any non-

compliance is the responsibility of the surgeon to take the most appropriate corrective action.

Useful Information to Avoid Risks Arising from Implantation: To decrease the risks of implantation it should be followed strictly: contraindications, instructions for use and all information contained in the "Instructions of Use" of the product.

Intrinsic Risk of the Implantation: Cemented Tibial Component - BM are manufactured with material of reorganized biological use, which follows the requirements of ASTM F757-07 rule Standard Specification for Cobalt-28 Chromium-6 Molybdenum Alloy Castings and Casting Alloy for Surgical Implants (UNS R30075). Contamination Risks: There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because Cemented Tibial Component - BM come in contact with body tissues and fluids. Explanted product should be treated as highly infective.

Sterilization: This product is supplied in sterile form. As the client request, it can be supplied with the following sterilization processes:

Process	Validity
Sterilization by ethylene oxide gas (ETO)	5 Years
Sterilization by Gamma Radiation (Gamma Ray)	5 Years

ETO process

Note: The differentiation of the sterilization type will be defined by the product code. Product sterilized by ETO gas have different codes of products sterilized by Gamma Radiation. If the product is with the sterilization expire date overdue and/or packaging damaged, it must be returned to the responsible distributor or directly for Biomechanics.

Sterilized Material Handling: When handling the material with sterile aseptic technique, it must obey some rules in order to keep it sterile:

- It is fundamental wash hands with soap and water before handling sterilized equipment;
- use material with intact packaging, dry, without spots, with identification (material type and sterilization date);
- working in front of the material;
- manipulate the material at or up the waist level;
- avoid coughing, sneezing, talking about the exposed material;
- not make moves on the sterile area;
- make sure the validity and the integrity of the material;
- work in a clean, quiet, dry and draft-free area;
- maintain a certain distance between the body and the material to be manipulated;
- comply with the other principles of asepsis.

The nursing technique advocated in sterilized material handling is: open it, starting with the opposite side of the hand; touch the exposed material with the sterile field that involves it; touch with hands only on the outside of the package;

do not keep as sterile material previously opened package;

Warranty: The warranty Cemented Tibial Component - BM family is related to component. It does not include information about the use of the product, its usage information, contraindications, warnings and precautions, warnings, possible side effects, packaging, sterilization, recommended cleaning and decontamination in this technical report and in instruction of use.

Product Disposal: Cemented Tibial Component - BM explanted of the patient must be properly disposed of by the hospital. It is the responsibility of the hospital full implant distortion preventing its re-use and method for distortion of the implant. Biomechanics recommends that the explanted implants are deformed mechanically with a hammer or treated with autoclave or ethylene oxide gas. The explantation process should be performed in a sterile environment. The patient should be informed of the risks associated with the use of this product. The patient should be informed of the risks associated with the use of this product.

Traceability: Along the packaging of implantable component follow 6 traceability labels containing the data of implant used. This label is the same label contained in the external labeling of the product. The information contained in this traceability label are listed in form "Presentation Form" of this technical report. One label must be strictly adhered to the patient clinical record and another in the report delivered to the patient. Other in the tax documents that creates billing, in AIH, in the case of patients treated by SUS, or in the bill of sale, if patient treated by complementary health system, other available to provider control (historical distribution record - RHD, other available for the control of responsible surgeon (principal) and the last for health insurance if any. We inform that patient's records must have detailed information to trace the implanted product. Collage of this label in the record allows traceability of the product used. Among the main information we highlight as essential the name of the manufacturer of the implant used, name of implant used, the implant code, batch number and product registration number at ANVISA. Such information are described in the traceability labels that follows with product and in the external labeling. Additional information should also be considered essential, as the surgery date, name of the patient receiving the implant, surgeon name, patient weight, patient age and other information requested in the patient record must also be fulfilled.

Marking

Cemented Tibial Component - BM contain the following information laser engraved, to allow traceability and monitoring the patient postoperative: Logo of Biomechanics - Batch Manufacturing Number; Acronym of the raw material used to manufacture the product (CA); Code; Dimension

After Sales (Customer Complaint): If there is need for a complaint of Cemented Tibial Component - BM family related to any adverse event that affects the safety of the user, as the product does not work, damage to the metallic implantable component, serious problems or death related to these components the surgeon must communicate this adverse event to the competent health organ and Biomechanics by email sac@biomecanica.com.br or by phone 0xx14 2104 7900. In cases of doubt the surgeon or health professional can make the communication of adverse event through Notification System for Sanitary Surveillance in ANVISA site: http://www.anvisa.gov.br/hotste/motivsa/index.htm

Submitting Material to Manufacturer Analyze: In case of sending implants to the manufacturer conduct analysis, this should be sanitized in the hospital using a broad-spectrum bactericidal and antifungal solution. Then must be disinfected or sterilized in a steam autoclave or ethylene oxide gas. The explanted product must be kept in intact packaging, identified with the method of cleaning, sterilization and product data.

Description of the efficacy and safety of medical products, in accordance with the regulations of ANVISA which addresses the Essential Efficacy and Safety Requirements of Medical Products.

The risks associated with Cemented Tibial Component - BM according to the Essential Safety and Efficacy Requirements of the Product according to RDC 56/01 for products covered by the Rule 8, risk class 3 and Risk Management of the Product and Process are conform to NBR ISO 14971-2009 - Risk Management for Health Products, as pursuant to procedure PRP003 - Risk Analysis Management.

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