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only once. Although it may appear undamaged, previous stress may create imperfections that may reduce the success of the implant. Improper selection of implant can cause unusual stress on the implant and may result in subsequent fracture of the implant. Because of their presentation and their characteristics, this product is not made to side effects; however, it is suggested to monitor the recovery of the patient through X-ray according to the International Safety Standards. Do not use MRI for the product manufactured in stainless steel and it may cause interference impairing the diagnostic test and may also cause displacement of the implant which may be problematic for the patient. The main factors that hinder the postoperative rehabilitation in ACL reconstructions (ACL) are:

Preoperative Factors
• Allow about three weeks between trauma and surgery as surgeries performed within the first two weeks post trauma usually evolve with greater functional limitation due to arthrofibrosis.
• ACL reconstruction is not recommended for patients with signs of inflammation such as severe pain, important swelling and / or effusion and especially with a MDA limitation greater than 10 degrees articulate. Some related injuries require early surgical solution (within two weeks), and grade III injury in the posterolateral or medial corner, because after this two week window it is more difficult to repair such lesions and the knee being locked leads to atrophied muscles.
• Physical and mental preparation of the patient is an important factor in preventing complications such as reflex sympathetic dystrophy and / or arthrobiosis. It is important to inform the patient that the time in appearance lost with fitness and / or emotional pre-op preparation increase the probability of success intraproprio factors:

• The major technical flaws that lead to a difficult postoperative rehabilitation of ACL reconstruction and even complete failure of the surgery are those related to intraoperative factors. Very anterior tibial tunnel leads to an important limitation in postoperative knee extension, with consequent loosening during physical activity, limiting the gain of extension capacity and/or even rupture of the graft in this situation. A Femoral tunnel brought forward, the most common technical error in the reconstruction leads to a higher friction of the graft against the femoral notch often causes limitations in the extension by formation of hypertrophic scars on the graft (cytopos) and / or loosening and breakage.

The tendency to widen the notch of the femur during surgery, for a better view of the correct femoral tunnel or a better position to accommodate the graft in this region can lead to formation of a hypertrophic scar tissue of fibrocartilaginous feature, which will end up pinching the graft leading to their loosening or breakage during the rehabilitation period.

Failures in graft fixation can be obtained from the use of proper screws as well as the use of unsuitable materials and fixings with poor fixation techniques performed will no doubt lead to complications in the postoperative period by early loosening or breakage of the graft.

• Extra care should be taken in removing the graft. Especially when using the patellar tendon (central 1/3) in ACL reconstructions, the lack of care in this donor area and especially its closure with high strain or with the knee in flexion, can lead to the development of pain or low patella postoperatively, which is difficult to resolve because of significant pain and limitation.

• Inadequate meniscus suture done with inappropriate instrument and equipment and without proper training of the surgeon may lead to the formation of intra articular areas of adhesion and pains which lead to very difficult recovery and in some cases may lead to further surgery to remove intra articular granulomas or darts fragments or arrows used in the procedure.

• Excessive repair of peripheral structures as well as of the PCL (Posterior Cruciate Ligament) can cause significant knee arthrofibrosis leading to great difficulty in postoperative recovery. Postoperative factors

• Use of rigid or fixed braces for more than three weeks lead to joint stiffness often difficult to solve.

• Accelerated rehabilitation protocol does not mean rushing back to the daily work or physical activities. The surgeon / physiotherapist combination is crucial for a good postoperative evolution of reconstructions.

• Reflex sympathetic dystrophy and infections are also causes of delay and difficulty in rehabilitation after reconstruction.
• The surgeon must be familiar with and hold sufficient knowledge of osteosynthesis and its limitations, including the pre and post-operative surgical technique used, precautions and potential risks. The surgeon must be sure the bone preparation for the introduction of ligament screws because the stability and success of fixation may be compromised. During handling of the implant scratches and nicks should be avoided, or anything that may damage or mark the implant, since these defects are stress concentrators and can be sites of crack nucleation and decrease the corrosion resistance and may result in implant fracture. To ensure proper deployment, Biomecánica equipment should be used, because they have been designed and manufactured specifically for the use of these implants. Using instruments from different manufacturers may compromise the surgery along with not being compliant with the product registration in the National Health Surveillance Agency (ANVISA). Do not use the product if it is damaged.

NON-STERILIZED PRODUCT - STERILIZE BEFORE USE. SINGLE USE PRODUCT - IT IS EXPRESSLY PROHIBITED TO RE-USE THIS PRODUCT.

CONTRAINDICACIONES: Store in cool, dry and away from light and away from the bad weather. Do not use the product if the package is damaged.
Note: Components from different manufacturers implants should not be used. We therefore recommend that the products have the same origin. Manufacturing date, expiry date and batch of the product: SEE LABEL

Product Conformation Limit: During handling of the implant scratches and nicks should be avoided, or anything that may damage or mark the implant, since these defects are stress concentrators and can be sites of crack nucleation and decrease the corrosion resistance and may result in implant fracture.

STERILIZATION: This product is supplied non-sterilized. It should be sterilized before use. We recommend steam sterilization in autoclaves at the hospital (ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices)

AUTOClave: It is an apparatus for sterilization by high pressure saturated steam.

OPERATING INSTRUCTIONS: To start the machine, check: if the circuit breaker is switched on; if the register of the water network is open; if the discharge register is closed. 1 - open the door of the equipment; 2 - accommodate the material to be sterilized properly; 3 - close the door equipment; 4 - select the desired cycle according to the material to be sterilized; 5 - turn on the main switch; 6 - the cycle will elapse, in sequence; 7 - when the "end cycle" beacon turns on, partially open the door for ten (10) minutes or so to cool the material.

It is recommended that the following physical parameters of sterilization in autoclaves (saturated steam) are applied:

Cycle	Temperature	Exposure Time (Minimum)	Drying Time (Maximum)
Convencional (1 atm de presión)	121°C (250ºF)	30 minutos	
High Vacuum	134°C (273ºF)	6 minutos	15 minutos

Adverse effects: Infection: According to the literature the incidence of infection is low, varying from 0.3 to 0.5%. Some known factors that can increase the risk of postoperative infection are: previous surgery, meniscus suture with posterior incision, lack of intraoperative antibiotic administration and immediate postoperative phase, prolonged duration of surgery and steroid injection at the end of the procedure. Thromboembolic disease: Patients undergoing ACL reconstruction have a lower risk, as they are usually are young and active. The chances of these patients to have predisposed local or systemic problems are lower. The discovery of risk factors and early stroll and movement are the main factors of prevention. Mechanical loosening can result from faulty fixation, stable reconstruction of the fractured bone or hidden infection. Sutures: metal sensitivity in patients have been reported, to be rare. Implantation of foreign material in tissues results in histological reactions involving various sizes of macrophages and fibroblasts. The chemical significance of this effect is uncertain, as similar changes can also occur as a precursor or during the healing process. The misplacement of screws resulting from technical error can cause many complications for the patient. The treatment consists of three risks: osteoarthritis, avascular necrosis and chondrolysis. Improper use of crutches after surgery while all the signs have not totally disappeared, can cause complications with the implant. Neurological complications, bleeding, reactions to foreign bodies, bone fractures, post-operative pain and pain with abnormal sensations are have also been reported.

Information to be Provided to the Patient: The wording of this statement in the items: Indications, Contraindications, Usage Information, Potential Adverse Effects, Warnings and Precautions. The patient should be informed about the importance of postoperative follow-up. The lack of monitoring prevents the detection of post-surgical problems. Not performing revision surgery when there is loosening or osteolysis of components can cause progressive loss of the joint and the lock. The patient should be informed about the type of implant used, brand, code, name, batch, surgeon's name and date of surgery.

Disposal of Product: The Metal Cannulated Screws Set for Ligament Fixation - BM explanted from patients should be properly disposed of by the hospital. It is the responsibility of the hospital to ensure the complete disposal of the implant preventing their reuse as well as the method used. Biomecánica recommends that explanted implants are deformed mechanically with the aid of hammer or press impact and should be identified by the words "Unfit for Use." As described in Resolution No. 2605 of 11/08/2006, implantable devices of any kind classified as single use, are prohibited from being reprocessed.

Traceability: Along with the packaging of implantable component are attached two labels of traceability containing the data of the implant used. These labels are the same label contained in the external labeling of the product. The information contained in

healing. Should the product be out of the correct position, or present any non-compliance it is the surgeon's responsibility to take the most appropriate corrective action.

Useful information to Avoid Risks Arising from Implantation: To lessen the risks of implantation should be followed strictly, contraindications, instruction for use and all information contained in the "Operating Manual" product.

Intrinsic Risk of Implantation: The Metal Cannulated Screws Set for Ligament Fixation - BM are manufactured with certified biomedical quality materials that follow the requirements of the ASTM F138 - Standard Specification For Wrought 18 Chromium-14 Nickel-2,5 Molybdenum Stainless Steel Bar And Wire For Surgical Implants (UNS S31673) and ASTM F136 - Standard Specification for Wrought Titanium-6 Aluminum-4 Vanadium EL(ExtraLow Interstitial) Alloy For Surgical Implant Applications (UNS R56401).

Risk of Contamination: There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the metal components of the Metal Cannulated Screws Set for Ligament Fixation - BM come in contact with body tissues and fluids. Explanted products should be treated as highly contaminating.

Decontamination, Cleaning and Sterilization: All models of Metal Cannulated Screws Set for Ligament Fixation - BM are non-sterilized, so the institution that will perform the surgery will be responsible for sterilization and should use methods that are regularly approved. The product is supplied in its original packaging previously cleaned. Before use, the packaging should be inspected to ensure the absence of breakage or damage. The product should only be unwrapped before sterilization in order to preserve the surface finish and configuration, handling it, if possible, as recommended by ISO 8828. To ensure the reduction of risk of infection, the implants (which were not used in surgery, but were placed on the operating table or that are no longer sterilized) and surgical instruments must be cleaned after surgery and before sterilization. Cleaning can be carried out by manual or mechanical methods.

Equipment for cleaning: - Ultrasonic Washer, - Sterilizing Washer, - Disinfecting Washer.
Manual Cleaning: it is the procedure in which the dirt is removed by physical action with the aid of detergent, water and supplies such as sponge and brush. Manual cleaning is most recommended for gentle treatment of implants and surgical instruments. For manual cleaning water should be used at room temperature, using nylon brushes, never steel or abrasive as it may harm the protective layer of the material. For decontamination a broad-spectrum bactericidal and antifungal solution should be used, avoiding aggressive cleaning agents, metal brushes so that the quality of metallic implants is not impaired. In case of use of cleaning materials, implants and surgical instruments should be placed in drawers so that there is no impact between them avoiding material damage. Those responsible for the cleaning of surgical instruments and implants should make sure the type of cleaning product used complies with the concentration recommended by the manufacturer. Implants and surgical instruments should be thoroughly rinsed until all residue is removed and then dried immediately after cleaning.

Types of cleaners: Enzymes are basically composed of enzymes, surfactants and solubilizers. The balanced combination of these elements causes product to remove the organic matter of the material in a short period of time. Enzyme solutions: have excellent cleaning action, but do not have bactericidal and bacteriostatic features. Enzymes: substances produced by living cells that govern the chemical reactions of the process. Once produced by the cells, an enzyme can be isolated and will maintain its catalytic properties if certain conditions are maintained during its manufacture. Enzymes are classified into three major functional groups depending on the type of substrate that they will affect: proteases, amylases and lipases that act on protein substrates, fats and carbohydrates, and which tend to solubilize and loosen up the surface of the material, the cleaning of enzymes of complex configuration is recommended to ensure cleanliness.

Disinfection: It is a process that destroys microorganisms, pathogenic or not, present on the articles, with the exception of bacterial spores, by physical or chemical means. Levels of disinfection: High level: destroys all microorganisms except the spores in high number => Glutaraldehído 2% - 20-30 minutos; Intermediate: preferred for hospital area. Medium: well-eliminate vegetative bacteria, most viruses, fungi and mycobacteria => 1% sodium hypochlorite - 30 minutes; Indication: for UBS, daycare, nursing homes, nursing home. Low level: eliminates most bacteria, some viruses and fungi, but does not eliminate mycobacteria => 0.025% Sodium Hypochlorite. Indication: nutrition.

Sterilization: This product is supplied non-sterilized. It should be sterilized before use. We recommend steam sterilization in autoclaves at the hospital (ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices)

AUTOClave: It is an apparatus for sterilization by high pressure saturated steam.

Cycle	Temperature	Exposure Time (Minimum)	Drying Time (Maximum)
Convencional (1 atm de presión)	121°C (250ºF)	30 minutos	
High Vacuum	134°C (273ºF)	6 minutos	15 minutos

Adverse effects: Infection: According to the literature the incidence of infection is low, varying from 0.3 to 0.5%. Some known factors that can increase the risk of postoperative infection are: previous surgery, meniscus suture with posterior incision, lack of intraoperative antibiotic administration and immediate postoperative phase, prolonged duration of surgery and steroid injection at the end of the procedure. Thromboembolic disease: Patients undergoing ACL reconstruction have a lower risk, as they are usually are young and active. The chances of these patients to have predisposed local or systemic problems are lower. The discovery of risk factors and early stroll and movement are the main factors of prevention. Mechanical loosening can result from faulty fixation, stable reconstruction of the fractured bone or hidden infection. Sutures: metal sensitivity in patients have been reported, to be rare. Implantation of foreign material in tissues results in histological reactions involving various sizes of macrophages and fibroblasts. The chemical significance of this effect is uncertain, as similar changes can also occur as a precursor or during the healing process. The misplacement of screws resulting from technical error can cause many complications for the patient. The treatment consists of three risks: osteoarthritis, avascular necrosis and chondrolysis. Improper use of crutches after surgery while all the signs have not totally disappeared, can cause complications with the implant. Neurological complications, bleeding, reactions to foreign bodies, bone fractures, post-operative pain and pain with abnormal sensations are have also been reported.

Information to be Provided to the Patient: The wording of this statement in the items: Indications, Contraindications, Usage Information, Potential Adverse Effects, Warnings and Precautions. The patient should be informed about the importance of postoperative follow-up. The lack of monitoring prevents the detection of post-surgical problems. Not performing revision surgery when there is loosening or osteolysis of components can cause progressive loss of the joint and the lock. The patient should be informed about the type of implant used, brand, code, name, batch, surgeon's name and date of surgery.

Disposal of Product: The Metal Cannulated Screws Set for Ligament Fixation - BM explanted from patients should be properly disposed of by the hospital. It is the responsibility of the hospital to ensure the complete disposal of the implant preventing their reuse as well as the method used. Biomecánica recommends that explanted implants are deformed mechanically with the aid of hammer or press impact and should be identified by the words "Unfit for Use." As described in Resolution No. 2605 of 11/08/2006, implantable devices of any kind classified as single use, are prohibited from being reprocessed.

Traceability: Along with the packaging of implantable component are attached two labels of traceability containing the data of the implant used. These labels are the same label contained in the external labeling of the product. The information contained in

this traceability label are listed under "Form Submission" in this operating manual. One label must be attached to the patient's records and one for the hospital records. The hospital must record this information in their files and pass this information to the patient. In addition, the following traceability labels with the product must have recorded information which allow the tracking of the implanted product. The attachment of the label to the medical record allows traceability of the product used. Among the essential information are the name of manufacturer of the implant used, name of implant used, the implant code, batch number and registration number of the product at ANVISA. This information is available in the following traceability labels with the product and the external labeling. Other information should also be considered essential, as the date of the surgery, name of the patient receiving the implant, surgeon's name, patient weight, patient age and other information requested from the patient.

Marking: The Metal Cannulated Screws Set for Ligament Fixation - BM all screws get marking on the screw head and the following information: Laser etched to enable traceability and post-operative monitoring of the patient: -Logof Biomecánica; - Manufacturing Batch Number; - Acronym of the raw material used to manufacture the product (S1 for stainless steel and T3 for titanium alloy); -Dimension

After Sales (Customer Complaint): If there is any complaint of Metal Cannulated Screws Set for Ligament Fixation - BM related to any adverse effect that affects user safety, to malfunctioning product, implantable metallic component damage, serious problems or deaths related to these components the surgeon should report this adverse event to the competent health care body and Biomecánica via email sac@biomecánica.com.br or by phone +55 11 214 7219. Biomecánica will assign a responsible professional responsible can report the adverse event through the Notifications in Health Surveillance System on the ANVISA website: http://www.anvisa.gov.br/hotsite/notisia/index.htm

Submitting Content for the manufacturer to Analyze: In case of shipping implants to the manufacturer for analysis, it must be sanitized at the hospital using a broad spectrum bactericidal and antifungal solution. Then it must be disinfected or sterilized in a steam autoclave or ethylene oxide. It should be forwarded to Biomecánica in an undamaged packaging, duly labeled with the method of cleaning, sterilization and product data indicated.

Description of the efficiency and safety of the medical device, in accordance with the ANVISA standards which address the Essential Requirements of Efficiency and Safety of Medical Products. Risks associated with the Metal Cannulated Screws Set for Ligament Fixation -BM comply with the Essential Requirements of Safety and Efficiency of the Product as defined by the RDC 56/01 for products covered by Rule 8, risk class 3 and Risk Management of Product and Process as defined by NBR ISO 14971 as provided in PRPD02 - Risk Analysis Management. The Metal T3 and T3 Screws Set for Ligament Fixation - BM are manufactured in Stainless Steel Bar and Wire as defined by ASTM F138 and titanium alloy as defined by ASTM F136 approved by laboratory test reports and of certified quality. These materials are acceptable for the manufacture of implants as defined in Annex A of the NBRISO21534 norm.

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