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
# FIXED ARTICULAR INSERT WITHOUT STABILIZER - BM

| Code        | Description  | Illustration |
|-------------|--|--------------|
| 2358-08-000 | BKS II - Large Tibial Insert with Preservation - 6.0 mm thick        |              |
| 2358-10-000 | BKS II - Large Tibial Insert with Preservation - 8.0 mm thick        |              |
| 2358-12-000 | BKS II - Large Tibial Insert with Preservation - 10.0mm thick        |              |
| 2358-15-000 | BKS II - Large Tibial Insert with Preservation - 13.0mm thick        |              |
| 2358-17-000 | BKS II - Large Tibial Insert with Preservation - 15.0mm thick        |              |
| 2358-19-000 | BKS II - Large Tibial Insert with Preservation - 17.0mm thick        |              |
| 2358-20-000 | BKS II - Large Tibial Insert with Preservation - 18.0mm thick        |              |
| 2358-24-000 | BKS II - Large Tibial Insert with Preservation - 22.0 mm thick       |              |
| 2359-08-000 | BKS III - Extra Large Tibial Insert with Preservation - 6.0 mm thick |              |
| 2359-10-000 | BKS III - Extra Large Tibial Insert with Preservation - 8.0mm thick  |              |
| 2359-12-000 | BKS III - Extra Large Tibial Insert with Preservation - 10.0mm thick |              |
| 2359-15-000 | BKS III - Extra Large Tibial Insert with Preservation - 13.0mm thick |              |
| 2359-17-000 | BKS III - Extra Large Tibial Insert with Preservation - 15.0mm thick |              |
| 2359-19-000 | BKS III - Extra Large Tibial Insert with Preservation - 17.0mm thick |              |
| 2359-20-000 | BKS III - Extra Large Tibial Insert with Preservation - 18.0mm thick |              |
| 2359-24-000 | BKS III - Extra Large Tibial Insert with Preservation - 22.0mm thick |              |

**Precautions of Handling and Transport of the Medical Product:** The storage of medical product should be clean, dry and properly lit in order to maintain the ideal conditions for storage and transport, as well as its physical and chemical integrity. The product should be stored and transported in a cool, dry place at room temperature (max. 35 °C), relative humidity around 30% and 70% and away from direct sunlight. The effects of vibration, shock, corrosion, defective seating during handling and transportation, improper stacking during storage should be avoided. Implants should be handled with care to avoid damage that might affect the quality of the material and the safety of the patient. Implants cannot be stored directly on the floor (minimum height = 20cm). Implants cannot be kept in very high places, near lamps, which can cause dryness or damage to the packaging label. Implants should not be stored in locations in which are stored contaminants such as, for example, cleaning materials, insecticides, pesticides, etc. The carrier should be informed about the content and date of delivery. Carry with care to avoid falling and friction in order not to damage the packaging and the device, thereby ensuring sterility of the material. The opening of the sealed blister packaging can be done manually, in an appropriate location (operating room), pulling the surgical grade paper from the opening. This opening should only be performed during surgery to prevent product contamination.

**The statement, purpose or use for which the medical product is intended, as indicated by the manufacturer.**

**Intended Use:** The Fixed Articular Insert without Stabilizer - BM aims for the replacement or restitution of the femorotibial knee joints when these are no longer able to maintain their functional capacity due to deformation, pain and blockage of movements. This procedure, together with the replacement of the patellofemoral joint, is known as "total knee arthroplasty", and seeks to create an artificial joint that has the same functionality of the original knee.

The total knee arthroplasty aims to eliminate the painful symptoms, correct deformities, allow functional breadth of motion and knee stabilization.

The Fixed Articular Insert without Stabilizer - BM is one of the components used in total knee arthroplasty and should therefore be used in conjunction with other implants. The Fixed Articular Insert without Stabilizer - BM has no locking system which preserves the posterior cruciate ligament. According to Villardi et al. (2008)2 and Villardi et al. (2005)1, there is no consensus on whether or not the preservation of the posterior cruciate ligament in the procedures for total knee arthroplasty, and that proceptive changes of the knee, patellar problems, joint line migration, improvement of breadth of motion, presence of rollback, greater efficiency in climbing stairs and maintenance of bone stock is insufficient to determine which is the best option between preserving or not the posterior cruciate ligament in total knee arthroplasty. But in their studies it became clear concerning the gait patterns and patellar height there was a slight difference giving advantage to preserving the posterior cruciate ligament. In studies from other authors quoted by Villardi et al. (2008)2 and Villardi et al. (2005)1 no difference was found between if there is preservation or not of the posterior cruciate ligament.

According to Weinstein et al. (1994) 4, as for elderly patients with advanced osteoarthritis, the best treatment is total knee replacement. In younger patients, total knee replacement is indicated only in cases of tricompartmental disease (medial and lateral tibiofemoral and patellofemoral), a future revision surgery will be necessary. In cases of rheumatoid arthritis total knee replacement is indicated for all ages because of the complications of the disease. As shown by the author, total knee replacement in cases of rheumatoid arthritis or osteoarthritis yields satisfactory results in about 90% of patients.

The Fixed Articular Insert without Stabilizer - BM is indicated in the total knee arthroplasty procedures, without preservation of the posterior cruciate ligament, due to:  
v Osteoarthritis: direct friction between the articular surfaces caused by wear of the articular cartilage that covers it, generating an extensive picture of pain, limited breadth of motion and deformity.  
v Rheumatoid Arthritis: a rheumatic disease that causes pain and deformity of the joints in general;  
v Consequences of comminuted fractures: many fragments of the knee where it is practically impossible to achieve a good coaptation.

## Bibliographical Reference

- Villardi A. M., Veiga L. T., Franco J. S., Cagy M., Silva P. J. G., Raptopoulos L., D'Angelo M. D., Gait analysis after total knee arthroplasty with and without preservation of the posterior cruciate ligament; 2005.
- Villardi A. M., Leite J. P. B., Santos B. F., Mandarino M., Avaliação comparativa de patelar height after total knee arthroplasty with and without preservation of the posterior cruciate ligament; 2008.
- Weinstein S. L., Buckwalter J. A., Turek T. R. Osteoporosis - Principles and Their Application; 2000.

## Contraindications

**Contraindications below do not limit the use of models Fixed Articular Insert without Stabilizer - BM, but should be analyzed by the surgeon:** Overall weakness of the quadriceps; infections or systemic diseases on near the joint; knee arthrosis after arthroplasty, revision surgery; prior patellofemoral arthritis; In younger patients, total knee replacement is indicated only in case of tricompartmental disease (medial and lateral tibiofemoral, and patellofemoral), a future revision surgery is required; obesity, fever or leukocytosis; intense neuromuscular, muscular or vascular deficiencies in the affected region; Severe lesions on bone structures generating an unacceptable risk in the implantation of a prosthesis; improper ligament system status for the selected component model; bone tumors near the component implantation location; muscle atrophy; bone immaturity; failure of the extensor mechanism of the knee; Failure of the structural soft tissue of the knee; Ankylosis of the knee; insufficient bone structure in the tibia, femur and patella; Patients unwilling or unable to abide by postoperative precautions and instructions.

**Precautions, restrictions, warnings, precautions and clarifications on the use of the medical product as well as its storage and transportation.**

## Precautions and Warnings

**SINGLE USE PRODUCT. DESTROY AFTER EXPLANATION. DO NOT REUSE THE PRODUCT. REPRESSING PROHIBITED.** The arthroplasty implant should be used on a single patient and only once. Although it may appear undamaged, previous stress may create imperfections that may reduce the success of the implant. Improper selection of implant can cause unusual stress on the implant and may result in subsequent fracture of the implant.

The surgeon must be familiar with and hold sufficient knowledge, including pre and postoperative surgical techniques used, precautions and potential risks.  
During handling of the implant should be avoided scratches, nicks, or anything that may damage or mark the implant, since these defects are stress concentrators and can be sites of crack nucleation and decrease the corrosion resistance and may result in implant fracture.

Check that the package containing the product is intact. Check the date of sterilization. Aseptically open the package after the date of sterilization is correct one. Validity of sterilization: 5 years (indicated in the inner and outer packaging of the product). To ensure proper implantation, only instruments manufactured by Biomecanica must be used, as these have been designed and manufactured specifically for the use of these implants. Using instruments from different manufacturers may compromise the surgery besides not being in compliance with the product registration in the National Health Surveillance Agency (ANVISA). The surgeon must be aware of the possibility of physical development of the patient.

Do not use the product if it is damaged.  
**SINGLE USE PRODUCT - IT IS EXPRESSLY PROHIBITED TO RE-USE THIS PRODUCT.**  
**Special storage conditions:** Store in a cool, dry place away from light and away from the action of weather.  
Do not use the product if the package is damaged.

**Note:** Components from different manufacturers implants should not be used. We therefore recommend that the products have the same origin.  
The date of manufacturing, expiry date and batch of the product: SEE LABEL.

**Adverse Effects:** Migration of the joint line; Amendment gait pattern; patellar tendon rupture; Deformation of patellar height; Braking; Suture dehiscence; progressive and asymmetric bone resorption (osteolysis) located around the component; Loosening, migration or fracture of the medical product; superficial or deep infection; vascular disorders including thrombosis and pulmonary embolism; shortening of limb on the basis of bone resorption; allergic reaction to the foreign body and may result in histological reactions involving various types of macrophages and osteoblasts. Pain, discomfort, or normal reactions depending on the presence or migration of the medical product; and also due to the surgical procedure itself: neural or neurological damage due to surgical trauma (injuries including paralysis and/or lesions of soft tissue); inability to perform normal daily physical activities.

**Warnings:** The Fixed Articular Insert without Stabilizer - BM aims to replace the damaged joint and not normal structures. The Fixed Articular Insert without Stabilizer - BM has no stabilization system which preserves the posterior cruciate ligament. The correct selection of the Fixed Articular Insert Stabilizer - BM model is extremely important and should be in accordance with the patient bone characteristics, analyzed by the surgeon. A wear off of the Fixed Articular Insert without Stabilizer - BM may occur due to poor positioning of the metal components, inappropriate ligament balance and aseptic loosening of the prosthesis. Infections generated are often related to comorbidity.

**Information to be Provided to the Patient:** The wording of this statement in the items: indications, Contraindications, Usage information, Potential Adverse Effects, Warnings and Precautions. The patient should be informed about the importance of postoperative follow-up. The lack of monitoring prevents the detection of post-surgical problems and possible loosening of

components or occurrence of osteolysis. The patient should be informed about the type of implant used, brand, code, name, batch, surgeon's name and date of surgery. The patient should be made aware that the bearing of weight, walking, excessive muscle activity , among others, before the recommended time for bone healing can cause complications.  
In the postoperative care of an orthopedic surgical procedure, the patient's willingness and ability to comply with instructions is extremely important. Children, elderly patients, patients with mental health problems or drug addictions, who are more likely to ignore the instructions and restrictions, can generate a high risk of failure of the procedure. The patient should also be informed that the component does not replace or has the same performance as normal bone. Consequently it can warp, break or loosen due to excessive efforts or activities and early loading.

**Disclaimer:** Biomecanica, as the manufacturer of this device, does not practice medicine or recommend this or any other surgical technique to be used on a particular patient. The surgeon performing any implant procedure is responsible for choosing and using the proper techniques to implant the prosthesis for each patient. Biomecanica is not responsible for choosing the appropriate surgical technique to be used on each patient. Surgical techniques depend on the technique used by the surgeon, this implies the choice of method, type and size of the products being placed.

**Assessment of the implemented Product:** After implantation, in the intraoperative phase the relevant personnel should perform radiological control to verify the correct positioning of the product. The relevant personnel should perform, clinical and radiological examinations after surgery at a frequency stipulated by the protocol to verify the status of the implant and the evolution of the bone healing, if the product is out of the correct position, or present any non-compliance it is the responsibility of the surgeon to take the most appropriate corrective action.

**Useful Information to Avoid Risks Arising from Implantation:** To lessen the risks arising from the implantation, the following should be strictly followed: contraindications, instruction for use and all information contained in the product "Operating Manual".

**Intrinsic Implantation Risk:** The Fixed Articular Insert without Stabilizer - BM is manufactured with materials appropriate for biomedical use, which follow the requirements of standard ISO 5834-2 version 08 - Implants for surgery - Ultra-high molecular weight polyethylene - Part 2: Products in the molded shape.

**Risk of Contamination:** There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the Fixed Articular Insert without Stabilizer - BM comes in contact with body tissues and fluids. Explanted products should be treated as highly contaminating.

**Sterilization**  
**This product is sterilized by Ethylene Oxide Gas (ETO) valid for:** Five years. If the product is valid for unsuccessful sterilization and/or damaged packaging, it must be returned to the distributor or directly to Biomecanica.

**Handling of Sterilized Material**  
**When handling the sterilized material through an aseptic technique, certain rules must be followed in order to keep it sterilized:**  
v washing your hands with soap and water before handling the sterilized equipment is critical;  
v use material that has a full packaging, dry, without spots, with identification labels (type of material and date of sterilization);  
v work facing the material;

v handle the material from the waist upwards;  
v avoid coughing, sneezing, talking above the exposed material;  
v avoid making moves above the sterilized area;  
v ensure the validity and integrity of the package;  
v work in a clean, quiet, dry and draft-free environment;  
v keep some distance between you and the material to be handled;  
v obey other principles of asepsis.

**The nursing technic recommended when handling the sterilized equipment is:**  
v open it, starting with the end opposite the handle;

v protect the material exposed to the sterilized area that around which is wrapped;  
v handle with your hands only the outside of the package;

v don't consider as sterilized equipment any package that has already been open.

**Warranty:** The warranty of the Fixed Articular Insert without Stabilizer - BM relies on the compliance with this use instruction. As in: indications and usage information, contraindications, warnings and precautions, warnings, possible adverse effects, packaging, sterilization, cleaning and decontamination recommended in this technical report and statement of use.

**Disposal of Product:** The Fixed Articular Insert without Stabilizer - BM explanted from patients should be properly disposed of by the hospital. It is the hospital's responsibility to ensure the complete disposal of the implant, preventing their re-use as well as the method used. Biomecanica recommends that explanted implants are deformed mechanically with a hammer or press impact, and should be identified with the words "Unfit for Use". As described in Resolution No. 2605 of 11/08/2006, implantable devices of any kind classified as single use product, are prohibited from being reprocessed.

**Traceability:** Along with the packaging of the implantable component come 6 traceability labels containing information about the implant used. This label is the same label contained on the external label of the product. The information contained in this traceability label are listed under "Presentation" of this technical report. One label must be attached to the patient clinical record and another in the report delivered to the patient, another with the tax documents that define the billing, in the case of AIH, of a patient treated by SUS, or in the bill of sale, if the patient has a complementary mutual fund, another label available for the supplier or control (historical record of distribution - RH-D), another available for the surgeon to control (principal) and a last one for health insurance if applicable. Please be advised that the mandatory patient record must have recorded information which allow the tracking of the implemented product. Attaching that label to the medical record allows traceability of the product used. Among the main information to be found on the label: the name of the manufacturer of the product, the name of the implant used, the implant code, batch number and product registration number at ANVISA. This information is described in the following traceability labels with the product and the external labeling. Other information may also be considered essential: the date of surgery, name of the patient receiving the implant, surgeon's name, patient weight, patient age and other information requested from the patient.

**Marking:** The Fixed Articular Insert without Stabilizer - BM contain the following information carved, by means of machining, to allow postoperative traceability and monitoring of the patient: - Identification of the manufacturer; Manufacture Batch; Dimension; Product Code.

**After Sales (Customer Complaint):** If there is need for a complaint about the Fixed Articular Insert without Stabilizer - BM related to some adverse effect that affects the security of the user, such as a product not working properly, implantable metallic component damage, serious injury or death related to these components the surgeon should report any adverse event to the competent health care body as well as Biomecanica via email sac@biomecanica.com.br or by phone +55 14 2104 7900. In cases of doubt the surgeon or health professional can report the adverse event through the Notifications in Health Surveillance System through ANVISA website: <http://www.anvisa.gov.br/hotsite/notvnsia/index.htm>

**Submitting Content to the Manufacturer for Analysis:** In case of sending implants for the manufacturer to perform analysis, it must be sanitized at the hospital using a broad spectrum bactericidal and antifungal solution. Then it must be disinfected or sterilized in a steam autoclave with ethylene oxide. It should be sent to Biomecanica in an undamaged packaging, clearly identified with the method of cleaning, sterilization and product data.

**Description of the efficiency and safety of the medical product in accordance with the rules of ANVISA which address the Essential Safety and Efficacy Requirements of Medical Products.**

The risks associated with the Fixed Articular Insert without Stabilizer - BM comply with the Essential Safety and Efficacy Requirements of the Product as defined by RDC 56/01 for products covered by Rule 8, risk class 3 and Product and Process Risk Management as defined by NBRISO14971 as provided in the PRPD03 - Risk Analysis Management.

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