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# PROSTHESIS FAMILY

**Special storage conditions:** Keep in a cool, dry place, away from light and away from the weather action.

Do not use the product if the packaging is damaged.

Note: It must not be used implant components from different manufacturers.

Therefore it is recommended that the products have the same origin.

Manufacturing date, expiration and product batch: SEE LABEL.

## Adverse Effects

**Allergies to contents of the implant:** The hypersensitive to the implants components may occur as a result of the body's reaction, primarily for nickel ions, and manifests as skin rash (skin redness), pruritus or urticaria. Although the nickel allergy in the general population lies between 10-15%, its manifestation after Hip Arthroplasty is rare and therefore it is believed that the body can develop tolerance to metal ions after implantation of orthopedic devices. However, the prevalence and significance of hypersensitivity reactions to metal implants are not well understood.

**Chondrolysis:** Although concern about the increase of the occurrence risk of tumoral processes in connection with orthopedic implants, there is no clinical evidence in the literature to support this assumption.

**Bone Cement Implantation Syndrome:** This rare syndrome, but potentially deadly, named as BCS - Bone Cement Implantation Syndrome - can occur intraoperatively in patients undergoing implantation of cemented orthopedic devices. There is no way to identify patients at risk for the occurrence of this uncommon condition, and thus no proper monitoring and support the patient during the cementing procedure is strongly recommended.

**Information to be Provided to the Patient**  
Talk to your orthopedic surgeon about any medical conditions that might affect the surgery. Total prostheses are successful surgeries in more than 90% of patients. When complications arise, most can be treated successfully. Among the complications that can arise, we have:

**Infection:** Infection may occur in the wound or can be deep (around the prosthesis), it can arise while the patient is in the hospital or at home. It can appear even years later. Minor infections in the wound are usually treated with antibiotics. Major or deep infections may require further surgery (for deep cleaning or even removal of the prosthesis. Any infection in your body (bladder, throat, teeth, ears, etc.) can lead germs through the blood to the prosthesis and cause infection.

**Thrombosis:** they are blood clots resulting from various factors, including your reduced mobility, and they made the blood flow more slowly in the veins, which may facilitate the appearance of thrombosis. Blood clots should be suspected if there is pain and edema in your thigh or calf. If it occurs, your orthopedic surgeon should do tests and exams to assess the risk of blood clots in the leg, and may require assessment and monitoring by vascular surgeon. Several measures should be used to reduce the possibility of blood clots:

- Medication for "thinning" the blood (anticoagulants)
- Compression stockings
- Exercise to increase circulation in the legs.

Plastic boots that inflate and squeeze the soles of the feet and calf, increasing the venous return.

Even with the use of these preventive measures, blood clots may still occur. If you notice edema (swelling), redness or pain in his calf after discharge from the hospital, you should contact your orthopedic surgeon.

**Loosening:** loosening of the prosthesis inside your body can occur after the surgery, it can cause pain if the loosening is significant. A surgical revision replacement prostheses may be required. New materials and new methods of attachment should minimize this problem.

**Dislocation:** possibly, after a hip prosthesis, the head of the prosthesis can leave behind the acetabulum that. This is what we call dislocation. Most of the cases the hip can be put back in place without the need for further surgery. To prevent dislocation is important that the muscles are strong (do exercises recommended by your orthopedist and does not flex the hip more than 90 degrees in the first months.

**Wear:** some kind of wear can be detected in any type of prosthesis. Excessive wear may contribute to loosening and may require further surgery (revision)

**Breakage of the Prostheses:** Once implanted, the bio-materials used in the prostheses necessarily come in contact with body fluids. These fluids, while seemingly harmless, can, over time, significantly degrading most materials with considerable "chemical inertia". Additionally, many of the orthopedic implants are subjected to static mechanical stress and/or cyclic, that may lead to material failure, i.e. breakage of the prosthesis.

With the materials currently used in prosthetic, implant fracture is very difficult. However, if this occurs, a new surgery to replace the fractured prosthesis is required.

**Soft tissue damage:** nerves in the vicinity of the prosthesis may be damaged during surgery, although this is very rare. This lesion is easier to occur when the patch surgery involves major deformity in the bowl or stretching of a very short leg due to wear. Over time these nerves injury can and generally improve complete recovery. Occasionally you can get for surgical exploration of the involved nerve.

**Reviews of the Implanted Product**  
After the implantation, interoperative responsible professional should perform radiological control to verify the correct positioning of the product. The professional responsible must make, eye to eye, the responsibility, the clinical and radiological examinations passed after the surgery procedure in frequency as it shall prescribe to check the status of the implant and the evolution of bone consolidation. If the product is outside of the proper placement, or shows any non-compliance, is the surgeon's responsibility take the most appropriate corrective action.

**Useful Information for Preventing Risks Arising from Implementation**  
To decrease the risks of implementation should be strictly follow: contraindications, instruction for use and all information contained in the product's "use instructions".

## Intrinsic Implementation Risks

The Hip Arthroplasty Prosthesis Family is manufactured with materials of recognized biomedical use. Being the Titanium alloy According to the ASTM F136 standard - Standard Specification For Wrought Titanium-6Aluminum-4Vanadium III (Extra Low Interstitial Alloy For Surgical Implant Applications.

## Contamination Risks

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis, because the modular stem comes into contact with tissue and body fluid. Explanted products should be treated as highly contaminating.

**Commercial Presentation of Product- Feature Package**  
The Hip Arthroplasty Prosthesis Family is provided in two types of sterilization, by gamma radiation, and Gamma Radiation For each type of sterilization, the package are composed by:

Type of Sterilization	Descrição das embalagens
Ethylene Oxide	Blisters de PVC or PET, cot., Tyvek, cartón and security seal
Gamma Radiation	Blisters PET non-PET, cot., Tyvek, cardboard box and PVC shrink film

For sterilization with ethylene oxide gas are used blister packaging made of PVC (Polyethylene Chloride) or PET (Polyethylene Terephthalate) and are used PET containers (polyethylene Terephthalate) for products sterilized by gamma radiation.

PET (Polyethylene terephthalate) and PVC (polyvinyl chloride) are raw materials suitable for the manufacture of blisters formed and cots term (vacuum - forming) in the pharmaceutical industry, food, toys, electronics, cosmetics and automobile. Both have the following main characteristics: good impact resistance; great dimensional stability; non-toxic, allowing contact with food and medical products; good moldability, good printability and self extinguishing. They are preferred for their mechanical properties, high transparency, with virtually no defects, good fluidity but not to present any danger to human health.

## Typk

This paper is an interlaced high-density polyethylene film (HDPE) to adhesive application. The Tyvek paper is suitable for the manufacture of packaging and Dental Hospital applicable for sterilization by ethylene oxide and Gamma Radiation. This paper is used in combination with the blister pack. Among its main features are:

- It supports high temperatures
- High tensile strength and punching
- Microbial barrier
- Limited use of high cost
- Disposable by incineration

## Security seal

Once packaged, the cardboard box that has a plastic finish cannot receive the security seal tampering, egshelling type with the following specification: adhesive label seal - measure 3xcm, length 9cm, printed 1 color.

**Presentations of the Medical Product**  
The Hip Arthroplasty Prosthesis Family after receive the double blister packaging is packed in a cardboard box. These boxes are customized with the company logo and sealed with the security seal, containing the labels of products and sterile product indication (red label), distributed uniformly, as shown in Figure below:

Figure - Form of presentation of distribution of the Hip Arthroplasty Prosthesis Family

**Presentations of the Medical Product**  
The external label is the product traceability tag that contains the following information:

Symbolism Traceability Label Label according to NBRISO15223:	
XXXX	Manufacturing Date
XXXXXX	Expiration
LOTXXXX	Product Batch Number
XXXXXX	Manufacturer's data
△	Non-sterile
REFXXXXXX	Code
△	"Single-Use Product"

The non-sterile product is identified with a green color label saying "Non-sterile"

Symbols of the outer packaging:	
△	Fragile, handle with care.
☀	Keep dry.
☀	Keep out of the sun.
⚠	Caution.
⊘	Do not use if package is damaged.

**Contained information in Cardboard Box (Outer Packaging)**  
The cardboard boxes (outer packaging) also have sayings of special conditions of storage, conservation and/or Product Handling with symbols. The table below describes the information in this cardboard box.

Table - Special conditions of storage, conservation and/or Product Handling contained on the outer packaging in the product	
☀	Open Packaging Aseptically
☀	Temperature Range 95°/32°C
☀	Store and transport in a cool dry place, at room temperature (Max. 95°/32°C) and relative humidity between 30% and 70%.
⊘	Do not store directly on the floor (minimum height = 20cm) neither in very high places, near areas which can cause dripping or damage to the label.
⊘	Do not store it where contaminants are stored such as, for example, cleaning supplies, insecticides, pesticides, etc.
⚠	Fragile, handle with care.
☀	Do not use if package is damaged
☀	Keep out of the sun
⚠	Attention: Instructions, warnings and precautions. See instructions for use.
☒	Check Use Instructions
☒	Do not use if package is damaged

The Hip Arthroplasty Prosthesis Family will be distributed as a unit according to the table below:

Table - Models of the Hip Arthroplasty Prosthesis Family		
Code	Description	Sterilization:
2518-06-000	Bionic® Stem 1 6 mm without support - not cemented	ETO
2518-07-000	Bionic® Stem 1 7 mm without support - not cemented	ETO
2518-08-000	Bionic® Stem 1 8 mm without support - not cemented	ETO
2518-09-000	Bionic® Stem 1 9 mm without support - not cemented	ETO
2518-10-000	Bionic® Stem 1 10 mm without support - not cemented	ETO
2518-11-000	Bionic® Stem 1 11,2 mm without support - not cemented	ETO
2518-12-000	Bionic® Stem 1 12,5 mm without support - not cemented	ETO
2518-13-000	Bionic® Stem 1 13,7 mm without support - not cemented	ETO
2517-06-000	Bionic® Stem 1 6 mm with support - not cemented	ETO
2517-07-000	Bionic® Stem 1 7 mm with support - not cemented	ETO
2517-08-000	Bionic® Stem 1 8 mm with support - not cemented	ETO
2517-09-000	Bionic® Stem 1 9 mm with support - not cemented	ETO
2517-10-000	Bionic® Stem 1 10 mm with support - not cemented	ETO
2517-11-000	Bionic® Stem 1 11,2 mm with support - not cemented	ETO
2517-12-000	Bionic® Stem 1 12,5 mm with support - not cemented	ETO
2517-13-000	Bionic® Stem 1 13,7 mm with support - not cemented	ETO
2618-06-000	Bionic® Stem 1 6 mm without support - not cemented	Gamma Ray
2618-07-000	Bionic® Stem 1 7 mm without support - not cemented	Gamma Ray
2618-08-000	Bionic® Stem 1 8 mm without support - not cemented	Gamma Ray
2618-09-000	Bionic® Stem 1 9 mm without support - not cemented	Gamma Ray
2618-10-000	Bionic® Stem 1 10 mm without support - not cemented	Gamma Ray
2618-11-000	Bionic® Stem 1 11,2 mm without support - not cemented	Gamma Ray
2618-12-000	Bionic® Stem 1 12,5 mm without support - not cemented	Gamma Ray
2618-13-000	Bionic® Stem 1 13,7 mm without support - not cemented	Gamma Ray

**Sterilization:**  
This product is provided Sterile. It is supplied with the following sterilization process:

Process	Expiration
Sterilization by Ethylene Oxide Gas (E.T.O.)	5 years
Sterilization by Gamma Radiation (Gamma Ray) dosage of 25 kGy	5 years

If the product is with sterilization expired and/or damaged packaging, it should be returned to the responsible distributor or directly for Biomechanica.

**Cares with Storage, Storage, Handling and Transportation of the Medical Product**  
The storage place of the medical product should be clean, dry and light in order to keep the ideal conditions for storage and transport, as well as their physical and chemical integrity. The product must be stored and transported in a cool dry place, at room temperature (Max. 95°F/35°C) and relative humidity between 30% and 70%.

The effects of vibration, shock, corrosion, defective laying during handling and transportation, improper stacking during storage should be avoided.

Do not store directly on the floor (minimum height = 20cm). Do not store in very high places, near lamps, which can cause dryness of the packaging or damage to the label.

Do not store it where contaminants are stored such as, for example, cleaning supplies, insecticides, pesticides, etc.

The carrier must be informed about the content and delivery. Transport carefully to avoid falling and friction in order to not damage the packaging and the device, thereby ensuring the sterility of the material.

The opening of the sealing blister pack can be done manually, in a suitable place (surgical environment), by pulling the Tyvek by the end around the blister. This opening should only be performed during surgery to prevent product contamination.

The product should not be used if the packaging is damaged.

Note: Should not be used implant components from different manufacturers; We therefore recommend that the products have the same origin. Manufacturing date, expiration and product batch: SEE LABEL.

**Handling of Sterilized Material**  
When handling the sterilized material with aseptic technique, one must obey certain rules in order to keep it sterile:

- It is essential to wash hands with soap and water before handling the sterilized material;
- Use the material packaging complete, dry, without spots, with ID type of material and sterilization dates);
- Working (keeping the material);
- Manipulate the material at the waist up;
- Avoid coughing, sneezing, talking over the exposed material;
- Do not make moves over the sterile field;
- Make sure of the validity and appropriateness of the package;
- Work in a clean, quiet, dry and without airflow place;
- Keep some distance between the body and the material to be manipulated;
- Comply with other principles of asepsis.

The naming technique recommended in sterilized material handling is:

- Open it, starting with the opposite end to the handle;
- Protect the exposed material to the sterile field that involves it;

- Touch with your hands only the outside of the package;
- Do not store as material sterilized an open package;

**Warranty**  
The warranty of the Hip Arthroplasty Prosthesis family is related to the compliance of the operating instructions. Where: indications and use information, contraindications, warnings and precautions, cautions, possible adverse effects, packaging, sterilization, cleaning and decontamination recommended in this instruction for use.

**Product Disposal**  
The implants that integrate the Hip Arthroplasty Prosthesis family explained from patients must properly discarded by the hospital. It is under the hospital's responsibility the complete mischaracterization of the implant preventing their reuse. It is under the hospital's responsibility the method used for the mischaracterization of the implant. The Biomechanica recommends that explanted implants are mechanically destroyed with hammer impact or press should not then identified with the mean "Unfit for Use".

In accordance with the Resolution No. 2605, of 11/08/2006, implantable devices of any kind classified as single use, are prohibited from reprocessing.

**Aftermarket (Customer Complaint)**  
If there is need for a complaint about Hip Arthroplasty Prosthesis Family related to some adverse effect that affects the safety of the user, as a product not working, damage to the implantable metallic component, serious problems or deaths related to these components the surgeon in charge must report this adverse reaction to the competent health agency and the Biomechanica by email tecnovigilancia@biomechanica.com.br or phone +55 14 2104 7900. In cases of doubt, the responsible surgeon or the health professional may make the communication of adverse events through the Notifications System Health Surveillance at the site of ANVISA: <http://www.anvisa.gov.br/notisae/notivisa/index.htm>

**Traceability**  
With the package of the implantable component follows 6 traceability labels containing information on the implant used. This label is the same label on the outside labeling of the product. The information in this traceability label are listed on the item "Delivery Form" of this technical report. One label must be affixed to the patient's record, another in the report delivered to the patient, another at the tax documents that generates the collection, on AIIH (authorization for hospital admission), in case of the patient being assisted by SUS, or on the bill of sale, in the case of a patient attended by complementary health system, and/or one made available to the control vendor (historical record distribution - INV), another one available for control of the attending surgeon (main) and the last to health insurance if one exists. We reported that necessary the patient's medical record must have recorded information which track the deployed product. Glue this label in the chart allows traceability of the product used. Among the main information highlighted as essential the implant manufacturer's name used, implant name used, highliged, lot number and product registration number at ANVISA. This information is described in the traceability labels that follow with the product and the external labeling. Other information must also be considered essential, as the date of surgery, patient's name who received the implant surgeon's name, patient weight, patient's age and other information requested in the patient must also be met.

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
**Material Submission for the Manufacturer Review**  
In case of sending implants for the manufacturer to perform analysis, it must be cleaned in the hospital using a bactericidal solution broad spectrum anti-fungal. Must then be disinfected or sterilized steam autoclave or ethylene oxide. Should be referred to Biomechanica in undamaged packaging, identified by the method of cleaning, sterilization and product data.

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**Regulations of ANVISA which deals with the Essential Requirements of Efficacy and Safety of Medical Products**  
The risks associated with the Hip Arthroplasty Prosthesis Family comply with the Essential Requirements of Product Safety and Efficacy as the RDC 560/14 for products covered by the rule, risk class 3 and the Management of Product Risk and Process (Appendix II) agree with NBR ISO 14971 as provided in the PRPD03 - Risk Analysis Management

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