FIXADOR EXTERNO BIOTUB (ETO) - BM

PORTUGUÊS

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ENGLISH

Product Description

The External Fixator has two basic functions: a) To join the fractured parts of the bone tissue, either diaphyseal or epiphyseal, creating compression or dynamic distraction or static stabilization between the fragments. b) To promote the stretching of the bone tissue, in case of their shortening. The External Fixator is a non-implantable medical product developed to aid bone fixation, which aims to: reduce, align, stabilize and fix various types of fractures of the lower and upper limbs of the human skeleton, according to the indication of use, and aims to transfer the load applied on the fractured bone to the apparatus, keeping the bone unloaded, while regenerating from the fracture. External fixators are provided with only pipe and connector structures, connected by means of wires and pins, with different designs and correct methods of application, but with general principles that guide their applications.

The Wires and Pins are ancillary components of the external fixators, however "They are not Object of this Register and are not part of this product", so they are registered in another process and sold separately.

Manufacturing Materials

The external fixators are manufactured with materials recognized in medical use according to ASTM standard for steels and NBR standard for aluminum. There are two types of raw material used, according to each component being Stainless Steel type 304 according to ASTM F899, Stainless Steel type 420 according to ASTM F899 and aluminum alloy 6065-T5 according to standard NBR 6834. Table 1 describes the codes, set components and raw material of each component of the fixator.

Table 1 - Relation of external fixators

| Code | Description | Assembly Components | Raw Material | Illustration of The Types of Assemblies |
|--------------|------------------------|-------------------------------|--|---|
| 50060-01-001 | Tubular Fixator TT-202 | 06 Bar-Pin Connector Assembly | Aluminum alloy 6065-T5 according to norm NBR 6834 and Stainless steel type 304 according to the standard ASTM F899 | |
| | | 03 Bar-Bar Connector Assembly | Aluminum alloy 6065-T5 according to norm NBR 6834 and Stainless steel type 304 according to the standard ASTM F899 | |
| | | 03 Tube Ø9,5 x 180mm | 304 stainless steel according to ASTM F899 | |
| | | 01 Tube Ø9,5 x 340mm | 304 stainless steel according to ASTM F899 | |

Symbol of the Traceability Label (Label) according to NBRISO15223

| | | | | | |
|--------------------|---------------------|------------|------------|---|----------------------|
| | | | | | |
| LOT xxxxx | xxxx | xxxx-xx | REF xxx-xx | STERILE EO | |
| Product Lot Number | Date of Manufacture | Shelf Life | Code | Identification for products sterilized in ETO Gas | "SINGLE USE PRODUCT" |

TECHNICAL INFORMATION

Indications, Precautions, Restrictions / Warnings, Adverse Effects and Contraindications: The external fixator is a product intended to be used on

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lower and upper limb bones, in corrective surgeries or for trauma treatment, maintaining bone fracture and stabilizing bone or adjacent structures, in order to provide bone fusion or bone repair. The surgeon should be familiar with and have sufficient knowledge of osteosynthesis and its limitations, including pre and post operative, surgical technique, precautions and potential risks. Preoperative and surgical procedures, including knowledge of the surgical technique, selection and appropriate placement, are important considerations for the successful use of temporary external fixation devices. Do not use the product if it is damaged. Do not use the product if the packaging is damaged. Possible adverse effects that may occur in the use of external fixator are: superficial or deep postoperative infection; pain, discomfort or abnormal sensations due to the presence or migration of the medical product, and also due to the surgical procedure and migration, limb edema, pseudoarthroses, movement limitation, non-union of the fracture. Due to the incompatibility of materials, designs and dimensions it is prohibited to use components from different manufacturers, the use of separate components can cause product problems and bring risks, eliminating the responsibility of the company to any warning that may occur. Components of different manufacturers should not be used. Therefore we recommend that the products have the same origin.

Sterilization: This product is supplied in sterile form. Sterilization is done by Ethylene Oxide. If the product has expired sterilization validity or damaged packaging, it must be returned to the responsible distributor or directly to Biomecanica.

Product storage and transport: The storage location must be clean, dry and illuminated in order to maintain the ideal storage and transport conditions, as well as its physical and chemical integrity. The product should be stored and transported in a dry and cool place, at room temperature (Max 35 ° C), relative humidity around 30% to 70% and protected from direct sunlight. The effects of vibration, shocks, corrosion, defective seating during handling and transportation, improper stacking during storage, should be avoided. The product should be handled with care to avoid damage that could impair the quality of the material and the safety of the patient. Can not be stored directly on the floor (minimum height = 20cm). They can not be placed in very high places, near bulbs, which could cause dryness of the packaging or damage to the label. They should not be stored in places where contaminants such as cleaning materials, insecticides, pesticides, etc. are stored. The carrier must be informed about the content and deadline of delivery. Carry carefully to avoid falling and friction in order not to damage the packaging and the device, thereby ensuring the sterility of the material. The opening of the sealed blister pack can be done manually, in a suitable place (surgical environment), by pulling the surgical grade paper by the leftover around the blister. This opening should be performed only in the surgical act to avoid contamination of the product.

Identification: The external fixators are labeled with the following information: Manufacturer's Data; Technical manager; Registration number in ANVISA; Technical Product Name; Product Commercial Name; Commercial Product Model; Raw Material; Number of items in the package; Lot number; Product code; Date of manufacture; Shelf life; Manufacturer's address; External fixators contain the following laser-engraved information to enable traceability: Biomecanica logo; Manufacturing Lot Number and Code. Manufacturer recommends Single Use.

Customer Service: If there is a need to make any complaint regarding the use of external fixators **related to any adverse effect that affects the safety of the user, the surgeon in charge should communicate this adverse event to the competent sanitary body and Biomecanica through the email sac@biomecanica.com. or by telephone 0xx14 2104 7900.** In case of doubts the responsible surgeon or the health professional can make the communication of the adverse event through the System of Notifications in Sanitary Surveillance in the ANVISA site: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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