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ENGLISH

USA

Detailed description of the medical product, including principles of operation and action, its content or composition, when applicable, and list of accessories intended to integrate the product.

Orthopedic Cement for Vertebroplasty - BM, Figure 1, is a self-curing acrylic compound, obtained **by mixing a powdered component**, based on poly (methyl methacrylate) and barium sulfate with a polymerization initiator, benzoyl peroxide; and liquid component, consisting of methyl methacrylate monomer with NN dimethyl-p-toluidine activator.

Figure 1 - Presentation of Orthopedic Cement for Vertebroplasty - BM, code 2010-00-000

Orthopedic Cement for Vertebroplasty - BM is specifically formulated to perform percutaneous vertebroplasty procedures, Figure 2. In this procedure, the collapsed vertebra is filled with bone cement by using a minimally invasive needle.

Fractured Vertebra

Bone Cement Injection

Figure 2 - Example of using the product Orthopedic Cement for Vertebroplasty - BM.

When the monomer and polymethyl *methacrylate* come into contact (mixed properly), *NN dimethyl-p-toluidine* activates the *benzoyl peroxide*, initiating the polymerization.

The process is progressive, and the liquid product in the initial moments becomes pasty until it reaches the solid state. The reaction is exothermic and with maximum temperatures that can vary depending on the volume and heat dissipation, reaching up to 90°C.

Presentation

- **01 liquid component (20ml):** colorless, flammable monomer with a characteristic odor, aseptically filled in a glass ampoule. Table 1 shows the composition of the liquid component.

Table 1 - Composition of the liquid component.

Component	Content
Methyl methacrylate	99% w/w
NN Dimethyl-p-toluidine*	1% w/w

Hydroquinone**	30ppm
<p>** Hydroquinone is added to prevent premature polymerization that can occur under certain conditions, such as high temperature and exposure to light (see pharmaceutical precautions).</p> <p>* The NN Dimethyl-p-toluidine is added to promote polymerization after mixing the liquid and powder components of the cement.</p> <p>*** Sterile liquid by microfiltration.</p>	

- **01 powder component (29g):** White powder, finely particulate. Table 2 shows the powder component composition.

Table 2 - Composition of the powder component.	
Component	Content
Polymethyl methacrylate	67,5% w/w
Barium Sulphate*	31,0% w/w
Benzoyl Peroxide	1,5% w/w
<p>*Barium Sulfate provides a radiopacity characteristic to the orthopedic cement.</p> <p>***Supplied sterile and in double blister.</p>	

Accessories

Orthopedic Cement for Vertebroplasty - BM has no accessory.

List of Instruments

To assist in the setting of Orthopedic Cement for Vertebroplasty - BM it is necessary to use instruments contained in the sterile kit: spatula and mixing container.

Ancillary Components

Orthopedic Cement for Vertebroplasty - BM does not have ancillary components.

Indication, purpose or use intended for the medical device, as indicated by the manufacturer.

Indication of Use

Vertebroplasty and / or percutaneous kyphoplasty is a minimally invasive procedure in which bone cement is injected into the vertebral body, i.e., into the bone of the spine. Vertebroplasty and / or kyphoplasty are indicated for patients with osteolytic lesions of the vertebral body resulting from metastases, myeloma, vertebral hemangioma and osteoporosis lesions, according to the following criteria:

- a) Hemangioma with aggressive clinical and / or radiological behavior.
- b) Severe and disabling pain, resistant to the use of analgesics.
- c) Progression of deformity and increase due to osteolytic lesions of the vertebral body.
- d) When the vertebral body is not completely destroyed, it must have at least 1/3 of its original height.

Liability Disclaimer

Biomechanica, as the manufacturer of this device, does not practice medicine or recommend this or any other surgical technique for use on a particular patient. The surgeon who performs any implant procedure is responsible for defining and using

the appropriate techniques to implant the prosthesis in each patient. BIOMECANICA is not responsible for choosing the appropriate surgical technique for use on each patient. Surgical Techniques depend on the technique adopted by the surgeon, who must choose the method, type and dimension of the products to be placed.

Preparation and Application

For the preparation, the product must be stored at a temperature of 23°C. If the product has been stored at a different temperature, it must be stored for at least 2 hours at 23°C to reach room temperature.

Orthopedic Cements are very sensitive to heat, and any increase in the room temperature or the cement components above 23°C may reduce preparation, handling, and polymerization times. On the other hand, lower temperatures will increase the time.

The packaging must be opened by a circulating nurse. Sterile content should be transferred aseptically to the sterile operating area. Break the ampoule containing the liquid with the aid of gauze, and pour ALL the liquid content into the powder in the blister. Use a sterile spatula that accompanies the kit to mix.

To prevent any possible contamination of the cement with glass fragments DO NOT BREAK THE AMPOULE ON THE CONTAINER.

ADD THE LIQUID TO THE POWDER, NOT THE POWDER TO THE LIQUID, IF YOU DO NOT USE A MIXER.

Mixing Instructions

The mixing of the Orthopedic Cement for Vertebroplasty - BM must be done for **1 minute (60sec)** after adding the liquid monomer to the powder component with the aid of a sterile spatula.

Homogenize at **15-second** intervals until achieving the working consistency.

Homogenization should be done with slow movements in order to avoid the incorporation of air bubbles into the product, which can impair its mechanical performance.

The correct working consistency of the product for application to the bone must be determined by the surgeon's experience and needs.

Digital Application

Select an appropriate method for applying the cement to the bone.

The Orthopedic Cement for Vertebroplasty - BM has an application work period of 03 to 09 minutes from the beginning of the mixing, a time of 06 minutes of application and curing time (polymerization time) between 20 - 30 minutes from the time mixing, considering the temperature of 23°C.

Handling and hardening characteristics of Orthopedic Cement for Vertebroplasty - BM may vary depending on the temperature and mixing technique.

Technical Information

(Varies with room and storage temperature. See Preparation and application).

Adhesion Time* (doughing time)	Polymerization Time* (setting time)
11 min - 15min	20 min - 30 min


Orthopedic Cement For Vertebroplasty - B, (23°C)

Adhesion and polymerization times are counted from the beginning of the mixing of the components at a temperature of 23°C.

Approximated Average Times And Polymerization Phases At Different Temperatures (°C) For Orthopedic Cement For Vertebroplasty - BM

The mass and hardening times of Orthopedic Cement for Vertebroplasty - BM vary with temperature, as shown in the graph below.

Precautions, restrictions, warnings, special care and, clarifications on the use of the medical product, as well as its storage and transportation.

Contraindications

The Procedure is contraindicated:

- A. In the presence of infection;
- B. In uncompensated coagulation disorders;
- C. In case of destruction of the posterior wall of the vertebral body, its use must be evaluated, since the potential for extravasation of the injection material into the space of the spinal canal can produce compression of nervous structures;
- D. When there is hypersensitivity to any of the components of the formula.

Interaction With Other Medicines

Unknown.

Adverse Effects

The most severe and frequent adverse effects, which can occur with the use of bone cement, are as follows. Surgeons must be aware of such reactions, and prepared to treat them if applicable.

- *Severe:*
 - 1) Myocardial infarction;
 - 2) Cerebrovascular accident;
 - 3) Cardiac arrest;
 - 4) Sudden death; and
 - 5) Pulmonary embolism.

- *More frequently:*
 - 1) Transient drop in blood pressure;

- 2) Thrombophlebitis;
 - 3) Hemorrhage and hematoma;
 - 4) Superficial or deep infection of the lesion;
 - 5) Short-term cardiac conduction irregularity.
- **Other Potential Adverse Events:**
Other potential associated adverse events include:
 - (1) allergic pyrexia;
 - (2) hematuria;
 - (3) dysuria;
 - (4) local neuropathy; and
 - (5) local vascular erosion and occlusion.

Use in Pregnancy and Breastfeeding

Unknown.

Warnings

Adverse cardiovascular reactions, including hypotension, hypoxemia, cardiac arrhythmia, bronchospasm, cardiac arrest, myocardial infection, pulmonary infarction, stroke and possible death: hypotensive reactions may occur between 10 and 165 seconds after application of PMMA bone cement, which can last from 30 seconds to more than 5 minutes. Some hypotensive reactions have progressed to cardiac arrest. Patients' blood pressure should be carefully monitored during and immediately after application of PMMA bone cement. In addition, under pressurization of PMMA bone cement during insertion should be avoided in order to minimize the occurrence of pulmonary embolism.

Experience and training of the surgeon: The surgeon must be fully familiar with the properties, handling characteristics and application of PMMA bone cement. As the handling and harde-

ning characteristics of the cement vary with temperature and application technique, they are better determined according to the surgeon's real experience.

Respiratory tract, eye and liver irritation: Care should be taken when mixing the liquid and powder components of the orthopedic cement to prevent excessive exposure to concentrated vapors of the liquid monomer, which can cause irritation in the respiratory tract, the eyes, and, possibly, the liver. People wearing contact lenses should not mix PMMA orthopedic cement.

Contact dermatitis: The liquid component is a powerful lipid solvent. It can cause contact dermatitis in susceptible individuals. The use of a second pair of surgical gloves and strict adherence to the mixing instructions can decrease the possibility of hypersensitivity reactions. The compound must not be allowed to come into direct contact with sensitive tissues, i.e., absorbed by the body. The liquid component must not be allowed to come into contact with rubber, including rubber gloves.

PRECAUTIONS

1. The use of Orthopedic Cement requires collaboration and consultation between surgeon and anesthetist.
2. The anesthetist must be informed, during the operation, of the moment of setting the cement.
3. The monomer can cause hypotension; the drop in blood pressure can cause cardiac arrhythmias or ischemic myocardium.
4. Hypotensive effects of methyl methacrylate can be enhanced if the patient suffers from hypovolemia.
5. The preparation of the bone marrow cavity results in marrow contents entering the bloodstream. By widening a hole in the medullary cavity, one can have similar effects on blood pres-

sure, such as the introduction of Orthopedic Cement for Vertebroplasty - BM.

6. The medullary cavities must be purged when the Orthopedic Cement for Vertebroplasty - BM is introduced digitally.

7. The complete irrigation of the medullary cavity during preparation reduces the risk of the medullary content being introduced into the vascular system while setting the Orthopedic Cement for Vertebroplasty - BM.

8. The implantation of a foreign body in the tissues increases the normal risk of infection after surgical interventions.

9. As for the application of cement, it is important to maintain the position of the prosthetic component until the completion of the polymerization process.

10. The use of an appropriate prophylactic antibiotic may be desirable in the clinician's prescription.

User Precautions

The liquid and powder components of Orthopedic Cement for Vertebroplasty - BM are produced for use in a single dose. Always add the entire liquid component to the powder component when mixing the Orthopedic Cement for Vertebroplasty - BM in the container, when using a mixer add the powder to the liquid. The liquid component is volatile. The operating room must be properly ventilated. Continuous inhalation of liquid vapor can have a soporific effect or also irritate the respiratory tract and the eyes, which should be avoided.

The liquid component is a powerful solvent that can cause contact dermatitis in sensitive people. The possibility of hypersensitive reactions is reduced by wearing a second pair of gloves and strictly following the preparation instructions.

If the liquid component has direct contact with the eyes, wash

them with plenty of water. Concentrated vapors from the liquid component may have an adverse reaction with soft contact lenses.

Storage

Store the product in a dry, ventilated place away from light, with a temperature of up to 25°C.

Otherwise, spontaneous polymerization of the liquid component may occur.

Sterilization

Product sterilized by Ethylene Oxide - ETO - Validity 3 years.

Sterilization by microfiltration.

This product must not be re-sterilized.

If the surgical packaging is damaged, it must be disposed of.

Note: If the product has expired sterilization and / or damaged packaging, it must be returned to the responsible distributor or directly to Biomecanica.

After Sales (Customer Complaint)

When required to make a complaint regarding the Orthopedic Cement for Vertebroplasty - BM related to any technical complaint or adverse event that affects the safety of the user, severe problems, or death related to these components, the responsible surgeon must report such adverse event to the competent health agency and Biomecanica through the email tecnovigilancia@biomecanica.com.br or by phone 0xx14 2104 7900. In cases of doubt, the surgeon in charge or the health professional may communicate the technical complaint or adverse event through the Health Surveillance Notification System on the AN-VISA website: <http://portal.anvisa.gov.br/notificacoes>.

Product Disposal


Eliminate the rest of the Orthopedic Cement for Vertebroplasty - BM and the partially used, expired or damaged content, following the rules and procedures applicable to this type of hospital waste.

Labeling

The external labeling is the product's traceability label that contains the following information:

	<p>PRODUCTS STERILIZED IN ETHYLENE OXIDE GAS Manufacturer Information: BIOMECANICA INDÚSTRIA E COMÉRCIO DE PRODUTOS ORTOPEDICOS LTDA. Address: Rua Luiz Pengo, 145; 1º Distrito Industrial; CEP: 17212-811 Jaú / SP Brazil - Phone: +55 (14) 2104-7900 CNPJ: 58.526.047 / 0001-73 / INSC EST: 401.042.207.113 Resp. Technician - Josiane da Silva Taletti - CRQ-SP nº 04272807 www.biomecanica.com.br / biomecanica@biomecanica.com.br Aut. Func. ANVISA / MS nº X353Y38XY66Y</p>
<p>Technical Person Responsible</p>	<p>Quantity of items per packaging</p>
<p>ANVISA registration</p>	<p>Quantity and Description of the Labels that come with the Product: - 6 adhesive labels for product identification and traceability</p>
<p>Technical Product Name</p>	
<p>Commercial Product Name</p>	
<p>Commercial Product Model</p>	
<p>Raw Material</p>	<p>In writing: "NO REPROCESSING"</p>

**Traceability Label (Label) Symbols
according to NBR ISO 15223**

	Product Lot Number.
	Date of Manufacturing.
	Use until.
	Code.
	Identification for products sterilized in ETO Gas.
	Identification for products sterilized through aseptic processing techniques.
	"NO REUSING".
	"NO RESTERILIZING".

Information in the Cardboard Box (External Packaging)

Cardboard boxes (external packaging) also have special phrases for Storage, Conservation and / or Product Handling with symbols. Table 3 describes the information contained in this cardboard box.

Table 3 - Special Conditions for Product Storage, Conservation and / or Handling contained in the external packaging

Aseptically Open the Packaging

Maximum Temperature Limit

Store and transport in a cool, dry place with a temperature lower at 25°C and away from light.

Do not store directly on the floor (minimum height = 20 cm) or in very high places, close to lamps, which could cause dryness of the packaging or damage to the label.

Do not store in places where contaminating substances such as cleaning materials, insecticides, pesticides etc. are stored.

Warning

Before using read the internal instructions.

Do not use if the packaging is damaged.

Keep dry.

	<p>Flammable liquid.</p>
	<p>Sterile Product Identification.</p>
	<p>Fragile. Handle with care.</p>
	<p>Keep out of sunlight.</p>
	<p>“SEE INSTRUCTIONS FOR USE”.</p>
	<p>Degree Celsius.</p>
	

Figure 3 illustrates the label model used externally and internally in the packaging to identify and maintain product traceability. Figure 3 illustrates the label of products sterilized in ETO.

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