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ENGLISH USA

DESCRIPTION

Bone Cement - Normal Viscosity, code 2005-00-000, is a self-curing acrylic compound obtained by mixing a powder component with a liquid component.

The powder component consists of Polymethyl Methacrylate and Barium Sulphate with a polymerization initiator, Benzoyl Peroxide. Barium sulphate is the agent that gives the cement its radiopaque properties.

The liquid component consists of the monomer Methyl Methacrylate with the activator N,N-Dimethyl-p-Toluidine and Hydroquinone. Hydroquinone is a stabilizer added to prevent premature polymerization that can occur under certain conditions, such as exposure to high temperatures and light.

When the powder and liquid components are brought into contact (mixed properly), the N,N-Dimethyl-p-Toluidine activates the Benzoyl Peroxide which starts the polymerization. The polymerization process is progressive, the mixture starts out liquid and becomes dough in the next few moments until it reaches a solid state. The mixture triggers an exothermic reaction with maximum temperatures that can vary depending on the volume and heat dissipation, reaching up to 90 °C.

FORM OF PRESENTATION

- 2005-00-000 - Bone Cement - Normal Viscosity

Each package of Bone Cement - Normal Viscosity consists of: 1 ampoule with 20 mL of the liquid component, 1 blister with 40 grams of the powder component and a spatula.

The set is supplied sterile in a double blister and packed in a cardboard box with a plasticized finish.

01 liquid component (20 mL): Colorless, flammable monomer with a characteristic odor, aseptically packaged in a glass ampoule. Sterile by microfiltration. Table 1 below shows the composition of the liquid component.

Table 1 - Composition of the liquid component.

	Component	Contents	Reference Standard
	Methyl Methacrylate (MMA) N,N-Dimethyl-p-Toluidine (DMPT)	99% m/m 1% m/m	ISO5833:2002 ASTM F451:21
	Hydroquinone (HQ)	30 - 50 ppm	

* N,N- Dimethyl-p-toluidine is added to promote polymerization after mixing the liquid and powder components of the cement.

** Hydroquinone is added to prevent premature polymerization which can occur under certain conditions, such as high temperature and exposure to light (see pharmaceutical precautions).

*** Sterile liquid by microfiltration.

01 powder component (40 g): White, finely particulate powder. Table 2 below shows the composition of the powder component.

Table 2 - Composition of the powder component.

	Component	Contents	Reference Standard
	Polymethylmethacrylate (PMMA)	88% m/m	ISO5833:2002 ASTM F451:21
	Barium Sulphate (BaSO ₄)	10% m/m	
Benzoyl Peroxide (BPO)	2% m/m		

*Barium sulphate provides the radiopacity characteristic of orthopedic cement.

**Supplied sterile and in a double blister pack.

01 Spatula: White, used only for mixing. Table 3 below shows the composition of the spatula.

Table 3 - Spatula composition.

	Component	Content
	Polypropylene	01 unit

Table 4 shows how the sterile set is presented.

Table 4 - Presentation of the sterile set.

	Component	Content
	Sterile Bone Cement Set - Normal Viscosity	01 unit

ACCESSORIES

Bone Cement - Normal Viscosity has no accessories.

LIST OF INSTRUMENTS

To assist in the implantation of Bone Cement - Normal Viscosity, it is necessary to use the spatula contained in the sterile set and the primary packaging (blister) of the powder component, which is used as a mixing container.

The Biomechanica Cementing Kit, which is not the subject of this registration, can be used to aid in the implantation of Bone Cement - Normal Viscosity.

ANCILLARY COMPONENTS

Bone Cement - Normal Viscosity has no ancillary components.

INDICATION, PURPOSE, OR USE FOR WHICH THE MEDICAL PRODUCT IS INTENDED, AS INDICATED BY THE MANUFACTURER.

INDICATIONS FOR USE

Bone Cement - Normal Viscosity is indicated for fixing prostheses to living bone in musculoskeletal orthopedic surgical interventions in cases of rheumatoid arthritis, osteoarthritis, traumatic arthritis, osteoporosis, avascular necrosis, collagen disease, severe secondary joint destruction following trauma or other conditions and in the revision of previous arthroplasty procedures.

PRECAUTIONS, CONTRAINDICATIONS, WARNINGS, SPECIAL CASES OF POLYMERIZATION OF THE LIQUID COMPONENT, PRECAUTIONS, AND ADVERSE EVENTS.

CONTRAINDICATIONS

The use of Bone Cement - Normal Viscosity is contraindicated in the presence of active or incomplete infection at the site where the bone cement is to be applied.

The use of Bone Cement - Normal Viscosity is contraindicated for patients who are allergic or hypersensitive to any of its components.

WARNINGS

Carefully monitor patients for any changes in blood pressure during and immediately after the application of bone cement. Adverse reactions affecting patients' cardiovascular systems have been associated with the use of bone cements. Hypotensive reactions occur between 10 and 165 seconds after application of the bone cement; they last from 30 seconds to 5 minutes or more. Some progress to cardiac arrest. Patients should be carefully monitored for any changes in blood pressure during and immediately after the application of bone cement.

Follow the handling instructions carefully.

- Care should be taken when mixing the two components to avoid excessive exposure to concentrated monomer vapors, which can cause irritation to the respiratory tract, eyes, and possibly the liver. People wearing contact lenses should not be near or involved in mixing this bone cement.
- The polymerization of bone cement is an exothermic reaction, which occurs while the cement is hardening *in situ*. The heat released can damage bone or other tissues around the implant.

- Inappropriate fixation or unforeseen post-operative events can affect the cement-bone interface and lead to micro-movement of the cement against the bone surface. A layer of fibrous tissue can develop between the cement and the bone and loosening of the prosthesis can occur, leading to implant failure. Regular long-term follow-up is recommended for all patients.
- Do not allow the liquid component to come into contact with rubber gloves. The liquid component is a powerful lipid solvent. If contact occurs, gloves may dissolve and tissue damage may occur. The use of a second pair of surgical gloves and strict adherence to the mixing instructions can reduce the possibility of hypersensitivity reactions. Once mixed, bone cement should not come into contact with the gloved hand until the cement has acquired the consistency of the mass. This usually occurs between one and two minutes after mixing the liquid and powder components.

Avoid excessive pressurization of the bone cement, as this can lead to extrusion of the bone cement beyond the intended application site and damage to the surrounding tissues.

The safety of bone cement in pregnant women or children has not been established. Bone cement can adversely affect bone growth and fetal health.

The hypotensive effects of methylmethacrylate can be potentiated if the patient suffers from hypovolemia.

Preparation of the bone marrow cavity results in marrow contents entering the bloodstream. Enlarging a hole in the marrow cavity can have similar effects on blood pressure as the introduction of Bone Cement - Normal Viscosity.

The medullary cavities must be purged when Bone Cement - Normal Viscosity is introduced digitally. Complete irrigation of the medullary cavity during preparation reduces the risk of medullary content being introduced into the vascular system during the placement of Bone Cement - Normal Viscosity.

The implantation of a foreign body in the tissues increases the normal risk of infection after surgical interventions.

When applying the cement, it is important to maintain the position of the implant until the polymerization process is complete.

PRECAUTIONS

This product is for single use only and must not be reprocessed.

Do not use the product if the packaging is damaged. If the sterile packaging is damaged, the product must be discarded. Do not reesterilize.

Follow the handling instructions to avoid contact dermatitis. Strict adherence to the instructions for mixing the liquid and powder components can reduce the incidence of this complication.

Ventilate the operating room properly to eliminate as much monomer vapor as possible. Liquid monomer is highly volatile and flammable. Ignition of monomer vapors caused by the use of electrocautery devices has been reported at surgical sites near recently implanted bone cements.

Do not use this product after the expiration date printed on the packaging. This device may not be safe or effective after the expiration date.

Do not mix components from different manufacturers. The use of components from different manufacturers introduces risks that have not been taken into account, eliminating any responsibility on the part of Biomecanica for adverse events that may occur.

The surgeon must, through specific training and experience, be completely familiar with the properties, handling characteristics, and application of bone cements. As the handling and curing characteristics of this bone cement range with temperature, humidity, and mixing technique, they are best determined by the surgeon's actual experience.

The use of bone cement requires collaboration and consultation between the surgeon and the anesthesiologist. The anesthesiologist must be informed during the operation when bone cement is to be implanted.

Biomecanica is not responsible for choosing the appropriate surgical technique for each patient. Surgical techniques depend on the technique adopted by the surgeon, who is responsible for choosing the method, type, and size of the products to be placed.

Special cases of polymerization of the liquid component: (1) in constant exposure to sunlight, directly or indirectly; and (2) in contact with or exposure to heat (consider additional precautions, as the product is flammable).

In view of the possible cases of polymerization, it is recommended that extra doses of Bone Cement - Normal Viscosity be taken to the surgeries.

The patient must be informed of the importance of post-surgical follow-up. Lack of follow-up prevents the detection of post-surgical problems.

The patient must be informed of the type of implant used, brand, code, name, batch, surgeon's name, and date of surgery.

ADVERSE EVENTS

Serious adverse events, some fatal, associated with the use of acrylic bone cements include myocardial infarction, cardiac arrest, stroke, and pulmonary embolism.


The most frequent adverse reactions reported with acrylic bone cements are a transient drop in blood pressure, thrombophlebitis, hemorrhage and hematoma, loosening or displacement of the prosthesis, superficial or deep wound infection, trochanteric bursitis, and short-term cardiac conduction irregularities. Other adverse reactions reported include heterotopic new bone formation and trochanteric separation.

Other adverse events reported for acrylic bone cements include pyrexia due to an allergy to the bone cement, hematuria, dysuria, bladder fistula, delayed sciatic nerve entrapment due to extrusion of the bone cement beyond the region of its intended application and adhesions and ileal stenosis due to the heat released during polymerization.

HANDLING

A. Preparation

Before starting to prepare Bone Cement - Normal Viscosity, the product must be at a temperature of 23 °C. If the product has been stored at a different temperature, it must be stored for at least 2 hours at 23 °C to balance to that temperature. Bone cements are very sensitive to heat and any increase in ambient temperature or in the temperature of their components above 23 °C can result in a reduction in the cement's adhesion and polymerization times and maximum temperature. In turn, lower temperatures increase these times.

- 1** Place your thumb on the indicated break point.
- 2** Apply force to the thumb while holding the ampoule completely.
- 3** With the aid of gauze, hold as indicated.
- 4** Apply force at the indicated point until it breaks.

To prepare the bone cement, open the carton pack. The sterile packaging is double-blistered. The first blister contains two blisters, one containing the powder component and the other containing the liquid component in a glass ampoule. This first blister must be opened by a circulating nurse. Both blister packs should be transferred aseptically to the sterile operating area. Open the blister containing the powder component and open the blister containing the glass ampoule. Proceed with breaking the ampoule using gauze, as shown in the figure below.

To prevent any possible contamination of the cement with fragments of glass, **DO NOT BREAK THE AMPOULE OVER THE OPEN BLISTER CONTAINING THE POWDER.**

Care must be taken during handling so that there is no contamination or insertion of other components during mixing that could modify the original characteristics of the product and change its properties.

B. Mixing

Empty ALL the liquid contents of the ampoule onto the blister containing the powder. Use the sterile spatula supplied with the product to mix. Mixing can be done in a blender.

If mixing by hand in the blister containing the powder, **ADD THE LIQUID ON TOP OF THE POWDER, NOT THE POWDER ON TOP OF THE LIQUID.**

If mixing in a blender, **ADD THE POWDER ON TOP OF THE LIQUID.**

Bone Cement - Normal Viscosity should be mixed by hand for **1 minute** after adding the liquid component to the powder component using a sterile spatula. Homogenize at **15-second** intervals. After **2 minutes**, when it is in a dough point, touch the surface of the mixture with a glove until the cement no longer adheres to the surface of the glove.

If the cement is applied to gloved hands too soon, it may initially stick to the gloves, making it necessary to knead (homogenize) the cement for longer.

Homogenization must be carried out with slow movements in order to avoid the incorporation of air bubbles into the product, which can impair its mechanical performance.

C. Digital Application

- Before applying Bone Cement - Normal Viscosity, the bone cavity must be thoroughly cleaned and dried to prevent blood and foreign bodies from mixing with the cement, which can impair mechanical performance.

- During the digital application of Bone Cement - Normal Viscosity, it is important to use a thin plastic probe to allow blood, air, and fat to escape. The probe must be removed before inserting the prosthesis.
- Bone Cement - Normal Viscosity has no adhesive properties but it does have a mechanical bond with the irregular surfaces of bone and prostheses.

Bone Cement - Normal Viscosity is ready for application in approximately **5 minutes**, counting from the time the liquid component is mixed with the powder. Digital application of the cement should take between **1 and 2 minutes** (after this time, the cement becomes too viscous to insert into the narrow cavity).

After digital application, the prostheses should be inserted and held firmly in position until the Bone Cement - Normal Viscosity finishes polymerizing, which takes between **8 and 11 minutes**.


The handling and hardening characteristics of cement can range depending on the temperature and mixing technique used.

Excess Bone Cement - Normal Viscosity must be removed before the cement has completely set.


TECHNICAL INFORMATION

Adhesion time <i>(doughing time)</i>	Polymerization time <i>(setting time)</i>
2' to 3'30"	8' to 11'

The adhesion and polymerization times are counted from the start of mixing the liquid and powder components at 23 °C. These times may vary with room temperature and storage temperatures.

The graph below shows the reaction times at temperatures between 20 °C and 24 °C for mixing and manual application; if using a digital applicator, consider the times for Phase I, Phase II, Phase III, and Phase IV.

TRACEABILITY

The packaging of the implantable component comes with three traceability labels containing the data of the implant used. The traceability label contains the following information: the product's commercial name and code, the manufacturer's company name, the batch number, and the ANVISA registration number. One label must be mandatorily affixed to the health record, another to the document to be given to the patient, and another to the tax document that generates the charge. Please note that the patient's health record must record information that allows the implanted product to be traced. Sticking this label on the patient's health record allows the product used to be traced. Among the main pieces of information, we highlight the following as essential: the name of the manufacturer of the implant used, the name of the implant used, the implant code, the batch number, and the product's ANVISA registration number. This information is described on the traceability labels that come with the product and on the external labeling. Other information should also be considered essential, such as the date of the surgery, the name of the patient who received the implant, the name of the surgeon, the patient's weight, the patient's age, and other information requested in the patient's health records.

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

Figure 1 - Traceability label model

STERILIZATION

Product sterilized with Ethylene Oxide (EtO) - Shelf life: Three years.

Sterilization of liquid by microfiltration.

This product must not be resterilized.

If the surgical packaging is damaged, it should be discarded.

Note: If the product has expired and/or damaged packaging, it must be returned to the distributor responsible or directly to Biomechanica.

HANDLING STERILE MATERIAL

When handling sterilized material with aseptic technique, certain guidelines must be followed in order to keep it sterile:

- It is essential to wash your hands with soap and water before handling sterilized material;
- Use material with intact, dry packaging, without stains, with identification (type of material and date of sterilization);
 - Work facing the material;
 - Handle the material from the waist up;
 - Avoid coughing, sneezing, or talking on exposed material;
 - Do not move over the sterilized area;
 - Make sure the packaging is valid and suitable;
 - Work in a clean, quiet, dry, and draught-free environment;
 - Keep a certain distance between the body and the material to be handled;
 - Obey the other principles of asepsis.

STORAGE AND TRANSPORTATION OF MEDICAL PRODUCTS

- Store the product in a cool, dry place away from light at a temperature of up to 25 °C. Otherwise, spontaneous polymerization of the liquid component may occur. Humidity should be between 30% RH and 70% RH and keep away from the elements.

- Do not store directly on the floor (minimum height = 20 cm) or in very high places near lamps, as this could cause the packaging to dry out or damage the label.

- Do not store in places where contaminating substances are stored, such as cleaning materials, insecticides, pesticides, etc.

- When transporting it, avoid bumps and improper stacking.

- Transport in a cool, dry place with an ambient temperature below 25 °C. Otherwise, spontaneous polymerization of the liquid component may occur. Storage humidity should be between 30% RH and 70% RH.

CUSTOMER COMPLAINT

If there is any need to make a complaint about Bone Cement - Normal Viscosity related to a technical complaint or adverse event that affects user safety, serious problems or death related to these components, the surgeon responsible must report this adverse event to the competent health agency and to Biomecanica by e-mailing tecnovigilancia@biomecanica.com.br or calling 0xx14 2104 7900. If in doubt, the responsible surgeon or healthcare professional can report the technical complaint or adverse event via the Health Surveillance Notification System on the ANVISA website: <http://portal.anvisa.gov.br/notificacoes>.

RECEIPT OF COMPLAINTS

Complaints can be received verbally (telephone/in person) or in writing (e-mail/letter). When the complaint is received, the person responsible must distinguish whether the complaint is related to an adverse event with an implantable product or instrument, a technical complaint that may lead to an adverse event and/or a technical complaint without an adverse event. Once the complaint has been recorded, it is forwarded to the Technovigilance Department, which in turn submits the complaint for evaluation by the Technical Chamber. This procedure is in harmony with that specified in the resolutions: RDC 67/2009 (amended by RDC 557/2021 and amended by RDC 702/2022).

Recall: If applicable, the Technical Manager will decide to recall the batch from the market and segregate the products still in stock, depending on the seriousness of the complaint.

If a manufacturing problem causes the affected batch(es) to be withdrawn from the market, this must be reported to the competent authorities, in accordance with procedure **PRTEC02 - Recall**.

DISPOSAL

According to Resolution RE No. 2605 of August 11, 2006, implantable materials of any kind classified as single-use are prohibited from being reprocessed.

If Bone Cement - Normal Viscosity has been removed from its internal packaging and inserted into the surgical environment, even if it has not been used or contaminated by other sources, and if it is not intended for further study and/or analysis, it must be properly disposed of by the hospital.

Expired products should also be discarded.

The rest of the Bone Cement - Normal Viscosity and the partially used, expired, or damaged contents must be disposed of following the rules and procedures applicable to this type of hospital waste.

The methods and procedures for disposing of the product used must ensure that it is completely uncharacterized, preventing any possibility of its reuse. Disposal of the medical product is the sole responsibility of the hospital, as are the disposal methods and procedures used. However, Biomecanica recommends that such products are clearly and visibly identified, indicating their status, i.e., unfit for use.


LABELING

Next to the secondary packaging of the implantable component, one label is attached containing the data of the implant used and one red label which indicates that this content is sterile. Three (03) traceability labels are also attached to the product packaging and one label is attached to the primary packaging containing the following information shown in Table 5.

Table 5 - Identification Labels

Symbolologies:	
	Catalog Number
 xxxxxx	Lot Number
 xxxx	Date of Manufacture
 xxxx	Use by the date
	Do not reuse
	Do not reesterilize
	Sterilized using aseptic processing techniques with double sterile barrier system
	Sterilized with ethylene oxide with double sterile barrier system
	Manufacturer

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Sterile Product Label

LABELING SECONDARY PACKAGING

The symbols on the outer packaging were designed in accordance with the requirements of the ANBT NBR ISO15223-1 standard, as shown below in Table 6.

Table 6 - Secondary packaging wording and symbols

Wording:	
<ul style="list-style-type: none"> - Open the package aseptically. - Store and transport in a cool, dry place, below 25 °C and away from light. - Do not store directly on the floor (minimum height = 20 cm) or in very high places near lamps, which could cause the packaging to dry out or damage the label. - Do not store in places where contaminating substances are stored, such as cleaning materials, insecticides, pesticides, etc. - Contains 01 Liquid Component (20 mL): Methylmethacrylate (MMA) 99% w/w; NN Dimethyl-p-toluidine (DMPT) 1% w/w; Hydroquinone (HQ) 30-50 ppm. Contains 01 Powder Component (40 g): Polymethylmethacrylate (PMMA) 88% w/w; Benzoyl Peroxide (BPO) 2% w/w; Barium Sulfate (BaSO₄) 10% w/w. Contains 01 polypropylene spatula. 	
Symbologies:	
	<p>Attention Read the internal instructions before use</p>
	Do not use if packaging is damaged
	Keep dry
	Flammable liquid
	Fragile. Handle with care
	Keep out of the sun
	See instruction and use
	Maximum Storage Temperature Limit
	Storage Humidity Limitation



	Manufacturer
<div style="border: 1px solid black; padding: 2px; display: inline-block;">STERILE</div>	Sterile with double sterile barrier system

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WARNING TO THE USER

The Instructions for Use for the correct use and handling of the product are available at www.biomecanica.com.br. They can be checked in the search field using the ANVISA registration number. The instructions for use provided will always be in accordance with the most current version. The revised Instructions for Use are identified on the product labeling. If the user is interested, the Instructions for Use can be supplied in printed form, free of charge. Requests for these instructions should be made to the manufacturer's Customer Service Channel, as indicated below:

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