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Blunt Non-Hinged Surgical Instruments - BM

RG.: 80128580118

**ENGLISH  
EUA**

**Table 1: Listing of Blunt Non-Hinged Surgical Instruments - BM**

Code	Description	Level (hip, knee, shoulder, etc.)	Raw Material	Indication for Use
6235-00-000	Femur Distal Support Intramedullary Rod	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in femur distal part.
6207-05-000	Extension Impactor	Trauma	Stainless Steel AISI 304 (NBR 13911)	Rod impaction support.
6206-01-000	Diastase Correction Adapter	Trauma	Stainless Steel AISI 304 (NBR 13911)	Fracture reduction adapter.
5168-10-000	Acetabular Measuring Handle	Hip	Stainless Steel AISI 304 (NBR 13911)	Acetabular measurement assistance handle
5103-00-000	Femoral Rod Impactor	Hip	Stainless Steel AISI 420 (NBR 13911)	Impact the rod.
6017-35-000	Plate Twister 3.5 (Pair)	Trauma	Stainless Steel AISI 440C (NBR 13911)	Model plate.
6017-05-000	Plate Twister 3.5	Trauma	Stainless Steel AISI 420 (NBR 13911)	Model plate.
6017-10-000	Plate Twister 3.5			
6006-00-000	Screw Tweezers	Trauma	Stainless Steel AISI 440C (NBR 13911)	Pick the screw.
6017-45-000	Plate Twister 4.5 (Pair)	Trauma	Stainless Steel AISI 304 (NBR 13911)	Model plate.
6234-00-000	Extension	Trauma	Stainless Steel AISI 420 (NBR 13911)	Extend the rod for impaction.
6234-05-000	Extension Impactor	Trauma	Stainless Steel AISI 420 (NBR 13911)	Rod impaction support.
6277-00-000	Measuring Ruler - Proximal Bilocking	Trauma	Stainless Steel AISI 304 (NBR 13911)	Measurement Instrumentation.
5116-00-000	Extractor Pin	Hip	Stainless Steel AISI 304 (NBR 13911)	Support to extract the rod.
5117-00-000	Insert Pressurizer	Hip	Stainless Steel AISI 304 (NBR 13911)	Pressurize insert.
	Handle		Aluminum (NBR ISO209:10)	
5392-00-000	Non-Cemented Femoral Rod Impactor	Hip	Stainless Steel AISI 420 (NBR 13911)	Impact the rod.

5147-00-000	Gigli Saw Handle (Pair)	Hip	Stainless Steel AISI 304 (NBR 13911)	Saw Handle.
8521-00-000	Twister	Column	Stainless Steel AISI 304 (NBR 13911)	Plate Modeler.
8540-00-000	Tightening Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Wrench Tightening
8512-00-00	Spondylolisthesis Screw Protection Guide	Column	Stainless Steel AISI 304 (NBR 13911)	Break Spondylolisthesis.
	Handle		Polyacetal	
8507-00-000	Countertorque Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Lock pedicle system.
	Handle		Polyacetal	
8523-00-000	Positioner Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Implant orientation and alignment.
8502-00-000	Rod Dropper Guide	Column	Stainless Steel AISI 304 (NBR 13911)	Implant alignment.
	Handle		Polyacetal	
8526-00-000	Rupture Wrench (Spondylolisthesis)	Column	Stainless Steel AISI 440C (NBR 13911)	Break Spondylolisthesis.
	Handle		Polyacetal	
8501-00-000	Hook Locking Screw Wrench	Column	Stainless Steel AISI 440C (NBR 13911)	Screw tightening/loosening wrench
	Handle		Polyacetal	
8530-00-000	Hexagon Regulating Torque Wrench Ø5.0 mm	Column	Stainless Steel AISI 440C (NBR 13911)	Screw tightening/loosening wrench
8534-00-000	Quick-Coupling Monoaxial Screw Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Screw tightening/loosening wrench
8535-00-000	Quick-Coupling Poliaxial Screw Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Screw tightening/loosening wrench
8536-00-000	Quick-Coupling Spondylolisthesis Screw Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Screw tightening/loosening wrench
8533-00-000	Quick-Coupling T Handle (Pedicle)	Column	Stainless Steel AISI 304 (NBR 13911)	Quick-coupling instrument assist handle.
	Handle		Polyacetal	
8533-10-000	Quick-Coupling T Handle (BR3)	Hip	Stainless Steel AISI 304 (NBR 13911)	Quick-coupling instrument assist handle.
	Handle		Polyacetal	
7172-00-000	Hammer with Slot	Hip	Stainless Steel AISI 304 (NBR 13911)	Impaction instrument assist
6210-00-000	Propeller Hammer	Trauma	Stainless Steel AISI 304 (NBR 13911)	Impaction instrument assist
5504-00-000	Base Wrench BR3	Hip	Stainless Steel AISI 304 (NBR 13911)	Rod tightening with base.
6229-00-000	Distal Locking Impactor	Trauma	Stainless Steel AISI 304 (NBR 13911)	Implant impaction
6278-00-000	Wire Guide	Trauma	Stainless Steel AISI 304 (NBR 13911)	Wire assist instrument.

6225-05-000	Side Guide	Trauma	Stainless Steel AISI 304 (NBR 13911)	Screw Proximal Guide.
8505-00-000	Hexagon handle wrench 5.0 mm	Column	Stainless Steel AISI 440C (NBR 13911)	Screw tightening/loosening wrench
8539-00-000	Distraction Rod Positioner Wrench	Column	Stainless Steel AISI 420B (NBR 13911)	Place the rod.
	Handle		Polyacetal	
8624-00-000	Hook Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Hook placement.
8624-00-001	Hook Wrench Impactor Guide	Column	Stainless Steel AISI 304 (NBR 13911)	Implant impaction
8625-00-000	Pedicle Hook Impactor	Column	Stainless Steel AISI 304 (NBR 13911)	Implant impaction
	Handle		Polyacetal	
6056-00-000	END IMPACTOR	Trauma	Stainless Steel AISI 304 (NBR 13911)	Implant impactor.
6058-00-000	T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6077-02-000	One-Step Rod	Trauma	Stainless Steel AISI 304 (NBR 13911)	Locks the implant at the time of introduction.
6077-01-000	One-Step Insertion Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Implant Introduction.
6245-00-000	Tibia Side Guide Intramedullary Rod	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in tibia lateral proximal part.
6237-00-000	Femur Anterior Support Intramedullary Rod	Trauma	Stainless Steel AISI 304 (NBR 13911)	Guide the Implant.
6272-01-000	Proximal Support Intramedullary Rod	Trauma	Stainless Steel AISI 440C (NBR 13911)	Implant Guide.
6272-02-000	Distal Support Intramedullary Rod - Proximal Biolocking	Trauma	Stainless Steel AISI 440C (NBR 13911)	Implant Guide.
6157-05-000	Quick-Coupling with Chuck Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Drilling assist.
6157-05-000	Quick-Coupling with Chuck Handle - 3.5 / 4.5	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument Assist Handle.
5607-00-000	Acetabular Measuring Handle	Hip	Stainless Steel AISI 304 (NBR 13911)	Acetabulum measurement assist cable.
	Handle		Aluminum (NBR ISO209:10)	
5605-00-000	Insert Positioner/Impactor	Hip	Stainless Steel AISI 304 (NBR 13911)	Impact the implant.
	Impactor Head		Aluminum (NBR ISO209:10)	
5608-01-000	Scrap Extension Cable (Chuck 1/4)	Hip	Stainless Steel AISI 304 (NBR 13911)	Assists acetabulum scraping.
	Ring		Polyacetal	
5608-00-000	Scrap Extension Cable (Chuck 1/4)		Stainless Steel AISI 304 (NBR 13911)	
	Ring		Polyacetal	
5602-10-001	Fixed Wrench 10/13	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument tightening.

8684-01-000	Impactor Positioner Rod          Spacer with Thread	Column	Stainless Steel AISI 304 (NBR 13911)	Locks the implant in wrench.
8587-00-000	Pin Fastening Wrench	Column	Stainless Steel AISI 420 (NBR 13911)	Instrument pin introduction.
	Handle		Polyacetal	
5321-00-000	Alignment Bar	Knee	Stainless Steel AISI 304 (NBR 13911)	Anatomic alignment.
5320-00-000	Alignment Bar with Coupler	Knee	Stainless Steel AISI 304 (NBR 13911)	Anatomic alignment.
5336-00-000	Tibial Drill Bit Guide	Knee	Stainless Steel AISI 304 (NBR 13911)	Bone Drilling Guide.
5341-00-000	Test Tibia Impactor	Knee	Stainless Steel AISI 304 (NBR 13911)	Test instrument impaction.
5340-00-000	BKS III - Test Tibia Impactor	Knee	Stainless Steel AISI 304 (NBR 13911)	Remove the test.
5319-00-000	BKS III - Angular Support	Knee	Stainless Steel AISI 304 (NBR 13911)	Mechanical alignment verification.
5326-00-000	BKS III - Tibial Cutting Guide Extensor - IM	Knee	Stainless Steel AISI 304 (NBR 13911)	Tibial alignment verification.
5326-75-000	BKS III - Tibial Cutting Guide Extensor - IM (75mm)			
5315-00-000	Anterior Femoral Cutting Guide	Knee	Stainless Steel AISI 304 (NBR 13911)	Cutting guide.
5359-00-000	BKS III - Universal Extractor	Knee	Stainless Steel AISI 304 (NBR 13911)	Extract the pins.
5367-00-000	BKS III - Instrumental Screw Extractor Wrench BKS III	Knee	Stainless Steel AISI 304 (NBR 13911)	Screw tightening/loosening wrench
	Handle		Aluminum (NBR ISO209:10)	
5307-00-000	BKS III - Quick-Coupling with Chuck Handle	Knee	Stainless Steel AISI 304 (NBR 13911)	Bone drilling assist.
5372-00-000	BKS III - Screw Pointer	Knee	Stainless Steel AISI 304 (NBR 13911)	Instrument Introduction.
8686-00-000	Lumbar Spacer Tuning Fork	Column	Stainless Steel AISI 304 (NBR 13911)	Impaction Instrument.
	Handle		Polyacetal	
6242-00-000	Tibia Distal Guide Intramedullary Rod	Trauma	Stainless Steel AISI 304 (NBR 13911)	Proximal part screw guide.
	Guide Base		Aluminum (NBR ISO209:10)	
5123-00-000	Head Impactor	Hip	Stainless Steel AISI 304 (NBR 13911)	Implant impaction
	Handle		Aluminum (NBR ISO209:10)	
	Impactor Head		Polyacetal	
5107-00-00P	Insert Initial Impactor - Small	Hip	Stainless Steel AISI 304 (NBR 13911)	Implant impaction

5107-00-00G	Insert Initial Impactor - Large			
5108-00-00P	Insert Final Impactor - Small	Hip	Stainless Steel AISI 304 (NBR 13911)	Implant impaction
5108-00-00G	Insert Final Impactor - Large			
5155-00-000	Acetabular Impactor	Hip	Stainless Steel AISI 304 (NBR 13911)	Implant impaction
6228-10-000	Combined Wrench 10/10	Trauma	Stainless Steel AISI 304 (NBR 13911)	Nut tightening and loosening.
6228-15-000	Combined Wrench 15/15			
6228-13-000	Combined Wrench 10/13			
6064-00-000	Tuning Fork	Trauma	Stainless Steel AISI 304 (NBR 13911)	Implant impaction assist.
6102-00-000	Ender Modeler (pair)	Trauma	Stainless Steel AISI 304 (NBR 13911)	Model implant.
5226-00-168	BKS III - Spring - 168 mm	Knee	Stainless Steel AISI 304 (NBR 13911)	Guide fastening.
5226-00-121	BKS III - Spring - 121 mm			
8512-00-000	Spondylolisthesis Screw Protection Guide	Column	Stainless Steel AISI 304 (NBR 13911)	Break Spondylolisthesis.
2319-00-000	Rod Ring Combined Wrench - GR	Hip	Stainless Steel AISI 304 (NBR 13911)	System assembly tightening wrench.
2318-00-000	Main Base Screw - GR	Hip	Stainless Steel AISI 304 (NBR 13911)	System assembly tightening wrench.
5188-00-000	Position/Extractor of Rod of Shoulder and Elbow	Elbow	Stainless Steel AISI 304 (NBR 13911)	Place implant impactor.
5189-00-000	Shoulder and Elbow Head Extractor	Elbow	Stainless Steel AISI 420 (NBR 13911)	Extract the implant.
7255-00-000	Anchor Screw Hammer	Arthroscopy	Stainless Steel AISI 304 (NBR 13911)	Impaction assist.
6310-00-000	Extension	Trauma	Stainless Steel AISI 440C (NBR 13911)	Implant impaction assist.
6322-00-000	Rod Impactor	Trauma	Stainless Steel AISI 440C (NBR 13911)	Implant impaction assist.
6319-00-000	Hammer	Trauma	Stainless Steel AISI 304 (NBR 13911)	Implant impaction assist.
6302-00-000	Distal Guide Intramedullary Rod	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in femur distal part
6302-10-000	Distal Guide Intramedullary Rod (Retrograde)	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in femur distal part
6305-05-000	Support Connector	Trauma	Aluminum (NBR ISO209:10)	Bone Drilling Guide.
6305-00-000	Anterior Locking Wrench Support	Trauma	Aluminum (NBR ISO209:10)	Bone Drilling Guide.
6307-00-000	Locking Wrench Handle	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument Assist Handle.

6345-00-000	Quick-Coupling	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument Assist Handle.
6350-00-000	Screw Tweezers	Trauma	Stainless Steel AISI 304 (NBR 13911)	Pick the screw.
6340-00-000	Intramedullary Guide Rod Proximal Drilling Guide	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in tibia proximal part.
6343-00-000	Intramedullary Guide Rod Distal Drilling Guide	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in tibia medial and distal part.
6342-00-000	Anterior Support Intramedullary Guide Rod	Trauma	Aluminum (NBR ISO209:10)	Guide screw drilling and placement in tibia distal part.
5126-00-000	Lock Tweezers	Hip	ASTM F138	Disassemble implant.
5339-00-000	BKS III - Tibial milling cutter handle	Knee	Stainless Steel AISI 304 (NBR 13911)	Tibial milling assist cable.
6157-00-000	T Handle with Quick-Coupling (3.5 / 4.5)	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument assist handle.
5135-00-000	Insert Positioner/Impactor	Hip	Stainless Steel AISI 304 (NBR 13911)	Impactor positioner.
5301-01-000	Biomec II Rod Positioner.	Hip	Stainless Steel AISI 304 (NBR 13911)	Impactor positioner.
5158-00-000	Prosthesis Positioner and Impactor.	Hip	Stainless Steel AISI 304 (NBR 13911)	Impactor positioner.
6207-00-000	Extension - Tibial/Femur/Proximal Bilocking	Trauma	Stainless Steel AISI 440C (NBR 13911)	Implant impaction assist.
5307-10-000	BKS III - Quick-Coupling (Chuck 1/4)	Knee	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)	Instrument handle
5608-10-000	Scrap Extension Cable (Modified Trinkle Coupling)	Hip	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)/Polyacetal	Femoral scraping assist instrument.
6210-02-000	Propeller Hammer (1Kg)	Hip	Stainless Steel AISI 304 (NBR 13911)	Impaction assist.
5115-10-000	Acetabular Positioner/Impactor CP3	Hip	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)	Place the implant.
5513-35-000	Torque Wrench Pointer (Hex.3.5)	Hip	Stainless Steel AISI 440C (NBR 13911)	Screw torque pointer.
8505-25-000	Hexagon handle wrench 2.5mm	Column	Stainless Steel AISI 440C (NBR 13911)	Screw tightening/loosening wrench
8542-00-000	Screw Positioner Wrench 3.5	Column	Stainless Steel AISI 440C (NBR 13911)	Screw tightening/loosening wrench
8690-35-000	Hexagon Wrench 3.5mm	Column	Stainless Steel AISI 440C (NBR 13911) / Aluminum (NBR ISO209:10) / Polymer TPE	Screw tightening/loosening wrench
8695-15-000	Guide with Inside Thread - Small	Column	Stainless Steel AISI 304 (NBR 13911)	Guide the screw.
8695-20-000	Guide with Inside Thread - Large	Column	Stainless Steel AISI 304 (NBR 13911)	Guide the screw.
8696-00-000	Hexagon T Wrench 10.0mm	Column	Stainless Steel AISI 304 (NBR 13911)	Nut tightening and loosening.
8521-01-000	Left Twister	Column	Stainless Steel AISI 304 (NBR 13911)	Plate Modeler

5508-50-000	Torque Wrench Pointer (Hex.5,0)	Column	Stainless Steel AISI 440C (NBR 13911)	Wrench Pointer for Tightening/Loosening Screw with Torquemeter
5393-10-000	Biomec III Rod Impactor Positioner	Hip	Stainless Steel AISI 304 (NBR 13911)/Polyacetal	Implant Placement.
5605-20-000	Biomec III Insert Positioner/Impactor (Constricted)	Hip	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)	Implant Placement.
5508-50-012	Torque Wrench - Pedicle 12N	Column	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)/Polyacetal	Wrench for Tightening/Loosening Screw with Torquemeter
5115-20-000	Constricted Acetabulum Impactor Positioner CP3	Hip	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)	Implant Placement.
5140-01-000	Charnley Acetabular Positioner	Hip	Stainless Steel AISI 304 (NBR 13911) / Aluminum (NBR ISO209:10)	Implant Placement.
5605-10-000	Insert Positioner/Impactor - Biomec III	Hip	Stainless Steel AISI 304 (NBR 13911) Stainless Steel AISI 420 (NBR 13911)	Implant Placement.
6342-10-000	Anterior Support Intramedullary Rod (Small)	Trauma	Stainless Steel AISI 304 (NBR 13911) / Aluminum (NBR ISO209:10)	Bone Drilling Distal Guide.
5612-00-000	Card Wrench with Spring	Hip	Stainless Steel AISI 304 (NBR 13911)	Screw tightening/loosening wrench
5140-00-000	Müller Acetabular Positioner	Hip	Stainless Steel AISI 304 (NBR 13911) / Aluminum (NBR ISO209:10)	Place the implant
5169-00-000	Moore Prosthesis Extractor	Hip	Stainless Steel AISI 304 (NBR 13911)	Remove the implant.
5314-00-000	Impactor	Knee	Stainless Steel AISI 304 (NBR 13911)	Rod Impaction.
8508-00-000	Reducer Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Implant alignment.
8528-00-000	Temporary Fastening Pin Wrench	Column	Stainless Steel AISI 440C (NBR 13911)	Pin Introducer.
	Handle		Polyacetal	
6500-03-011	Cleaning wire 2.5 x 270	Trauma	Stainless Steel ASTM F138	Instrument used to clean cylindrical milling cutter.
6500-08-008	Cleaning wire 1.5	Trauma	Stainless Steel ASTM F138	Instrument used to clean cylindrical milling cutter.
6500-09-017	Cleaning wire 2.5	Trauma	Stainless Steel ASTM F-138	Instrument used to clean cannulated instrument.
6500-17-022	Moldable Bar (template)	Column	Aluminum (NBR ISO209:10)	Distraction rod modeling test instrument.
6500-18-002	Drill Bit Handle	Trauma	Aluminum (NBR ISO209:10)	Drill bit use assist.
6423-00-000	Proximal SP2 extension	Trauma	Stainless Steel AISI 420B (NBR 13911)	Instrument to connect proximal SP2 implant with rod guide

6409-00-000	Impaction base	Trauma	Stainless Steel AISI 630 (NBR 13911)	SP2 proximal guide assist item, that enables using the hammer to femoral canal implant introduction.
6408-00-000	Extractor impactor hammer	Trauma	Stainless Steel AISI 304, 630 (NBR 13911)	Hammer to be used against impaction base, assisting implant introduction in femoral canal.
6416-00-000	Fracture reduction wrench (pin without guide)	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument used as proximal femur fracture reduction element. It must be connected to sliding screw.
	Handle		Polypropylene and Polychloroprene	
6417-00-000	Fracture reduction wrench (pin with guide)	Trauma	Stainless Steel AISI 304 (NBR 13911); Stainless Steel ASTM F138	Instrument used as proximal femur fracture reduction element. It must be connected to sliding screw with guide.
	Handle		Polypropylene and Polychloroprene	
6425-00-000	Rod Extractor	Trauma	Stainless Steel AISI 420B (NBR 13911)	Instrument to proximal SP2 rod extraction.
6460-00-000	Humerus SP2 extension	Trauma	Stainless Steel AISI 420B (NBR 13911)	Instrument to connect humerus SP2 implant with rod guide
6461-00-000	Impaction base	Trauma	Stainless Steel AISI 630 (NBR 13911)	SP2 humerus guide assist item, that enables using the hammer to intramedullary canal implant introduction.
6462-00-000	Extension wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Extension tightening wrench.
6463-00-000	Hammer	Trauma	Stainless Steel AISI 304; 630 (NBR 13911)	Hammer to be used against impaction base, assisting implant introduction in intramedullary canal.
6474-00-000	Guide wire meter	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument used to guide wire measurement.
6451-00-000	Wire Guide	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument to place guide wires.
6451-01-000	Wire Guide	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument to place guide wires.
	Body		Polypropylene and Polychloroprene	
6500-20-006	Hammer	Hip	Stainless Steel AISI 420 (NBR 13911)	Impaction instrument, assisting hip cement removal.
6500-20-011	Femoral prosthesis extractor	Hip	Stainless Steel AISI 420 (NBR 13911)	Instrument to assist femoral prosthesis removal.
6500-01-017	Plate Modeler	Trauma	Stainless Steel AISI 420 (NBR 13911)	Implant anatomic modeling instrument
6500-01-022	Screw Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-01-023	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-02-015	Plate Modeler (large)	Trauma	Stainless Steel AISI 420 (NBR 13911)	Implant anatomic modeling instrument
6500-02-019	Screw Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers

6500-02-020	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-03-001	Compression wrench	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument used to fracture compression in DHS/DCS kit.
6500-03-008	Wire meter	Trauma	Stainless Steel AISI 304 (NBR 13911)	Measurement instrument
6500-03-012	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-03-013	Quick-Coupling Straight Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-03-021	Screw Tweezers	Trauma	Stainless Steel AISI 304 (NBR 13911)	Screw Tweezers
6500-03-031	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-04-016	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-04-017	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-04-019	Screw Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-05-009	Screw Holder Tweezers	Pelvis	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-05-011	Screw Tweezers	Pelvis	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-06-015	Cable with Torquemeter 1.5N.M	Trauma	Stainless Steel AISI 304 (NBR 13911); (NBR ISO209:10)	Torque restriction instrument handle
6500-06-019	Screw tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-06-020	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-06-023	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-06-030	L Wrench (for drill bit stop)	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6500-07-004	L Wrench for drill bit guide 4.5	Trauma	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-07-010	Cable with Torquemeter 3.0N.M	Trauma	Stainless Steel AISI 304; 440C (NBR 13911); 440C; (NBR ISO209:10)	Torque restriction instrument handle
6500-07-012	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-07-015	Screw Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-07-016	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw Tweezers
6500-07-033	L Wrench (drill bit stop)	Trauma	Stainless Steel AISI 304 (NBR 13911)	Drill bit stop tightening/loosening wrench
6500-08-009	Screw Tweezers	Trauma	Stainless Steel AISI 304 (NBR 13911)	Screw Tweezers
6500-09-012	Screw Holder Tweezers	Trauma	Stainless Steel AISI 420 (NBR 13911)	Screw picking tweezers
6500-09-013	Screw Tweezers	Trauma	Stainless Steel AISI 304 (NBR 13911)	Screw picking tweezers
6500-09-018	Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench

6500-10-000	Guide - part 01	Trauma	Stainless Steel AISI 630 (NBR 13911)	Drilling guide to introduce intramedullary rod and relevant locking screws.
6500-10-001	Guide - part 02	Trauma	Stainless Steel AISI 304 (NBR 13911), (NBR ISO209:10)	
6500-10-002	Guide - part 03	Trauma	Stainless Steel AISI 304 (NBR 13911), (NBR ISO209:10)	
6500-10-003	Guide - part 04	Trauma	Stainless Steel AISI 304 (NBR 13911), (NBR ISO209:10)	
6500-10-004	Guide - part 05	Trauma	Stainless Steel AISI 304 (NBR 13911), (NBR ISO209:10)	
6500-10-005	Guide - part 06	Trauma	Stainless Steel AISI 304 (NBR 13911), (NBR ISO209:10)	
6500-10-007	Hexagon L Wrench SW3.0 /SW5.0	Trauma	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-10-013	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-10-017	Small extractor adapter	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument to intramedullary rod extraction.
6500-10-018	Sliding hammer	Trauma	Stainless Steel AISI 420 (NBR 13911)	Impaction hammer.
6500-10-022	Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6500-11-000	Guide - part 01	Trauma	Stainless Steel AISI 630 (NBR 13911)	Drilling guide to introduce intramedullary rod and relevant locking screws.
6500-11-001	Guide - part 02	Trauma	Aluminum (NBR ISO209:10)	
6500-11-002	Guide - part 03	Trauma	Stainless Steel AISI 440C (NBR 13911)	
6500-11-003	Guide - part 04	Trauma	Stainless Steel AISI 440C (NBR 13911)	
6500-11-005	Rod introducer adapter	Trauma	Stainless Steel AISI 420 (NBR 13911)	Intramedullary rod introduction adapter
6500-11-020	Blade extractor	Trauma	Stainless Steel AISI 420 (NBR 13911)	Lock blade extraction instrument.
6500-11-021	Hammer	Trauma	Stainless Steel AISI 304, 630 (NBR 13911)	Impaction hammer.
6500-11-031	Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6500-11-032	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-12-000	Guide - part 01	Trauma	Stainless Steel AISI 304, 630 (NBR 13911)	Drilling guide to introduce intramedullary rod and relevant locking screws.
6500-12-001	Guide - part 02	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-12-002	Guide - part 03	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-12-003	Guide - part 04	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	

6500-12-004	Guide - part 05	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-12-005	Guide - part 06	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-12-007	Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6500-12-008	Hexagon L Wrench	Trauma	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-12-010	Unlocking screw extractor	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument to remove damaged screw.
6500-12-011	Sliding screw impaction wrench	Trauma	Stainless Steel AISI 304, 420 (NBR 13911)	Sliding screw placement wrench.
6500-12-012	Extractor adapter	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument to intramedullary rod extraction.
6500-12-013	Sliding hammer	Trauma	Stainless Steel AISI 420 (NBR 13911)	Impaction hammer
6500-12-014	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-13-000	Guide - part 01	Trauma	Stainless Steel AISI 304, 630 (NBR 13911)	Drilling guide to introduce intramedullary rod and relevant locking screws.
6500-13-001	Guide - part 02	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-13-002	Guide - part 03	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-13-003	Guide - part 04	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-13-004	Guide - part 05	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-13-005	Guide - part 06	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-13-007	Hexagon L Wrench	Trauma	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-13-019	Extractor adapter	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument to intramedullary rod extraction.
6500-13-020	Sliding hammer	Trauma	Stainless Steel AISI 420 (NBR 13911)	Impaction hammer
6500-13-025	Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6500-13-027	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Milling cutter instrument handle
6500-14-000	Guide - part 01	Trauma	Stainless Steel AISI 304, 630 (NBR 13911)	Drilling guide to introduce intramedullary rod and relevant locking screws.
6500-14-001	Guide - part 02	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-14-002	Guide - part 03	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-14-003	Guide - part 04	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	


6500-14-004	Guide - part 05	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-14-005	Guide - part 06	Trauma	Aluminum AISI ISO209:10), Stainless Steel AISI 440C (NBR 13911)	
6500-14-007	Hexagon L Wrench	Trauma	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-14-008	Hexagon L Wrench	Trauma	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-14-019	Large extractor adapter	Trauma	Stainless Steel AISI 420 (NBR 13911)	Instrument to intramedullary rod extraction.
6500-14-020	Sliding hammer	Trauma	Stainless Steel AISI 420 (NBR 13911)	Impaction hammer
6500-14-024	Wrench	Trauma	Stainless Steel AISI 304 (NBR 13911)	Tightening/Loosening Wrench
6500-14-027	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-17-004	Countertorque	Column	Stainless Steel AISI 304 (NBR 13911)	Countertorque Wrench
6500-17-006	Rotation Wrench	Column	Stainless Steel AISI 420 (NBR 13911)	Instrument used to distraction rod placement.
6500-17-007	Reduction Wrench	Column	Stainless Steel AISI 420 (NBR 13911)	Instrument used to distraction rod screw reduction.
6500-17-008	Wrench with tear	Column	Stainless Steel AISI 440C (NBR 13911)	Tightening/Loosening Wrench
6500-17-009	Quick-Coupling T Handle	Column	Stainless Steel AISI 304 (NBR 13911)	Instrument handle
6500-17-010	Wrench	Column	Stainless Steel AISI 304, 420 (NBR 13911)	Tightening/Loosening Wrench
6500-17-024	T Wrench	Column	Stainless Steel AISI 304, 440C (NBR 13911)	Tightening/Loosening Wrench
5117-00-001	Insert Pressurizer	Column	Stainless Steel AISI 304 (NBR 13911)	Pressurize insert in acetabulum
	Handle		Polypropylene and Polychloroprene	
8585-00-001	Impaction positioning wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Implant positioning and impacting wrench
	Handle		Polypropylene and Polychloroprene	
8686-00-001	Tuning Fork	Column	Stainless Steel AISI 420 (NBR 13911)	Vertebra milling assist wrench
	Handle		Polypropylene and Polychloroprene	
8507-00-001	Countertorque Wrench	Column	Stainless Steel AISI 304 (NBR 13911)	Countertorque Wrench
	Handle		Polypropylene and Polychloroprene	
5508-35-008	Hexagon Wrench 3.5mm with Torquemeter (BR3 - 8N)	Hip	Aluminum AISI ISO209:10), Stainless Steel AISI 304, 440C (NBR 13911)	Screw tightening screw with torque limit.
	Handle		Polyphenylsulfone	

5348-00-000	BKS III - Femoral Impactor	Knee	Stainless Steel AISI 304 (NBR 13911)	Instrument to final impaction of femoral component BKS3.
	Impactor		Polyacetal	
5349-00-000	BKS III - Tibial Impactor	Knee	Stainless Steel AISI 304 (NBR 13911)	Instrument to final impaction of tibial component BKS3.
	Impactor		Polyacetal	
6078-00-000	Plate Impactor	Trauma	Stainless Steel AISI 304 (NBR 13911)	Instrument to DHS/DCS type plate impaction.
	Impactor		Polyacetal	
6455-00-000	Quick-Coupling T Handle	Trauma	Stainless Steel AISI 304, 302 (NBR 13911)	Connection cable with various instruments.
	Handle		Polypropylene and Polychloroprene	
6455-01-000	Quick-Coupling Straight Handle	Trauma	Stainless Steel AISI 304, 302 (NBR 13911)	Connection cable with various instruments.
	Handle		Polypropylene and Polychloroprene	
5418-00-000	Wrench 13.8 mm for counternut	Knee	Stainless Steel 304 (NBR 13911)	Tighten nut in test rod.
5435-00-000	Extramedullary Alignment Bar Extension	Knee	Stainless Steel 420 (NBR 13911)	Extend alignment bar.
5436-00-000	Tibial Wedge Cutting Guide Angular Support.	Knee	Stainless Steel 420 and 304 (NBR 13911)	Intramedullary support
5437-00-000	Initial Tibial Cutting Angular Support	Knee	Stainless Steel 420 and 304 (NBR 13911)	Intramedullary support
5438-00-000	Tibial Positioner Bushing DIAM. 10.0 mm (Eccentric)	Knee	Stainless Steel 420 (NBR 13911)	Guide drill bit drilling Ø 10.0 mm.
5440-00-000	Tibial Positioner Bushing DIAM. 10.0 mm	Knee	Stainless Steel 420 (NBR 13911)	Guide drill bit drilling Ø 10.0 mm.
5444-00-000	Tibial Drill Bit Guide DIAM. 16.8 mm	Knee	Stainless Steel 420 (NBR 13911)	Guide drill bit drilling Ø 16.8 mm.
5450-00-000	T Handle with Adapter	Knee	Stainless Steel 304 (NBR 13911)	Handle to be coupled in drill bits and milling cutters, for instrument manual use.
5458-00-000	Fastening pin guide	Knee	Stainless Steel 420 (NBR 13911)	Drilling Guide
5459-00-000	Plateau Connection Screw	Knee	Stainless Steel 304 (NBR 13911)	Fasten Test Plateau in Test Tibia.
5476-00-000	Posterior Coupler for Femoral Drilling Guide Box	Knee	Stainless Steel 304 and 420 (NBR 13911)	Femur posterior cutting alignment.
5482-35-000	Torque Wrench Pointer (Hex.3.5mm)	Knee	Stainless Steel 440C (NBR 13911)	Tighten the screw with torque wrench.
5482-35-001	Torque Wrench Pointer (Hex.3.5mm)	Knee	Stainless Steel 420 (NBR 13911)	Tighten the screw with torque wrench.
5483-00-000	Femoral Revision Guide Elevator 6°	Knee	Stainless Steel 420 (NBR 13911)	Connection with coupler instrument for distal femoral cutting
5486-35-000	Torque Wrench 10N (Revision Knee)	Knee	Stainless Steel 420, 440C and 304 (NBR 13911) / Aluminum (NBR ISO209:10)	Controls the torque during screw tightening

	Handle		Polyacetal	
8505-35-000	Hexagon Wrench 3.5 mm	Knee	Stainless Steel 440C (NBR 13911)	Screw picking wrench and start the tightening.
6500-18-001	Handle Tightening Wrench	Trauma	Stainless Steel 420 (NBR 13911)	Tighten screw to fasten drill bit in handle
6500-18-004	Piston	Trauma	Stainless Steel 420 (NBR 13911)	Cannula Cleaner
	Handle		Polyacetal	
6500-18-014	Handle	Trauma	Stainless Steel 304 (NBR 13911)	Tightening assist handle.
6500-03-011	Cleaning wire 2.5 x 270	Trauma	Stainless Steel ASTM F138	Instrument used to clean cylindrical milling cutter.
6500-08-008	Cleaning wire 1.5	Trauma	Stainless Steel ASTM F138	Instrument used to clean cylindrical milling cutter.
6500-09-017	Cleaning wire 2.5	Trauma	Stainless Steel ASTM F-138	Instrument used to clean cannulated instrument.
6410-00-000	Guide For Proximal Locking Screw	Trauma	Stainless Steel 304	Guide for insertion of proximal locking screw.
6472-01-008	Flexible Shaft For Reamer 420mm	Trauma	Stainless Steel 304	Support for milling the medullary canal.
6472-01-016	Flexible Shaft For Reamer 500mm	Trauma	Stainless Steel 304	Support for milling the medullary canal.
6472-01-024	Flexible Shaft For Reamer 550mm	Trauma	Stainless Steel 304	Support for milling the medullary canal.

### Traceability Label Symbols (Label) as per NBR ISO 15223

#### TECHNICAL INFORMATION

**Product Description:** Surgical instruments can have several designs. Some have been developed specifically to comply with specific surgeon preferences. Several bear the name of the surgeon to whom they have been developed. New instruments have been developed to improve several surgical procedures or even enable new techniques. Even though innovation means new designs and even new instrument categories, in general instruments can be grouped by type or applications. Surgery techniques vary according to surgeon knowledge and choice. The surgeon is responsible for the final choice of method, type and dimension of instruments and implants to be used, as well as after-surgery result evaluation criteria. Biomecânica brand instruments should be used during the surgical procedure, according to purpose each instrument has been developed. Such instrument used is indicated to be used by professionals only. Surgeons that manage such product use need to know perfectly implant processes, and implant instrument and component handling. Sharp components can only be used for their specific purposes. Misuse can cause blade and edge cutting loss, and impact their performance within surgical procedure. Instruments that fall on the floor are invariably damaged. Whenever that happens, the instrument should be separate, inspected carefully and submitted to repair. Instrument counting is critical. At first, for the patient protection, and secondly, to prevent the possibility of being accidentally submitted to laundry along with surgical drapes, where they can be lost or damage, in case they are not noticed.

**Manufacturing materials:** Most surgical instruments are made of stainless steels, as provided by NBR 13911:1997 or ASTM F899:07; that characterize their chemical composition, mechanical properties and final behavior of stainless steel, to produce surgical instruments. Some components are made of aluminum, polyphenylsulfone, polypropylene and polychloroprene.

**Indications, Precautions, Restrictions/Warnings, Adverse Effects and Contraindications:** Blunt Non-Hinged Surgical Instruments - BM are indicated as critical assistance to prepare the bone segment to receive the implant. Biomecânica brand instrument can only be used to assist placing Biomecânica brand implants, as they have been manufactured and designed to be used as a set. Using different manufacturer instruments may compromise the surgery. Through time, instruments tend to wear naturally due to regular use. Instruments should only be used as specified. As Blunt Non-Hinged Surgical Instruments - BM are reused in another patient, it is required that

it is cleaned, sanitized and then re-sterilized. This process has to be performed every time in case of use in new surgery. The reuse without undergoing sanitization and re-sterilization process is expressly forbidden. Worn and broken instruments can cause significant problems if not detected during inspection and test. They can damage other instruments during cleaning/sterilization processes. If they remain in instrument set, from traumas in organic tissues to critical moment operation failure within surgical procedure can be caused. Identify as soon as possible instruments with problems. Separate them from the other and place labels indicating they have to be replaced or repaired. Contraindication is connected to implant to be used, as provided in Instructions For Use supplied with implants. Instruments should not be used for different purposes than intended.

**Cleaning:** Appropriate care with surgical instruments starts with appropriate cleaning. Appropriate care also means instrument regular maintenance, including sharpening and settings. There is not standard schedule; maintenance should be defined according to frequency of usage. Apart from instrument cleaning care, these instructions mention several surgical instrument enemies, including: blood, tissue in general, surgical residues (that are pitting corrosion point preliminary causes), stains and discoloration of instruments. Water and dampness also have damaging effects, enabling that such substances dry or soak into your instruments, causing unwanted stains. Another enemies used in instrument washing with inadequate solutions include: soap, discoloring agents, disinfectants and other solutions not recommended. For instrument correct conservation, recommended cleaning methods should be used, as well as understanding unwanted effect causes, like stains. Stains appear with orange or brown coloration. The idea is assuring appropriate care, to limit or extinguish them. Blood, pus and other surgical secretions contain chloride ions, that lead to corrosion, appearing more often with orange-brown color. If the instrument remains for a long period (one to four hours) in contact with residues, marks and stains will appear in it, especially if such residues dry along with the instrument. Instruments should be thoroughly cleaned and dried upon use. Only sterilize a clean instrument. Autoclave high temperature will cause chemical reactions that can cause permanent stains in instruments and/or their yellowing.

**Instrument Sterilization:** Instruments are supplied non-sterile. They have to be sterilized before use.

**We recommend autoclave steam sterilization at the hospital** (*ISO 17665-1: 2006 Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*)

**AUTOCLAVE:** It is a piece of equipment for sterilization through saturated steam under pressure.

**OPERATING INSTRUCTIONS:** For equipment start-up, check: if the circuit breaker is on; if the water supply valve is open; if the discharge valve is closed.

1 - open the equipment door; 2 - store the material to be sterilized properly 3 - close the equipment door; 4 - choose the cycle wanted according to the material to be sterilized; 5 - turn on the general switch; 6 - the cycle will be executed automatically, in sequence; 7 - as 'end of cycle' lamp is on, open the door partially for ten (10) minutes, approximately, for material chilling.

It is recommended that the following physical sterilization parameters are applied in autoclaves (saturated steam):

**Table - Autoclave physical sterilization parameters**

<b>Cycle</b>	<b>Temperature</b>	<b>Exposure Time (minimum)</b>	<b>Drying Time (maximum)</b>
Standard (1atm pressure)	121°C (250°F)	30 minutes	-
High Vacuum	134°C (273°F)	6 minutes	15 minutes

Remark: Tim must be written down whenever sterilization chamber heat achieves temperature wanted.

**For further clarification, refer to instruction manual supplied with autoclave.**

**Another sterilization method:** that can be used apart from autoclave is **Ethylene Oxide Sterilization** (E.T.O.) - parameters and procedures provided in validation protocol and ISO 11135-1 - Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

**NOTE:** The hospital institution is responsible for sterilization method, equipment, controls, and sterilization instructions used. Sterilized article cautions.

Sterilized article storage conditions: - environment: must be clean, ventilated, dry, restricted to sector personnel - article: upon sterilization process, do not place it on cold surface, (stone or stainless steel), use hollow baskets or containers until they chill - casing (raw cotton textile, non-woven textile, surgical grade paper, crepe paper, paper with film, tyvec or drilled metal cases) must remain intact and not handled much, to prevent torn packages or seal break - be stored in closed cabinets with shelves - labeled shelves, in order to enable material removal - the material must be stored according to sterilization valid date, for easier distribution and preventing expired material in stock - store non-sterile separately to decrease external contaminant level.

**Storage:** Non-sterile surgical instruments must be stored so that the packaging is not impacted. Packages cannot be compressed, crushed, drilled or exposed to possible damaged caused by water. Storage area cannot have impurities, dust or pollutants of any kind, and it cannot be exposed to extreme temperature or dampness. Area traffic should be reduced to a minimum, Instruments must be stored in trolleys or shelves, away from floor, ceiling or walls. Keep the in ventilated, dry place, away from light and severe weather action. As they are transported, inadequate bumps and stacking must be avoided. Store and transport in dry and ventilated place, with room temperature (max. 35°C) and relative humidity around 30% to 70%. Do not store directly on the floor (minimum height = 20 cm) or in very high places, close to lamps, as the packaging

might get dried or the label might be damaged. Do not store in places where contaminant substances are stored, including, for instance, cleaning materials, insecticides, pesticides, etc.

**Identification:** Instruments are provided with labels bearing the following information: Manufacturer Name; Product Trade Name; Technical Name; Batch Number; ANVISA Registration Number; Product Code; Quantity; Product Description in the package (model); Manufacturing Date; Valid Date; Description of Raw Material used to manufacture the product (Composition); Manufacturer Address; Technical Responsible Name; Statements: 'Non-Sterile Product' ;'Before using, see instructions for use; Warnings/Precautions/Special Cautions/Symbols: See Instructions For Use.

**Consumer Service:** In case it is required to perform any complaint related to instrument use related to any adverse effect that impacts the user safety, the responsible surgeon must report such adverse event to the competent health surveillance agency and Biomecânica through email [sac@biomecanica.com.br](mailto:sac@biomecanica.com.br) or telephone 0xx14 2104 7900. In case of doubt, the responsible surgeon or health care professional can report the adverse event through Health Surveillance Notification System on ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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